
# Barrier-Free<sup>®</sup> Exam Chair

[English](#)  
[Deutsch](#)  
[Français](#)  
[Español](#)  
[Dansk](#)  
[Suomi](#)  
[Italiano](#)  
[Svenska](#)  
[Norsk](#)  
[Nederlands](#)  
[Polski](#)  
[Türkçe](#)

For Models:

627 (-011)

## User Guide

0 0 3 - 1 0 0 2 8 - 9 9

Style U

003-10028-99 Rev. AB4 (07/03/24)

TP202 20-42-FO-00014 Rev A1 C2169

# Product Information


**Dealer:**

**Date of Purchase:**

**Model / Serial Number:**

**Midmark Authorized Service  
Company:**

**Model / Serial Number Label Located  
on the Back Mounting Frame**  
(exact location may vary)

MA100411

# Table of Contents

## **Important Information**

<u>Safety Symbols</u> .....	4
<u>Product Registration</u> .....	5
<u>Transportation / Storage Conditions</u> .....	5
<u>Operating Conditions</u> .....	5
<u>Disposal of Equipment</u> .....	5
<u>Authorized Representatives</u> .....	5
<u>Intended Use</u> .....	7
<u>Electrical Requirements</u> .....	7
<u>Electromagnetic Interference</u> .....	7
<u>Safety Instructions</u> .....	9

## **Operation**

<u>Active Sensing Technology®</u> .....	11
<u>Base Up / Down Function</u> .....	12
<u>Back Up / Down Function</u> .....	13
<u>Paper Roll and Tear Strap</u> .....	14
<u>Stirrups</u> .....	15
<u>Footrest and Treatment Pan</u> .....	16
<u>Pelvic Tilt</u> .....	17
<u>Clean Assist Roller System</u> .....	18

## **Maintenance**

<u>Calling for Service</u> .....	19
<u>Cleaning</u> .....	19
<u>Periodic Maintenance</u> .....	19

## **Specifications**

<u>Specifications Chart</u> .....	20
<u>Model Identification / Compliance Chart</u> .....	21
<u>Range of Motion</u> .....	22
<u>Dimensions</u> .....	23

## **Warranty Information**

<u>Limited Warranty</u> .....	24
-------------------------------	----

# Important Information - Safety Symbols

## **WARNING**

Indicates a potentially hazardous situation which could result in serious injury.

## **Caution**

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices

## **Equipment Alert**

Indicates a situation which could result in equipment damage.

## **Note**

Amplifies a procedure, practice or condition.

## Symbol Glossary

Duty Cycle (motor run time)  
30 seconds ON, 5 minutes OFF

Serial Number

Proper shipping orientation

Keep dry


Fragile

Maximum stacking height  
(palletted units)

Type B, Applied Part  
(Upholstery)

Fuse rating specification

Protective earth ground

Dangerous voltage / shock hazard

Patient Weight Limit

Refer to manual

Manufacturer

Catalogue Number

Operating Instructions

Medical Device

# Product Registration

## Transportation / Storage Conditions

Ambient Temperature Range:.....-30°C to +60°C (-22°F to +140°F)  
Relative Humidity..... 10% to 90% (non-condensing)

## Operating Conditions

Ambient Temperature Range:.....+10°C to +40°C (+50°F to +104°F)  
Relative Humidity.....30% to 75% (non-condensing)  
Altitude.....3000m or less

## Disposal of Equipment

At the end of product life, the chair, accessories, and other consumable goods may have become contaminated due to normal medical use. Consult local codes and ordinances for proper disposal of equipment, accessories and other consumable goods.

## Authorized Representatives

### WARNING

*In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.*

Customers in the EU should direct all questions, incidents and complaints to the Midmark Authorized Representative listed below.

CEpartner4U  
Esdoornlaan 13  
3951 DB Maarn, The Netherlands  
Phone: +31 343 442 524  
Fax: +31 343 442 162

Customers in the UK should direct all questions, incidents, and complaints to the Midmark UK Responsible Person listed below.

UKCApartner4U Ltd.,  
7 CAMPION WAY, BINGHAM, NOTTINGHAM NG13 8TR, UK.  
www.UKCApartner4U.co.uk  
Phone : +44 7905 384429

Customers in Australia should direct all questions, incidents, and complaints to the Midmark Sponsor listed below.

ICONA Pty Ltd.  
1A/2A Westall Road  
Clayton, Victoria 3168  
Phone: 1300 442 662

## ***Authorized Representatives - continued***

Customers in India should direct all questions, incidents and complaints to the Midmark Authorized Representative listed below.

Midmark India  
Kalpataru Point, Unit No. 12, 1st Floor Kamani road, Sion (E)  
Mumbai, 400022  
Tel: +91 22 4915 3000  
Fax: +91 22 4915 3100

Customers in Saudi Arabia should direct all questions, incidents and complaints to the Midmark Authorized Representative listed below.

Attieh Medico Ltd  
Abdullah M. Al Khalifi Street  
Al Nakheel Dist. II  
P.O. Box 116105  
Jeddah (East)  
21391  
Saudi Arabia  
Tel: + 966 2 286 4707  
Fax: + 966 2 286 4744

Customers in the United Arab Emirates should direct all questions, incidents and complaints to the Midmark Authorized Representative listed below.

GULF DRUG LLC.  
Al Barsha,  
Dubai, UAE  
Phone: +971 4 501 4000  
Fax: +971 4 501 4100

Customers in Hong Kong should direct all questions, incidents and complaints to the Midmark Authorized Representative listed below.

Associated Medical Supplies Company LTD. 沒有中文版使用說明  
Room 1201-1202, 12/F, Fo Tan Industrial Center  
26-28 Au Pui Wan Street, Fo Tan, Shatin  
New Territories, Hong Kong  
Phone: +852 2604 9389  
fax: +852 2694 0866

Customers in Qatar should direct all questions, incidents and complaints to the Midmark Authorized Representative listed below.

Tadmur Trading  
Midmac Roundabout  
Opposite Al Mana Petrol Station  
Doha, Qatar  
Tel: +974.44337000  
Fax: +974.44337100

Customers in Israel should direct all questions, incidents and complaints to the Midmark Authorized Representative listed below.

Medtechnica Ltd.  
HaTnufa St 7  
Petah Tikva, Israel  
Tel: +972 3.925.4040  
Fax: +972 3.924.9977

## ***Intended Use***

The examination chair is intended to be used as a chair / table to provide positioning and support of patients during general examination procedures conducted by medical professionals.

## ***Electrical Requirements***

Before connecting the Midmark device power cord into an electrical power source outlet, it is recommended to contact a local licensed electrician or a local Authority Having Jurisdiction over the internal building power distribution to verify that the device when powered up will comply with all local electrical codes. Use the electrical power rating marked on the device when determining the appropriate electrical branch circuit and outlet requirements. For 115V rated device models with a current rating at or above 12 Amps, a dedicated branch circuit may be required if the branch circuit protection is provided by a 15 Amp breaker.

### ***Caution***

*To completely isolate the chair from electrical mains supply, power cord must be unplugged.*

### ***Caution***

*Prior to use of high frequency surgical devices or endocardial catheters, consult the operating instructions for those devices. Failure to comply may result in electric shock or burns to the patient.*

## ***Electromagnetic Interference***

This product is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this product:

- Remove interfering device from room
- Plug chair into isolated circuit
- Increase separation between chair and interfering device
- Contact Midmark if interference persists

### ***Equipment Alert***

*Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Exam Chair including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.*

# EMC - Manufacturer's Declaration and Guidance

## Note

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals. The exam chair is not intended to be used in a residential environment.

Immunity Test	Immunity Test Level
Electrostatic Discharge	± 8 kV Contact, ±15 kV Air
Radiated RF EM Fields	3V/m 80MHz-2.7 GHz 80% AM at 1kHz
Rated Power Frequency Magnetic Fields	30 A/m
Electrical Fast Transients	±2 kV at 100kHz repetition frequency
Surges	Line to Line: ±1 kV Line to Ground: ±2 kV
Conducted RF	3V: 0.15MHz -80MHz 6V in ISM Bands between 0.15MHz - 80MHz 80% AM at 1kHz
Voltage Dips	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°
Voltage Interruptions	0 % UT; 250/300 cycle

Radiated Test	Radiated Test Level
Conducted and Radiated RF Emissions	CISPR 11 Group 1 Class A
Harmonic Distortion	Tested per IEC 61000-3-2
Voltage Fluctuations and Flicker	Tested per IEC 61000-3-3

## Safety Instructions

### **WARNING**

*No modification of this equipment is allowed.*

### **WARNING**

*Connecting equipment to the multiple socket outlet effectively leads to creating a Medical Electrical System and the result can be a reduced level of safety.*

### **WARNING**

*To avoid a fire hazard or risk of electric shock, this equipment must only be connected to a supply mains with protective earth ground.*

### **WARNING**

*Ensure that the surrounding area is clear of obstructions prior to each use. Otherwise improper operation and serious injury could occur.*

### **WARNING**

*Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.*

### **WARNING**

*Equipment is not suitable for use in the presence of a flammable anesthetic mixture.*

*Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.*

## Safety Instructions - continued

### **WARNING**

*Be sure that all personnel and equipment are clear of the chair before activating any function. Keep a minimum 6 in. (15.25 cm) clearance surrounding the chair. Failure to do so could result in personal injury.*

At the chair's lowest height [18 in. (46 cm)], the clearance between the bottom of the drawer and the floor is approximately 5 in. (12.7 cm). Make sure patient (esp. feet) and all objects are clear of this area before lowering chair.


### **Equipment Alert**

*Chair must be positioned in normal use such that the power cords can be easily accessed to unplug.*

If the chair malfunctions, immediately release controller, unplug the chair power cord, and assist the patient from the chair. If chair continues to malfunction, call for service.

Failure to comply with these instructions may result in injury to patient and/or damage to equipment.

# Operation


## Active Sensing Technology®

### **WARNING**

*Be sure that all personnel and equipment are clear of the chair before activating any function. Failure to do so could result in personal injury.*

*The Active Sensing Technology prevents the chair from moving Down if the bottom of the drawer contacts another object. If the Active Sensing Technology is activated, it will sound a series of “Beeps,” disable the Base Down function and move chair in the Base Up direction. To resume Base Down operation, remove object from under the bottom of the drawer and press the Base Up / Down buttons on controller.*

# Operation - continued

## Base Up / Down Function

### Attention!

**The Active Sensing Technology® will stop the Base Down function and retract Up if the bottom of the drawer contacts another object.**

*(You will hear a series of "Beeps" if the Active Sensing Technology is activated.)*

To resume Base Down operation:


Remove object from under the bottom of the drawer.

Press Base Down button.

### Note

*The chair automatically stops at its maximum and minimum height.*

**To move the chair Up or Down...**  
*Press and hold the appropriate button.*

ArtHalf

# Operation - continued

## Back Up / Down Function

### Note


The back section automatically stops at the top and bottom of its travel.

### To position the back section...

Press and hold the appropriate button.

ArtHalf


MA10045i

# Paper Roll And Tear Strap

Paper Roll Size (maximum):..... 24 in. wide x 3.5/8 in. diameter  
(61 cm x 9.2 cm)


## Note

One additional paper roll may be stored under the back section.

## WARNING

To avoid injury, do not put head or extremities under the paper tear strap.

**To install paper tear strap...**

Attach strap to snaps on each side of chair.

MA10065I

# Stirrups


## **Caution**

*Be sure the stirrups are locked in place before using.  
The stirrups will not support the patient's entire weight.*

## **To position the stirrups...**

- A) *Pull the stirrup out, then unfold.*
- B) *Lift the stirrup slightly, then move it left or right as desired.*
- C) *Release stirrup to lock in position.*

MA10067i

## Footrest and Treatment Pan


### **WARNING**

*Do not use footrest to support patient's entire weight.*

### **To extend footrest...**

*Pull footrest out to desired position.*

### **To access the treatment pan...**

*Pull treatment pan out to desired position.*

## Pelvic Tilt


### **WARNING**

To avoid risk of injury do not put extremities under the seat.

### **To elevate seat section...**

Lift seat until pelvic support drops into place.


**Pelvic Support**

MA10050i

### **To return seat section to flat position...**

- A) Lift seat section slightly.
- B) Pivot pelvic support handle upward.
- C) Lower seat section.

**Pelvic Support**

MA10051i

# Clean Assist Roller System

Control levers located under the back section on each side, activate and deactivate the **Clean Assist Roller System**. The rollers are located on each corner of the chair to allow positioning of chair.

## Caution

To prevent injury, never activate / deactivate Clean Assist Roller System or move chair with patient aboard. Use caution when moving chair.

## Equipment Alert

The intended use of the Clean Assist Roller System is for chair positioning within the exam room only.

## Equipment Alert

Be sure to activate both control levers.


## WARNING

To avoid risk of fall, the Clean Assist Roller System must be disengaged while a patient ingresses or egresses the chair

## To activate Clean Assist Roller System...

- Raise back section 30 to 45 degrees.
- Ensure chair is at least two inches higher than its lowest position.
- Move control levers on each side of chair to the Clean Assist Roller System engagement position.
- Move chair to the lowest position to activate rollers.


## Clean Assist Roller System Engagement Position

MA10098i

Clean Assist Roller System

## Clean Assist Roller System Disengagement Position

MA10099j

## Equipment Alert

After Clean Assist Roller System has been activated, disconnect power cord and foot control and stow away from chair to prevent damage to cords.

## To deactivate Clean Assist Roller System...

- Move chair to desired location.
- Reconnect power cord.
- Raise chair at least 2" to deactivate Clean Assist Roller System.

Note: An audible click will occur to confirm deactivation.

## Note

Control levers will automatically reset to the Clean Assist Roller System disengagement position.

# Maintenance

## Calling for Service

### Note

Model / serial number information is required when calling for service.

If service is required, contact Midmark directly:  
US 844.856.1230 / Canada 937.526.8585  
Monday-Thursday 8:00 AM until 6:00 PM  
Friday 8:00 AM until 5:00 PM (ET)

## Cleaning

### WARNING

The upholstery should be inspected before each use and replaced if punctured or otherwise damaged.

### Equipment Alert

The upholstery is resistant to most medicinal-type stains, but may be damaged by solvents and dyes. Immediately remove any fluids spilled on the upholstery.

## Upholstery

Wash your upholstery weekly with a mild liquid soap and water mixture, rinse with clear water and dry completely to remove disinfection cleaner build-up.

Disinfect your upholstery with a solution of standard bleach and water mixed 1 in 10 (10%) or chlorine based cleaners. Follow this with a clear water rinse and thorough drying of material. See current CDC Guideline for Disinfection and Sterilization in Healthcare Facilities.

To minimize damage caused by disinfectant cleaner residue build-up, do not allow disinfectants to pool on the upholstery surface. Once the approved contact time has been obtained, remove and dry excess liquid remaining on the surface.

Detailed care and maintenance instructions are included with your product. This information is also available on [midmark.com](http://midmark.com) in the Technical library under the User Information tab for your product.

## Painted Metal / Plastic Surfaces

Clean the painted metal and plastic surfaces weekly using a clean soft cloth and mild cleaner.

## Periodic Maintenance

Periodically inspect the following areas:

- Power cord(s) should be free of cuts or other visible damage.
- All fasteners should be in place and tightened securely.
- All mechanical functions should operate properly.

Periodically lubricate the back hinge to maintain quiet, smooth, operation (use light machine oil).

Have a qualified service technician inspect your equipment every twelve months.

### WARNING

Should damage and or loose parts be identified, the examination chair should be removed from use and serviced. Failure to do so may result in injury.

# Specifications



<b>Specifications Chart</b>	
<b>Patient Weight</b> ( <i>maximum</i> )	500 lbs (227 kg)
<b>Paper Roll</b> ( <i>maximum size</i> ):	24 in. long x 3.5/8 in. diameter (61 cm x 9.2 cm)
<b>Range of Motion &amp; Dimensions:</b>	[See Range of Motion & Dimensions pages]
<b>Weight of Chair:</b> w/upholstery w/packaging & skid ( <i>no uph.</i> ) Uph. w/packaging ( <i>shipped separately</i> )	450 lbs (204.1 kg) 460 lbs (208.7 kg) 40 lbs (18.1 kg)
<b>Power Cord Length:</b>	8 ft. (244 cm) long
<b>Electrical Requirements:</b>	[See Model Identification / Compliance Chart]
<b>Foot / Hand Control Voltage:</b>	3.3 VDC, SELV (Safety Extra Low Voltage)
<b>Fuse Ratings:</b> <i>IEC inlet fuses</i> F1 ( <i>on Main PC board</i> ) F2 ( <i>on Main PC board</i> )	T8AH, 250V, 5 x 20 mm T6.3AH, 250V, 5 x 20 mm T250mAH, 250V, 5 x 20 mm
<b>Duty Cycle</b> <b>(Motor Run Time):</b>	Intermittent Operation [30 seconds ON / 5 minutes OFF]
<b>Classifications:</b>	Class I, Type B Applied Part, Intermittent Operation [30 seconds ON / 5 minutes OFF]
<b>Certifications:</b>	ES 60601-1:2012, CSA C22.2 No. 60601-1:14; (Safety Standards) IEC 60601-1-2:2014 (EMC Standards)
<b>Protection against ingress of fluids:</b>	IPX0 Chair IPX1 Foot Control

# Model Identification / Compliance Chart

Model	Description	Complies To:				Electrical Ratings:			
		ANSI/AAMI ES60601-1	CAN/CSA-C22.2 No. 60601-1:14	IEC 60601-1-2	NFPA 99	CE	VAC +/-10%	Amps	Cycles (Hz)
627-011	Midmark Chair - Export w/ Pelvic Tilt	•	•	•	•	•	230	5	50/60




Fire Code Ratings: All upholstery complies with California Bureau of Home Furnishing Technical Bulletin 117 and California Code of Regulations, Sect. 93120-93120.12, Title 17. Optional upholstery is available that complies with California Bureau of Home Furnishing Technical Bulletin 133.

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© Midmark Corporation 2024
**Midmark Corporation**