# 204

## Manual Examination Table

# Ritter®

[English](#)

[Deutsch](#)

[Français](#)

[Español](#)

[Dansk](#)

[Suomi](#)

[Italiano](#)

[Svenska](#)

[Norsk](#)

[Nederlands](#)

[Polski](#)

[Türkçe](#)

ArtFull

CE

# User Guide

Style U

003-10363-99 Rev AA5 (6/26/23)

TP202 20-42-FO-00014 Rev A1 C2169

# Product Information

(The information below is required when calling for service)

**Dealer** [name / phone]:

**Date of Purchase:**

**Model / Serial Number:**

**Midmark Authorized Service Company** [name / phone]:

**Model / Serial  
Number Location**

# Table of Contents

## **Important Information**

<i>Safety Symbols</i> .....	4
<i>Warnings</i> .....	5
<i>Intended Use</i> .....	5
<i>Electrical Requirements</i> .....	5
<i>Electromagnetic Interference</i> .....	6
<i>Disposal of Equipment</i> .....	6
<i>Transportation / Storage Conditions</i> .....	6
<i>Authorized Representatives</i> .....	6

## **Operation**

<i>Back Section Positioning</i> .....	8
<i>Stirrups</i> .....	9
<i>Footrest &amp; Treatment Pan</i> .....	10
<i>Footstep</i> .....	10
<i>Paper Roll &amp; Tear Strap</i> .....	11
<i>Pelvic Tilt (optional)</i> .....	12
<i>Receptacle (optional)</i> .....	13

## **Maintenance**

<i>Cleaning</i> .....	14
<i>Periodic Maintenance</i> .....	14
<i>Calling For Service</i> .....	14

## **Specifications**

<i>Weights &amp; Dimensions</i> .....	15
<i>Upholstery Dimensions</i> .....	16
<i>Voltage, Current, Frequency Ratings &amp; Compliance Chart</i> .....	17

## **Warranty Information**

<i>Scope of Warranty</i> .....	18
--------------------------------	----

# Important Information

## Safety Symbols

### **WARNING**

Indicates a potentially hazardous situation which could result in serious injury.

### **Caution**

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

### **Equipment Alert**

Indicates a potentially hazardous situation which could result in equipment damage.

### **Note**

Amplifies a procedure, practice or condition.

Proper shipping orientation

Keep dry


Fragile

Maximum stacking height  
(palletted units)

Type B, Applied Part  
(Upholstery)

Fuse rating specification

Protective earth ground

Dangerous voltage / shock hazard

Patient Weight Limit

Refer to manual

Manufacturer

Catalogue Number

Serial Number

Medical Device

## Warnings

### **WARNING**

*To avoid injury, do not put head or extremities under the paper tear strap.*

### **WARNING**

*Connecting equipment to the multiple socket outlet effectively leads to creating a Medical Electrical System and the result can be a reduced level of safety.*

### **WARNING**

*Table must be positioned in normal use such that the power cords can be easily accessed to unplug.*

### **WARNING**

*To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.*

### **WARNING**

*This product can expose you to chemicals including Antimony Oxide (Antimony Trioxide), which is known to the state of California to cause cancer. For more information go to [www.p65warnings.ca.gov](http://www.p65warnings.ca.gov)*

### **WARNING**

*Equipment is not suitable for use in the presence of a flammable anesthetic mixture.*

*Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.*

### **WARNING**

*No modification of this equipment is allowed.*

## Intended Use

This product is intended to be used as a table to provide positioning and support of patients during examination procedures conducted by qualified medical professionals in an office environment.

## Electrical Requirements

### **Caution**

*To completely isolate the table from electrical mains supply, both power cords must be unplugged.*

*When using high frequency surgical devices or endocardial catheters:*

- Use non-conductive material to insulate patient from metal portions of the table.*
- Consult operating instructions for the device before using in conjunction with drawer or upholstery heater.*

*Failure to comply may result in electric shock or burns to patient.*

# Electromagnetic Interference

This product is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this product:

- Remove interfering device from room.
- Plug table into an isolated circuit.
- Increase separation between table and interfering device.
- Contact Midmark if interference persists.

## Disposal of Equipment

At the end of product life, the table, accessories, and other consumable goods may become contaminated from normal use. Consult local codes and ordinances for proper disposal of equipment, accessories and other consumable goods.

## Transportation / Storage Conditions

Ambient Temperature Range:.....-5°C to +38°C (+23°F to 100°F)  
Relative Humidity.....10% to 90% (non-condensing)  
Atmospheric Pressure .....500hPa to 1060hPa (0.49atm to 1.05atm)

## Authorized Representatives

### WARNING

*In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.*

Customers in the EU should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

CEpartner4U  
Esdoornlaan 13  
3951 DB Maarn, The Netherlands  
Phone: +31 343 442 524  
Fax: +31 343 442 162

Customers in the UK should direct all questions, incidents, and complaints to Midmark's UK Responsible Person listed below.

UKCApartner4U Ltd.,  
7 CAMPION WAY, BINGHAM, NOTTINGHAM NG13 8TR, UK.  
www.UKCApartner4U.co.uk  
Phone : +44 7905 384429

Customers in Australia should direct all questions, incidents, and complaints to Midmark's Sponsor listed below.

ICONA Pty Ltd.  
1A/2A Westall Road  
Clayton, Victoria 3168  
Phone: 1300 442 662

## ***Authorized Representatives - continued***

Customers in India should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Midmark India  
Kalpataru Point, Unit No. 12, 1st Floor Kamani road, Sion (E)  
Mumbai, 400022  
Tel:+91 22 4915 3000  
Fax:+91 22 4915 3100

Customers in Saudi Arabia should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Attieh Medico Ltd  
Abdullah M. Al Khalifi Street  
Al Nakheel Dist. II  
P.O. Box 116105  
Jeddah (East)  
21391  
Saudi Arabia  
Tel: + 966 2 286 4707  
Fax: + 966 2 286 4744

Customers in the United Arab Emirates should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

GULF DRUG LLC.  
Al Barsha,  
Dubai, UAE  
Phone: +971 4 501 4000  
Fax: +971 4 501 4100

Customers in Hong Kong should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Associated Medical Supplies Company LTD. 沒有中文版使用說明  
Room 1201-1202, 12/F, Fo Tan Industrial Center  
26-28 Au Pui Wan Street, Fo Tan, Shatin  
New Territories, Hong Kong  
Phone: +852 2604 9389  
fax: +852 2694 0866

Customers in Qatar should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Tadmur Trading  
Midmac Roundabout  
Opposite Al Mana Petrol Station  
Doha, Qatar  
Tel: +974-44337000  
Fax: +974-44337100

Customers in Israel should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Medtechnica Ltd.  
HaTnufa St 7  
Petah Tikva, Israel  
Tel: +972 3-925-4040  
Fax: +972 3-924-9977

# Operation


## Back Section Positioning

### To position the back section...

- A) While holding the back section in place, squeeze the Back Release Handle (either side).
- B) Move the back section to desired position.
- C) Release handle to lock back section in position.

**Back Release Handle**

MA697500i

# Stirrups


## Caution

*Be sure the stirrups are locked in place before using.  
The stirrups will not support the patient's entire weight.*

## To position the stirrups...

- A) Pull the stirrup out, then unfold.
- B) Lift the stirrup slightly, then move it left or right as desired.
- C) Release stirrup to lock in position.

MA6963001

## Footrest & Treatment Pan

### **Caution**

*Do not use footrest to support patient's entire weight.*

### **To access the treatment pan...**

*Lift footrest pad slightly, then slide pad back into table.*

### **To extend footrest...**

*Pull footrest out to desired position.*

## Footstep


### **Caution**

*Ensure patients place their feet fully on the flat surface of the footstep when getting aboard and exiting the table.*

### **To access the footstep...**

*Lift footstep slightly, then pull outward.*

## Paper Roll & Tear Strap

### Note

Paper Roll Size (max.):

Seamless Upholstery.....21 in. wide x 3,5 in. diameter  
(53 cm x 9 cm)

Soft Touch Upholstery .....21 in. wide x 4.75 in. diameter  
(53 cm x 12 cm)

**To install the paper roll...**

*Insert rod into sockets under head end of upholstery.*

MA9895p

**To install the paper tear strap...**

*Attach strap to snaps on each side of table.*

ArtF

## Pelvic Tilt (optional)

### **WARNING**



To avoid risk of injury do not put extremities under the seat.

### **To elevate seat section...**

Lift seat until pelvic support drops into place.

### **To return seat section to flat position...**

- A) Lift seat section slightly.
- B) Pivot trigger handle upward.
- C) Lower seat section.



## Receptacle (optional)

The duplex receptacle on the side of the table, provides power for accessories used during medical procedures. There are two circuit breakers located under the seat section. If the receptacle's maximum load is exceeded, the circuit breakers interrupt power to the receptacle.

**Maximum Load** ..... 5 amps

**To reset circuit breaker(s)...**  
Push the circuit breaker toggle switch.

MA9894p

# Maintenance

## Cleaning

### Upholstery

#### **Equipment Alert**

*The upholstery is resistant to most medicinal-type stains, but may be damaged by solvents and dyes. Immediately remove any fluids spilled on the upholstery.*

Wash your upholstery weekly with a mild liquid soap and water mixture, rinse with clear water and dry completely to remove disinfection cleaner build-up.

Disinfect your upholstery with a solution of standard bleach and water mixed 1 in 10 (10%) or chlorine based cleaners. Follow this with a clear water rinse and thorough drying of material. See current CDC Guideline for Disinfection & Sterilization in Healthcare Facilities.

To minimize damage caused by disinfectant cleaner residue build-up, do not allow disinfectants to pool on the upholstery surface. Once the approved contact time has been obtained, remove and dry excess liquid remaining on the surface.

Detailed care and maintenance instructions are included with your product. This information is also available on [midmark.com](http://midmark.com) in the Technical library under the User Information tab for your product.

### **Painted Metal / Plastic Surfaces**

Clean the painted metal and plastic surfaces weekly using a clean soft cloth, and mild cleaner.

## Preventative Maintenance

Periodically inspect the following areas:

- Power cord should be free of cuts or other visible damage.
- All fasteners should be present and tightened securely.
- All mechanical functions should operate properly.

Have a qualified service technician inspect your equipment every twelve months.

## Calling For Service

If service is required, contact your authorized Midmark dealer.

### **Note**

*Model / Serial Number information is required when calling for service.*


If service is required, contact Midmark directly:

US 844.856.1230 / Canada 937.526.8585  
Monday-Thursday 8:00 AM until 6:00 PM  
Friday 8:00 AM until 5:00 PM (ET)  
[midmark.com](http://midmark.com)

# Specifications

## Weights & Dimensions




<b>Maximum Patient Weight:</b>	500 lbs (226 kg)
<b>Paper Roll</b> ( <i>maximum size</i> ):	[See Paper Roll & Tear Strap page]
<b>Range of Motion &amp; Dimensions:</b>	(See illustration below)
<b>Weight of Table:</b>	220 lbs (99.8 kg)
<b>Power Cord:</b>	6 ft. (1.8 m) long

MA8376i

# Specifications (continued)

## Upholstery Dimensions

MA83771

# Specifications *(continued)*

## Voltage, Current, Frequency Ratings & Compliance Chart

Model	Description	Compliance			Electrical Ratings:		
		ES 60601-1	CAN / CSA C22.2, 60601-1	CE	VAC +/- 10%	Amps	Cycles (Hz)
204-005 204-013	Ritter® Exam Table - Export <i>(Pass-thru Drawers) w/ Pelvic Tilt</i>	•	•	•	-	-	-
204-006 204-014	Ritter® Exam Table - Export <i>(Pass-thru Drawers) w/ Receptacle &amp; Pelvic Tilt</i>	•	•	•	230	5	50/60

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**Midmark Corporation**