Midmark® Mobile Dental Delivery System

[English](#)

[Español](#)

[Français](#)

For Models:

AH9000-001

AH9000-002

AH9000-003

AH9000-004

AH9000-005

ArtHalf

User Guide

Style U

TP202 20-42-FO-00014 Rev A3 ECO088522
003-10510-99 Rev. AA8 (6/6/24)

Product Information


(The information below is required when calling for service.)

Dealer:

Date of Purchase:

Model / Serial Number:

**Midmark Authorized Service
Company:**

Model / Serial Number Location

ArtFull

Table of Contents

Important Information

Safety Symbols.....	4
Intended Use	5
Electromagnetic Interference.....	5
Disposal of Equipment.....	5
Transportation / Storage / Operating Conditions	5

Transport

Transporting Dental System	6
----------------------------------	---

Component Location

Component Location and Description	7
Access Procedures.....	8

Operation

Electrical System	10
Water System	11
Water Bottle Level Indicator	12
Air System	13
Delivery System.....	14
Low-Speed Handpiece Operation	16
High-Speed Handpiece Operation.....	17
Syringe Tip Replacement and Handpiece Installation	18
Drive Air Adjustment and Locking Caster	19
Air Pressure Regulator Adjustment and Water Flow Adjustment	20
Coolant Air Adjustment	21
Locking Casters and Storage Area.....	22

Cleaning / Maintenance

Midmark Delivery System Cleaning/Maintenance	23
Midmark High-Speed Handpiece Cleaning/Maintenance	33
Midmark Low-Speed Handpiece Cleaning/Maintenance.....	39
Midmark Scaler Handpiece Cleaning/Maintenance.....	40
Dental Suction Accessory Cleaning/Maintenance	42
Calling For Service	43
Model Identification / Compliance Chart.....	43

Specifications

Weights / Dimensions	44
Electrical	45
Electromagnetic Compatibility	45
Emissions	46

Warranty

Warranty Information & Registration	47
---	----

For installation information go to: [Technical Library at Midmark.com](https://www.midmark.com/technical-library)

Important Information

Safety Symbols

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage if not avoided.

Note

Amplifies a procedure, practice, or condition.

Proper Shipping Orientation

Handle With Care

Temperature Limit

Maximum stacking height
(Do not stack)

Humidity Limit

Type B, Applied Part

Fragile

Fuse Rating Specifications

Type BF, Applied Part

Keep Dry

Fuse Rating Specifications

Protective Earth Ground

Refer to Instruction Manual/Booklet

Pressure Limit

TUV - Domestic Models Only

Do Not Tumble

WARNING

This product can expose you to chemicals including Nickel (Metallic), which is known to the state of California to cause cancer. For more information go to www.p65warnings.ca.gov

Intended Use

The Midmark Mobile Dental Delivery System is intended to give qualified veterinarians and/or technicians a self-contained, compact, portable system for performing veterinary dental procedures.

Disposal of Equipment

At the end of product life, dental delivery(s), accessories and other consumable goods may be contaminated from normal use. Consult local codes and ordinances for proper disposal of equipment and other consumable goods.

Authorized Representatives

Warning
In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.

Customers in the EU should direct all questions, incidents and complaints to Midmark's Authorized Representative listed below.

CEpartner4U
Esdoornlaan 13
3951 DB Maarn, The Netherlands
Phone: +31 343 442 524
Fax: +31 343 442 162

Intended Operating Environment

Temperature:20°C to 24.5°C (68°F to 76°F)
Humidity:30% to 65% (non-condensing)
Atmospheric Pressure:70 kPa (20 in Hg) to 106 kPa (31 in Hg)

Safe Operating Environment

Temperature:10°C to 40°C (50°F to 104°F)
Humidity:30% to 75% (non-condensing)
Atmospheric Pressure:70 kPa (20 in Hg) to 106 kPa (31 in Hg)

Transportation / Storage Environment

Temperature:-30°C to 60°C (-22°F to 140°F)
Humidity:10% to 90% (non-condensing)
Atmospheric Pressure:50 kPa (14 in Hg) to 106 kPa (31 in Hg)

Transport


WARNING

Exercise care when transporting over thresholds and other low-lying obstacles. Failure to do so can lead to injury.


When transporting over large thresholds or steps.

- 1) Use the recessed handle pockets as shown to pull the dental delivery system up and over onto level surface.

Mobile Dental Delivery

Component Location and Description


ArtFull

Lid Access Procedure

To access under lid.

- 1) Remove mat from top of lid.
- 2) Place hand onto lever and push or pull towards back of system.
- 3) Pull on front of lid to open.


To close lid.

- 1) Rotate lid forward onto inside latches.
- 2) Push down firmly on the tray until click is heard.

Door Access Procedure

To access door.

- 1) Remove water bottles by pushing up slightly and turning.
- 2) When bottle is lined up with pin pull down.
- 3) Open access door by moving water sensors out of way and pulling on access door handle.

ArtSU5

Operation


Electrical System

On/Off Switch

On/Off Switch is located on the back of the dental delivery. It controls the main power source and must be ON to operate the delivery system.

Circuit Breakers

There are two circuit breakers located on the lower back panel of the dental delivery. If the dental delivery's maximum load is exceeded, the circuit breakers interrupt power to the dental delivery. Maximum load is 115 VAC, 6 amps (domestic) 230 VAC, 3 amps (export).

Mains Disconnect

Power can be terminated by unplugging the mains disconnect or turning switch off. Main disconnect is meant to fit loosely in the socket for quick disconnect to disengage power.

Operation

Water System


Caution

Failure to properly seal water bottles can lead to user injury.

Water Bottles

The dental delivery has two one-liter bottles. The water bottles supply water to the handpieces and air/water syringe. Air flows through the water regulator and the toggle switch to the water bottles. The water system is regulated to 40 psi. Use the water knobs to adjust water flow for each instrument. The toggle switch between the water bottles controls which bottle will supply the water. When the toggle is in the center neither water bottle is pressurized. Move the toggle toward the left to pressurize the left bottle and depressurize the right bottle. Toggle to the right to pressurize the right bottle and depressurize the left bottle.

Art-Half

To fill water bottles.

- 1) Move toggle to center.
- 2) Remove water bottle and fill. (**Distilled water is recommended**)
- 3) Replace water bottle and move toggle toward filled water bottle.


Note: Red LED indicator light will illuminate when filling is necessary.

Operation

Water Bottle Level Indicator

Water Level Sensor

The dental delivery has indicator lights that display the level of water in each bottle. When illuminated in blue, there is a sufficient supply of water. When an indicator light turns red during a procedure, 25% of the water remains. When an indicator light turns red or the water is depleted, toggle to the other bottle. Once a bottle is refilled light will return to blue.

EQUIPMENT ALERT
Lights are not working.

To Reset.

- 1) Turn delivery power switch off.
- 2) Wait approximately 1 minute.
- 3) Turn delivery power switch on.

To adjust water level light intensity.


Use push button under the lid to toggle the intensity of light brightness.

Operation

Air System

Compressor

The air system is an air compressor with a one gallon air tank and regulator. It supplies air for the handpieces and air/water syringe. The compressor will cycle between 70 to 100 PSI. The regulator is factory set at 65 PSI.

ArtHalf

Operation

Delivery System

Caution

Handpieces and Syringe Tip should be changed or disinfected and sterilized between each patient.


Handpieces

Depending on model, the dental delivery powers three or four handpieces and an air/water syringe. The handpiece pressure is controlled by handpiece drive air adjustment valves. The gauge on front of the delivery system indicates the handpiece pressure. Handpiece pressures are initially set at the factory. The scaler tool has an adjustment knob for the intensity level.

Equipment Alert

Refer to handpiece manufacturer for suggested psi settings. Only check one handpiece at a time.


Low-Speed Handpiece

When using the low-speed handpiece for polishing, turn the water knob to the off position. This handpiece position can also accommodate a high-speed handpiece if the pressure is adjusted to the correct psi setting.

Operation

Air/Water Syringe

There are air and water buttons on top of the syringe. When both buttons are depressed you will get a mist of both water and air. Air and water adjustments are located under the top cover. Insert flat head screwdriver into adjustment screw(s). Counter-Clockwise to increase flow and clockwise to decrease flow. Turn screw until air/water flow is satisfactory.

ArtSU5

Low-Speed Handpiece

Operating Recommendations

Air Pressure: Dental delivery air pressure for low-speed handpiece should be between 35 and 50 psi.

Air Supply: Moist or dirty air from the compressor can damage Midmark low-speed handpieces. If dirt or moisture is present, bleed system (see page 29) and replace filter (see page 32).

Empty Operation: Low-speed handpiece will not operate unless the chucking system is fully engaged. This minimizes risk of damage to the chuck. (see next page)


Operation

Low-Speed Handpiece Operation

Installation Notes

- This product is not sterile. Clean and sterilize before its first use. Follow manufacturer's sterilization and maintenance instructions.
- Check air line filter in dental delivery system. (see page 31) If dirt or moisture is present, bleed system (see page 29) and replace filter (see page 32).
- Preset low-speed handpiece air pressure between 35 and 50 psi.
- Align drive air and exhaust tubes with outlets in the supply tube and screw the delivery hose nut onto the Low-speed handpiece, making sure contact between handpiece and gasket is flush.
- Secure handpiece bur or prophy angle attachment, if applicable, to the low-speed handpiece, using directions see below.

Note: Midmark low-speed handpiece will operate only when the chuck ring is in the locked position.

ArtSU3

Placing and Removing Handpiece Burs or Prophy Angle Attachments and Heads

- 1) Hold handpiece securely, rotate the chuck ring clockwise (while looking at the nose cone) until it clicks into place. The arrow on the handpiece should now align with the unlock symbol on the chuck ring.
- 2) Align prophy angle slot with the pin on nose cone
- 3) Slide the prophy angle or prophy attachment fully onto the nose cone up to the shoulder. Rotate chuck ring into lock position.
- 4) Double check to be sure the attachment is secured by gently pulling on it.
- 5) To remove the prophy angle, prophy attachment or cutting bur, repeat step 1 to unlock the chuck and pull the attachment straight out.
- 6) Detach handpiece from the delivery tube at the end of each day and prior to other extended periods of non-use.

Note: See low-speed handpiece users manual for conditioning instructions.





Operation

High-Speed Handpiece Operation

To Change Bur in High-Speed Handpiece.

- 1) Push button on back of handpiece and pull the bur out.
- 2) Push bur into handpiece. (The shaft will only go in half-way)
- 3) Push button on back of handpiece and push the bur in until shaft is fully seated.
- 4) Pull lightly on the bur to ensure it is seated properly and does not come out.

Note: Friction grip burs (FG Style) are designed so that the entire straight portion of the shaft should be inside the chuck of the handpiece. You should only see the tapered portion of the bur if it is installed properly. Carbide burs are for one-time use. Dispose after procedure. Diamond burs can be re-used. Sterilize before each use.


ArtFull

Syringe Tip Replacement

To Change Syringe Tip.

- 1) Press down on large collar.
- 2) When a soft “click” is felt, pull tip straight out.
- 3) Hold collar down and insert new tip. Be sure to press it all the way in, then release collar.
- 4) Give a firm pull on tip to assure it is inserted all the way and locked in place.

Handpiece Installation

Install handpieces.

Line up components properly to prevent damage resulting in hand piece failure.

Low-Speed
Handpiece
Connection

High-Speed
Handpiece
Connection


High-Speed with
Swivel Coupler
Connection

High Speed

Low Speed

Lightly spray the exposed o-rings
with one-step conditioner.

Tighten coupler firmly with wrench.


ArtHalf

Drive Air Adjustment

To adjust drive air.

- 1) Open lid, place hand onto lever and push or pull towards back of system and lift lid.
- 2) Adjust by turning knob until proper air pressure displays on the pressure gauge.

ArtSU5

Operation

Air Pressure Regulator Adjustment


EQUIPMENT ALERT

Do not adjust regulator to more than 65 psi.

To Adjust Regulator...

- 1) Pull adjustment knob upward to unlock.
- 2) Rotate adjustment knob clockwise to increase pressure, counterclockwise to decrease pressure.
- 3) Push adjustment knob downward to lock.


ArtSU3

Water Flow Adjustment

To adjust water flow.

Rotate water adjustment knob for corresponding handpiece until desired flow is achieved.

ArtSU5


Operation

Coolant Air Adjustment


To adjust coolant air.

Insert flat head screwdriver into the corresponding screw to adjust a handpiece's coolant air. Counter-Clockwise to increase flow and Clockwise to decrease flow. Adjust until desired air pressure is reached.

Note: Drive air adjustment screws are located under lid.

Locking Caster


ArtSU5

Storage Areas:

To store objects.

- 1) Store foot control in the rear of the dental delivery system. The handle can be used to store the foot control tubing and power cord.
- 2) Storage drawer will hold prophylaxis angles, gauze, handpieces and other items as needed.

ArtSU5

Mobile Dental Delivery Cleaning/Maintenance

Note

Do not perform maintenance on machine while dental delivery system is in use.

- The following products are considered acceptable for cleaning and disinfecting exterior surfaces and handpiece tubing. They must be used according to the manufacturer's instruction unless otherwise specified:

Active Ingredient	Name
Sodium Hypochlorite	Clorox Healthcare Bleach Germicidal
	Diluted Bleach Solution (Details below)
Potassium Peroxy-monosulfate	Trifectant Broad Spectrum Disinfectant
Hydrogen Peroxide	Clorox Healthcare Hydrogen Peroxide Cleaner Disinfectant Wipes
	Rescue Concentrate
Quaternary ammonium	Clorox Disinfecting Wipes
	KennelSol Germicidal Detergent and Deodorant

- Bleach solution of 1 part household bleach (5.5% or 6% hypochlorite concentrate) and 9 parts water. Wet surface for 3 minutes, then wipe off with a clean damp cloth.
 - To properly care for these surfaces, it is extremely important to adhere to the contact time limits and thoroughly wipe these products off as directed.
 - After cleaning, visually inspect the product for deterioration of covers. If excessive discoloration, cracking or other signs of wear are noticeable, refer to the section (Calling for Service).

For external surfaces

American Veterinary Medical Association

<https://www.avma.org>

American Animal Hospital Association

<https://www.aaha.org>

For dental delivery plumbing

Organization for Safety and Asepsis Procedures

<http://www.osap.org>

Dept. of Health and Human Resources Centers for Disease Control and Prevention (CDC)

<http://www.cdc.gov>


Mobile Dental Delivery Cleaning/Maintenance

continued...

Waterline Treatment and Maintenance

Dental delivery waterline treatment is important for the health and safety of the patient and all clinic personnel. Due to the long, narrow waterline tubing in dental delivery systems, standing water can lead to microbial colonization and biofilm growth. A biofilm is a community of bacterial cells and other microbes that adhere to surfaces and form a layer of slime. If left untreated these microbes can break off, transferring bacteria to the patient or even be aerosolized into the clinic environment.

All waterlines must be flushed.

Prior to first time use of your dental delivery system, perform initial treatment to shock your waterlines and review current clinic routine waterline maintenance protocol.

If no maintenance protocol exists, one should be created that includes testing, shocking and regular germicidal treatment.

Midmark recommends testing your waterlines monthly. If a failure occurs, shock all waterlines.

Waterline shocking and maintenance treatments can come in many forms. Currently, the most popular methods on the market are shock tablets and straw/cartridge-based systems for regular germicidal treatment. Midmark has tested and recommends the use of the following HuFriedyGroup products for waterline shocking and regular germicidal treatment:

- VistaTab™ Dental Waterline Cleaner Tablets (part number: IMS-1451)
- DentaPure™ Independent water cartridge (part number: DP40B).
 - o When your system needs to be shocked, remove the cartridge(s) and attach DentaPure Cartridge Shock Accessory Kit (part number: DPBKIT)

For details on using these products, refer to the manufacturer's instructions for use. Contact HuFriedy Customer Care at 1.800.483.7433 for questions specific to waterline maintenance treatment protocols.

For customers outside of the United States, use comparable products.

General Use for Routine Maintenance.

When using a straw/cartridge-based system, install it on the primary bottle or purchase a straw/cartridge for each bottle. When shocking your system, replace with a shocking accessory. Do not shock through a functional straw or cartridge product.

If you are using IMS-1451 for routine maintenance, ensure the product is completely flushed out of the system before patient use.

Refer to the QR code located on the resource label for waterline treatment and maintenance details.


Mobile Dental Delivery Cleaning/Maintenance

continued...

After Procedure:

- Remove high-speed, low-speed and scaler handpieces.
 - Remove bur(s), disposable prophyl angle, and tip.
 - Discard used carbide bur(s) and disposable prophyl angles.
- Flush waterlines. (See next page)
- Perform handpiece maintenance.
 - Clean, condition and sterilize.
- Sterilize diamond bur(s) and scaler tip(s).
- Disinfect exterior and air/water syringe.
- Fill water bottles as necessary.

Note: The resource label provides a summary of tasks to perform after each procedure. Refer to manufacturer's instructions for details.

Equipment Alert

HANDPIECE STERILIZATION

Between all procedures complete manufacturer's instructions for "Infection Control / Sterilization" found in the documentation included with the handpieces. Failure to do this will cause premature failure of handpieces.

Equipment Alert

AUTOCLAVABLE SYRINGE TIP STERILIZATION

The autoclavable syringe tips supplied with the delivery system must be sterilized prior to use with each patient, including initial use. Be sure to thoroughly rinse and clean syringe tips prior to sterilization. Any debris may reduce the effectiveness of the sterilization process. The recommended sterilization process is steam autoclave. Recommended parameters are 125°C (250°F) and 106 kPa (15 psi) for 40 minutes at temperature and pressure.

Mobile Dental Delivery Cleaning/Maintenance

continued...

At the beginning of each work day.


- Fill water bottles with fresh distilled water.
- Flush handpiece waterlines.

Note

Distilled water is highly recommended to help prevent the development of biofilms and corrosion.

To Flush Delivery System.

- 1) Remove all handpieces, including low-speed.
- 2) Open water adjustment knob.
- 3) Hold end of tubing over a container and depress foot control.
- 4) Allow water to run for 20 to 30 seconds.
- 5) Install newly disinfected handpieces and syringe tip.

ArtSU5

Mobile Dental Delivery Cleaning/Maintenance


continued...

End of Day Maintenance.

- Remove and empty water bottles.
- Wipe entire water level sensors (2) from top to bottom with a dry lint-free cloth.
- Drain water system. (see next page for complete details)
- Drain air regulator. (see page 29)
- Drain condensation from air tank. (see page 29)
- Remove and sterilize handpieces.
- Disinfect exterior.

Ensure waterline maintenance is being followed.

Note: The resource label provides a summary of tasks required after each procedure.
Refer to handpiece manufacturer's instructions for details.

Cleaning Mat.


- Wash mat with soapy water and wipe dry. Follow-up with disinfectant.

Mobile Dental Delivery Cleaning/Maintenance

continued...

Drain Water System.

- 1) Move water bottle toggle switch to center.
- 2) Remove and empty water bottles.
- 3) Install water bottles.
- 4) Select either water bottle
- 5) Remove all handpiece tubing from holders and step on foot control to blow air through waterlines.
- 6) Replace handpiece tubing in holders.
- 7) Hold air water syringe and press water button to blow air through waterline.


ArtSU5

Mobile Dental Delivery Cleaning/Maintenance

continued...

To Drain Air Regulator.

Place a cloth under the regulator and depress the button until no water comes out.

ArtHalf

To Drain Condensation from Air Tank.


- 1) Turn power switch off.
- 2) Hold end of tube with cloth.
- 3) Turn drain valve slowly to open.
- 4) When water stops coming out, close valve.

Mobile Dental Delivery Cleaning/Maintenance

continued...

Monthly Maintenance cleaning of air/oil separator.

- 1) Turn power switch OFF.
- 2) Remove container by squeezing sides in while pulling out.
- 3) Remove fluid and saturated gauze.
- 4) Clean container and mounting cap.
- 5) Install clean gauze and reinstall the container.


ArtHalf

Mobile Dental Delivery Cleaning/Maintenance

continued...

Monthly Maintenance cleaning of regulator filter.

- 1) Turn power switch OFF.
- 2) Open door on dental delivery.
- 3) Remove regulator cover and regulator filter.
- 4) Clean and inspect regulator cover and filter.
- 5) If damaged, replace. Install regulator filter and cover.


ArtHalf

Mobile Dental Delivery Cleaning/Maintenance

continued...

Replacement intake filter.

- 1) Turn power switch OFF and unplug power cord.
- 2) Remove screws (4) and compressor assembly from back of unit.
- 3) Remove screws (4) per sound shield cap (2).
- 4) Remove old filter. Replace new filter into cap.
- 5) Install cap (2) onto compressor with screws (4) per cap.
- 6) Reinstall compressor assembly with screws (4).
- 7) Plug power cord into back of unit.

ArtHalf

Midmark High-Speed Handpiece Cleaning/Maintenance

WARNING

Handpiece should be sterilized between patients to prevent cross-contamination which could result in serious illness or death from infectious diseases.

Equipment Alert


All Midmark Handpieces must be sterilized and conditioned after every patient. Conditioning will extend life, increase performance and reduce debris within the handpiece. Conditioning can be performed between handpiece cleanings if the handpiece appears to be underperforming. If equipped, the swivel coupler can be sterilized in a bag at the same temperatures as a handpiece.

- 1) Do not use Flash sterilization techniques. This can damage components and will void the warranty.
- 2) Never submerge components in any cleaning or disinfecting solution.
- 3) Do not use Ultrasonic Cleaners on handpieces or couplers.
- 4) Do not Autoclave or Chemiclave for extended / Unusual periods of time (example: overnight).
- 5) Do not leave handpiece components in sterilizer after cycle is completed.


To clean the water port(s)...

- A) Use water spray tool to periodically clean spray ports as needed.

Clearing The Single Port Spray

Clearing The Quad Port Spray

midmark®

ArtHalf

Midmark High-Speed Handpiece Cleaning/Maintenance *continued...*



To clean the handpiece.

- 1) Detach handpiece from hose coupler (if equipped with a swivel quick disconnect remove handpiece then remove swivel from coupler with wrench).
- 2) Remove bur from handpiece.
- 3) Keeping the head of handpiece upright, clean external surface of handpiece thoroughly to remove saliva, blood and other organic material. Scrub handpiece with a soft brush and warm water to remove debris. Rinse and dry.

Equipment Alert

Align components properly to prevent damage resulting in the handpiece failure.


ArtHalf

Midmark High-Speed Handpiece Cleaning/Maintenance *continued...*

To clean/condition the handpiece drive air ports.

(Applicable on Non-Swivel High-Speed Handpieces Only).

- 1) Assemble spray tube onto spray can head by firmly pushing tube into conditioner can spray tip as far as possible.
- 2) Insert tube into drive air port of handpiece.
- 3) Direct handpiece away from your body and towards an absorbent towel. Spray until conditioner flows from both ends of handpiece head (1 to 2 seconds).
- 4) If overspray on towel appears discolored, repeat spray for an additional two seconds. Wipe excess conditioner from handpiece with a clean dry towel.


ArtHalf

Midmark High-Speed Handpiece Cleaning/Maintenance *continued...*

To clean/condition the handpiece swivel coupler cavity.

(Applicable on Swivel High-Speed Handpieces Only).

- 1) Screw swivel nozzle onto spout of a well-shaken can of Midmark One-Step Handpiece Cleaner/Conditioner.
- 2) Remove swivel coupler from handpiece and Insert nozzle tip into handpiece where swivel coupler was removed.
- 3) Direct handpiece away from your body and towards an absorbent towel. Spray until conditioner flows from both ends of handpiece head (1 to 2 seconds).
- 4) If over spray on towel appears discolored repeat spray for an additional two seconds. Wipe excess conditioner from handpiece with a clean dry towel.
- 5) Remove swivel nozzle from cleaner/conditioner can and lightly spray o-rings of swivel coupler. Re-install swivel coupler onto handpiece.


ArtHalf

Midmark High-Speed Handpiece Cleaning/Maintenance *continued...*

To clean/condition the handpiece push button and chuck.

(Applicable on ALL High-Speed Handpieces).

- 1) Screw push button chuck conditioner nozzle onto spout of a well-shaken can of Midmark One-Step Handpiece Cleaner/Conditioner. Orient curved tip of nozzle above can spout.
- 2) Insert nozzle tip into handpiece spindle bore (Where the Bur is inserted and removed) until it stops.
- 3) Hold very slight pressure on handpiece towards nozzle tip to keep it seated in position. Direct handpiece and nozzle away from your body and towards an absorbent towel. Spray for a minimum of two full seconds.
- 4) If overspray on the towel appears discolored repeat spray for an additional two seconds. Wipe excess conditioner from handpiece with a clean dry towel.
- 5) Re-attach handpiece and install a bur or bur blank. Run without water for 45 seconds. Wipe off excess with a dry towel.

ArtHalf

Midmark High-Speed Handpiece Cleaning/Maintenance *continued...*


To sterilize a Midmark High-Speed Handpiece.

- 1) Remove bur. (If equipped, the swivel coupler must be detached and autoclaved separately at the same temperatures as a handpiece.)
- 2) Place handpiece in autoclave bag. The use of autoclave bags dramatically reduces fiber optic light output deterioration and handpiece cosmetic damage.
- 3) Load autoclave bag containing handpiece into autoclave or chemiclave. Load autoclave bags for maximum penetration of steam or chemical vapor.
- 4) Cycle the handpiece according to the autoclave / chemiclave manufacturer's instructions. Dynamic air sterilizers autoclave the handpieces at 270°F (132° C) for at least 4 minutes. (Do not exceed 275°F (135° C).
- 5) Immediately remove handpiece from autoclave or chemiclave. Allow to cool sufficiently prior to handling.
- 6) As soon as handpiece is cool enough to handle, wipe fiber optic light transmitting surfaces clean on both ends of handpiece with isopropyl alcohol using a cotton swab. This step is particularly critical to maintaining light output (Fiber Optic High-Speed Handpieces Only).
- 7) If cleaning fiber optics of a swivel high-speed handpiece, remember to clean surfaces at both ends of the swivel coupler.
- 8) If equipped, reattach swivel coupler to handpiece.

Equipment Alert

For all other handpieces, refer to handpiece manufacturer's instructions for suggested cleaning and maintenance procedures.

Midmark Low-Speed Handpiece Cleaning/Maintenance


WARNING

Handpiece should be sterilized between patients to prevent cross-contamination which could result in serious illness or death from infectious diseases.

Maintenance - Low-Speed Handpiece

The Midmark low speed handpiece is designed to operate without the need for repeated lubrication. However, a build-up of debris inside the motor may affect its performance over time. Should that happen, apply Midmark one-step handpiece cleaner/conditioner to the drive air inlet tube and run the handpiece for 15 seconds to flush the debris out of the motor. Performance should be restored. Wipe the Midmark low speed handpiece with a clean, dry cloth or gauze to remove excess conditioner. Lubricant should be applied to the joint between the nose cone and the motor.

To Clean Handpieces

When autoclaved according to the following instructions, the Midmark 0-5k Straight low speed handpieces will withstand up to 1200 autoclave cycles before requiring service.

- 1) Prepare the Midmark low speed handpieces for autoclave by carefully scrubbing the external surface with a soft brush and warm tap water $40^{\circ}\text{C} \pm 5^{\circ}\text{C}$ ($104^{\circ}\text{F} \pm 10^{\circ}\text{F}$) to remove visible debris and particles. Do not run water into back end tubes. Dry with a towel or dry compressed air.
- 2) Seal handpieces in separate autoclave bags. Make sure to remove dental attachments from handpiece before sealing.
- 3) Process in dynamic air sterilizer as follows:
 - Place handpiece(s) in autoclave tray(s).
 - Autoclave at a temperature of 270°F (132°C) for at least 4 minutes.
 - NOTE: Temperature should not exceed 275°F (135°C).
- 4) After autoclave cycle is completed, allow to cool for 30 minutes prior to handling.

Midmark Scaler Handpiece Cleaning/Maintenance


WARNING

Handpiece should be sterilized between patients to prevent cross-contamination which could result in serious illness or death from infectious diseases.

To Clean Scaler Handpiece

- 1) Remove scaler handpiece from delivery system.
- 2) Refer to manufacturer's instructions for cleaning and maintenance.
- 3) Reinstall sterilized scaler for new patient when needed.


Art-Half

Midmark Scaler Adjustment Knob Cleaning/Maintenance

To Clean Under.

- 1) Remove knob cap.
- 2) Loosen knob clamp screw with screwdriver.
- 3) Remove knob and clear indicator disk from adjuster shaft. Clean disk and surface behind disk.
- 4) Install clear indicator disk and knob. Be sure the three locating prongs on the disk align with the corresponding holes in the knob.
- 5) Lightly tighten the clamp screw. The knob should be tight enough to turn the adjuster shaft but loose enough to slip on the shaft when it reaches either travel limit. Once the shaft reaches its travel limit, whether to the right or left, continue to rotate the knob in the same direction until the indicator line aligns with the corresponding end of the color dial. Tighten the clamp screw and verify the indicator stops at each end of travel, from "1" to "10".
- 6) Install knob cap.


ArtHalf

Dental Suction Accessory - Cleaning/Maintenance

To clean vacuum tubing and HVE instrument.

- 1) Empty the collection container and reinstall.
- 2) Mix chosen treatment system, per manufacturer's instructions, into a clean bucket.
- 3) Turn unit ON and suction the treatment solution through the HVE instrument.
- 4) Empty solution from the collection container and reinstall.
- 5) Fill bucket with 16 oz. of water.
- 6) Turn unit ON and suction water from bucket through the HVE instrument.
- 7) Empty the collection container and reinstall.

ArtHalf

Calling For Service

Direct all service inquiries to your authorized Midmark dealer. When calling for service, you must provide the following information:

- Model / serial number
- Date of purchase
- Symptom(s) of malfunction

Service and Parts Manual is available at www.midmark.com.

Model Identification / Compliance Chart

Model	Description	Electrical Ratings:		
		VAC +/- 10%	Amps	Cycles (Hz)
AH9000-001	Midmark Mobile Dental Delivery System, Integrated 4 Position	115	6	50/60
AH9000-002	Midmark Mobile Dental Delivery System, High-Speed Swivel, Integrated 4 Position	115	6	50/60
AH9000-003	Midmark Mobile Dental Delivery System, High-Speed Swivel, Integrated 5 Position	115	6	50/60
AH9000-004	Midmark Mobile Dental Delivery System, Integrated 4 Position (Export)	230	3	50/60
AH9000-005	Midmark Mobile Dental Delivery System, Integrated 4 Position, Canada	115	6	50/60

Specifications

Cart	
Width	44.2 cm (17.4 in)
Length	48.26 cm (19 in)
Height	85.47 cm (33.8 in)
Weight (Mass)	58.97 kg (130 lb)
Shipping Weight	68.04 kg (150 lb)
Top of Cart Weight Restriction	11.33 kg (20 lb)
Compressor	
Motor	Oil Free 1/2 HP
Working Pressure	80-100 psi
Flow Rate	2.4 scfm @ 60 psi
Air Tank	1 gal
Water Module	
Reservoirs	Two - 1 L - Toggled
Pressurized System	30-40 psi
Aseptic	Non-Retracting HPC Coolant
Handpiece Flush System	Manual
Other	
Duty Cycle	Continuous Operation Note: Exception, Units with LED Scaler, 10 min ON, 5 min OFF
Classifications	Class 1, cart treated as Type B Applied Part, except as noted for accessories, Ordinary Equipment [IPX0]
Accessories:	
Acteon Satelec Scaler	Type BF Applied Part
High-Speed Handpiece with Fiber Optics	Type B Applied Part
Low-Speed Handpiece	Type B Applied Part
Air Water Syringe	Type B Applied Part

Specifications *continued...*

Electrical (See Regulatory Compliance Chart - Next Page)	
Fuse (Power Box) 115VAC	800mA 250V, SLO BLOW
Fuse (Power Box) 230VAC	400mA 250V, SLO BLOW
Fuse (Main Board)	250V 3.5, FA
Circuit Breakers (2) 115VAC	7A 250V, Push to Reset
Circuit Breakers (2) 250VAC	4A 250V, Push to Reset
Power Cord (115VAC)	15A/125VAC x 10 ft. Hospital Grade
Power Cord (230VAC)	10A/250VAC x 2.5 meters

Electromagnetic Compatibility (EMC) Declaration and Guidance		
Immunity		
Test	Standard	Level
Electrostatic Discharge	IEC 61000-4-2	± 8 kV Contact, ± 15 kV Air
Radiated RF EM Fields	IEC 61000-4-3	3 V/m (80 MHz -2.7 GHz)
Rated Power Frequency Magnetic Fields	IEC 61000-4-8	30 A/m, 60 Hz
Proximity Magnetic Fields	IEC 61000-4-39	134.2 kHz, Pulse Modulation 2.1kHz 65A/m13.56MHz, Pulse Modulation 50kHz, 7.5A/m
Electrical Fast Transients / Bursts	IEC 61000-4-4	± 2kV @ 100kHz Rep. (AC Mains)± 1kV @ 100kHz Rep. (I/O Ports)
Surges, Line to Line	IEC 61000-4-5	± 0.5 kV, ± 1 kV L-L
Surges, Line to Ground	IEC 61000-4-5	± 0.5 kV, ± 1 kV, & ± 2 kV L-PE
Conducted RF	IEC 61000-4-6	3 V m 0.15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage Dips	IEC 61000-4-11	0% UT for 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°.
Voltage Interruptions	IEC 61000-4-11	0% UT for 1 cycle 70% UT for 25/30 cycles Single phase: at 0° Notes: 30 cycles for 60 Hz 0% UT for 250/300 cycles Notes: 30 cycles for 60 Hz

Emissions		
Test	Standard	Level
Conducted and Radiated RF Emissions	CISPR11	Class A, Group 1
Harmonic Distortion	IEC 61000-3-2	Up to 16A per phase
Voltage Fluctuations and Flicker	IEC 61000-3-3	Up to 16A per phase

WARNING

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) from any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING

Equipment not suitable for use in the presence of a flammable anesthetic mixture. Isoflurane and Sevoflurane are recommended for use.

WARNING

Do not modify this equipment without authorization from the manufacturer.

WARNING

To avoid the risk of electric shock, this equipment must ONLY be connected to an outlet with a ground socket.

WARNING

The Acteon Newtron® LED scaler complies with current electromagnetic compatibility standards. However, it must not be used by or on those with cardiac stimulators or other active implants unless adequate electromagnetic interference shielding is present. Check with a cardiologist to verify the EMI shielding capabilities of the implant. See enclosed Newtron® LED manual for further advisories.

WARNING

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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Midmark Corporation