Smart Air Compressor


English
Español
Français

Applies to Models:

P21-050 - 115 Volts

P22-050, P32-050 - 230 Volts

P52-050, P72-050 - 230 Volts

User Guide

AAF2824

Style P

003-10538-99 Rev. AA5 (9/23/24)

TP202 20-42-FO-00014 Rev A3 ECO088522


Product Information

Dealer:

Date of Purchase:

Model / Serial Number:

Midmark Authorized Service Company:

**Model / Serial
Number Location**

Table of Contents

Product Information

Model / Serial Number Location 2

Important Information

Safety Symbols 4

Symbol Glossary 4

Product Registration 5

Warranty Information 5

Intended Use 5

Electromagnetic Interference 5

Disposal of Equipment 5

Transportation / Storage Conditions 5

Electromagnetic Compatibility / Emissions 6

Compressor Main Component Locations 7

Site Requirements

Electrical 8

Plumbing 8

Connectivity 8

Environmental 8

Installation Kit 9

Specifications

Weights / Dimensions 10

Electrical Specifications 10

Maintenance

Maintenance Components Location 11

Weekly 12

Monthly 12

Annually 12

Bi-Annually 12

Maintenance Kits / Parts 12

Calling for Service 12

Connectivity Module Information

Intended Use 13

Safety Instructions 13

Electromagnetic Interference 13

Transportation / Storage Conditions 13

Operating Conditions 13

IC Compliance Statement (Interference): 14

IC Compliance Statement (Antennas): 14

FCC Compliance Statement: 14

Classifications 14

Wi-Fi and Ethernet technology to transmit data 15

Cybersecurity: 16

Important Information

Safety Symbols

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage.

Note

Amplifies a procedure, practice, or condition.


Symbol Glossary

Keep dry

Maximum stacking height
(palletted units)

Fuse rating specification

Caution hot surface

IPXO

Ordinary Equipment

FCC

Refer to instruction manual

Proper shipping orientation

Fragile

Protective earth ground

Manufacturer

Serial Number

Catalogue Number

This product has been evaluated with respect to electrical shock, fire and mechanical hazards only, in accordance with UL60601-1 and CAN/CSA C22.2 NO. 601.1.

DENTAL AIR COMPRESSOR

UL 60601-1, IEC60601-1,
CAN/CSA C22.2 No.601.1

Product Registration

The following link and QR code are for access to the Smart Air Compressor warranty registration.
<https://www.midmark.com/service-support/terms-conditions/warranty-registration>

Warranty Information

The following link and QR code are for access to the Smart Air Compressor warranty.
<https://www.midmark.com/service-support/terms-conditions>

Intended Use

To provide compressed air during general examinations and procedures conducted by qualified dental professionals. Install compressor in a room separate from patients or professionals.

WARNING
*This product can expose you to chemicals including nickel (metallic), which is know to the state of California to cause cancer.
For more information go to [www.p65 warnings.ca.gov](http://www.p65warnings.ca.gov).*

Electromagnetic Interference

This Midmark Smart Air Compressor is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and these units: Remove interfering device from room.

- Remove interfering device from room
- Plug interfering device into an isolated circuit
- Increase separation between unit and interfering device
- Contact Midmark if interference persists

Disposal of Equipment

At the end of product life, the unit(s), accessory and other consumable goods may become contaminated from normal use. Consult local codes and ordinances for proper disposal of equipment, accessories and other consumable goods.


Transportation / Storage Conditions

Ambient Temperature Range:.....-30°C to 60°C (-22°F to 140°F)
Relative Humidity.....10% to 90% (non-condensing)
Atmospheric Pressure500hPa to 1060hPa (0.49atm to 1.05atm)

Electromagnetic Compatibility (EMC) Declaration and Guidance		
Immunity		
Test	Standard	Level
Electrostatic Discharge	IEC 61000-4-2	± 8 kV Contact, ± 15 kV Air
Radiated RF EM Fields	IEC 61000-4-3	3 V/m (80 MHz -2.7 GHz)
Rated Power Frequency Magnetic Fields	IEC 61000-4-8	30 A/m, 60 Hz
Proximity Magnetic Fields	IEC 61000-4-39	134.2 kHz, Pulse Modulation 2.1kHz 65A/m13.56MHz, Pulse Modulation 50kHz, 7.5A/m
Electrical Fast Transients / Bursts	IEC 61000-4-4	± 2kV @ 100kHz Rep. (AC Mains)± 1kV @ 100kHz Rep. (I/O Ports)
Surges, Line to Line	IEC 61000-4-5	± 0.5 kV, ± 1 kV L-L
Surges, Line to Ground	IEC 61000-4-5	± 0.5 kV, ± 1 kV, & ± 2 kV L-PE
Conducted RF	IEC 61000-4-6	3 V m 0.15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage Dips	IEC 61000-4-11	0% UT for 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°.
Voltage Interruptions	IEC 61000-4-11	0% UT for 1 cycle 70% UT for 25/30 cycles Single phase: at 0° Notes: 30 cycles for 60 Hz 0% UT for 250/300 cycles Notes: 30 cycles for 60 Hz

Emissions			
Test	Standard	Level	
		P22, P32	P21, P52, P72
Conducted and Radiated RF Emissions	CISPR11	Group 1, Class A	Group 1, Class A
Harmonic Distortion	IEC 61000-3-2	Up to 16A per phase	N/A
Voltage Fluctuations and Flicker	IEC 61000-3-3	Up to 16A per phase	N/A

Main Component Locations

MAF12228

P21-050 / P22-050 Model Shown
Components Shown are Typical on all Smart Air Compressor Models

Site Requirements

Note

Compressor should be installed in a room separate from patients or professionals.

Electrical

Supply	P21-050	P22-050	P32-050	P52-050	P72-050
Voltage +/-10%	115	230	230	230	230
Single Phase					
Disconnected Switch Box - Fused - User Supplied - Box(es) must be located within 3 ft. of compressor	20 Amps	20 Amps	20 Amps	30 Amps	40 Amps

Plumbing

Fresh Air Intake Line	
Type	PVC
Line Sizing	1" Minimum (Provided by Contractor)
Compressor Connections	1" MPT or 1 in. PVC Slip Connectors
Drain	
Type	Poly Tube (Provided with Install Kit)
Delivery Line	
Main Trunk Line	3/8" to 3/4" Copper, as site layout dictates. Terminate in mechanical room with 1/2" FPT. Flex hose provided for service connection. Trunk line to slope 1/4" / 10ft toward source compressor.
Operatory Branch Lines	3/8" to 3/4" Copper - Same size as trunk line. Terminate 1/2" or 5/8" FPT at Operatory junction box location.

Connectivity

Connections (all models)	1x RJ45 cable, 1x RJ45 outlet port
--------------------------	------------------------------------

Environmental

Equipment Room Ambient Temperature - Operational	40° - 104° Fahrenheit / 4° - 40° Celsius
Relative Humidity	10% to 90% (non-condensing)
Atmospheric Pressure	700hPa to 1060hPa (0.69atm to 1.05atm)


Site Requirements - continued

Installation Kit Provided with Compressors

**Exhaust
Silencer**


**Fresh Air Intake Hose
Length - 10'**


**Elbow for
Operatory Line
Connection**

Feet with Hardware

**Operatory Air Line
Length - 6'**

Specifications

Classifications: Class 1, No Applied Part

Weights / Dimensions

Model	# Users	CMF @ 80 PSI	Tank Capacity	Sound Level (dBA)	Sound Level (dBA) With Quiet Cover	Dimensions H x W x D (IN.)	Dimensions With Sound Cover H x W x D (IN.)	Product Weight (LBS.)	Fresh Air Intake Pipe Size
P21-050	1 - 3	5.2	10	65	60	27 x 26 x 24	32 x 26 x 25	225	2"
P22-050	1 - 3	5.2	10	65	60	27 x 26 x 24	32 x 26 x 25	225	2"
P32-050	3 - 5	7.8	20	67	62	29 x 33 x 26	33 x 33 x 28	295	2"
P52-050	5 - 7	10.6	20	70	65	31 x 33 x 26	33 x 33 x 28	345	2"
P72-050	7 - 10	15.9	32	73	68	33 x 40 x 27	39 x 40 x 30	425	2"

Electrical Specifications

Model	Total HP	Voltage +/-10%	Total Amps	Breaker Size	Motor Fuses	PCB Fuse	Transformer Fuse
P21-050	1.5	115	16.0	20	15A, 250V, Slo-Blo	1/2A, 250V, Slo-Blo	400mA, 250V, Slo-Blo
P22-050	1.5	230	7.0	20	15A, 250V, Slo-Blo	1/2A, 250V, Slo-Blo	200mA, 250V, Slo-Blo
P32-050	2.25	230	11.0	20	15A, 250V, Slo-Blo	1/2A, 250V, Slo-Blo	200mA, 250V, Slo-Blo
P52-050	3.2	230	20.0	30	15A, 250V, Slo-Blo	1/2A, 250V, Slo-Blo	200mA, 250V, Slo-Blo
P72-050	4.8	230	29.0	40	15A, 250V, Slo-Blo	1/2A, 250V, Slo-Blo	200mA, 250V, Slo-Blo

Maintenance


Equipment Alert

Never lubricate this oil-less compressor.
The motor / crankshaft bearings and Teflon
compressor rings are permanently lubricated.

Note

Upon initial startup, the compressor
will have a 15-20 second delay before
heads begin to run.

P21-050 / P22-050 Model Shown

Duty cycle: When properly sized, the compressor should run
a maximum of 35% of the time.

Maintenance - continued

Weekly

Check Moisture Indicator

- If indicator is *BLUE*, no action is required.
- If indicator is *PINK*, marginal life or compressor is overcycled.
- If indicator is *WHITE*, dryer tank needs replacing.

Monthly

Check Pressurization Time

Watch the Storage Air gauge as the compressor runs. Record the time it takes to go from 80 PSI to 100 PSI Max (+/- 5psi). It should take approximately 15 to 45 seconds. Note: For P52 and P72, Listen to hear fan running. When motors are cycled off, fan should be heard running.

[Contact your dealer if significantly out of range]

Clean / Dust Compressor

Wipe with a clean cloth to remove dust & debris.

Annually

Replace Air Intake Filters

The number of air intake filters varies by model. You must replace all of these filters annually.

Replace Coalescing Filter Element and O-ring

All models have only one coalescing filter element.

Bi-Annually

Replace Exhaust Muffler

Maintenance Kits / Parts

P21-050 - P22-050 - P32-050

(Order Kit #: 002-10807-00)

Kit Includes:

- Air Intake Filters (3)
- Coalescing Filter Element. (1)
- O-Ring (2)

P52-050 (Order Kit #: 002-10805-00)

Kit Includes:

- Air Intake Filters (4)
- Coalescing Filter Element. (1)
- O-Ring (2)

P72-050 (Order #: 002-10806-00)

Kit Includes:

- Air Intake Filters (6)
- Coalescing Filter Element. (1)
- O-Ring (2)

Calling for Service

Note

Model / serial number information is required when calling for service.

If service is required, contact your authorized Midmark dealer.

To contact Midmark directly:

1.800.MIDMARK (1.800.643.6275)
Monday-Thursday 8:00 AM until 6:00 PM
Friday 8:00 AM until 5:00 PM (ET)
midmark.com

Connectivity Module Information

Intended Use

The Connectivity Module is intended to be integrated within a Midmark device that is used in a healthcare professional environment in order to transmit, store, and convert non-PHI data to and from a remote server. This device is not intended to interact or come in contact with patients.

Safety Instructions

WARNING
No modification of this equipment is allowed.

Electromagnetic Interference

WARNING
Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Equipment Alert
Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The ambient environmental conditions inside the enclosure in which the Connectivity Module is installed must be within the ranges of or must not exceed:

Transportation / Storage Conditions

Ambient Temperature Range:.....-30°C to 60°C (-22°F to 140°F)
Relative Humidity.....10% to 90% (non-condensing)
Atmospheric Pressure500hPa to 1060hPa (0.49atm to 1.05atm)
Altitude.....-380 meter to 5,570 meter

Operating Conditions

Ambient Temperature Range:.....4°C to 65°C (40°F to 149°F)
Relative Humidity.....10% to 90% (non-condensing)
Atmospheric Pressure700hPa to 1060hPa (0.69atm to 1.05atm)
Altitude.....-380 meter to 3,000 meter

Connectivity Module Information - continued

IC Compliance Statement (Interference):

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device. Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes: (1) appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

IC Compliance Statement (Antennas):

Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.

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FCC Compliance Statement:

This equipment has been tested and found to comply with the limits for a class A digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their own expense. Any modifications to this product may violate the rules of the Federal Communications Commission and make operation of the product unlawful.

Classifications

Protection against electrical shock: Class II.

Operating mode: Continuous operation.

- i. Only products certified to one of the following may be connected to the Ethernet port
 - 1. IEC 60950-1 and complying with SELV or TNV-1 limits
 - 2. IEC 62368-1 and complying with ES 0 or ES 1 limits
- ii. Connecting the product/connectivity module to an IT-network could result in previously unidentified risks to patients, operators or third parties and that the responsible organization should identify, analyze, evaluate and control these risks.

Connectivity Module Information - continued

The Smart Air Compressor has a remote connection with Midmark Phone Application. The Smart Air Compressor includes an Ethernet connection and a Wi-Fi module that allows the device to connect to the local network. In rare instances, the Smart Air Compressor may be unable to connect to a wireless network at powerup due to operating location/environment frequency interference on certain wireless channels. Moving the Smart Air Compressor slightly may resolve this issue. Note: Interference conditions may change. Connecting the Smart Air Compressor via Ethernet is another option. Please consult the IT department/contact to verify the proper Wi-Fi coverage in the area where the device will be used.

The Smart Air Compressor supports Wi-Fi and Ethernet technology to transmit data

Objects	Description	Value/Comment
<i>Application and Operating System</i>		
Operating System/Version		
Network Protocols		Link Layer: IEEE 802.3, IEEE 802.11 Internet Layer: IP Transport: TCP, UDP App Layer: TLS, MQTT, NTP, HTTP, etc.
Network data	LAN Physical Connection	RJ-45
	LAN Cabling	Ethernet (min. CAT 5)
	LAN Max Speeds	10/100 Mbps
	LAN Protocols	IEEE 802.3u
	WiFi Connection	The END PRODUCT supports connectivity to Wi-Fi using the following standards: • 802.11a • 802.11b • 802.11g • 802.11n • 802.11ac
	Operational Frequency	Connects to a 2.4 GHz or 5.0 GHz
	WiFi Security Standards	When connecting to Wi-Fi, the END PRODUCT uses the following security standards: • WPA2-PSK • WPA2 PEAP/MSCHAPv2 • WPA2 EAP-TLS
	Network interface card parameters	Not Applicable

Specifications Chart

Certifications Midmark is an ISO 13485 and ISO 9001 Certified Company	IEC 60601-1 Edition 3.1 (IEC 60601-1:2005+A1:2021) - Certification for Connectivity Module IEEE 802.11 IEEE 802.3
---	--

Connectivity Module Information - continued

Cybersecurity:


Connecting the Smart Air Compressor to a local network could result in previously unidentified cybersecurity risks that organization responsible for the network should identify, analyze, evaluate, and control. The product captures relevant system and security logs to identify date, time, and source of user activity and a software bill of material (SBOM) is available upon request. Contact Midmark Technical Service for assistance if a cybersecurity vulnerability or incident is detected while using the Smart Air Compressor or for any other cybersecurity needs.

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Midmark Corporation