Midmark Dental Chair


Model: CH-1000

User Guide

*English
Español
Français*

Model / Serial Number Location

**Model / Serial
Number Location**

ArtFull

Table of Contents

Model / Serial Number Location

Model / Serial Number Location	2
--------------------------------------	---

Important Information

Symbols Glossary.....	4
Intended Use / Indication for Use.....	5
Transportation and Storage Conditions.....	5
Operating Conditions.....	5
Disposal of Equipment	5
FCC Compliance Statement.....	5
Midmark's Online Help	5
Calling for Service.....	5
Electromagnetic Interference.....	6
Wifi and Ethernet Connections	6
Unregistering a Chair.....	7
Cybersecurity.....	7
Software Update.....	7
Warnings	9-10
Authorized Representatives.....	11

Component Location

Component Location and Description.....	12
---	----

Product Features

POWER ON Light.....	13
Collision Protection System.....	14

Product Operation

Foot Control.....	15
Home Functions	16
Stop Functions.....	16
Back Functions	17
Base Functions	18
Programming Chair Positions	19
Chair Lockout	20
Hand Control	21
Heat and Massage	22
Chair Back.....	22
Heat and Massage Hand Control	22
Hand Control Details	23
Massage Mode Selections.....	23
Heat Mode Selections.....	24
Adjusting the Headrest.....	25
Positioning Headrest for Wheelchairs	26

Product Operation continued...

Adjusting Back Height	27
Rotating the Chair.....	28

Barriers, Cleaning, and Disinfecting

Clinical Contact Surfaces.....	29
Non-Clinical Contact Surfaces	29
Barriers Protection	29
General Purpose Cleaning.....	29
Cleaning and Disinfecting Assistance	30
Cleaning and Disinfecting	30

Specifications

Specifications.....	31
Certification and Compliance.....	31

Accessories

Cuspidor.....	31
Left-Right Monitor.....	31

Wi-Fi, Ethernet, Cybersecurity

Wi-Fi Connection Conditions	32
Cybersecurity	32
Connectivity Options Chart.....	33
Network-Related Hazardous Situations.....	34

EMC - Manufacturer's Declaration and Guidance

Electromagnetic Emissions.....	35
Recommended Separation Distances	35
Electromagnetic Immunity.....	36

Warranty Information

Warranty Information QR Code.....	38
Warranty Registration QR Code.....	38

Important Information

Symbols Glossary

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage.

Note

Amplifies a procedure, practice, or condition.

Refer to manual

Pressure Limit

Manufacturer

Consult IFU

Temperature Limit

Catalogue Number

Type B, Applied Part

Humidity Limit

Duty Cycle

Protective earth ground

Maximum stacking height
(palletted units)

TUV

Proper shipping orientation


Keep dry

Patient Weight

Fragile

Transportation
Orientation

The heated and massage upholstery is a Type BF, Applied Part. Do not penetrate upholstery with sharp objects, such as needles, knives, etc

Important Information

Intended Use / Indication for Use

The Next Generation Dental Chair is intended for use by professional dental practitioners ranging in stature, while providing treatment to dental patients in a dental operatory environment. It is designed to properly position the patient for said treatments.

Transportation and Storage Conditions

Transportation and Storage Temperature -30°C to 60°C (-22°F to 140°F)
Relative Humidity 10% to 90% (non-condensing)
Atmospheric Pressure 50 kPa to 106 kPa (7.3 psig to 15.3 psig)

Operating Conditions

Ambient Operating Temperature Range 15°C to 35°C (59°F to 95°F)
Relative Humidity 10% to 90% (non-condensing)
Atmospheric Pressure 69 kPa to 106 kPa (10 psig to 15.3 psig)

Disposal of Equipment

At the end of product life, the unit, accessories and other consumable goods may become contaminated from normal use. Consult local codes and ordinances for proper disposal of this equipment and other consumable goods.

FCC Compliance Statement

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense. Any modifications to this product may violate the rules of the Federal Communications Commission and make operation of the product unlawful.

Midmark's Online Help

All service, user, setup, and accessories documents for Midmark products are available in the Midmark Technical Library at midmark.com.

Calling for Service

If service is required, contact your authorized Midmark dealer. When calling for service you must provide the following information:

- Model / Serial Number
- Date of purchase
- Symptom(s) of malfunction

Important Information

Electromagnetic Interference

This dental chair is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this product:

- Remove interfering device from room
- Plug dental chair into isolated circuit
- Increase separation between dental chair and interfering device
- Contact Midmark if interference persists

Wifi and Ethernet Connections

- A remote connection with the myMidmark app is available for 30 minutes after powering on the dental chair.
- Disconnect and reconnect power to the dental chair if a connection to the myMidmark app is needed to modify network settings or to re-register the device.
- This dental chair includes an Ethernet connection and a Wi-Fi module that allows the device to connect to the local network.
- In rare instances, the dental chair may be unable to connect to a wireless network at power-up due to operating location/environment frequency interference on certain wireless channels. Moving the dental chair slightly may resolve this issue. Note: Interference conditions may change. Connecting the dental chair via Ethernet is another option.
- Please consult your IT department/contact to verify the proper Wi-Fi coverage in the area where the device will be used.
- Connecting the dental chair to an IT-network could result in previously unidentified risks to patients, operators or third parties and the system integrator should identify, analyze, evaluate and control these risks.

Important Information

Unregistering a Chair

The Midmark chair can be unregistered by calling Midmark.

- An unregistered device will be inactivated in the Midmark cloud services. Past activity and data will not be available to any new customer associated to that device. Future activity or data is not associated with a prior customer account.
- Device activity and data is retained in the cloud for a prior customer so that historical data may be retrieved by that customer.

Cybersecurity

Connecting the dental chair to a local network could result in previously unidentified cybersecurity risks that the organization responsible for the network should identify, analyze, evaluate, and control. The product captures relevant system and security logs to identify date, time, and source of user activity. A software bill of material (SBOM) is available upon request.

It is the responsibility of the user's organization to:

- Identify, analyze, evaluate, and control IT-network risks associated with this connection.
- Ensure the network provides the required characteristics and security measures defined in the device documentation.
- Review any network configuration changes, updates, or the addition/removal of connected equipment, as these may introduce new risks requiring re-assessment.
- Maintain access control, monitoring, and patching practices aligned with the facility's cybersecurity policy.

If a cybersecurity vulnerability or incident is detected while using the dental chair, or for any other cybersecurity needs, contact Midmark Technical Service for assistance (see Calling for Service).

Software Update

Within myMidmark.com, a dental chair software update may be available, and the update may be scheduled to be applied to the dental chair either immediately or for a specific time.

- If the dental chair is being used during the scheduled time, the software update will be canceled and will need to be rescheduled.
- The dental chair is considered in use if the chair back is reclined at least 45 degrees, a connected dental light is illuminated, or delivery system master switch is on.

Important Information

Caution

The Midmark Dental Chair has not been tested in all environments that may contain unique medical emitters such as electrocautery, MRI, electrosurgical units, and diathermy devices.

Federal law restricts this device to sale by or on the order of a licensed dental practitioner.

Note

Service performed on fuses and / or power cords shall be performed by a Midmark recommended Service Technician.

Note

Service and maintenance should not be performed by the user while in use with a patient.

Note

By connecting the chair into the 029-11724-00 power supply it creates a medical electrical system.

Note

Assembly of medical electrical systems and modifications during the service life of this equipment require evaluation to the requirements of the 60601 series of medical electrical equipment safety standards.

Important Information

WARNING

No modification of this equipment is allowed.

When using dental systems or patient chairs with other movable equipment, ensure all equipment is stationary before operation to avoid collisions or injury. Maintain a safe distance between devices, secure movable equipment to prevent unintended movement and regularly inspect for proper functionality. Train all personnel in safe operation procedures to minimize risks. Failure to follow these precautions may result in injury or equipment damage.

Use of the dental chair adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the dental chair and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by Midmark of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of the Dental Chair, including cables specified by the manufacturer. Otherwise, degradation on the performance of this equipment could result.

Equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen, air, or nitrous oxide.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

During transportation, the chair must be at its lowest height and all attachments secured in their lowest and most central position possible.

Failure to comply may result in injury and / or damage to equipment.

This product can expose you to chemicals including Nickel (Metallic), which is known to the state of California to cause cancer.

For more information go to www.p65warnings.ca.gov

Important Information

WARNING

Only connect approved Midmark dental chair, delivery systems, lights, monitor and power supply when assembling the dental operatory (medical electrical system).

Midmark-approved products listed below:

CH-1000	Midmark Dental Chair
DS-1100-001	Midmark Left-Right Chair-Mounted Doctor's Dental Delivery System
DS-1200-001	Midmark 12 O'clock Column/Cabinet-Mounted Doctor's Dental Delivery System
ADS-1100-001	Midmark Chair Mounted Assistant's Dental Delivery System
ADS-1200-001	Midmark Cabinet Mounted Assistant's Dental Delivery System
LT-1100-001	Midmark Left-Right Chair-Mounted Dental LED Light
LT-1200-001	Midmark Cabinet-Mounted Dental LED Light
LT-1300-001	Midmark Wall-Mounted Dental LED Light
LT-1400-XXX	Midmark Ceiling-Mounted Dental LED Light
LT-1500-XXX	Midmark Track-Mounted Dental LED Light with Monitor
LT-1600-XXX	Midmark Track-Mounted Dental LED Light
9A738002	Midmark Left-Right Chair-Mounted Dental Monitor
029-11724-00	Midmark Dental Power Supply

Other combinations of equipment have not been evaluated for safety.

Important Information

Authorized Representatives

WARNING

In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.

Customers in Colombia should direct all questions, incidents and complaints to the Midmark Authorized Representative:

Emergo Colombia
Calle 76 No. 4-28 of. 803
Bogotá - Colombia 110221
Tel: (57) 313 4419785
CC: LST.AUS.ProgramCommColombia@ul.com

Customers in Costa Rica should direct all questions, incidents and complaints to the Midmark Authorized Representative:

Emergo Costa Rica S.R.L.
Edificio Ribleh 2 Piso, Diagonal a la Biblioteca Municipal
Heredia, Belén
La Ribera, Costa Rica
Tel: +(506) 2201-14 53
LAvigilance@ul.com

Customers in Panama should direct all questions, incidents and complaints to the Midmark Authorized Representative:


TECNIDENTAL PANAMA S.A.
VILLA MAGNA, VILLA ZAITA
OFFICE B8
PANAMA
Tel: 507 6416-2235
+507 6492-0884
tecnidentalpanama@gmail.com
<https://tecnidentalpa.com/> [tecnidentalpa.com]

Customers in Mexico should direct all questions, incidents and complaints to the Midmark Authorized Representative:

(Exclusive Dental Dealer)
Operadora Nac e Int (ONIPO)
Av. Tepeyac 4437
Col. Prados Tepeyac
Zapopan, Jalisco 45050
México
Tel: 52.3331240542
Cell: +52 33 3570 4299
<https://onipo.dental/> [onipo.dental]

Product Overview


Component Location and Description

Product Features

POWER ON Light

The chair has a POWER ON light that illuminates when electrical power is supplied to the chair.


**POWER ON
Light**

Red	Warning immediate action required
Yellow / Orange	Maintenance required auto restart of chair
Green	Ready for use

Product Features

Collision Protection System

Midmark Dental Chairs have a Collision Protection System, with Active Sensing Technology™ that is designed to prevent damage to personnel or equipment. If a chair cover hits an object as it is lowering, it stops and raises slightly to disengage contact; the chair also beeps until the control button is released. Normal operation can then be resumed.

ArtFull

Product Operations

WARNING

Ensure all dental equipment, patient chairs, and movable devices are stationary before use to prevent collisions or injury. Keep devices at a safe distance. Failure to comply may cause injury or damage.


Foot Control

Controls can be used to quickly and accurately position the dental chair.

- Control devices are membrane-covered keypads.
- Connects to the rear of chair base by a cord.

Foot Control Functions

- Home button is programmed to position chair back to original position.
- Stop buttons stops motion of chair.
- Buttons to activate Up and Down functions for both Chair and Back.
- Numbered buttons (1-3) are programmable for positioning patient in common used positions.

ArtFull

Operation


WARNING

Ensure all dental equipment, patient chairs, and movable devices are stationary before use to prevent collisions or injury. Keep devices at a safe distance. Failure to comply may cause injury or damage.


Home Functions

Pushing the Home button takes the dental chair back position to all the way up and chair seat slightly above parallel. If the Collision Protection System, with Active Sensing Technology™ is press during Home function movement it will stop the Home function movement and the base will retract up slightly. This may also occur if a chair cover hits an object as it is lowering.

Stop Button Functions

Pushing Stop Button make dental chair stop any motion it is functioning immediately.

Operation


WARNING

Ensure all dental equipment, patient chairs, and movable devices are stationary before use to prevent collisions or injury. Keep devices at a safe distance. Failure to comply may cause injury or damage.

Back Functions

Pushing the Back buttons will make adjustments to the chair back. Press and hold the appropriate button on the foot control until the desired position is reached. If the Collision Protection System, with Active Sensing Technology™ is press during back function movement it will stop the back function movement and the base will retract up slightly. This may also occur if a chair cover hits an object as it is lowering; the chair also beeps until the control button is released.

- All control devices are membrane-covered keypads.
- All control devices invoke the same chair movements.

Back Up / Down

To Raise or Lower the Back.

Press and hold the desired Back Up or Back Down button on any control device. The back stops moving when the button is released or it reaches its range of motion limit.

Operation


WARNING

Ensure all dental equipment, patient chairs, and movable devices are stationary before use to prevent collisions or injury. Keep devices at a safe distance. Failure to comply may cause injury or damage.

Base Functions

Pushing the Base buttons will make adjustments to the chair base. Press and hold the appropriate button on the foot control until the desired position is reached. If the Collision Protection System, with Active Sensing Technology™ is press during base function movement it will stop the movement and the base will retract up slightly. This may also occur if a chair cover hits an object as it is lowering; the chair also beeps until the control button is released.

- All control devices are membrane-covered keypads.
- All control devices invoke the same chair movements.

Chair Up / Down

To Raise or Lower the Chair.

Press and hold the desired Chair Up or Chair Down button on any control device. The chair stops moving when the button is released or it reaches its range of motion limit.

Programming Chair Positions

WARNING

Ensure all dental equipment, patient chairs, and movable devices are stationary before use to prevent collisions or injury. Keep devices at a safe distance. Failure to comply may cause injury or damage.

To program a button to a specified chair position.


- 1) Use the arrow buttons to move the chair to desired position.
- 2) Hold the desired Program Button (1, 2, 3) down until 3 beeps are heard (around 3 seconds). Chair is now programmed for that Program Button to the assigned position.

Note

The chair can be stopped any time during a programmed positioning sequence by pressing any other button on the control device.

Note

The maximum base height allowed for saving a preset position is less than 66 cm (26 inches.)

Chair Lockout

WARNING


Ensure all dental equipment, patient chairs, and movable devices are stationary before use to prevent collisions or injury. Keep devices at a safe distance. Failure to comply may cause injury or damage.

Lockout

To lock movement of chair.

Press and hold the Stop and Back Chair Down buttons simultaneously for 3 seconds a beep will sound when chair is locked. This will prevent controls and movement of chair.

- Control functions are locked when an operator presses any control button and a beeping sound is heard with no movement of chair.

Unlock

To unlock and allow movement of chair.

Press and hold the Stop and Back Chair Down buttons simultaneously for 3 seconds a continuous beeping will occur until chair is unlocked. Press any of chair controls to test movement of chair.



Product Operations

Heat and Massage for Patient's Control

- A patient can use hand control for heat or massage which can be used separately or together for patient comfort.

Hand Control

- Connects to the Chair by a cord.
- Buttons to activate the Heat and Massage functions.

Heat and Massage Control rests on magnetic mount on chair base when not in use

ArtQuart

ArtFull


Heat and Massage (Optional)

Caution

- Ensure that the surface of the back upholstery is free from mechanical damage prior to each application.
- Do not penetrate back upholstery with sharp objects, such as needles, knives etc.
- The heated back upholstery has not been evaluated for safety when used in combination with other heat sources such as heated blankets. No other heat sources should be used.

Chair Back

ArtQuart

WARNING

- Do not use the heat and massage feature in conjunction with high frequency electrosurgical equipment, such as hyfrecators, etc. Failure to comply may result in electric shock or burns to the patient.
- Warming transdermal medications (patches) can increase drug delivery, resulting in possible harm.

The back upholstery heats from 24 °C (75 °F) to 34 °C (93 °F) in approximately 25 - 28 minutes. To reach this temperature, the heater must operate continuously over this duration. The Heat and Massage components are embedded in the back upholstery and are operated by a patient with the hand control. Heat or Massage can be used separately or together at the same time.

Heat and Massage Hand Control

ArtQuart

All buttons on the hand control activate with touch and release operation. Just touch and release, do not hold down a button. Heat and massage functions each have timers that turn off automatically after 20 minutes without user input.

Hand Control Details

Note

Indicator lights beside each button illuminate when the function is ON.

Massage Mode Selections

PULSE Button

- Press the PULSE button to select the Pulse Massage mode.
- Press the PULSE button once for the LOW intensity pulse massage, one light will turn on.
- Press the PULSE button twice for the MEDIUM intensity pulse massage, two lights will turn on.
- Press the PULSE button three times for the HIGH intensity pulse massage, three lights will turn on.
- Press the PULSE button a fourth time for the pulse massage to turn on.

Message Indicator Lights

WAVE Button

- Press the WAVE button to select wave massage mode.
- Press the WAVE button once for the LOW intensity wave massage, one light will turn on.
- Press the WAVE button twice for the MEDIUM intensity wave massage, two lights will turn on.
- Press the WAVE button three times for the HIGH intensity wave massage, three lights will turn on.
- Press the WAVE button a fourth time for the wave massage to turn off.
- Use SPD arrows to increase or decrease how fast massage motors turn ON and OFF.

Power Button

- Press the PWR button once, the HEAT LOW setting will turn on.
- Press the PWR button again, all functions will turn OFF.

Note

After approximately 20 minutes with no user interaction, massage functions turn OFF automatically.

Hand Control Details

Note

Indicator lights beside each button illuminate when the function is ON.

Heat Mode Selections

PULSE Button

- With the PWR button pressed once, the HEAT LOW light will be on. Leave the heat on or press the HEAT button twice to turn it off completely.


WAVE Button

- With the PWR button pressed once, the HEAT LOW light will be on. Leave the heat on or press the HEAT button twice to turn it off completely.

Heat Indicator Lights

HEAT Button

- With the PWR button pressed once, the HEAT LOW light will turn on.
- Press the HEAT button again for HIGH heat, two lights will turn on.
- Press the HEAT button again to turn the heat completely off.

ArtQuart

Power Button

- Press the PWR button once, the HEAT LOW setting will turn on.
- Press the PWR button again, all functions will turn OFF.


Note

After approximately 20 minutes with no user interaction, heat functions turn OFF automatically.

Adjusting the Headrest

To adjust the Headrest.

- 1) Loosen (turn counterclockwise) the clamp knob and rotate the headset at the pivot points to achieve the desired position.
- 2) Tighten (turn clockwise) the clamp knob to lock the headrest in the desired position.
- 3) Pull up or push down on the headrest, sliding the mounting bar up or down to the desired height.




ArtFull

Positioning for Wheelchairs

The Double-Articulating Gliding Headrest can be used to accommodate patients in wheelchairs.

To position the headrest to accommodate a wheelchair:

- 1) Slide the headrest mounting bar up until it is free from the Dental chair.
- 2) Rotate the headrest 180 degrees and slide the mounting bar into the backrest. Ensure that the mounting bar is pushed completely into the backrest.
- 3) Position headrest and Dental chair to its full back-up position.
- 4) The patient should be positioned so the wheelchair and the Dental chair are back-to-back. Adjust the height of the chair by moving the chair up or down for comfort of patient.
- 5) Ensure that the wheelchair wheels are locked.

Operation

Adjusting Back Height

The height of the back is only adjustable on dental chairs with the Adjustable Back feature.


Back Raise / Lower

To Raise or Lower the Back.

A) While holding the back section in place, squeeze the Back Release Handle.

B) Raise or lower the back section into position.




C) Release handle to lock.

Rotating the Chair

To rotate the Chair.

- Push the brake handle to unlock symbol and rotate the chair at the pivot hub to achieve the desired position.
- Push the brake handle to lock symbol to Lock.

Rotation limits differ for different Chair Models		
Model	° of rotation from center line	
Chair Only	30° Left and 30° Right	shown
L/R Chair	30° Left and 30° Right	not shown

Barriers, Cleaning, and Disinfecting

Clinical Contact Surfaces

Clinical Contact Surfaces are most likely to become contaminated with pathogens during dental treatment through touch or transfer surfaces.

(Source: Centers for Disease Control and Prevention. Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care. Centers for Disease Control and Prevention, U.S. Dept of Health and Human Services; October 2016.)

Non-Clinical Contact Surfaces

Non-Clinical Surfaces carry less opportunity of cross contamination than clinical contact surfaces and can be cleaned with soap and water or cleaned and disinfected if visibly contaminated with blood.

(Source: Centers for Disease Control and Prevention. Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care. Centers for Disease Control and Prevention, U.S. Dept of Health and Human Services; October 2016.)

Barrier Protection

Midmark recommends applying legally-marketed disposable barriers to all dental chairs and headrest. Only use barrier material that is intended for use with dental chairs and headrest. Follow barrier manufacturer instructions for proper use of these products. The use of uncontaminated barriers significantly reduces the need for cleaning and disinfecting equipment in-between patients, thereby extending the life of the equipment.)

General Purpose Cleaning of Non-Clinical Surfaces

Midmark recommends using cleaners that are appropriate for the situation, such as warm water and mild detergents for housekeeping surfaces. The clinical contact surface instructions should be followed if there is visible contamination present.

Barriers, Cleaning, and Disinfecting

Cleaning and Disinfecting Assistance

Local, State and federal guidelines should always be followed and considered when evaluating cleaning and disinfecting protocols. Below are suggested resources.

Association for Dental Safety (ADS): https://www.myads.org	Dept. of Health and Human Resources Centers for Disease Control and Prevention (CDC): https://www.cdc.gov
American Dental Association: https://www.ada.org	

Cleaning and Disinfecting

In addition to the use of barriers, Midmark recommends using an EPA registered cleaner/disinfectant to be used on all clinical contact surfaces that may be contaminated during dental procedures. Follow cleaner/disinfectant manufacturer instructions for proper use of the product. Proper Personal Protection Equipment (PPE) should be evaluated before use, and extreme caution when handling blood or bodily fluid.

Equipment Alert

The upholstery is resistant to most medicinal-type stains, but may be damaged by solvents and dyes. Immediately remove any fluids spilled on the upholstery. To clean and disinfect the upholstery follow the material manufacturers recommendations. Do not use Metrex CaviWipes™ or Cavicide™ to clean upholstery surfaces.

Midmark recommends the use of the below cleaner/disinfectants, which are low alcohol-based:

- Hufriedy Group AdvantaClear™
- Metrex CaviWipes™ or CaviCide™

All product names, trademarks, and registered trademarks are the property of their respective owners. Use of these names does not imply any affiliation with or endorsement by them.

- Surfaces should be cleaned and disinfected between patients and at the end of the day.
- Care should be taken to avoid excessive application and pooling if liquid spray, although Midmark recommends the use of wipes over liquid spray.
- After each cleaning/disinfecting application, complete a visual inspection to confirm that no additional contamination remains, and repeat as necessary.
- Cleaner/Disinfectants can leave a film residue on equipment, so all treated areas should be thoroughly cleaned with mild soap and warm water at least once daily. This routine washing helps minimize the accumulation of chemical residues that could harm the equipment over time.

Equipment Alert

Cleaners used outside of recommendations could cause serious degradation of equipment, like peeling paint, cracked plastic, etc.

Specifications

Specification	
Base Weight of Chair:	138.4 kg (305 lbs)
Shipping Weight:	190.5 kg (420 lbs)
Patient Weight (max):	226.8 kg (500 lbs)
Power Cord Length:	1.52 or 3.66 meters (5 or 12 feet)
Electrical Rating:	120 V~ 60Hz, 8A
Fuses:	F1 & F2 (Line Input) F 8A H 250 V
	F3 (Motors) T 6.3A H 250 V
Mode of Operation:	Intermittent Operation
Duty Cycle:	30 seconds ON / 5 minutes OFF, Base and Back Functions
Classifications:	
Protection Against Electric Shock:	Class I Type B Applied Part Type BF Applied Part, Heated Back Upholstery only
Ingress Protection (IEC 60529):	IPX0, Chair IPX1, Foot Control IPX2, Heated Back Upholstery

Certification and Compliance	
Certifications:	ANSI/AAMI ES60601-1:2005 & A1:2012 & A2:2021 CAN/CSA-C22.2 No. 60601-1:14 (R2022) IEC 60601-1-6:2010+A1:2013+A2:2020 CAN/CSA-C22.2 NO. 60601-1-6:11/A2:21 (R2021) IEC 80601-2-60:2019 CSA C22.2 No. 80601-2-60:23
Compliance:	IEC 60601-1-2:2014+A1:2020 IEC TS 60601-4-2:2024

Accessories

Accessory Name:	Part Number:
Cuspidor	9A738001
Left-Right Monitor	9A738002

Wi-Fi, Ethernet, Cybersecurity

WiFi Connection Conditions

The Dental Chair has a remote connection with Midmark Phone Application.

The Dental Chair includes an Ethernet connection and a Wi-Fi module that allows the device to be connected to the local network.

In rare instances, the Dental Chair may be unable to connect to a wireless network at power-up due to operating location/environment frequency interference on certain wireless channels.

Moving the Dental Chair slightly may resolve this issue. Note: Interference conditions may change.

Connecting the Dental Chair via Ethernet is another option.

Please consult your IT department/contact to verify the proper Wi-Fi coverage in the area where the service will be used.

Cybersecurity

Connecting the Dental Chair to a local network could result in previously unidentified cybersecurity risks that organization responsible for the network should identify, analyze, evaluate, and control. The product captures relevant system and security logs to identify date, time, and source of user activity. A software bill of material (SBOM) is available upon request. The Dental Chair communicates over a local area network (LAN) or Wi-Fi to enable device monitoring, configuration, and software maintenance functions.

Characteristics of the Dental Chair data include:

- Primary flow: Dental Chair ↔ Midmark Connectivity Module ↔ Midmark Cloud Services (when enabled).
- Directionality: All telemetry and configuration data originate from the device; remote commands are limited to authenticated maintenance functions.
- Data type: Equipment status, diagnostic logs, software update, and non-patient operational information only.
- No patient-identifiable information is transmitted. The device continues normal treatment operation even when the network or cloud connection is unavailable.

Contact Midmark Technical Service for assistance if a cybersecurity vulnerability or incident is detected or for any other cybersecurity needs.

Wi-Fi, Ethernet, Cybersecurity

Connectivity Options Chart		
Network Protocols		Link Layer: IEEE 802.3, IEEE 802.11 Internet Layer: IP Transport: TCP, UDP App Layer: TLS, MQTT, NTP, HTTP, etc.
Network Data	LAN Physical Connection	RJ-45
	LAN Cabling	Ethernet (min. CAT 5)
	LAN Max Speeds	10/100 Mbps
	LAN Protocols	IEEE 802.3u
	Wi-Fi Connection	The device supports connectivity to Wi-Fi using the following standards: <ul style="list-style-type: none"> - 802.11a - 802.11b - 802.11g - 802.11n - 802.11ac
	Operational Frequency	Connects to a 2.4 GHz or 5.0 GHz
	Wi-Fi Security Standards	When connecting to Wi-Fi, the device uses the following security standards: <ul style="list-style-type: none"> - WPA2-PSK - WPA2 PEAP/MSCHAPv2 - WPA2 EAP-TLS
	Network Interface Card Parameters	Not Applicable

* Any network administration would be dependent on each site. The fields are not predefined and are intended to be customized as needed.

Wi-Fi, Ethernet, Cybersecurity

<i>Network-Related Hazardous Situations and Mitigations</i>			
Network Deviation	Possible Effect	Hazardous Situation	Mitigation / Residual Risk
Network outage or disconnection.	Remote diagnostics or updates unavailable.	Delay in maintenance actions.	Device operates locally without safety or performance degradation.
Excessive latency or packet loss.	Time-sensitive telemetry delayed.	Delayed fault reporting.	Local indicators and logs remain active; safety and essential performance unaffected.
Incorrect time source (unsynchronized clock).	Event timestamps misaligned.	Inaccurate maintenance records.	Internal real-time clock maintains correct operational timing; does not affect therapy functions.
Unauthorized network access attempt.	Possible exposure of service interface.	Potential cybersecurity incident.	WPA2/WPA2-Enterprise, TLS 1.2+, and mutual authentication reduce likelihood; logs available for investigation.
<p>NOTE: In all cases, the safety of the Dental Chair is maintained. Network functions are supportive but not safety-critical.</p>			

EMC - Manufacturer's Declaration and Guidance

Guidance and manufacturer's declaration – electromagnetic emissions

The Midmark Dental Chair is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Dental Chair should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11 edition 6.2 (2019)	Group 1	The Midmark Dental Chair uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11 edition 6.2 (2019)	Class A	The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
Harmonic emissions IEC 61000-3-2	NA	Not applicable. The Midmark Dental Chair is 120 VAC 60 Hz and does not exceed 220 VAC.
Voltage fluctuations/flicker emissions IEC 61000-3-3	NA	

Recommended separation distances between portable and mobile RF communications equipment and the Midmark Dental Chair

The Midmark Dental Chair is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Midmark Dental Chair can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the Midmark Dental Chair as recommended below, according to the maximum output power of the communications equipment.

Radiated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \times \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.69	3.69	7.38
100	11.67	11.67	23.34

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity

The Midmark Dental Chair is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Dental Chair should assure that it is used in such an environment.

Immunity Test	Compliance Level	Electromagnetic environmental - guidance
Electrostatic Discharge (ESD): IEC 61000-4-2 edition 2.0 (2008)	± 8 kV Contact; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM Fields: IEC 61000-4-3 edition 3.2 (2010)	3 V/m (80 MHz -2.7 GHz) 80% AM at 1kHz	
Proximity Fields from RF Wireless Communications Equipment: IEC 61000-4-3 edition 3.2 (2010)	IEC 60601-1-2, Table 9	
Electrical Fast Transient / Burst: IEC 61000-4-4 edition 3.0 (2012)	± 2 kV @ 100kHz Rep. (AC Mains) ± 1 kV @ 100kHz Rep. (I/O Ports)	Mains power quality should be that of a typical commercial or hospital environment.
Surge: IEC 61000-4-5 edition 3.1 (2017)	± 0.5 kV, ± 1 kV L-L ± 0.5 kV, ± 1 kV, ± 2 kV L-PE	Mains power quality should be that of a typical commercial or hospital environment.
Conducted RF: IEC 61000-4-6 edition 4.0 (2013)	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	
Power frequency (50/60 Hz) magnetic field: IEC 61000-4-8 edition 2.0 (2009)	30 A/m, 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines: IEC 61000-4-11 edition 2.1 (2017)	0% UT for 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°. 0% UT for 1 cycle 70% UT for 30 cycles 0% UT for 300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Midmark Dental Chair requires continued operation during power mains interruptions, it is recommended that the Midmark Dental Chair be powered from an uninterruptible power supply or a battery.
Proximity Magnetic Fields: IEC 61000-4-39 edition 1.0 (2017)	IEC 60601-1-2, Table 11	

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The Midmark Dental Chair is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Dental Chair should assure that it is used in such an environment.

Immunity Test	Compliance Level	Electromagnetic environmental - guidance
---------------	------------------	--

Portable and mobile RF communications equipment should be used no closer to any part of the Midmark Dental Chair including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Recommended separation distance:

Conducted RF IEC 61000-4-6 edition 4.0 (2013)	3 V	$d = 1.2 \times \sqrt{P}$
Radiated RF IEC 61000-4-3 edition 3.2 (2010)	3 V/m	$d = 1.2 \times \sqrt{P}$ 80 MHz to 800 MHz
		$d = 2.3 \times \sqrt{P}$ 800 MHz to 2.5 GHz

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation in meters (m).

Field strength from fixed RF transmitters, as determined by the electromagnetic site survey. a) Should be less than the compliance level in each frequency range. b) Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location is which the Midmark Dental Chair is used exceeds the applicable RF compliance level above, the Midmark Dental Chair should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Midmark Dental Chair.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Equipment Alert

Use of accessories, transducers and cables other than those specified or provided could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warranty Information

Scan the following QR code for access to the Dental Chair warranty.

Warranty Registration

Scan the following QR code for access to the Dental Chair warranty registration.

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