Cuspidor

User Guide


*English
Español
Français*

For Models: 9A738001

Serial Number Location

(The information below is required when calling for service.)

Serial Number Label

Table of Contents

Important Information

Symbols Glossary.....	4
Intended Use	5
Transportation / Storage	5
Operating Conditions	5
Disposal of Equipment	5
Midmark's Online Help.....	5
Calling for Service.....	5
Electromagnetic Interference.....	6
Authorized Representatives.....	7

Operation

Item Location	8
Bowl rinse and adjustment	9
Cup fill and adjustment.....	10
Flow time adjustment	11
Bowl Rotation Drag Adjustment.....	12

Maintenance

Clinical Surfaces	13
Non-clinical Surfaces	13
Barriers.....	13
Cleaning and Disinfecting of Clinical Surfaces.....	14
General Purpose Cleaning of Non-clinical Surfaces.....	14
Visual Inspection.....	14
Cleaning and Disinfecting Assistance	14
Cleaning Solid Collector	15
Cuspidor Daily Cleaning and Maintenance.....	16
Specifications.....	17
Replacement Part.....	17

Warranty

Warranty Information.....	18
Warranty Registration.....	18

Important Information

Symbols Glossary

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage.

Note

Amplifies a procedure, practice, or condition.

Refer to manual

Pressure Limit

Manufacturer

Keep dry

Temperature Limit

Catalogue Number

Protective earth ground

Humidity Limit

Serial Number

Proper shipping orientation

Maximum stacking height
(n = palletted units)

Fragile

Important Information

Intended Use

The Midmark cuspidor is intended to provide water for oral rinsing, and a drain to dispose of waste water and oral debris during dental examinations.

WARNING

Do not pour water into the Cuspidor from external sources. It is designed to safely eliminate only the amount of water supplied from the Cuspidor fill and flush operations.

Transportation / Storage

Transportation / Storage Temperature:

-30°C to 60°C (-22°F to 140°F)

Relative Humidity:

10% to 90% (non-condensing)

Atmospheric Pressure:

50 kPa to 106 kPa (7.3 psi to 15.3 psi)

Operating Conditions

Operating Temperature Range:

15°C to 35°C (59°F to 95°F)

Relative Humidity:

10 to 90% (non-condensing)

Atmospheric Pressure:

69 kPa to 106 kPa (10 psi to 15.3 psi)

Disposal of Equipment

At the end of this product's life, the unit, accessories and other consumable goods may be contaminated from normal use. Consult local codes and ordinances for proper disposal of this equipment and other consumable goods.

Midmark's Online Help

All service, user, setup, and accessories documents for Midmark products are available in the Midmark Technical Library at midmark.com.

Calling for Service

If service is required, contact your authorized Midmark dealer. When calling for service, you must provide the following information:

- Model/Serial Number
- Date of purchase
- Symptom(s) of malfunction

Important Information

Electromagnetic Interference

Midmark dental cuspidor components are designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this cuspidor, remove the interfering device from room and / or plug product into an isolated circuit.

- Remove interfering device from room
- Plug operatory into isolated circuit
- Contact Midmark, if interference persists

Note

The Midmark Dental Cuspidor is MR Unsafe and is not recommended for use in a Magnetic Resonance (MR) or MRI environment.

Safety Risk

WARNING

Do not modify this equipment without authorization of the manufacturer.

Should damaged or loose parts be identified, the equipment should be removed from use and serviced. Failure to do so may result in injury.

Equipment is not suitable for use in the presence of a flammable anesthetic mixture.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

This product can expose you to chemicals including Nickel (Metallic), which is known to the state of California to cause cancer.

For more information go to www.p65warnings.ca.gov

Caution

Federal Law restricts this device to sale by or on the order of a licensed dental practitioner.

Important Information

Authorized Representatives

WARNING

In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.

Customers in Colombia should direct all questions, incidents and complaints to the Midmark Authorized Representative:

Emergo Colombia
Calle 76 No. 4-28 of. 803
Bogotá - Colombia 110221
Tel: (57) 313 4419785
CC: LST.AUS.ProgramCommColombia@ul.com

Customers in Costa Rica should direct all questions, incidents and complaints to the Midmark Authorized Representative:

Emergo Costa Rica S.R.L.
Edificio Ribleh 2 Piso, Diagonal a la Biblioteca Municipal
Heredia, Belén
La Ribera, Costa Rica
Tel: +(506) 2201-14 53
LAvigilance@ul.com

Customers in Panama should direct all questions, incidents and complaints to the Midmark Authorized Representative:


TECNIDENTAL PANAMA S.A.
VILLA MAGNA, VILLA ZAITA
OFFICE B8
PANAMA
Tel: 507 6416-2235
+507 6492-0884
tecnidentalpanama@gmail.com
<https://tecnidentalpa.com/> [tecnidentalpa.com]

Customers in Mexico should direct all questions, incidents and complaints to the Midmark Authorized Representative:

(Exclusive Dental Dealer)
Operadora Nac e Int (ONIPO)
Av. Tepeyac 4437
Col. Prados Tepeyac
Zapopan, Jalisco 45050
México
Tel: 52.3331240542
Cell: +52 33 3570 4299
<https://onipo.dental/> [onipo.dental]

Operation

Item Location



Operation

Bowl rinse and adjustment

To rinse the cuspidor bowl...

Press Bowl Rinse Button.

Note: After releasing the button, water will run for approximately 5 seconds before shutting off automatically.

To increase the flow volume...

Rotate bowl rinse flow adjustment counter-clockwise as necessary.

To decrease the flow volume...

Rotate bowl rinse flow adjustment clockwise as necessary.


Operation

Cup fill and adjustment


To fill a cup with water...

Press Cup Fill Button.

Note: After releasing the button, water will run for approximately 5 seconds before shutting off automatically.

ArtSU3

ArtSU3

To increase the flow volume...


Rotate cup fill flow adjustment counter-clockwise as necessary.

To decrease the flow volume...

Rotate cup fill flow adjustment clockwise as necessary.

Operation

Flow time adjustment


Operation

Bowl Rotation Drag Adjustment

Adjusting bowl rotation drag.

Use an Allen wrench to turn set screw clockwise to increase drag and counter-clockwise to decrease drag on bowl rotation.

ArtFull

Maintenance

Clinical Contact Surfaces

Clinical Contact Surfaces are most likely to become contaminated with pathogens during dental treatment through touch or transfer surfaces.

(source: Centers for Disease Control and Prevention. Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care. Centers for Disease Control and Prevention, U.S. Dept of Health and Human Services; October 2016.)

Non-Clinical Contact Surfaces

Non-Clinical Surfaces carry less opportunity of cross contamination than clinical contact surfaces and can be cleaned with soap and water or cleaned and disinfected if visibly contaminated with blood.

(source: Centers for Disease Control and Prevention. Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care. Centers for Disease Control and Prevention, U.S. Dept of Health and Human Services; October 2016.)

Barrier Protection

Midmark recommends applying legally-marketed disposable barriers to all clinician contact surfaces. Only use barrier material that is intended for use with dental equipment. Follow barrier manufacturer instructions for proper use of these products. The use of uncontaminated barriers significantly reduces the need for cleaning and disinfecting equipment in-between patients, thereby extending the life of the equipment.

Maintenance

Cleaning and Disinfecting Assistance

Local, State and federal guidelines should always be followed and considered when evaluating cleaning and disinfecting protocols. Below are suggested resources.	
Association for Dental Safety (ADS): https://www.myads.org	Dept. of Health and Human Resources Centers for Disease Control and Prevention (CDC): https://www.cdc.gov
American Dental Association: https://www.ada.org	
American Dental Association: https://www.ada.org	European Dental Association: https://www.eda-eu.org

Cleaning and Disinfecting of Clinical Surfaces

In addition to the use of barriers, Midmark recommends using an EPA registered cleaner/disinfectant to be used on all clinician controls, or surfaces that may come in contact with dental instruments during dental procedures. Follow cleaner/disinfectant manufacturer instructions for proper use of the product. Care should be taken to avoid excessive application and pooling of liquids. Midmark recommends the use of wipes over liquid spray. Disinfectants can leave a film residue on equipment, so all treated areas should be thoroughly cleaned with mild soap and warm water at least once daily. This routine washing helps minimize the accumulation of chemical residues that could harm the equipment over time.

Midmark recommends the use of the below cleaner/disinfectants, which are low alcohol-based:

- Hufriedy Group AdvantaClear™
- Metrex CaviWipes™ or CaviCide™

All product names, trademarks, and registered trademarks are the property of their respective owners. Use of these names does not imply any affiliation with or endorsement by them.

Surfaces should be cleaned and disinfected between patients and at the end of the day.

- Care should be taken to avoid excessive application and pooling if liquid spray, although Midmark recommends the use of wipes over liquid spray.
- After each cleaning/disinfecting application, complete a visual inspection to confirm that no additional contamination remains, and repeat as necessary.
- Cleaner/Disinfectants can leave a film residue on equipment, so all treated areas should be thoroughly cleaned with mild soap and warm water at least once daily. This routine washing helps minimize the accumulation of chemical residues that could harm the equipment over time.

Equipment Alert

Cleaners used outside of recommendations could cause serious degradation of equipment (e.g., peeling paint, cracked plastic).

General Purpose Cleaning of Non-Clinical Surfaces


Midmark recommends using cleaners that are appropriate for the situation, such as warm water and mild detergents for housekeeping surfaces. The clinical contact surface instructions should be followed if there is visible contamination present

Maintenance

Cleaning Solid Collector

To clean solid collector.

- A) Remove solid collector basket by pulling it straight out.
- B) Clean basket housing.
- C) Install solid collector basket.

Caution
Always dispose of biohazardous debris
according to local regulations.


Maintenance

Cuspidor Daily Cleaning and Maintenance

WARNING

To reduce the risk of cross-contamination between patients, it's essential to clean and disinfect after every dental procedure. Always follow manufacturer's guidelines and recommendations for the proper disinfection process.


Between patients.
Clean and disinfect solid collector and clinical surfaces (See Cleaning and Disinfecting).

End of each day.
A) Clean and disinfect solid collector and clinical surfaces (See Cleaning and Disinfecting).
B) Clean all non-clinical surfaces (See General Purpose Cleaning).
C) Visually inspect cuspidor system (See Visual Inspection).

ArtHalf

Maintenance

Specifications	
Width:	20.6 cm (8 1/8 in)
Length:	45.72 cm (18 in)
Height:	74.17 cm (29 3/16 in)
Weight (Mass):	10.9 kg (24 lbs)
Air Supply Requirements:	Maximum 80 - 100 PSI (5.52 - 6.9 bar)
Water Supply Requirements:	Maximum 30 - 50 PSI (2.1 - 3.4 bar)

Replacement Part	
Solids Collector Basket Part Number 053-1082-00	 <small>ArtSU3</small>

Warranty

Warranty Information

Scan the following QR code for access to the Dental Cuspidor warranty.

Warranty Registration

Scan the following QR code for access to the Dental Cuspidor warranty registration.

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(source : Centers for Disease Control and Prevention. Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care. Centers for Disease Control and Prevention, US Dept of Health and Human Services ; octobre 2016).

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Association for Dental Safety (ADS) : https://www.myads.org	Dept. of Health and Human Resources Centers for Disease Control and Prevention (CDC) : https://www.cdc.gov
American Dental Association : https://www.ada.org	
American Dental Association : https://www.ada.org	European Dental Association : https://www.eda-eu.org

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