
Procenter Delivery Systems

[Includes: Delivery Units and Assistant's Units (All Styles)]

[English](#)
[Español](#)
[Français](#)

Serial Number Locations

Note

Work surface shapes may vary per cabinet, your exact shape may not be shown. Control knobs and adjustments remain the same, or similar, for all units.

Compliance Information and Symbols:

Intended Use	2
Electromagnetic Interference	2
Disposal of Equipment	2
Transportation / Storage /	
Operating Conditions	2
Calling for Service	2
Symbols	3

Product Features:

Work Surface Weight Limit	4
Self-Contained Water System	4
Manual Shut-Off Valves for	
Air and Water Supply	5
Pressure Regulator Valves	5
Master ON / OFF Switch	6
Water Selector Switch	7
Flex Arm (Lock / Unlock)	8
Quick-Connect Ports	9
HVE / Saliva Ejectors	10

Operation and Adjustments:

Operating Handpieces	11
Handpiece Controls	12
Syringe Operation & Adjustments	13
Adjustment Knob Location	14 - 16

Cleaning and Maintenance:

Cleaning & Maintenance	17 - 23
Specifications	24
Replacement Parts	24

EMC Information:

EMC - Manufacturer's Declaration and Guidance	25
---	----

Warranty:

Warranty Information	28
--	----

For installation information go to:
[Technical Library at Midmark.com](#)

Compliance Information and Symbols

Intended Use

Midmark delivery systems are intended to provide dental professionals with air, water and suction to operate dental handpieces, syringes, and Midmark authorized accessories during dental examinations and procedures.

Electromagnetic Interference

Midmark dental operatory components are designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this operatory, remove the interfering device from room and / or plug product into an isolated circuit.

Transportation / Storage / Operating Conditions

Transportation / Storage Temperature:.....23°F to 100°F (-5°C to 38°C)

Relative Humidity:..... 10% to 90% (non-condensing)

Atmospheric Pressure: 7.2 PSI to 15.3 PSI (50 kPa to 106 kPa)

Operating Temperature Range:59°F to 95°F (15°C to 35°C)

Calling For Service

Direct all service inquiries to your authorized Midmark dealer. When calling for service, you must provide the following information:

Model / serial number

Date of purchase

Symptom(s) of malfunction

Disposal of Equipment

At the end of this product's life, the unit, accessories and other consumable goods may be contaminated from normal use. Consult local codes and ordinances for proper disposal of this equipment and other consumable goods.

WARNING

Do not modify this equipment without authorization of the manufacturer.

Caution

Federal law restricts this device to sale by or on the order of a licensed dental practitioner.

Symbols

These symbols may appear on your equipment and/or in the manuals. Warnings and cautions are provided in the manuals where applicable.

WARNING

Indicates a potentially hazardous situation which could result in serious injury if not avoided.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury if not avoided. It may also be used to alert against unsafe practices

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage if not avoided.

NOTE

Amplifies a procedure, practice, or condition.

Type B,
Applied Part

Type BF,
Applied Part

Protective
Earth Ground

Proper Shipping
Orientation

Fragile

Catalogue Number


Serial Number


Manufacturer


Refer to instruction
manual/booklet

Pressure
Limit

Temperature
Limit

Humidity
Limit

Maximum stacking
height (Refer to “n”
number on package.)

Keep Dry

Product Features


Work Surface Weight Limit


Caution

Exceeding the design weight limit may result in personal injury or damage the equipment.

Maximum load on Flex Arm mounted units is 7 lbs (3.2 KG)

DA2075

DA24701i

Maximum load for cabinet mounted arms is 40 lbs (18 KG)

Self-Contained Water System

The Self-Contained Water System gives you control of the quality of water used in your delivery system. It also provides a means for cleaning the system's internal tubing.

Caution

Residual air pressure may cause the water bottle to pop off and spill or spray the contents if not released.

Before you begin filling the water system...

- A) Turn master switch OFF..
- B) Press syringe air button and hold until all residual air pressure is expelled. (You can hear the air being expelled.)


NOTE

Distilled water is not required. Water used must meet all local requirements for drinking water.

To fill the Self-Contained Water System bottle...

- A) Unscrew water bottle.
- B) Fill bottle with water.
- C) Be sure Gasket is in place (in Mounting Cap).
- D) Screw bottle into Mounting Cap.


Note: If Gasket is improperly sealed, you will hear a hissing noise when the system is turned ON.

DA101201i

Manual Shut-Off Valves for Water and Air Supply


Manual shut-off valves allow you to stop the air and/or water supply at the point of input to the operatory. This is recommended during extended periods of non-use (ex. vacation, holidays, etc.), or in the event of an equipment malfunction.

Pressure Regulator Valves

Pressure regulator valves allow you to control the air and water pressure supplied to the instruments of the delivery.

Shut-off Valves
Rotate clockwise 1/4
turn to Shut off.


To Adjust the Pressure Regulators...

- Pull up knob and turn to adjust.
- Watch regulator gauge as you turn knob to achieve desired setting
- Push knob DOWN to lock.

Turn Clockwise
for more pressure

Turn Counter-Clockwise
for less pressure

Recommended Settings:

City Water Regulator Gauge Setting	30 psi
Air Pressure Regulator Gauge Setting	80 psi

Equipment Alert



Delivery components were designed to operate at the recommended settings. Poor performance or damage to the equipment may result if recommended settings are not maintained.

DA3820i

Master ON / OFF Switch

The Master ON / OFF Switch must be ON to operate the delivery system.

To turn the Delivery System ON or OFF...
Move the Master Switch to the desired position.





Water Selector Switch

The Water Selector Switch allows you to choose the water source for the delivery system. You may select either City Water (tap water), or water from the Self-Contained System's water bottle.

To select a water source....

Toggle the switch to select the desired water source.



DA229401i

Flex Arm (Lock / Unlock)

A lock mechanism exists on the Flex Arm to prevent unexpected downward movement, which could result in injury to the patient or staff. The Master Switch must be ON for the air-powered Flex Arm lock to work.


The Flex Arm can still be raised upwards when locked.

Procenter Units

To set the Flex Arm lock to a new height...




- A) Press and hold either Lock/Unlock Button.
- B) Raise flex arm slightly (to disengage the lock).
- C) Lower to desired height.
- D) Release the Lock/Unlock Button to engage the lock.

Quick-Connect Ports

Quick-Connect Ports (air [yellow] and water [blue]) allow you to quickly connect accessories to the delivery system.

Each Quick-Connect Port has a hinged clip that securely affixes the connecting line to a raised ring found on each port opening.


DA229602i

HVE and Saliva Ejectors

Before you begin...

- A) Turn facility vacuum system ON.
- B) Move lever to towards the OPEN position to activate and moderate suction for the HVE / Saliva Ejector.


Operation and Adjustments

Operating Handpieces

Handpieces can deliver water and air dependant on the setting of the Water Selector footswitch.

To Operate handpiece...

- A) Turn Master Switch ON
- B) Move Water Selector Switch to desired setting.
- C) Remove handpiece from holder.
- D) Press foot pedal to deliver the selected drive air/water to the handpiece.


Da100202i

Handpiece Controls


Each handpiece has its own Drive Air control knob and its own Coolant Water control knob.

There is also a Coolant Air control knob that adjusts the coolant air **for all handpieces**.

To adjust the Coolant Water...


Rotate the Coolant Water control knob for the desired handpiece, until the desired flow is achieved.

To adjust the Drive Air pressure...

Rotate the Drive Air control knob for the desired handpiece, until the desired air pressure registers on the pressure gauge.

Note: Refer to the handpiece manufacturer's literature for the recommended pressure setting.

DA2531

To adjust the air/water Spray Pattern...

Rotate the Coolant Air control knob until the desired pattern is achieved.



Note: If your handpiece has an internal coolant air connection this adjustment will have no effect.

Syringe Operation and Adjustments

This page shows operation and adjustment for a standard 3-way syringe.

Before you begin...

- A) Turn the Master Switch ON.
- B) Position the Water Selector Switch to Water or Air.

To Adjust the Air Spray pressure and Water Spray pattern...

Rotate the syringe Air control knob as desired.

To Adjust the Water stream pressure and Water Spray pattern...


Rotate the syringe Water control knob as desired.

DA25951

Adjustment Knob Location


ProCenter Delivery Units

Adjustment Knob Location - continued...

Assistant's Units


The Assistant's units have screw-type control knobs to adjust the air and water flow to the syringe.

DA2303

Adjustment Knob Location - continued...

Cabinet Mounted Units

Ka998503

Cleaning, Disinfecting and Maintenance

Barriers


Midmark recommends the use of disposable barriers on all clinician controls that may be in contact with clinician hands and fingers during dental procedures. The use of barriers significantly reduces the need for chemical cleaners, thus prolonging the life of the equipment.

Only use barrier material that is intended for use with dental equipment. Midmark recommends the use of an FDA market-cleared barrier such as Pinnacle Cover-all™. Follow barrier manufacturer instructions for proper use of these products.

Cleaning and Disinfecting

In addition to the use of barriers Midmark recommends the use of an EPA registered and FDA market-cleared cleaner/disinfectant such as Cavicide™ to be used on all clinician controls or surfaces that may come in contact with dental instruments during dental procedures.

Follow cleaner/disinfectant manufacturer instructions for proper use of the product. Care should be taken to avoid excessive application and pooling of liquids.

Handpiece Accessories

Only use dental handpiece accessories with the delivery system that are FDA market-cleared and refer to manufacturer's instructions for proper cleaning and disinfecting. Either an autoclavable syringe tip or a single use disposable syringe tip may be used.

Equipment Alert

AUTOCLAVABLE SYRINGE TIP STERILIZATION

The autoclavable syringe tips supplied with the delivery system must be sterilized prior to use with each patient, including initial use. Be sure to thoroughly rinse and clean syringe tips prior to sterilization, any debris may reduce the effectiveness of the sterilization. Recommended sterilization process is steam autoclave. Recommended parameters are 125°C (250°F) and 106 kPa (15 PSI) for 40 minutes at temperature and pressure.

General Purpose Cleaning

For general purpose cleaning, use cleaners that are appropriate for the situation, such as warm water and mild detergents, or a 10% solution of bleach with water.

Visual Inspection

After cleaning, visually inspect the product for deterioration of covers and touch pads. Do not use the delivery system if excessive discoloration, cracking, or other signs of wear are noticeable (see Calling for Service instructions).

Waterline Maintenance

Waterline maintenance is necessary to keep the count of heterotrophic bacteria from rising higher than desired levels. The desired level for a specific location should be determined by any local or regional guidelines. For example, The United States Centers for Disease Control and Prevention (CDC) guideline for heterotrophic bacteria is less than or equal to 500 CFU/mL (colony forming units per milliliter). Midmark recommends keeping this level under 200 CFU/mL.

Treatment can come in many forms. The most popular methods on the market currently are tablets and straw/cartridge based systems. Midmark recommends the use of a straw/cartridge based system that keeps the bacteria levels in check.

Regular monitoring should also take place to ensure that heterotrophic bacteria is not exceeding the desired limit. If the level is higher than desired, a shock treatment of the waterlines should be performed. When performing a shock treatment, be sure to check with the manufacturer of the regular treatment regimen being used to ensure chemical compatibility. Monitoring frequency should be established by your practice. As a suggestion, Midmark would recommend that you begin by monitoring on a monthly basis, and make adjustments to the frequency based on test results.


Per the CDC, routine flushing of the waterlines should be performed between every patient. Extra flushing maybe needed within Midmark equipment when tablets are used. Undissolved tablet particles can gather over time in places within the waterlines, obstructing the line and causing water flow to slow. By flushing the waterlines, water flow is maximized and should push any undissolved particles through.

Cleaning and Disinfecting Assistance

For assistance with cleaning and disinfecting instructions contact the Midmark Technical Service Department at 1-800-Midmark; it is helpful to provide the delivery system model number and serial number when asking for assistance.

Additional information is available from the organizations listed below:	
Organization for Safety & Asepsis Procedures: http://www.osap.org	Dept. of Health & Human Resources, Centers for Disease Control & Prevention (CDC): http://www.cdc.gov
American Dental Association: http://www.ada.org	European Dental Association: http://www.eda-eu.org

Cleaning the Delivery System

DA101201i

Note:

Water must be safe for drinking. Distilled water is not required.

For each new patient...


Replace disposable tips, instruments, etc. and perform a Purging Procedure.

Note: See Purging Procedure described in this manual.

At the beginning of each day...


Fill the Self contained water bottle with fresh water and perform a Purging Procedure.

Note: See Purging Procedure described in this manual.

DA2536i


Procenter Units

DA253501i

Cleaning the Delivery System - At the end of each day...

A) Remove disposable tips, instruments, etc.


DA2536i

B) Clean and disinfect the delivery system (see Cleaning and Disinfecting instructions).

C) Fill the water bottle with fresh water and perform a purging procedure.

Note: See Purging Procedure described in this manual.


DA253501i

DA2399_a

D) Turn Master Switch OFF.
Press and hold the foot control pedal until all pressure is released.

DA2538i



Purging Procedure for the Delivery

NOTE:

The purging procedure removes debris from the tubing to the handpieces and syringe. Performing this procedure frequently may help reduce the accumulation of biofilm on your instruments.

To begin the purging procedure...

- A) Turn Master Switch ON.
- B) Turn Water Selector Switch to bottle setting
- C) Move the foot control switch to the water setting.
- D) Disconnect handpiece from tubing.

Water
Selector
Switch


Caution

Hold the tubing and syringe over a container or drain while flushing.


Flush the tubing to the handpieces ...

- A) Press and hold the foot control pedal for 30 seconds.
- B) Press and hold the flush button for 30 seconds.
- C) Repeat for all tubing to the handpieces.

Flush Buttons


Water

Foot Control Switch

Flush the Syringe Tubing...

- A) Press and hold both syringe buttons (air and water) for 30 seconds.


Air/Oil Separator - Cleaning and Maintaining

Periodically check the fluid level in the air / oil separator container.
When the container is approximately 2/3 full, clean the air/oil separator as shown below.

To clean the air/oil separator...


- A) Turn master ON/Off switch Off.
- B) Remove (unscrew) air/oil separator container.
- C) Dispose of the fluid and saturated gauze.
- D) Clean container and mounting cap.
- E) Install clean gauze and reinstall the container.

DA2629i

Maintaining and Replacing the Regulator Filters


Step 1: Shut off manual valves to turn off water and air supplies to the operatory.


Step 3: Unscrew retainer nut and remove filter.

Retainer Nut

Filter

Filter Cap

Step 2: Unscrew filter cap.

Note: Use 9/16" wrench.

Step 4: Install new filter.
Secure with retainer nut.

Note: Install the filter & retainer nut with the ridged side up (as shown).

Step 5: Reinstall filter cap.

DA3821i


Assistant's Units - Cleaning

To clean the facility vacuum system...

Refer to the instructions provided by the vacuum system's manufacturer.

Note:

Every dental practice is different, and no single cleaner is the best choice for every facility. See the list of organizations that may assist you in choosing the best cleaners available for your practice, in the front of this manual.

To clean the solids collector...

- A) Turn facility vacuum OFF.
- B) Remove lid and basket.
- C) Clean basket and housing.
- D) Reinstall basket and lid.






Caution

Always dispose of biohazardous debris according to local regulations.

DA147102i

Specifications	
Electrical Rating for Elevance and Midmark Ultra Chair:	See power supply output (153941 ...)
Complies to the applicable requirements of:	ES/IEC/EN 60601-1, CAN/CSA C22.2 No. 60601-1 (Safety Standards), EN/IEC 60601-1-2 (EMC Standards)
Electrical Rating for Knight Chair:	24 VAC, 3.2 amp 13 VAC, 2 amp
Complies to the applicable requirements of:	UL/IEC/EN 60601-1, CAN/CSA C22.2 No. 601.1, IEC 60601-1-1
Air Supply Requirements:	Maximum 80-100 PSI (5.52 - 6.9 bar)
Water Supply Requirements:	Maximum 30 - 50 PSI (2.1 -3.4 bar)
Classifications:	Class I, Type B Applied Part, Ordinary Equipment (IPX0), Continuous Operation
Certifications:	Midmark Corporation.....ISO-9001
Optional Accessories:	
Acteon Satelec Scaler	Type BF Applied Part
Acteon Satelec Curing Light	Type B Applied Part
Acteon Sopro Intra-oral Camera	Type BF Applied Part
Bien Air Electric Motor	Type B Applied Part
Dentsply Cavitron Scaler	Type B Applied Part
Handpiece tubing and connectors intended to be used with ISO 7785-1 or ISO 7785-2 compliant air-driven handpieces	Type B Applied Part

Replacement Parts			
These replacement parts can be purchased from your authorized Midmark dealer.			
Saliva Ejector Valve Part Number 029-3975-00		HVE Handpiece Valve Part Number 029-3974-00	
Solids Collector Basket Large Part Number (1) 053-1106-00 (50) 002-10258-00		Syringe Tip Part Number 120514	
Solids Collector Basket Small Part Number (1) 118613 (50) 002-10259-00			
Saliva Ejector Tip Part Number 016-1116-00		Regulator Valve Filter Part Number 053-1166-00	
Saliva Ejector Lever Part Number 016-1113-00		HVE Lever Part Number 016-1109-00	 DA2304-2

WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen, air, or nitrous oxide.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

EMC - Manufacturer's Declaration and Guidance

Guidance and manufacturer's declaration – electromagnetic emissions

The Midmark Procenter Unit dental device is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Procenter Unit dental device should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The Midmark Procenter Unit dental device uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The Midmark Procenter Unit dental device is suitable for use in all establishments, including domestic establishment and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communications equipment and the Midmark Procenter Unit dental device

The Midmark Procenter Unit dental device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Midmark Procenter Unit dental device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the Midmark Procenter Unit dental device as recommended below, according to the maximum output power of the communications equipment.

Radiated maximum output power of transmitter <i>W</i>	Separation distance according to frequency of transmitter <i>m</i>		
	150 kHz to 80 MHz $d = 1.2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \times \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.69	3.69	7.38
100	11.67	11.67	23.34

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity

The Midmark Procenter Unit dental device is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Procenter Unit dental device should assure that it is used in such an environment.

<i>Immunity Test</i>	<i>IEC 60601 Test Level</i>	<i>Compliance Level</i>	<i>Electromagnetic environmental - guidance</i>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for AC and DC power lines I/O lines not tested, all less than 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line(s) to line(s) ±2kV lines(s) to earth	±1kV line(s) to line(s) ±2kV lines(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T for 0.5 cycle) 40% U_T (60% dip in U_T for 5 cycles) 70% U_T (30% dip in U_T for 25 cycles) <5% U_T (>95% dip in U_T for 5 s)	V Dip >30% of U_T for 500ms V Dip < 60% of U_T for 100ms V Dip > 95% of U_T for 5000ms and 10ms	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Midmark Procenter Unit dental device requires continued operation during power mains interruptions, it is recommended that the Midmark Procenter Unit dental device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The Midmark Procenter Unit dental device is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Procenter Unit dental device should assure that it is used in such an environment.

<i>Immunity Test</i>	<i>IEC 60601 Test Level</i>	<i>Compliance Level</i>	<i>Electromagnetic environmental - guidance</i>
			Portable and mobile RF communications equipment should be used no closer to any part of the Midmark Procenter Unit dental device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3V rms 150kHz to 80MHz	3V	$d = 1.2 \times \sqrt{P}$
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	$d = 1.2 \times \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \times \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation in meters (m). Field strength from fixed RF transmitters, as determined by the electromagnetic site survey, ^a Should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location is which the Midmark Procenter Unit dental device is used exceeds the applicable RF compliance level above, the Midmark Procenter Unit dental device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Midmark Procenter Unit dental device.

^b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.

Warranty Information

Midmark Limited Warranty - Dental Products

SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will at Midmark's option repair or replace components of the dental products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this limited warranty is limited to the repair or replacement of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the first retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows: Effective March 1, 2018 these applicable warranty periods are measured from the date of invoice to the original user, shall be as follows:

1. OPERATORY PRODUCTS

- a. Five (5) years for all products (except for the items in (b) through (e)).
- b. Two (2) years for upholstery (chairs and stools).
- c. "KINK-VALVE" module carries a ten (10) year warranty.
- d. The original light bulb on a new light carries a one (1) year warranty.
- e. Accessories not manufactured by Midmark are excluded including but not limited to Bien Air handpiece systems, Dentsply Cavitron scaler, Satelec scaler and curing light, and Sopro cameras.

2. ORAL SURGERY PRODUCTS are warranted for a period of one (1) year.

3. STERILIZER PRODUCTS are warranted for a period of one (1) year.

4. ULTRASONIC CLEANER PRODUCTS are warranted for a period of two (2) years.

5. AIR AND VACUUM PRODUCTS

- a. PowerAir® oil-less compressors – Five (5) years or 3,500 hours of use, whichever occurs first.
- b. PowerVac® and PowerVac® G dry vacuums – Five (5) years or 10,000 hours of use, whichever occurs first (except that the vacuum pump warranty term is ten (10) years or 20,000 hours of use, whichever occurs first).
- c. Classic Series® wet-ring vacuums – Five (5) years or 10,000 hours of use, whichever occurs first.
- d. PowerMax surgical suction – Two (2) years.
- e. Hg5 Series Amalgam Separator - One (1) year. (f) Midmark manufactured accessories - One (1) year.

6. SYNTHESIS™ DENTAL CASEWORK AND ARTIZAN® EXPRESSIONS PRODUCT

- a. Five (5) years for all products and components including door and drawer fronts, casters and slides, except for the items in (b), (c) and (d).
- b. Three (3) years for electrical components such as task lights/LED lights, cords, controls and accessories.
- c. Two (2) years for sliding track monitor mount and components and upholstery. (d) One (1) year for countertops and resin, including accessories.

7. IMAGING PRODUCTS are warranted for a period of two (2) years except for the ClearVision CR reader which is warranted for a period of one (1) year.

8. MIDMARK Replacement Parts and Accessories carry a ninety (90) day warranty

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for the following:

1. defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
2. products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environment conditions and electrical power requirements;
3. products considered to be of a consumable or sterile nature;
4. accessories or parts not manufactured by Midmark;
5. charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark
6. costs and expenses of routine maintenance and cleaning;
7. representations and warranties made by any person or entity other than Midmark;
8. matching of color, grain or texture except to commercially acceptable standards;
9. changes in color caused by natural or artificial light;
10. custom manufactured products;
11. alterations or modifications to the product by any person or entity other than Midmark; and
12. Products that would otherwise be covered under Sections 1 and 2 of this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical, animal health or dental market, as the case may be, in which purchaser intends to use the product.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

WARRANTY DISCLAIMER

THIS LIMITED WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, or implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

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