# M3 Steam Sterilizer

[English](#)

[Español](#)

[Français](#)

For Models:

M3 (-001 thru -004)  
(-040 thru -042)

# User Guide

# Product Information

## Attention Canadian Users - Action Required!


Thank you for purchasing an M3 Steam Sterilizer.

To comply with changes in Canadian regulatory requirements you must now provide documentation, a Manufacturer's Data Report (MDR), to your appropriate provincial regulatory agency when purchasing and using a steam sterilizer.

To obtain the MDR from Midmark, record the sterilizer Serial Number and the ASME National Board Number of your sterilizer. See the illustration below for the location of these numbers. Then, contact Midmark customer service at 1.800.MIDMARK and provide the recorded numbers.

Midmark will then complete the MDR and send it to you. Once you have the MDR, you must submit it to the appropriate provincial regulatory agency.

We apologize for any inconvenience and appreciate your assistance in helping maintain compliance to Canadian regulations for sterilizers. If you have any questions or concerns, please contact 1.800.MIDMARK.

# Table of Contents

## **Important Information**

Component Location .....	4
Safety Instructions .....	5
Intended Use .....	5
Sterilization Technology .....	5
Electromagnetic Interference .....	5
Safety Symbols .....	6
Transportation / Storage Conditions .....	7
Accessory Chart .....	7
Sterilization Monitoring Guidelines .....	8

## **Installation**

Location Requirements .....	9
Operating Environment .....	10
Electrical Ratings / Requirements .....	10
Power Cord Connection .....	11
External Condensing Tank Connections .....	12
External Condensing Tank Draining Procedures .....	14

## **Operation**

Quick Reference .....	15
Power Switch .....	16
Filling the Reservoir .....	16
Qualification Testing .....	17
Loading the Tray .....	17
Cycle Selection .....	21
Post Sterilization Processing .....	23
'Additional Heat' Cycle .....	23
Adjusting the Dry Time .....	24

## **Maintenance**

Maintenance Messages .....	24
Daily Care .....	25
Periodic Maintenance .....	26
Monthly Maintenance .....	28
Extended Use Maintenance .....	29

## **Printer (optional)**

Printer Harness Connection .....	30
Printer Readout .....	31
Paper Roll Installation .....	32
Printer Ribbon Cartridge Replacement .....	33

## **Troubleshooting**

Error Codes .....	34
User Diagnostic Mode .....	35
Calling for Service .....	36

## **Specifications**


Specification Information .....	37
Water Purity Specifications .....	38

## **Warranty Information**

Scope of Warranty .....	39
-------------------------	----

# Important Information

## Component Location

SA1795-1i

## ***Safety Instructions***

The primary concern of Midmark is that this equipment is operated and maintained with the safety of the patient and staff in mind. To assure safe and reliable operation:

- Read and understand this manual before attempting to install or operate the sterilizer.
- Assure that the appropriate personnel are informed on the contents of this manual. (This is the responsibility of the purchaser).
- Assure that this manual is located near the sterilizer, or if possible, permanently affixed to the sterilizer.

## ***Intended Use***

The M3 Steam Sterilizer can be used in medical, dental and veterinary offices, hospitals, clinics, nursing homes, laboratories and other facilities to sterilize heat and moisture stable reusable items (including dental handpieces) that are compatible with steam sterilization.

Refer to 'Loading the Tray' and 'Cycle Parameters' in this manual for detailed information.

## ***Sterilization Technology***

The Sterilization Technology used in the M3 utilizes a dynamic air removal system called Steam Flush Pressure Pulse to remove air from the chamber.

## ***Electromagnetic Interference***

The M3 is designed and built to minimize electromagnetic interference with other devices.

However, if interference is noticed between another device and this sterilizer:

- Remove interfering device from room
- Plug sterilizer into a dedicated circuit
- Increase separation between sterilizer and interfering device
- Contact Midmark if interference persists

## Safety Symbols

### **WARNING**

*Indicates a potentially hazardous situation which could result in serious injury.*

### **Caution**

*Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices*

### **Equipment Alert**

*Indicates a potentially hazardous situation which could result in equipment damage.*

### **Note**

*Amplifies a procedure, practice, or condition.*


## Shipping Symbols

Caution Shipping Damage

Keep dry

Proper shipping orientation

Fragile

2 Person Lift

## Transportation / Storage Conditions

### **Equipment Alert**

*All water must be removed from the reservoir before transporting or storing at +32°F (0°C) or below.*

Ambient Temperature Range:..... -22°F to +140°F (-30°C to +60°C)

Relative Humidity:..... 10% to 90% (non-condensing)

Atmospheric Pressure: ..... 49.6 kPa to 106.4 kPa (7.2 psi to 15.4 psi)

## Accessories

Accessories	Part Number	Model
Printer	9A401001	M3-00x
Printer	9A401002	M3-04x
Door Tray	9A402001	M3-00x
Door Tray	9A686001	M3-04x

Common Service Parts/Tools	Part Number
M3 Internal Filter Screen	016-1095-00
2 1/2" Paper Roll	060-0008-00
Printer Cartridge	053-0505-00
Printer Refill Kit	002-0371-00
VistaCool™ Direct-To-Drain Thermal Reduction System	(Single) 9A586001
VistaCool™ Direct-To-Drain Thermal Reduction System	(Twin) 9A586002

VistaCool Direct-to-Drain System for Autoclave Wastewater is a trademark of Crosstex International, Inc., a Cantel Medical Company, Hauppauge, NY.

# Sterilization Monitoring Guidelines

## Note

*The information below is provided for reference only. Contact appropriate state / local agencies for specific sterilization guidelines for your office. Additional information on infection control is available from the Centers for Disease Control and Prevention (CDC), Organization for Safety and Asepsis Procedures (OSAP) and the American Dental Association (ADA).*

## Physical Monitors

Temperature and pressure measuring devices can help detect sterilizer malfunctions. The sterilizer's control system aborts the cycle and displays a message if physical conditions go outside established limits. The optional printer can be used to create a record of each load's actual cycle time, temperature and pressure.

## Note

*Use only FDA cleared chemical and biological indicators designed for steam sterilization that are compatible with the particular sterilization cycle temperature and exposure time being monitored. Use sterility monitors with each sterilization load. If a sterilization cycle is terminated prematurely, reprocess instruments to ensure sterility of the load. Process the load according to your regular practice, placing indicators near the handle side of tray. Follow manufacturer's instructions for proper disposal of used indicators.*

## Chemical Indicators

Chemical indicators are designed to verify that conditions in the sterilizer chamber were adequate to achieve sterilization. They do not validate that a processed item is sterile. If a chemical indicator shows a failure, items in that load are considered non-sterile. Potential causes for sterilization failures include: improper cleaning, packing, loading, or a sterilizer malfunction. Determine the cause of any sterilization failure and remedy the situation before running the next cycle. Only FDA cleared chemical indicators labeled for use with the nontraditional steam sterilization cycle parameters, e.g. temperature and exposure time, of the M3 Sterilizer should be used for monitoring the three M3 cycles. Follow the chemical indicator's instructions for proper storage, use, interpretation and disposal.

## Biological Indicators


Biological indicators are microbiological devices designed to accompany items being sterilized to monitor adequacy of the sterilization process. If a biological indicator shows a failure, items in that load are considered non-sterile. Potential causes for sterilization failures include: improper cleaning, packing, loading, or a sterilizer malfunction. Determine the cause of any sterilization failure and remedy the situation before running the next cycle. Only FDA cleared biological indicators labeled for use with the nontraditional steam sterilization cycle parameters, e.g. temperature and exposure time, of the M3 Sterilizer should be used for monitoring the three M3 cycles. Follow the biological indicator's instructions for proper storage, use, interpretation and disposal.

# Installation

## Location Requirements

(Allow clearance on both sides)

SA1796-1i

### Support Surface

- Material should be water-resistant material  
(Ex. laminate, stainless steel, stone, etc.)
- Surface must be level to ensure proper operation.
- Surface should meet minimum dimensions listed below:

#### Dimensions

Depth (front to back).....24 in. (61 cm)  
Width (side to side).....22 in. (56 cm)

### Clearance Requirements

To ensure proper air circulation and to allow access to the reservoir fill port and drain coupling, adhere to the minimum clearance requirements listed below.

### Clearance Requirements

Back of Unit - Back Wall ..... 4 in. (10 cm)  
Front Sterilizer Feet - Front of Support Surface ..... 4 in. (10 cm)  
Side of Unit - Side Wall..... 2 in. (5 cm) each side  
Distance Above Unit\* ..... 2 in (5 cm)\*

\* The minimum clearance for proper air circulation is listed.  
However, be sure to allow access to the reservoir fill port located on top of the sterilizer.

## Operating Environment

### Equipment Alert

*Unit should be allowed to reach room temperature before operating.  
Failure to do so could result in damage.*

Ambient Temperature Range:..... +68°F to 104°F (+20°C to 40°C)  
Relative Humidity:..... less than 80% (non-condensing)  
*(Pollution Degree 2, in accordance to IEC664)*  
Normal Operating Altitude: ..... less than 6000 ft. (1829 m) above sea level

- Approved for indoor use only
- Environment should be relatively dust-free

## Electrical Ratings / Requirements

### Note

*To ensure unit is properly grounded, it must be connected to a matching grounded, dedicated, correctly polarized receptacle.*

### WARNING

***Use 104-127 VAC, 50/60 HZ alternating current only for 115 VAC rated models and 207-253, 50/60 HZ alternating current only for 230 VAC rated models. Failure to do so could result in electrical shock to personnel and will result in damage to sterilizer.***

M3 (115V model): .....115 VAC, 50/60 Hz, 12 amp  
Max. Power Consumption: .....1400 Watts  
Requires\*:.....Dedicated supply circuit rated at 120 VAC, 50/60 Hz, 12 amp

M3 (230V model): .....230 VAC, 50/60 Hz, 6 amp  
Max. Power Consumption: .....1400 Watts  
Requires\*:.....Dedicated supply circuit rated at 230 VAC, 50/60 Hz, 6 amp

*\*Power source must have over voltage limits less than 1500 watts from mains to ground.  
(Installation Category II in accordance to IEC 664)*

## Power Cord Connection

### **WARNING**

*Equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen, air, or nitrous oxide.*

*Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.*

### **WARNING**

*Check the serial number label on back panel of sterilizer to verify voltage rating for the unit. Failure to connect sterilizer to an appropriate power supply could result in damage to the unit and electrical shock to personnel.*


### **Equipment Alert**

*For optimal performance, allow sterilizer to reach room temperature before operating.*

### **To connect the power cord...**

- A) *Plug power cord into receptacle on back of sterilizer.*
- B) *Plug power cord into a properly polarized and grounded receptacle rated for a minimum of 15 amps. A dedicated circuit only used for the sterilizer is recommended.*

SA1816i

# External Condensing Tank Connections


## Note

External Condensing Tank is not needed if M3 is connected to the VistaCool™ Direct-To-Drain Thermal Reduction System.

## Equipment Alert

The tank hose must be properly connected (and not kinked). Improper connection will cause water / steam leaks or a sterilizer malfunction if water flow to the tank is restricted.


## To connect the condensing tank...

- If connected, disconnect sterilizer power cord.
- Push tank hose into compression fitting on lid. Tighten nut. (If step B is difficult - see NOTE)
- Connect tank hose to back of sterilizer.
- Fill condensing tank with tap water to minimum water level indicator line.
- Secure the lid on the condensing tank by aligning the tabs and rotate lid clockwise to secure.


## Note

### If tubing is difficult to install...

- Remove compression nut and two ferrules.
- Install nut and ferrules onto tank hose as shown. (Large and small ferrules must be positioned as shown)
- Insert hose into tank fitting, then tighten nut.


Compression Fitting

Min. Water Level Indicator Line

SA1826i

## External Condensing Tank Connections - continued

### **To connect the condensing tank - continued...**

- F) Place tank on a level surface, preferably below the sterilizer but in no case should it be higher than the sterilizer support surface.
- G) Connect sensor plug.
- H) Connect sterilizer power cord.

### **Note**

#### **Clearance...**

Maintain a minimum of 6 inch clearance above the condensing tank for proper steam ventilation. If enclosed in a cabinet the support surface and surrounding surfaces, should be protected with a water resistant material (e.g. plastic, laminate, stainless steel, etc.). If enclosed in a cabinet, it is recommended that the door be vented to avoid heat, moisture build up and potential damage to the inside of cabinet.

SA1827i

# External Condensing Tank Draining Procedure


## Caution

Water that is discharged to external condensing tank can be **VERY HOT**; person emptying pitcher should allow the temperature to cool. Always use carrying handle and use caution when emptying.


### To remove lid...

- A) Rotate lid counterclockwise.
- B) Lift lid assembly out of condensing tank pitcher.

### To drain pitcher...

- A) Grasp pitcher by handle.
- B) Pour cool water into an approved drain.

### To replace lid...

Note: Lid can be reinstalled in one of four locations to facilitate safe water disposal.

- A) Insert lid assembly into pitcher.
- B) Align the four tabs for correct positioning.
- C) Rotate lid clockwise to secure.


# Operation

## Quick Reference

(Detailed instructions for each step are outlined in the following pages of the Operation section).

### Basic Operation...

- A) Turn power switch ON (I).
- B) Fill reservoir.
- C) Load tray.
- D) Press desired cycle button.
- E) Press <Start> button.

### Note


There is a 10 minute (dry) heat cycle that can be used to pre-heat the chamber, or for additional drying time. Refer to: '**Additional Heat**' Cycle in this section.

Cycle  
Buttons

## Power Switch

The power cord must be connected and the power switch must be ON (I) for the sterilizer to operate.



SA1800i

## Filling the Reservoir

### To fill reservoir...

Pour one (1) gallon of distilled water into fill port. Do not fill above lower lip of fill port.

### Equipment Alert

Use distilled water or water that meets the referenced water purity specifications. Failure to comply may result in sterilizer malfunction and/or premature failure due to excessive corrosion.

SA1801i


## Qualification Testing

Your sterilizer should be tested after sterilizer installation, malfunctions, relocation, major repairs and after sterilization process failure. Qualification testing should be performed prior to placing the sterilizer in service. If multiple cycle types are used, e.g. "Pouches" and "Low Temp" each cycle type should be qualified. Qualification testing should include at least one Biological Indicator (BI) (sometimes referred to as Spore Tests) and one Chemical Indicator (CI). The test pack should be placed near the front of the tray and performed with items routinely processed and considered to be the most difficult to sterilize. Additional items should be placed in the chamber along with the Biological Indicator and Chemical Indicator so that chamber is fully loaded (don't exceed the maximum capacities listed in the tables under "Guidelines for Loading" in this manual). Three consecutive test runs, for each cycle type tested, with negative results from the BIs and the appropriate readings from all physical monitors and chemical indicators demonstrating complete sterilization, provide verification that the sterilizer has been properly installed (or reinstalled after relocation) or repaired to the manufacturer's specifications and that it will function effectively in the facility in which it is installed. All items processed during qualification testing should be quarantined until the results of the biological testing for all three test runs are available.

## Loading the Tray

### Types of Items

Before placing any instrument in the M3, check with the instrument manufacturer to be sure the materials are compatible with steam sterilization and to verify the acceptability of sterilization parameters.

SA1802i

The M3 is designed to sterilize the following:

- High and low speed handpieces
- Metal instruments
- Rubber / plastic devices (ex. suction cannulas, impression trays, etc.)
- Wrapping / bundling materials (ex. CSR wrap, instrument pouches, etc.)
- Cassettes (Hu-Friedy Signa-Stat [6.5" x 10.5" x 1.25"] or smaller)
- Surgical instruments (ex. ophthalmologic instruments)

### Equipment Alert

Do **not** sterilize items composed of any of the following materials in the M3:

- Corrosion sensitive metal (ex. carbon steel, iron, etc.)
- Fragile items susceptible to breaking under pressure / high temperature
- Liquids
- Biomedical waste
- Textiles (including towels, gauze, etc.)
- Plastics that may break down or produce residue when exposed to steam / high temperatures.

#### Examples

Polyethylene

PVC

PPO (Noryl™)

Styrene

Textiles

Latex

Cellulosics

Acrylic (Plexiglass™)

Neoprene

ABS

## Loading the Tray - continued

### Immediate Use Sterilization


The M3 is capable of immediate use sterilization - sterilizing unwrapped instruments for immediate use. Please consider the following when choosing whether or not to sterilize your instruments unwrapped:

- The sterility of unwrapped instruments is compromised upon exposure to a non-sterile environment. Follow CDC guidelines for using unwrapped, sterilized instruments.
- Due to the sensitive nature of some types of surgery (including, but not limited to ophthalmological), instruments used in such procedures must be wrapped or pouched in order to reduce their exposure to sterilization process residues. The water reservoir should also be drained and refilled with fresh distilled water on a daily basis when processing instruments for these procedures on a routine basis.

### Pouching and Wrapping Items

The M3 is capable of sterilizing pouched or wrapped items.

- When pouching or wrapping items, use only sterilizer pouches and wraps that have been cleared by the FDA and labeled for use with the nontraditional steam sterilization cycle parameters, e.g. temperature and exposure time, of the M3 sterilizer. Follow the manufacturer's instructions for use.
- When using cassettes in the M3 follow the manufacturer's instructions for use.
- Pouched items to be sterilized should be placed lengthwise with plastic side of pouch facing up in the M3 tray. The pouches may overlap slightly, but items must not be layered. Refer to diagram below.

#### **WARNING**

*Failure to comply with these guideline may adversely affect sterilization and/or drying.*

## Loading the Tray - continued

### Load Size

The M3 Steam Sterilizer can accommodate loads weighing up to **2.4 lbs (1.1 kg)**.

[Note: This is the weight of the contents in the tray (ex. instruments, cassettes, pouches, etc.). The weight of the tray itself has already been accounted for].

### WARNING

**Do not overload the sterilizer tray. Failure to allow adequate space around items for circulation will compromise sterilization and drying.**

Use the table below as a general guideline for weights of commonly used items. Consult manufacturer's specifications for the exact weight of any particular instrument.

Item Description	Weight*	
	lbs.	kg
Scissors	0.066	0.030
Dental Scalars	0.044	0.020
Forceps	0.033	0.015
Dental Handpiece	0.121	0.055
Suction Cannula	0.022	0.010
Plastic Mouth Mirror	0.018	0.008
Impression Tray	0.033	0.015
Plastic X-Ray Positioning Ring	0.044	0.020
Hu-Friedy Signa-Stat Cassette	1.500	0.680

### Packing the Tray

(\*actual weights may vary)

### Equipment Alert

**Do not use towels or packaging which contains chlorine bleach residue. Chamber and/or tray rusting or discoloration may occur. The life of the sterilizer may be shortened significantly.**

### WARNING

**Clean and dry instruments thoroughly before placing them into tray. Improper cleaning may result in non-sterile instruments or damage to the unit. Follow instrument manufacturer's guidelines and CDC recommendations for handling and cleaning instruments prior to sterilization.**

In addition to total load weight outlined above, all items must be processed in accordance with Centers for Disease Control and Prevention (CDC), 'Guidelines for Infection Control in Dental Healthcare Settings' - 2003, MMWR 2003; 52 (no. RR-17), which states:

"Items to be sterilized should be arranged to permit free circulation of the sterilizing agent (e.g., steam, chemical vapor, or dry heat); manufacturer's instructions for loading the sterilizer should be followed."


- All items must fit in M3 tray.
- Loaded tray must slide into chamber opening without scraping.
- Items must not touch one another.
- Pouched items should be loosely packed.
- Pouches may overlap slightly, but items must not be layered.

### WARNING


**Failure to comply with these guideline may adversely affect sterilization and/or drying.**

## Loading the Tray - continued

# Cycle Selection

(The parameters for each cycle are outlined on the following page.)

## Cycle Selection - continued

### Cycle Parameters

(Before sterilizing any items in the M3, refer to **Loading the Tray** in this section.)

CYCLE	STERILIZATION PARAMETERS	DRY TIME	ITEMS TO BE STERILIZED
 <b>Unwrapped</b>	Temperature: 270°F (132°C)  Pressure: 27.1 psi (186 kPa)  Time: 3:30 Minutes	Time: 25 Minutes	<ul style="list-style-type: none"> <li>Dental instruments / handpieces loose on a tray.</li> <li>Items manufacturers recommend for exposure at 270°F (132°C), loose on tray.</li> </ul>
 <b>Pouches</b>	Temperature: 270°F (132°C)  Pressure: 27.1 psi (186 kPa)  Time: 5:30 Minutes	Time: 30 Minutes	<ul style="list-style-type: none"> <li>Dental instruments / handpieces in pouches, wrapped, or in a wrapped cassette.</li> <li>Items manufacturers recommend for exposure at 270°F (132°C), in pouches, wrapped, or in a wrapped cassette.</li> </ul>
 <b>Low Temp</b>	Temperature: 250°F (121°C)  Pressure: 15.0 psi (104 kPa)  Time: 20:00 Minutes	Time: 50 Minutes	<ul style="list-style-type: none"> <li>Rubber or plastic dental devices, dental instruments / handpieces loose on a tray, in pouches, wrapped, or in a wrapped or unwrapped cassette.</li> <li>Items manufacturers recommend for exposure at 250°F (121°C), loose on tray, in pouches, wrapped, or in a wrapped or unwrapped cassette.</li> </ul>

## Post-Sterilization Processing


After sterilization is complete, all items must be handled in accordance with accepted and documented standards, such as the Centers for Disease Control and Prevention (CDC) document, 'Guidelines for Infection Control in Dental Healthcare Settings' - 2003, MMWR 2003; 52 (no. RR-17), as well as any local requirements that may apply.

Qualified personnel responsible for infection control should prepare a protocol for handling sterilized items. This protocol should be followed by all personnel responsible for handling sterilized items.

### 'Additional Heat' Cycle

The Additional Heat Cycle activates the dry heaters for ten minutes.

This cycle can be used to pre-heat the chamber at the beginning of the workday, or for extended drying time at the end of a cycle.

SA1804i

#### **To pre-heat chamber prior to running a cycle...**


- A) Press <Start> button when 'SELECT CYCLE' appears on display.
- B) During the ten minute pre-heat mode, 'ADDITIONAL HEAT' will flash on the display.
- C) When 'ADDITIONAL HEAT' stops flashing, press desired cycle button, then press <Start>.

#### **For extended drying time at the end of a cycle...**

- A) Press <Start> button when 'SELECT CYCLE' appears on display.
- B) During the ten minute drying mode, 'ADDITIONAL HEAT' will flash on the display.

## Adjusting the Drying Time

The M3 allows the operator to adjust the drying time from 20 - 60 minutes using 1 minute increments for the three pre-programmed cycles.

SA1805I

### ***To adjust the drying time for a pre-programmed cycle...***

*A) After pressing desired cycle button, press <P> button.*

*[Display will show current setting. (ex. DRY TIME: 30 MINUTES)]*

*B) Press the <+> / <-> buttons to increase / decrease the drying time.*

*C) Press the <P> button to save your changes.*

*(Pressing the <Stop> button cancels the changes and returns to last saved setting.)*

## Maintenance

### Maintenance Messages


To assure correct operation and maximum sterilizer life, the M3 provides the operator with reminders when it's time to perform operator maintenance. After the M3 is powered ON for 7, 14 and 21 days, a message "Perform Periodic Maintenance" will be displayed. After 28 days, a "Perform Monthly Maintenance" message will be displayed. Refer to the appropriate maintenance instructions in this manual. The maintenance reminders are removed from the display when a cycle is started. If power is turned OFF, the timer will reset, initiating a new cycle of messages.

# Daily Care


- **Clean External Surfaces / Tray and Chamber**

- A. Wash unit according to your facility's procedure for clinical contact surfaces noting the following:  
(Use only quaternary disinfectants to disinfect unit. Staining, pitting, discoloration, or softening could occur if phenolic, iodophor, or glutaraldehyde-based disinfectant is used on plastic surfaces of the unit. Also, use of alcohol or aerosol spray cleaner / disinfectant containing substantial amounts of alcohol in the formula can damage the faceplate.)
- B. Wring excess solution from cloth.
- C. Using soft cloth, wipe all external surfaces.
- D. Follow the instructions provided with the cleaner / disinfectant used regarding rinsing and drying of external surfaces.

### **WARNING**

Check to assure the chamber filter is not blocked by debris and that it is properly seated; in the bottom hole in the back of the chamber. Failure to make sure the filter is in place and clear of debris could result in serious personal injury and or equipment damage.

- **Clean Chamber Filter**


Check that the chamber filter is free of debris and properly positioned in the hole at the back of the sterilizer chamber. If the filter is blocked by debris follow the procedures under "Monthly Maintenance" to remove and clean the filter.

- **Clean Door Gasket / Mating Surface**

- A. Wash with a damp cloth.
- B. Inspect gasket for damage.
- C. Replace gasket if necessary.

SA1806i

## Periodic Maintenance

- **Drain / Refill Reservoir (As Needed)**

- Remove drain hose from storage location.
- Place open-end of drain hose into container or sink.
- Connect adapter-end of drain hose to coupling as shown.
- Once water has drained, press release lever and remove hose.
- Return drain hose to storage location.
- Refill reservoir with distilled water.

### **Equipment Alert**

Use distilled water or water that meets the referenced water purity specifications. Failure to comply may result in sterilizer malfunction and/or premature failure due to excessive corrosion.


## Periodic Maintenance - continued

### • Empty / Clean External Condensing Tank (As Needed)

#### Note

This process is not necessary if the sterilizer is connected to a direct-to-drain thermal reduction system.

- A. Empty water from tank. (Do not reuse water!)
- B. Clean tank with diluted bleach solution (1/4 cup bleach : 1 gallon water) and a brush.
- C. Rinse tank thoroughly.
- D. Refill tank to minimum water level indicator line.

# Monthly Maintenance

- **Remove and Clean Filter**


Wash with mild soap solution to remove debris. Rinse with distilled water.  
(Use a stiff brush to scrub, or place in ultrasonic cleaner if necessary.)

- **Clean Condensing Tank Level Sensors**

**Note**


This process is not necessary if the sterilizer is connected to a direct-to-drain thermal reduction system.

Clean two sensors with mild soap solution, then wipe dry.

SA1809I


## Monthly Maintenance - continued

### To remove filter...

- A) Slide notch in filter tool (provided) over filter.
- B) Pull out to remove.

### Tapered End

### To install filter...

- A) Position filter in notch of tool as shown.
- B) Align filter with hole in back of chamber\* and press in gently.
- C) Flip tool over, then gently press filter in to secure.

\*Tip: Slide tool along bottom of chamber to align filter with hole.

SA1810i

## Extended Use Maintenance

The M3 is designed and tested to provide exceptional reliability throughout its service life. However, like all electro-mechanical devices it is subject to wear and degradation with use.


To ensure the integrity, performance and safety of all major components it is the responsibility of the user to have the sterilizer performance / operation verified by a Midmark Authorized Service Provider after 8 years or 20,000 cycles of use, whichever comes first. After 8 years or 20,000 cycles of use an annual inspection by a Midmark Authorized Service Provider is recommended.

# Printer (optional)

## Printer Harness Connection

### To connect printer harness...

- A) Turn sterilizer power switch OFF (O).
- B) Connect printer harness as shown.
- C) Turn sterilizer power switch ON (I).

SA1813i

# Printer Readout

```

BEGIN
UNWRAPPED
270 deg F 3.5 MINUTES

MM - DD - YYYY
__ - __ - __
HH : MM
__ : __
TOTAL CYCLES 140

HEATING

MM:SS degF PSI
0:00 151.5 0.1
0:30 193.1 0.2
1:00 218.1 3.8
1:30 243.2 11.9
2:00 244.6 13.5
2:30 268.5 27.8

STERILIZING

MM:SS degF PSI
0:00 271.0 28.8
0:30 274.8 31.6
1:00 272.9 29.2
1:30 275.1 25.8
2:00 273.7 30.8
2:30 273.8 30.0
3:00 273.5 29.3
3:30 275.2 25.8

MIN TEMP 271.0 degF
MAX TEMP 275.7 degF
MIN PRESS 23.2 PSI
MAX PRESS 32.1 PSI

VENTING CHAMBER

DRYING

MM:SS
0:00

TOTAL PROCESS TIME
9:50

COMPLETE
    
```

Selected Cycle

Set Temperature and Time

Date

Total # of cycles run on unit

During Heat-Up phase, the printer records chamber temperature and pressure in 30-second increments.


During Sterilization phase, the printer records chamber temperature and pressure in 30-second increments.

Summary of Sterilization phase

Indicates Venting phase initiated

Indicates Drying phase initiated (Printer records drying time in 5 minute increments).

Paper Feed Button

## If printer readout becomes light ...


- A) Separate lid / printer from housing.
- B) Press Paper Feed button for three seconds. (This advances the ribbon inside the cartridge).
- C) Repeat step B as necessary\*.

\*Note: If pressing the Paper Feed button fails to correct the problem, replace the printer cartridge.

## Paper Roll Installation

### To install paper roll...

- A) Loosen thumbscrew, then separate lid / printer from housing.
- B) Slide paper into slot as shown\*.
- C) Press paper feed button until approx. 2 in. (5 cm) of paper is fed thru printer.
- D) Pull paper thru slot in lid.
- E) Position paper roll / spindle in slots.
- F) Reassemble lid / printer / housing.




# Printer Ribbon Cartridge Replacement

## To remove old printer ribbon cartridge...

- A) Turn sterilizer power switch OFF (O).
- B) Separate lid from printer.
- C) Press down on side of printer cartridge labeled "EJECT".

## To install new printer ribbon cartridge...

- A) Install new cartridge as shown\* (it will "snap" into place).

\*Note: The printer paper must be between the cartridge ribbon and the metal plate.

- B) Turn Adjustment Knob clockwise until ribbon is tight.

# Troubleshooting

## Error Codes

If a malfunction is detected during a cycle, a numeric error code will appear on the display panel. Use the chart below to diagnose and correct the most common, maintenance-related error codes. If you encounter an error code not identified below, follow the instructions on the display panel. If error code persists, contact your authorized service provider.

**Example:**

SA18041

Error Code	Probable Cause	Corrective Action
<b>C010</b>	Sterilizer lost power during cycle.	Press <STOP> button, then restart cycle.
<b>C100-series (all C101, C102, etc.)</b>	<STOP> button was pressed during the cycle.	Press <STOP> button, then restart cycle.
<b>C231, C232</b>	Not enough water in reservoir to complete the cycle.	Fill reservoir with distilled water, or water that meets the referenced water purity specifications.
<b>C441, C442</b>	External condensing tank is full.	Empty external condensing tank.
<b>C533, C633</b>	Water pump needs primed.	Put the sterilizer in the User Diagnostic Mode and prime the pump. (see following page)

# Troubleshooting - continued

## User Diagnostic Mode

The User Diagnostic Mode is used:

- To set “English” or “Metric” units on the display panel,
- To retrieve the last five (5) error codes stored in the unit memory.
- To prime the water pump if needed.

### **To activate User Diagnostic Mode...**

- A) Turn Power Switch OFF (O).
- B) Press and hold the <START> button.
- C) Turn Power Switch ON (I).
- D) Press the <START> button when the display shows “USER DIAGNOSTIC”.

### **To change Display Units...**

- A) Put the unit in User Diagnostic Mode.
- B) Press the <P> button to select units.
- C) Press the <+> button to change the temperature and pressure display from English-to-Metric or Metric-to-English units. (Factory default setting is English).
- D) Press the <START> button to continue.
- E) Turn the power switch OFF (O) to exit User Diagnostic Mode.

### **To retrieve the five (5) most recent error codes...**

- A) Put the unit in User Diagnostic Mode.
- B) Press the <STOP> button to recall errors.
- C) The last five (5) error codes will be displayed.
- D) Press the <START> button to return to the User Diagnostic Mode display.
- E) Turn the power switch OFF (O) to exit User Diagnostic Mode.

### **To prime the sterilizer pump...**

- A) Put the unit in User Diagnostic Mode.
- B) Press the <START> button to start the pump priming progress.  
The unit will automatically cycle through a pre programmed priming cycle...
  - Closing the sterilizer door.
  - Heating the boiler.
  - Cycling the pump ON and OFF until the pump is primed.When finished the 2nd line of the display will show “PRIMING COMPLETE”.
- C) Press the <START> button to return to the User Diagnostic Menu.
- D) Turn the power switch OFF (O) to exit User Diagnostic Mode.

# Calling for Service

**Note**

Please mark down any displayed Code(s) and be sure to relay this information to the service technician.

Contact your **Midmark Authorized Dealer**, or log on to [www.midmark.com/technical-library](http://www.midmark.com/technical-library).  
**Model and serial number information will be required** when calling for service.

To contact Midmark directly:

1.800.MIDMARK (1.800.643.6275) or 937.526.3662  
8 am to 5 pm EST (Monday thru Friday)  
*[excluding standard US holidays]*

# Specifications

## Fuse Ratings:

### 115 VAC Unit

F1 ..... 15 Amp, 250 V, Fast Acting, 1/4" x 1 1/4"  
F2 ..... 0.25 Amp, 250 V, Slo-blo, 1/4" x 1 1/4"

### 230 VAC Unit

F1 ..... 8 Amp, 250 V, Fast Acting, 5 x 20 mm  
F2 ..... 0.125 Amp, 250 V, Slo-blo, 5 x 20 mm

## Certifications:

ASME Boiler and Pressure Vessel Code, Section VIII, Division 1  
Canadian Registration Number Available  
UL 61010-1, 2nd Edition  
IEC 61010-2-040  
CAN/CSA-C22.2 No. 61010-1, 2nd Edition  
FCC Part 15, Sub-part B

## Physical Dimensions:

Overall Length: ..... 21 in. (53.3 cm)  
Overall Width: ..... 17.8 in. (45.2 cm)  
Overall Height: ..... 7.1 in. (18 cm)  
Shipping Carton Length: ..... 25 in. (63.5 cm)  
Shipping Carton Width: ..... 22 in. (55.9 cm)  
Shipping Carton Height: ..... 16.6 in. (42.2 cm)  
Counter Area: ..... 24 in. (61 cm) deep x 22 in. (55.9 cm) wide  
Chamber Volume: ..... 0.49 gal (1.8 liter)

## Weight:

Empty Reservoir: ..... 71 lbs. (32.2 kg)  
Full Reservoir: ..... 80 lbs. (36.3 kg)  
With Shipping Carton: ..... 80 lbs. (36.3 kg)

**Water Reservoir Capacity:** ..... 1.20 gal (4.5 liter)

**Pressure Relief Valve Setting:** ..... 40 PSI (275.8 kPa)

## Chamber Pressure:

@ 270°F (132°C) ..... 27.1 psi. (186.2 kPa)

## Water Purity Specifications

<b>Water Purity Specification Chart</b>	
	<b>AAMI ST-46 (ref.)</b>
Evaporate Residue	≤ 15 mg/l
Silica	≤ 2 mg/l
Iron	≤ 0.2 mg/l
Cadmium	≤ 0.005 mg/l
Lead	≤ 0.05 mg/l
Rest of heavy metals, excluding Iron, Cadmium and Lead <sup>1</sup>	≤ 0.1 mg/l
Chloride	≤ 3 mg/l
Phosphate	≤ 0.5 mg/l
Conductivity (at 20°C) <sup>2</sup>	≤ 50 μs/cm
pH	6.5 to 8
Appearance	Colorless, clean, w/o sediment
Hardness	≤ 0.1 mmol/l (.585 gr/gal)

# Warranty Information

**SCOPE OF WARRANTY** Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

Midmark warrants to the original retail purchaser that during the applicable warranty period it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications.

**APPLICABLE WARRANTY PERIOD** The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product and shall be one (1) year for all warranted products and components.

**OBTAINING WARRANTY SERVICE** Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at [www.midmark.com](http://www.midmark.com); by phone at 1.800.MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380. It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

**EXCLUSIONS:** This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Products considered to be of a consumable or sterile nature;
- (4) Accessories or parts not manufactured by Midmark;
- (5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) Costs and expenses of routine maintenance and cleaning; and
- (7) Representations and warranties made by any person or entity other than Midmark.
- (8) Matching of color, grain or texture except to commercially acceptable standards;
- (9) Changes in color caused by natural or artificial light;
- (10) Custom manufactured products;
- (11) Alterations or modifications to the product by any person or entity other than Midmark; and
- (12) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market.

**SOFTWARE; WITH RESPECT TO SOFTWARE THAT IS A PRODUCT OR COMPONENT THEREOF, MIDMARK DOES NOT WARRANT THAT THE SOFTWARE: (1) IS ERROR FREE; (2) CAN BE USED WITHOUT PROBLEMS OR INTERRUPTIONS; OR (3) IS FREE FROM VULNERABILITY TO INTRUSION OR ATTACK BY VIRUSES OR OTHER METHODS.**

**EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER** MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN. WARRANTY DISCLAIMER THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

**STATUTE OF LIMITATIONS** No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.


**NO AUTHORIZATION** No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

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