Midmark® 1000 Mobile Dental Delivery System

[English](#)

[Español](#)

[Français](#)

For Models:

8000-001
8000-002
8000-003
8000-004
8000-005
8000-006
8000-006S
8000-007
8000-007S
8000-008
8000-008S

User's Guide

Product Information


(The information below is required when calling for service.)

Dealer:

Date of Purchase:

Model / Serial Number:

Midmark Authorized Service Company:

Model / Serial Number Location

Table of Contents

Important Information

Safety Symbols.....	4
Intended Use	5
Electromagnetic Interference.....	5
Disposal of Equipment.....	5
Transportation / Storage / Operating Conditions	5

Transport

Transporting Dental System	6
----------------------------------	---

Operation

Electrical System	7
Water System	8
Air System	9
Delivery System.....	10
Low-Speed Handpiece Operation	11
High-Speed Handpiece Operation.....	12
Syringe Tip Replacement and Handpiece Installation	13
Drive Air Adjustment and Locking Caster	14
Air Pressure Regulator Adjustment and Water Flow Adjustment	15
Coolant Air Adjustment	16

Cleaning / Maintenance

Midmark 1000 Cleaning/Maintenance	17
Midmark High-Speed Handpiece Cleaning/Maintenance.....	23
Midmark Low-Speed Handpiece Cleaning/Maintenance.....	29
Calling For Service	30

Specifications

Drawings.....	31
Weights / Dimensions.....	32
Electrical	33
Model Identification / Compliance Chart	34

Warranty

Warranty Information	35
----------------------------	----

For installation information go to:

Technical Library at Midmark.com

Important Information

Safety Symbols

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage if not avoided.

Note

Amplifies a procedure, practice, or condition.

Proper Shipping Orientation

Handle With Care

Type B, Applied Part

Maximum stacking height
(Do not stack)


Humidity Limit

Type BF, Applied Part

Fragile

Fuse Rating
Specifications


Protective Earth Ground

Keep Dry

Pressure Limit

LED Scaler Handpiece
Duty Cycle


Refer to Instruction
Manual/Booklet

TUV - Domestic
Models Only

Do Not Tumble

Temperature Limit

Intended Use

The Midmark 1000 is a mobile dental delivery system intended to give qualified veterinarians and/or technicians a self-contained, compact, portable unit for performing veterinary dental procedures.

Electromagnetic Interference

This product is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this cart:

- Remove interfering device from room
- Plug cart into an isolated circuit
- Increase separation between Midmark product and interfering device
- Contact Midmark if interference persists

Disposal of Equipment

At the end of product life, the unit(s), accessories and other consumable goods may be contaminated from normal use. Consult local codes and ordinances for proper disposal of equipment and other consumable goods.

Authorized Representatives

EC **REP**

Warning

In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.

Customers in the EU should direct all questions, incidents and complaints to Midmark's Authorized Representative listed below.

CEpartner4U
Esdoornlaan 13
3951 DB Maarn, The Netherlands
Phone: +31 343 442 524
Fax: +31 343 442 162

Transportation / Storage / Operating Conditions

Transportation / Storage Temperature Range:	23°F to 100°F (-5°C to 38°C)
Operating Temperature Range:	50°F to 104°F (+10°C to +40°C)
Transportation / Storage Relative Humidity:	10% to 90% (non-condensing)
Operating Relative Humidity:	30% to 75% (non-condensing)
Transportation / Storage Atmospheric Pressure:	50 kPA (7.2 psi) to 106 kPA (15.3 psi)
Operating Atmospheric Pressure:	.700 hPa (20 in Hg) to 1060 hPa (31 in Hg)


Transport

WARNING

Exercise care when transporting over thresholds and other low lying objects.
Failure to do so can lead to injury.

When transporting over large thresholds or steps...
A) Pull cart up and over onto level surface.

VA1626

Operation


Electrical System

On/Off Switch

On/Off Switch is located on the back of the cart. It controls the main power source and must be ON to operate the delivery system.

Circuit Breakers

There are two circuit breakers located on the lower back panel of the cart. If the unit's maximum load is exceeded, the circuit breakers interrupt power to the unit. Maximum load is 115 VAC, 6 amps (domestic) 230 VAC, 3 amps (export).

Mains Disconnect

Power can be terminated by unplugging the mains disconnect or turning switch off.

Operation

Water System

Water Bottles

The unit has two 1 liter bottles. The water bottles supply water to the handpieces, air/water syringe and a quick disconnect on the back of the cart. Air flows through the water regulator and the toggle switch to the water bottles. The water system is regulated to 40 PSI. Use the water needle valve to adjust water flow for each instrument. The toggle switch between the water bottles that controls which bottle will supply the water. When the toggle is in the center neither water bottle is pressurized. Move the toggle toward the left to pressurize the left bottle and depressurize the right bottle. Toggle to the right to pressurize the right bottle and depressurize the left bottle.


Caution

Failure to properly seal water bottles can lead to patient injury.

Water Quick Disconnect

The water quick disconnect has a blue ring. You can attach accessories to the cart that will utilize the water from the cart. You may need to add a regulator if attaching equipment that can not handle 40 PSI.

To fill water bottles...

- A) Move toggle to center.
- B) Remove water bottle and fill. (Distilled Water Is Recommended)
- C) Replace water bottle and move toggle toward filled water bottle.


Note: Water bottles should be checked in between each patient and filled if necessary.

Operation

Air System

Compressor

The air system is a 1/3 HP compressor with a 2 liter air tank and regulator. It supplies air for the handpieces, air/water syringe and a quick disconnect on the back of the cart. The compressor will cycle between 70 to 100 PSI. The regulator is factory set at 65 PSI.

Air Quick Disconnect



The air quick disconnect has a red ring. You can attach accessories to the cart that will utilize the air from the cart. You may need to add a regulator if attaching equipment that can not handle 100 PSI.

Operation

Delivery System


Handpieces

The delivery system controls up to four handpieces. The handpiece pressure is controlled by handpiece drive air adjustment valves. The gauge on front of the delivery system indicates the handpiece pressure. Handpiece pressure should not exceed 60 PSI under normal operation. The scaler tool (optional) has an adjustment knob for the intensity level.

Air/Water Syringe

There are air and water buttons on top of the syringe. When both buttons are depressed you will get a mist of both water and air. Air and water adjustments are located under the syringe holder. Insert flat head screwdriver into adjustment screw(s). Counter-Clockwise to increase flow and Clockwise to decrease flow. Turn screw until desired air/water flow is coming through the instrument.

Operation

Low-Speed Handpiece Operation

Operating Recommendations

Air Pressure: Dental unit air pressure for low-speed handpiece should be between 35 and 50 psi.


Air Supply: Moist or dirty air from the compressor can damage Midmark low-speed handpieces. Check filters often. If dirt or moisture is present, replace filter and bleed system.

Empty Operation: Low-speed handpiece will not operate unless the chucking system is fully engaged. This minimizes risk of damage to the chuck.

Installation Notes

- This product is not sterile. Clean and sterilize before its first use. Follow sterilization and maintenance instructions.
- Check air line filter in dental machine. If dirt or moisture is present, bleed system and replace filter.
- Preset low-speed handpiece air pressure between 35 and 50 psi.
- Align drive air and exhaust tubes with outlets in the supply tube and screw the Low-speed handpiece to the dental unit hose nut, making sure contact between handpiece and gasket is flush.
- Secure handpiece bur or doriot attachment, if applicable, to the low-speed handpiece, using directions found below.

Note: Midmark low-speed handpiece will operate only when the chuck ring is in the locked position.

ArtHalf

Placing & Removing Handpiece Burs or Doriot Attachments & Heads

- 1) Hold handpiece securely, rotate the chuck ring to the left (clockwise while looking down at the nosecone) until it clicks into place. The arrow on the handpiece should now align with the unlock symbol on the chuck ring.
- 2) Align prophy angle or Doriot slot with the pin on nose cone
- 3) Slide the prophy angle or doriot attachment fully onto the nose cone up to the shoulder.
- 4) Hold handpiece securely, rotate the chuck ring to the right (counter clockwise while looking down at the nosecone) until the arrow on the handpiece aligns with the lock symbol on the chuck ring.
- 5) Double check to be sure the attachment is secured in place by gently pulling on it.
- 6) To remove the prophy angle, Doriot attachment, or cutting bur, repeat step 1 to unlock the chuck and pull straight out.
- 7) Detach handpiece from the delivery tube at the end of each day and prior to other extended periods of non-use.




Operation

High-Speed Handpiece Operation


To Change Bur in High-Speed Handpiece...

- A) Push button on back of handpiece and pull the bur out.
- B) Push bur into handpiece. (The shank will only go in half-way)
- C) Then push button on back of handpiece and push the bur in until shank is all the way in.

Note: Friction grip burs (FG Style) are designed so that the entire straight portion of the shaft should be inside the chuck of the handpiece. You should only see the tapered portion of the bur if installed properly. **New bur is recommended for each procedure.**

Syringe Tip Replacement


To Change Syringe Tip...

- A) Press down on large collar.
- B) When a soft "click" is felt, pull tip straight out.
- C) Hold collar down and insert new tip. Be sure to press it all the way in, then release collar.
- D) Test the installation by giving a firm pull on tip, to assure it is inserted all the way and locked in place.

VA1807i

Handpiece Installation


Install handpieces.

Drive Air Adjustment and Locking Caster

To adjust drive air...

Insert flat head screwdriver into knob for Drive Air Adjustment (for desired handpiece). Turn knob until proper air pressure displays on the pressure gauge.

Operation


Air Pressure Regulator Adjustment

EQUIPMENT ALERT
Do not adjust regulator over
65 PSI.

To Adjust Regulator...

- Pull adjustment knob upward to unlock.
- Rotate adjustment knob clockwise to increase pressure, counterclockwise to decrease pressure.
- Push adjustment knob downward to lock.



Water Flow Adjustment

To adjust water flow...

Rotate Water Needle Valve (for desired handpiece) until desired flow is achieved.

Water Needle Valve


Operation

Coolant Air Adjustment

To adjust coolant air...

Insert flat head screwdriver into knob for Coolant Air Adjustment (for desired handpiece). Counter-Clockwise to increase flow and Clockwise to decrease flow. Turn screw until desired air pressure is coming through the instrument.

Note: Drive air adjustment screws are located below tool holders for each instrument.

Midmark 1000 Cleaning/Maintenance

- Use cleaners that are appropriate for the situation, such as warm water and mild detergents.
- Clean the exterior surfaces with a clean, lint free cloth, dampened with disinfectant. Do not drip any liquid into open vents, plugs or connectors. Dry the surfaces with a clean cloth or paper towel.
- After cleaning, visually inspect the product for deterioration of covers. Do not use the delivery system if excessive discoloration, cracking, or other signs of wear are noticeable (See Calling for Service instructions).

Note

Dental unit waterline treatment is important for the health and safety of the patient and all clinic personnel. Due to the long, narrow waterline tubing in dental units, standing water can lead to microbial colonization and biofilm growth. Biofilm is a community of bacterial cells and other microbes that adhere to surfaces and form a protective slime layer. If left untreated these microbes can break off, transferring bacteria to the patient, or even become aerosolized into the clinic environment.

Waterline treatment is necessary to keep the count of heterotrophic bacteria from rising higher than desired levels. The desired level for a specific location should be determined by any local or regional guidelines. Regular monitoring should also take place to ensure that heterotrophic bacterial is not exceeding the desired limit. Monitoring frequency should be established by your facility. As a suggestion, Midmark would recommend that you begin by monitoring on a monthly basis and make adjustments to the frequency based on test results.

Treatment can come in many forms. The most popular methods on the market currently are tablets and straw/cartridge based systems. Midmark recommends the use of DentaPure by Crosstex, part number DP365B. Refer to the manufacturer's instructions for use.

Extra flushing may be needed within Midmark equipment when tablets are used. Undissolved tablet particles can gather over time in places within the waterlines, obstructing the line and causing water flow to slow. By flushing the waterlines, water flow is maximized and should push any undissolved particles through.


Every veterinary practice setting is different and no one disinfectant is the best choice for every facility. There are several organizations educated in disinfection procedures. Refer to web sites of accredited sources that can assist veterinary personnel in choosing what is best for their practice. Listed below are a few of the sites.

Organization for Safety and Asepsis Procedures
<http://www.osap.org>

Dept. of Health and Human Resources Centers for Disease Control and Prevention (CDC)
<http://www.cdc.gov>

When using disinfectants...

- Carefully read the product label and directions for use.
- Do not exceed the dilution rate.

DA148100i

Midmark 1000 Cleaning/Maintenance *continued...*

At the beginning of each work day...


- Fill water bottles with fresh distilled water.
- Flush handpieces.

Note

Distilled water is not mandatory. However, water must meet requirements for drinking.

To Flush Delivery System...

- Remove all handpieces that use water.
- Open water needle valve.
- Hold end of tubing over a container and depress foot control.
- Allow water to run 20 to 30 seconds.
- Install newly disinfected handpieces and syringe tip.
- Refill water bottle.

VA1226

Between each patient...

- Remove disposable tips, instruments, etc.
- Flush delivery system.
- Disinfect handpieces and scaler per manufacturers instructions.
- Disinfect syringe tip.
- Fill water bottles.

Equipment Alert

HANDPIECE STERILIZATION

Between each procedure complete manufacturer's instructions for "Infection Control / Sterilization found in the documentation included with the handpieces. Failure to do this will cause premature failure of handpieces.

Equipment Alert



AUTOCLAVABLE SYRINGE TIP STERILIZATION

The autoclavable syringe tips supplied with the delivery system must be sterilized prior to use with each patient, including initial use. Be sure to thoroughly rinse and clean syringe tips prior to sterilization, any debris may reduce the effectiveness of the sterilization. Recommended sterilization process is steam autoclave. Recommended parameters are 125°C (250°F) and 106 kPa (15 PSI) for 40 minutes at temperature and pressure.

Midmark 1000 Cleaning/Maintenance *continued...*

At the end of each day...


- Clean and disinfect the delivery system (see *Cleaning and Disinfecting instructions*).
- Drain air regulator.
- Drain condensation from air tank.

Midmark 1000 Cleaning/Maintenance *continued...*

To Drain Air Regulator...


Place a cloth under the regulator and depress the button until there is no water coming out.

To Drain Condensation from Air Tank...

- A) Turn power switch off.
- B) Position end of red tube so it drains into a container.
- C) Turn drain valve slowly to open.
- D) When no water is coming out, slowly close valve.


Midmark 1000 Cleaning/Maintenance *continued...*

Midmark 1000 Cleaning/Maintenance *continued...*

To clean the air/oil separator...

- A) Turn master switch OFF.
- B) Remove container by squeezing sides in while pulling out.
- C) Remove fluid and saturated gauze.
- D) Clean container and mounting cap.
- E) Install clean gauze and reinstall the container.

Midmark High-Speed Handpiece Cleaning/Maintenance

WARNING

Handpiece should be sterilized between patients to prevent cross-contamination which could result in serious illness or death from infectious diseases.

Equipment Alert


All Midmark Handpieces must be sterilized and lubricated after every patient. Lubrication will elongate life, increase performance and eliminate debris from within handpiece. Lubrication can be performed between handpiece cleanings if the handpiece appears to be under performing. If equipped, the swivel coupler can be sterilized in a bag at the same temperatures as a handpiece.

- 1) Do not use Flash sterilization techniques or damage to components will occur, voiding warranty.
- 2) Never submerge components in any cleaning or disinfecting solution.
- 3) Do not use Ultrasonic Cleaners on handpieces or couplers.
- 4) Do not Autoclave or Chemiclude for extended / Unusual periods of time (such as overnight).
- 5) Do not leave handpiece components in sterilizer after cycle is completed.


To clean the water port(s)...

- A) Use water spray tool to periodically clean spray ports as needed.

Clearing The Single Port Spray

Clearing The Quad Port Spray

midmark®

ArtHalf

Midmark High-Speed Handpiece

Cleaning/Maintenance *continued...*


To clean the handpiece...

- A) Detach handpiece from hose coupler (if equipped with a swivel quick disconnect remove handpiece then remove swivel from coupler with wrench).
- B) Remove bur from handpiece.
- C) Keeping the head of handpiece upright, clean external surface of handpiece thoroughly to remove saliva, blood and other organic soil. Scrub handpiece with a soft brush and warm water to remove debris. Rinse and dry.


Equipment Alert

Line up components properly to prevent damage resulting in hand piece failure.

Bur Removal

ArtFull



Midmark High-Speed Handpiece

Cleaning/Maintenance *continued...*

To clean/condition the handpiece drive air ports...

(Applicable on Non-Swivel High-Speed Handpieces Only).

- A) Assemble spray tube onto spray can head by firmly pushing tube into conditioner can spray tip as far as possible.
- B) Insert tube into drive air port of handpiece.
- C) Direct handpiece away from your body and towards an absorbent towel. Depress handpiece conditioner can top to spray until conditioner flows from both ends of handpiece head (1 to 2 seconds).
- D) If overspray on towel appears discolored repeat spray for an additional two seconds. Remove handpiece and wipe excess conditioner from handpiece with a clean dry towel.

ArtFull


Midmark High-Speed Handpiece

Cleaning/Maintenance *continued...*

To clean/condition the handpiece swivel coupler cavity...

(Applicable on Swivel High-Speed Handpieces Only).

- A) Screw swivel nozzle onto spout of a well-shaken can of Midmark One Step Handpiece Cleaner/Conditioner.
- B) Remove swivel coupler from handpiece and Insert nozzle tip into handpiece where swivel coupler was removed
- C) Direct handpiece away from your body and towards an absorbent towel. Depress handpiece conditioner can top to spray until conditioner flows from both ends of handpiece head (1 to 2 seconds).
- D) If overspray on towel appears discolored repeat spray for an additional two seconds. Remove handpiece and wipe excess conditioner from handpiece with a clean dry towel.
- E) Remove swivel nozzle from cleaner/conditioner can and lightly spray o-rings of swivel coupler. Re-install swivel coupler to handpiece.

Midmark High-Speed Handpiece


Cleaning/Maintenance *continued...*

To clean/condition the handpiece push button and chuck...

(Applicable on ALL High-Speed Handpieces).

- A) Screw push button chuck conditioner nozzle onto spout of a well-shaken can of Midmark One Step Handpiece Cleaner/Conditioner using knurled base. Orient curved tip of nozzle above can spout.
- B) Insert nozzle tip into handpiece spindle bore/ (Where the Bur is inserted and removed.) completely until it stops.
- C) Hold very slight pressure on handpiece towards nozzle tip to keep it seated in position. Direct handpiece and nozzle away from your body and towards an absorbent towel. Depress handpiece conditioner can top to spray for a minimum of two full seconds.
- D) If overspray on the towel appears discolored repeat spray for an additional two seconds. Remove handpiece from nozzle tip and wipe excess conditioner from handpiece with a clean dry towel.

Re-attach handpiece and install a bur or bur blank. Run WITHOUT water for 45 seconds. Wipe off excess conditioner with dry towel.

ArtFull

Midmark High-Speed Handpiece Cleaning/Maintenance *continued...*


To sterilize a Midmark high-speed handpiece...

- A) Remove bur. (If equipped, the swivel coupler must be detached and autoclaved separately at the same temperatures as a handpiece.)
- B) Place handpiece in autoclave bag. The use of autoclave bags dramatically reduces fiber optic light output deterioration and handpiece cosmetic damage.
- C) Load autoclave bag containing handpiece into autoclave or chemiclave. Load autoclave bags for maximum penetration of steam or chemical vapor.
- D) Cycle the handpiece according to the autoclave / chemiclave manufacturer's instructions. Do not exceed 275°F (135° C).
- E) Immediately remove handpiece from autoclave or chemiclave. Allow to cool sufficiently prior to handling.
- F) As soon as handpiece is cool enough to handle, wipe fiber optic light transmitting surfaces clean in all locations (both ends of handpiece) with isopropyl alcohol and cotton swab. This step is particularly critical to maintaining light output (Fiber Optic High-Speed Handpieces Only).
- G) If cleaning fiber optics of a swivel high-speed handpiece, remember to clean surfaces in both locations on the handpiece and swivel coupler.
- H) If equipped, reattach swivel coupler to handpiece.

Equipment Alert

Refer to handpiece manufacturer on all other handpieces for suggested cleaning and maintenance procedures.

Midmark Low-Speed Handpiece Cleaning/Maintenance


WARNING

Handpiece should be sterilized between patients to prevent cross-contamination which could result in serious illness or death from infectious diseases.

Maintenance - Low-Speed Handpiece

The Midmark low speed handpiece is designed to operate without the need for repeated lubrication. However, over time a build-up of debris inside the motor may affect its performance. Should that happen, apply Midmark one-step handpiece cleaner/conditioner to the drive air inlet tube and run the handpiece for 15 seconds to flush the debris out of the motor. Like new performance should then be restored. Wipe the Midmark low speed handpiece with a clean, dry cloth or gauze to remove excess conditioner. For the nose cone swivel, the lubricant should be applied to the joint between the nose cone and the motor.

AerEcol

To Clean Handpieces

When autoclaved according to the following instructions, the Midmark 0-5k Straight low speed handpieces will withstand up to 1200 autoclave cycles prior to requiring service.

- 1) Prepare the Midmark low speed handpieces for autoclave by carefully scrubbing the external surface with a soft brush and warm tap water $40^{\circ}\text{C} \pm 5^{\circ}\text{C}$ ($104^{\circ}\text{F} \pm 10^{\circ}\text{F}$) to remove visible debris and particles. Be careful not to run water into backend tubes. Dry with a towel or dry compressed air.
- 2) Seal handpieces in separate autoclave bags. Make sure to remove dental attachments from handpiece prior to sealing.
- 3) Process in gravity steam autoclave as follows:
 - Place handpiece(s) in autoclave tray(s).
 - Autoclave at a temperature of 270° F (132° C) for at least 15 minutes.
 - NOTE: Temperature should not exceed 275° F (135° C).
- 4) After autoclave cycle is completed, allow to cool down 30 minutes prior to handling.


Calling For Service

Direct all service inquiries to your authorized Midmark dealer. When calling for service, you must provide the following information:

- Model / serial number
- Date of purchase
- Symptom(s) of malfunction

Service and Parts Manual is available at www.midmark.com.

Midmark 1000

Specifications

Cart	
Width	23 inches (58.42 cm)
Length	18 1/2 inches (46.99 cm)
Height	36 1/2 inches (92.71 cm)
Weight	130 lbs (58.97 kg)
Shipping Weight	150 lbs (68.04 kg)
Top of Cart Weight Restriction	25 lbs (11.33 kg)
Compressor	
Motor	Oil Free 1/3 HP
Working Pressure	70-100 PSI
Flow Rate	2.5 scfm @ 40 PSI
Sound Level (full load)	40 dB
Air Tank	1 Gallon
Water Module	
Reservoirs	Two - 1 Liter - Toggled
Pressurized System	30-40 PSI
Aseptic	Non-Retracting HPC Coolant
Handpiece Flush System	Manual
Air and Water Connections	External Quick Connects
Other	
Duty Cycle	Continuous Operation Note: Exception, Units with LED Scaler, 10 min. ON, 5 min. OFF
Classifications	Class 1, Type B Applied Part, except as noted for optional accessories, Ordinary Equipment [IPXO]
Optional Accessories:	
Acteon Satelec Scaler	Type BF Applied Part
High-Speed Handpiece with and without out Fiber Optics	Type B Applied Part
Low-Speed Handpiece	Type B Applied Part
Air Water Syringe	Type B Applied Part

Specifications *continued...*

Electrical (See Regulatory Compliance Chart - Next Page)	
Fuse (Scaler Transformer)	T0.8AL, 250V, 3AG
Circuit Breakers (2) 115VAC	6A 250V, Push to Reset
Circuit Breakers (2) 250VAC	3A 250V, Push to Reset
Power Cord (115VAC)	15A/125VAC x 10 ft. Hospital Grade
Power Cord (230VAC)	10A/250VAC x 2.5 meters
Handpiece	
Connection Type	Midwest 4-Hole Type Handpiece

WARNING

Equipment not suitable for use in the presence of a flammable anesthetic mixture.

WARNING

*Do not modify this equipment without authorization
of the manufacturer.*

WARNING

*To avoid the risk of electric shock, this equipment must ONLY
be connected to supply mains with protective earth.*

WARNING

*The Acteon Newtron® LED scaler complies with current
electromagnetic compatibility standards. However, it must not
be used by or on those with cardiac stimulators or other active implants
unless adequate electromagnetic interference shielding is present.
Check with a cardiologist to verify the EMI shielding capabilities of the
implant. See enclosed Newtron® LED manual for further advisories.*

Model Identification / Compliance Chart

Model	Description	Complies To:				Electrical Ratings:		
		UL 60601-1, Ed. 1	CAN / CSA 22.2, #601.1-M90	IEC 60601-1, Ed. 2	CE	VAC +/- 10%	Amps	Cycles (Hz)
8000-001	Midmark 1000 w/ Scaler	•	•	•	•	115	6	60
8000-002	Midmark 1000 w/o Scaler	•	•	•	•	115	6	60
8000-003	Midmark 1000 w/ Scaler and Fiber Optics	•	•	•	•	115	6	60
8000-004	Midmark 1000 w/o Scaler and Fiber Optics	•	•	•	•	115	6	60
8000-005	Midmark 1000 w/ Scaler, Export				•	230	3	50
8000-006	Midmark 1000 w/ LED Scaler and Fiber Optics. Export				•	230	3	50
8000-006S	Midmark 1000 w/ LED Scaler and Swivel Fiber Optics. Export				•	230	3	50
8000-007	Midmark 1000 w/ LED Scaler and 2 Position Fiber Optics	•	•	•	•	115	6	60
8000-007S	Midmark 1000 w/ LED Scaler and 2 Position Swivel Fiber Optics	•	•	•	•	115	6	60
8000-008	Midmark 1000 w/ Scaler and Fiber Optics	•	•	•	•	115	6	60
8000-008S	Midmark 1000 w/ Scaler and Swivel Fiber Optics	•	•	•	•	115	6	60

Warranty

SCOPE OF WARRANTY Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will repair or replace components of the animal health products manufactured by Midmark (except for products and components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable products or components. This limited warranty shall only apply to defects that: (i) are reported to Midmark within the applicable warranty period; and (ii) are determined to exist upon examination by Midmark. This limited warranty extends only to the original retail purchaser of a product and is not transferable or assignable.

APPLICABLE WARRANTY PERIOD The applicable warranty period, measured from the date of delivery to the original user, shall be as follows: **The VetPro® 1000 unit is warranted against defect in material and workmanship for a period of two years from the time of delivery.**

PRODUCTS are warranted against defect in material and workmanship for the period indicated below from the time of delivery:

- a) 2 years for all products and components (except for the periods in b and c).
- b) Accessories not manufactured by Midmark are excluded (i.e. Fiber optic systems, Scalers, Handpieces, etc.).
- c) Replacement parts and accessories carry a 90 day warranty

OBTAINING WARRANTY SERVICE Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at www.midmark.com, by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380, or by phone at: 1.800.MIDMARK. It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident (including animal acts of any kind), freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) matching of color, grain or texture except to commercially acceptable standards;
- (3) changes in color caused by natural or artificial light;
- (4) products which are not installed, used and properly cleaned and maintained as required in the installation and operation manuals for the applicable product (imaging product must be installed by a certified Midmark installer);
- (5) products considered to be of a consumable nature;
- (6) accessories or parts not manufactured by Midmark;
- (7) specially manufactured products;
- (8) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (9) costs and expenses of routine maintenance and cleaning;
- (10) all sinks, faucets and plumbing accessories;
- (11) representations and warranties made by any person or entity other than Midmark; and
- (12) with respect to software that is a product or a component thereof, that the software will be error free, can be used without problems or interruptions, or will be free from vulnerability to intrusion or attack by viruses or other methods.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER Midmark's only obligation under this limited warranty is the repair or replacement of defective parts. Midmark shall not be liable for and hereby disclaims any direct, special, indirect, incidental, exemplary or consequential damages or delays including, but not limited to, damages for loss of profits or income, loss of use, downtime, cover and employee or independent contractor wages, payments and benefits. This disclaimer shall survive any failure or asserted failure of the essential purpose of this limited warranty or its remedies specified herein.

NO AUTHORIZATION No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with Midmark products.

WARRANTY DISCLAIMER: THIS LIMITED WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

STATUTE OF LIMITATIONS No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the warranty period. In the event multiple warranty periods exist with respect to a product, the ninety (90) day period provided for herein shall begin to run from expiration of the warranty period for the component to which the claim relates.

SEVERABILITY In the event any provision of this limited warranty is determined to be invalid or otherwise unenforceable: (i) the provision shall be enforced in a manner that closest holds to its intent, while at the same time curing the invalidity or unenforceability and (ii) the balance of this limited warranty, being severable, shall not be affected in the event the provision that is invalid or unenforceable cannot be enforced in any respect.

Midmark Corporation
60 Vista Drive
Versailles, OH 45380 USA
1-800-643-6275
1-937-526-3662
www.midmark.com

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