

Warranty Information



SCAN ME

Warranty Registration



SCAN ME

General Information

Attention: Wired Foot Controls

Wired foot controls are shipped standard with all models. A wired foot control is required to perform the Association Procedure. This procedure may be required in the event of a malfunction, do not discard the wired foot control.

(Refer to: **Association Procedure** in this User's Guide)

Intended Use

This product is intended to control patient positioning while providing better patient access by eliminating the cord from the foot control.

Calling for Service

If service is required, contact Midmark directly:

US 844.856.1230 / Canada 937.526.8585
Monday-Thursday 8:00 AM until 6:00 PM
Friday 8:00 AM until 5:00 PM (ET)

Compliance Information

FCC Compliance Statement

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Any modifications to this product may violate the rules of the Federal Communications Commission and make operation of the product unlawful.

Midmark Corporation
60 Vista Drive
Versailles, OH 45380-0286
Phone: 937-526-3662
Fax: 937-526-5542
midmark.com



MA7942revB

This guide outlines basic operation of the wireless foot control accessory. If service is required, contact your authorized Midmark dealer.

Dealer (name / phone):

Date of Purchase:

Association Procedure

The wireless controller will only function with the chair it has been “associated” with. This prevents unwanted movement in locations with more than one chair.

Each wireless controller’s association is preset at the factory. This procedure is only required in the event of a malfunction, if a new controller is purchased, or if you are attempting to use the controller with a chair other than the one it was originally associated with.

Association Procedure:

A) Using the wired foot control, press & hold the **Stop & QuickChair** buttons simultaneously.

(After two seconds, you will hear a single “beep”)

B*) Using the wireless foot control, press & hold the **Stop & QuickChair** buttons until you hear three “beeps”.

(One long “beep” indicates the association failed. Unplug unit for one minute and repeat Steps A & B.)

* Note: Step B must be performed within 3 seconds of Step A.


C) Check for proper operation of wireless foot control.



MA7942revB

Specifications

Range of Operation..... 20 ft. (6 m) radius
Batteries..... AA (qty. 2)



MA7935i

Troubleshooting

- When any control button is pressed, that controller sounds a single “beep”. (All functions operate)

Cause: Low battery

Solution: Replace batteries

- When any control button is pressed, nothing happens. (No “beeps”, no movement, etc.)

Cause: 1. Batteries completely dead
2. Wireless controller has lost “association” to chair

Solution: 1. Replace batteries
2. Perform **Association Procedure***

* The wireless control will only function with the chair it has been “associated” with. This prevents unwanted movement in locations with more than one chair.

If further service is required, contact your authorized Midmark dealer.