Elevance® Chair Setup Guide

Applies to Models:
153810 (all)

Special Tools:
none

Note

Chair is shipped on a box covered skid. Hardware and the components ordered are in boxes film-strapped to the chair. Upholstered seat and back ship in a separate box. Head rest comes in another box.

Unpacking the Chair

WARNING


During transportation, the chair must be at its lowest height and all attachments secured in their lowest and most central positions possible. Failure to comply may result in injury and/or damage to equipment.


Equipment Alert

We recommend anchoring your Elevance Chair to the floor for optimal stability when mounting a delivery unit.

Anchor holes (choose one)

DA3085i

Step 1: Remove all six skid clamps.

Step 2: Remove boxes strapped to chair seat.

Note: You may have one or two boxes dependant on the components ordered.

Step 3: Remove and discard two desiccant packets (not shown).

Step 4: Remove chair from skid.

Step 5: Remove shipping strap from chair.

Note: Strap holds chair to the base, preventing air from entering the drive cylinder. Only remove it after chair is positioned on the floor.

DA2322

Installing Upholstery and Connecting Chair to Power Supply

Step 1: Install the upholstery. Reference the Elevance Upholstery Set Up Guide packaged with the upholstery.

Step 2: Connect chair to the proper voltage, grounded Power Supply. The Power ON indicator light illuminates when power is supplied to chair.

Step 3: Tap any one of the lift arm covers 3 times to raise the chair to the pre-programmed “home” position, a convenient height for installation or patient access.


Note: The chair beeps as it moves and stops beeping when it reaches the home position.

Step 4: Remove the shipping pin from the chair pivot.

Note: The shipping pin prevents the chair from pivoting during shipping.

Note

A chair ordered alone comes with a standard 3-prong electrical plug. A chair ordered with delivery units or lights comes with an IEC plug that connects to a Power Supply unit.

Step 5: Test the base and back movements by pressing buttons on each control device ordered.

Special Programming Features - Button 4



The number 4 Position Button on the control devices is capable of being set to one of three operation modes. Which operating mode the number 4 button invokes depends on the dipswitch settings of the PROG SW on the PC Board. The operation modes and their respective dipswitch settings are shown below.

Equipment Alert

Operation Mode Switches should be set before the chair is powered on. Changes in Operation Mode Switch settings do not “register” until the power is cycled off and back on.

Push in with a screw driver to release clips and lift Upper Arm Lift Cover up to free front clips from shaft. Lift Cover must be removed to access Chair PC Board and cover.


Cuspidor Position Mode Setting

Press button 4 to raise the chair (patient) to the cuspidor position, a position convenient for cuspidor use. The cuspidor solenoid output also activates on the PC board.

The cuspidor solenoid output (J5) can be connected to an optional automatic bowl flush solenoid, initiating an automatic bowl flush function as the chair raises to the cuspidor position.

Press button 4 again to return chair to last position used before the cuspidor position.


Return to the Last Position Mode Setting

Press button 4 to move the chair back and forth between two “remembered” positions.

The chair remembers any position it sits motionless in for 2 seconds or more.

For example: Move chair to desired position, leave it sit still for 2 seconds. Use arrows to move chair to a 2nd desired position and leave it sit still for 2 seconds. Now pressing button 4 moves the chair back and forth between these two desired, remembered positions.

Normal Mode Setting - [This is the factory default setting for button 4.](#)

Press button 4 to move the chair to a programmed position.

Program button 4 the same as any other position button.

Program position 4 button by moving chair to desired position, press and hold button 4 (2 seconds) until you hear 3 beeps, indicating the position is programmed in memory.

DA2318

Chair Calibration Procedure

Note


Calibration is only required if the chair is not operating at all, or if it is not raising/lowering to the proper extents. Power to the chair must also be cycled after moving the #4 dipswitch to complete the procedure. Calibration should always be performed after replacing a chair's PC Board, a tilt sensor or any lift mechanism.

Movement of the chair is monitored by sensory devices that keep the chair within the range of motion accommodated by the drive mechanisms controlling the chair. **The Calibration Mode** is a programmed routine of chair movements that locates the end of travel points, calculates working values for these points, stores these values in memory and then tests the results of the routine. Follow these steps to calibrate a chair.


Step 1: Unplug the chair from the power supply.

Step 2: Remove the Upper Lift Arm Cover and PC Board cover (2 screws). Locate S1 and S2 switch components on the Chair PC board.

Push in with a screw driver to release clips and lift Upper Arm Lift Cover up to free front clips from shaft. Lift Cover must be removed to access Chair PC Board and cover.

Step 3: Move the #4 dipswitch on the S2 program switch to the Calibration Mode position (ON).

DA2319

Step 6: Install PC Board Cover and Upper Lift Arm Cover. Plug chair back into power supply. Activate control buttons to test chair operation.

Step 5: Unplug the Chair and move dipswitch #4 on S2 program switch out of Calibration Mode (OFF position).

Equipment Alert

Control buttons do not work normally during a Calibration Mode. Pressing any button (including the S1 Test pushbutton) will stop the chair movement and end the Calibration Mode.

Step 4: Plug chair into power supply and press the TEST button on switch S1.

As soon as the TEST button is pressed, Calibration begins. Step back and allow the chair to complete all 3 up/down cycles of the calibration routine.

Note: The chair emits a short beep every 2 seconds while calibrating. When finished calibrating the chair emits 3 long, confirmation beeps.

The Calibration Mode has failed if the chair emits 1 long beep at the end of the calibration mode.

Manual Calibration - use to set your own high and low extents for base and back travel.

Note

Manual Calibration is only required to modify one or more of the limit of travel extents the base and back are set to originally.
Power to the chair must be cycled after dipswitch #4 is returned to OFF position to complete the procedure.

Base and back movements are monitored by sensory devices that keep them within a controlled range of extents. **Manual Calibration** is a modified calibration routine that allows you to set your own end of travel points, stores these values and then tests the results of the routine. Below is the Manual Chair Calibration procedure.

Manual Calibration...

- A) Unplug chair from power supply.
- B) Remove upper lift arm cover and PC Board cover (two screws). Locate S1 and S2 switches on the PC board.
- C) Move dipswitch #4 on the S2 Prog SW to ON position.
- D) Plug chair into power supply, press and hold the TEST button on S1 switch (for five seconds). Listen for a beep to indicate it is in Manual Calibration mode. Beeps continue, one every 5 seconds, as long as you are in manual calibration.
- E) Extents can be manually set or skipped for each of these base and back travel points using the same keypad selections.

Extents must always be set in this order:	To set an extent:	
	Move base or back to desired extent, then press and hold button 1 for two seconds to set the extent.	Two beeps indicate extent is set.
	OR, to skip setting current extent:	
	Press and hold button 4 for two seconds.	Two beeps indicate extent remains at chair limit.
Base Low		
Back High		
Base High		
Back Low		


Note: You can abort the manual calibration procedure any time by pressing the TEST button on S1 switch.

- F) As soon as the final (Back High) extent is set the calibration routine begins. Step back and allow the chair to complete the necessary number of up/down cycles, which takes a few minutes.

Note: The chair emits a short beep every 2 seconds while calibrating. When calibration is successful, the chair emits three long, confirmation beeps.

If calibration failed, the chair emits one long beep at the end of the routine. Repeat entire Manual Calibration.


Push in with a screw driver to release clips and lift Upper Arm Lift Cover up to free front clips from shaft. Lift Cover must be removed to access Chair PC Board and cover.


Manual Calibration continued...

- G) Unplug chair from power supply and move dipswitch #4 on S2 program switch out to OFF position.
- H) Install covers (see B).
- I) Plug chair back into to power supply.
- J) Activate keypad and remote buttons to test operation.


Chair Dimensions

Chair Dimensions continued...

Set Up Checklist

DA2344

Test	First with an empty chair, then with a person seated in the chair, verify the following:	✓
1	Power On light illuminates when chair is plugged in.	
2	Tapping any Lift Arm cover 3 times moves chair to home position.	
3	Foot Control works (optional, may or may not have one).	
4	Remote Control works (optional, may or may not have one).	
5	Brake release pedals work (depress pedal and rotate chair, release pedal to stop rotation).	
6	Chair Lockout works - depress brake release pedal 5 times to lockout all chair movement and rotation - depress pedal 5 times to release the chair lockout.	
7	Headrest moves with ease and holds in place.	

Elevance Chair

Base Weight of Chair:	426 lbs / 193 kg
Shipping Weight:	506 lbs / 230 kg
Patient Weight (max):	450 lbs / 204 kg
Power Cord Length:	4 or 8 feet / 1.22 or 2.44 meters
Electrical Rating:	115 V~
	60Hz, 9A
	230 V~
	50 Hz, 6A
Duty Cycle:	Intermittent operation
	(Pump Motor and Brake: 30 seconds ON / 5 minutes OFF)
Fuses - 115 VAC Models:	F1 & F2 (line input) 250 V~, T6.3AL
	F3 (transformer) 250 V~, T0.16AL
Fuses - 230 VAC Models:	F1 & F2 (line input) 250 V~, T3.15AL
	F3 (transformer) 250 V~, T0.080AL
Classifications:	Class I, Type B Applied Part, except as noted for optional heated back upholstery, Ordinary Equipment, Intermittent Operation
Certifications:	ES/IEC/EN 60601-1, CAN/CSA C22.2 No. 60601-1 (Safety Standards) EN/IEC 60601-1-2 (EMC Standards)
	Optional heated back upholstery, Type BF Applied Part, complies with IEC 60601-2-35
	IPX0 Dental Chair IPX1 Foot Control IPX2 Heated Back Upholstery

This page is intentionally left blank

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