

# Warranty Information



SCAN ME

# Warranty Registration



SCAN ME

## General Information

### Intended Use

*This product is intended to control patient positioning while providing better patient access by eliminating the cord from the controller(s).*

### Calling for Service

*If service is required, contact Midmark directly:*

US 844.856.1230 / Canada 937.526.8585  
Monday-Thursday 8:00 AM until 6:00 PM  
Friday 8:00 AM until 5:00 PM (ET)

## Compliance Information

### FCC Compliance Statement

*This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their own expense.*

*Any modifications to this product may violate the rules of the Federal Communications Commission and make operation of the product unlawful.*

*This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.*

Midmark Corporation  
60 Vista Drive  
Versailles, OH 45380-0286  
Phone: 937-526-3662  
Fax: 937-526-5542  
midmark.com



## User's Guide 625 Wireless Controls 9A429



*This guide outlines basic operation of the wireless hand / foot control and the associated base station. If service is required, contact your authorized Midmark dealer.*

**Dealer (name / phone):** \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

**Date of Purchase:** \_\_\_\_\_

# Compliance Information - continued...

## IC Compliance Statement

### Interference

“This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.”

“Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes: (1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.”

### Antennas

“Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.”

“Conformément à la réglementation d'Industrie Canada, le présent émetteur radio peut fonctionner avec une antenne d'un type et d'un gain maximal (ou inférieur) approuvé pour l'émetteur par Industrie Canada. Dans le but de réduire les risques de brouillage radioélectrique à l'intention des autres utilisateurs, il faut choisir le type d'antenne et son gain de sorte que la puissance isotrope rayonnée équivalente (p.i.r.e.) ne dépasse pas l'intensité nécessaire à l'établissement d'une communication satisfaisante.”

## Association Procedure

The wireless controller will only function with the chair it has been “associated” with. This prevents unwanted movement in locations with more than one chair.

Each wireless controller's association is preset at the factory. This procedure is only required in the event of a malfunction, if a new controller is purchased, or if you are attempting to use the controller with a chair other than the one it was originally associated with.

### Association Procedure:


A) Disconnect power to the table for at least 3 seconds.

B\*) Reconnect power, wait for the power up beep sequence from the table, then press & hold the **Stop & Home** buttons on the **wireless** controller until you hear three “beeps”.

\* Note: Step B must be performed within 10 seconds of the power up beep sequence from the table.

C) Check for proper operation of wireless controller.

D) If associating a hand **and** foot control, repeat procedure using the other wireless controller.




# Specifications

Range of Operation..... 20 ft. (6 m) radius

### Batteries:

Hand Control.....AA (qty. 2)

Foot Control.....AA (qty. 2)



# Troubleshooting

When any control button is pressed, that controller sounds a single “beep”. (All functions operate)

Cause: Low battery

Solution: Replace batteries

When any control button is pressed, nothing happens. (No “beeps”, no movement, etc.)

Cause: 1. Batteries completely dead  
2. Wireless controller has lost “association” to chair

Solution: 1. Replace batteries  
2. Perform **Association Procedure\***

The wireless controllers will only function with the chair they have been “associated” with. This prevents unwanted movement in locations with more than one chair.

If further service is required, contact your authorized Midmark dealer.