Barrier-Free® Examination Table

[English](#)
[Español](#)
[Français](#)

For Models:
625 (-001 / -003)
625 (-005 / -006)

User's Guide

Product Information

Dealer :

Date of Purchase:

Model / Serial Number:

**Midmark Authorized Service
Company:**

**Model / Serial Number Label
Located on the Back Mounting Frame**
(exact location may vary)


Table of Contents

Important Information

| | |
|--|---|
| <u>Safety Symbols</u> | 4 |
| <u>Transportation / Storage Conditions</u> | 4 |
| <u>Compliance Information</u> | 5 |
| <u>Intended Use</u> | 6 |
| <u>Electrical Requirements</u> | 6 |
| <u>Electromagnetic Interference</u> | 6 |
| <u>Disposal of Equipment</u> | 6 |
| <u>Safety Instructions</u> | 7 |

Operation

| | |
|--|----|
| <u>Base / Back Functions</u> | 8 |
| <u>Home Function</u> | 9 |
| <u>Quick Exam® Function</u> | 10 |
| <u>Exam Assistant® Drawer System</u> | 11 |
| <u>Paper Roll & Tear Strap</u> | 12 |
| <u>Stirrups</u> | 13 |
| <u>Footrest & Treatment Pan</u> | 14 |
| <u>Drawer Heater (optional)</u> | 14 |
| <u>Pelvic Tilt</u> | 15 |
| <u>Duplex Receptacle</u> | 16 |

Maintenance

| | |
|---------------------------------------|----|
| <u>Calling for Service</u> | 17 |
| <u>Cleaning</u> | 17 |
| <u>Preventative Maintenance</u> | 17 |

Specifications

| | |
|---|----|
| <u>Weights, Dimensions, Electrical Specifications</u> | 18 |
| <u>Model Identification / Compliance Chart</u> | 19 |
| <u>Range of Motion & Dimensions</u> | 20 |

Warranty Information

| | |
|-------------------------------|----|
| <u>Limited Warranty</u> | 22 |
|-------------------------------|----|

Important Information

Safety Symbols

DANGER

Indicates an imminently hazardous situation which will result in serious or fatal injury. This symbol is used only in the most extreme conditions.

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage.

Note


Amplifies a procedure, practice, or condition.

Consult User Guide for important information

Protective earth ground

Duty Cycle (motor run time):
30 seconds ON, 5 minutes OFF

Fuse rating specification

Keep dry

Type B, Applied Part

Proper Shipping Orientation

Maximum stacking height
(palletted units)

Fragile

Emergency Stop

Patient weight limit

Transportation / Storage Conditions

Ambient Temperature Range:.....41°F to 104°F (5°C to 40°C)
Relative Humidity.....10% to 90% (non-condensing)
Atmospheric Pressure500hPa (.49 atm) to 1060hPa (1.04 atm)

Compliance Information

FCC Compliance Statement

This equipment has been tested and found to comply with the limits for a class A digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their own expense.

Any modifications to this product may violate the rules of the Federal Communications Commission and make operation of the product unlawful.

IC Compliance Statement

Interference

“This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.”

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Antennas

“Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.”

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Intended Use

This exam table is intended to be used as a chair / table to provide positioning and support of patients during general examination procedures conducted by medical professionals.

Electrical Requirements

Caution

When using high frequency surgical devices or endocardial catheters:

- *Use non-conductive material to insulate patient from metal portions of the chair.*
- *Consult operating instructions for the device before using in conjunction with drawer or upholstery heaters.*

Failure to comply may result in electric shock or burns to patient.

Electromagnetic Interference

This product is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this product:

- Remove interfering device from room
- Plug chair into isolated circuit
- Increase separation between chair and interfering device
- Contact Midmark if interference persists

Disposal of Equipment

At the end of product life, the chair, accessories, and other consumable goods may become contaminated from normal use. Consult local codes and ordinances for proper disposal of equipment, accessories and other consumable goods.

Safety Instructions



Caution

Be sure that all personnel and equipment are clear of the table before activating any function. Failure to do so could result in personal injury.

At the table's lowest height [18 in. (45.7 cm)], the clearance between the bottom of the drawer and the floor is approximately 2.1 in. (5.4 cm). Make sure patient (esp. feet) and all objects are clear of this area before lowering table.

If the table malfunctions, immediately release foot pedal, unplug the table power cord, and assist the patient from the table. If table continues to malfunction, call for service.

Failure to comply with these instructions may result in injury to patient and/or damage to equipment.

MA8022i

Operation

Base / Back Functions

To operate the Base / Back functions...

Press & hold the Up or Down arrow button for the desired function.

Note: The table automatically stops at its maximum and minimum height.
In the event of a malfunction, press the Stop button to terminate movement.

Equipment Alert

If the table's weight limit is exceeded, you will hear a "beep" and the function will be disabled. If the problem occurs with a normal load, contact your service technician.

'Home' Function

The Home function automatically lowers the table to approximately 20" (51 cm).

To activate the Home function...
Press & release the Home button.

Note: In the event of a malfunction, press the Stop button to terminate movement.

Quick Exam® Function

The Quick Exam® function automatically commands the table seat height to a programmed position.

To program the Quick Exam® function...

- Press the Base Up button to move the table to the desired position.
- Press & hold the Quick Exam® & Base Up buttons simultaneously for three seconds.

Note: You will hear a brief chirp when the new Quick Exam® position is saved.

To activate the Quick Exam® function...

Press & release the Quick Exam® button.

Note: In the event of a malfunction, press the Stop button to terminate movement.

Equipment Alert

If the table's weight limit is exceeded, you will hear a "beep" and the function will be disabled. If the problem occurs with a normal load, contact your service technician.

Exam Assistant® Drawer System


To use the drawer front as a shelf...

Pull drawer front down to release the magnetic latch.

Equipment Alert

The maximum weight the shelf will support is 10 lbs (4.5 kg).


If drawer front folds past horizontal...

Pull drawer front upward until it snaps back into the horizontal position.

Equipment Alert

If drawer front folds past horizontal, do not operate table.

MA7996i

Paper Roll & Tear Strap

Upholstery sets are available in two widths: 28" (standard) & 32" (optional).
Paper rolls of common size may be stored under the back section as follows:

| <u>Uph. Width</u> | <u>Paper Roll Size</u> | <u>Max. Qty.</u> |
|-------------------|--|------------------|
| 28 in. | 21 in. long x 3.0 in. diameter (53.3 cm x 7.6 cm)..... | 4 |
| 28 in. | 21 in. long x 3.5 in. diameter (53.3 cm x 8.9 cm)..... | 2 |
| 32 in. | 24 in. long x 3.0 in. diameter (70 cm x 7.6 cm)..... | 2 |
| 32 in. | 24 in. long x 3.5 in. diameter (70 cm x 8.9 cm)..... | 1 |


Warning

To avoid injury, do not put head or extremities under the paper tear strap.


To install paper roll...

Insert rod ends into sockets
under head end of upholstery.

To install paper tear strap...

Attach strap to snaps on each side of table.

MA8018i

Stirrups


Caution

*Be sure the stirrups are locked in place before using.
The stirrups will not support the patient's entire weight.*

To position the stirrups...

- A) *Pull the stirrup out, then unfold.*
- B) *Lift the stirrup slightly, then move it left or right as desired.*
- C) *Release stirrup to lock in position.*


MA7995i

Footrest & Treatment Pan

Caution

Do not use footrest to support patient's entire weight.

To access the treatment pan...

Lift footrest pad slightly, then slide pad back into table.

To extend footrest...

Pull footrest out to desired position.


MA7994i

Drawer Heater

The drawer heater warms the contents of the drawer to approximately 98°F (37°C).

To turn the drawer heater ON / OFF...

*Press and release the drawer heater switch.
(Switch will illuminate when heater is ON)*


MA7990i

**Drawer Heater
Switch**

Pelvic Tilt


To elevate seat section...

Lift seat until pelvic support drops into place.

To return seat section to flat position...

- A) *Lift seat section slightly.*
- B) *Pivot pelvic support handle upward.*
- C) *Lower seat section.*


MA7991i

Duplex Receptacles

The duplex receptacles on each side of table provide power for accessories used during medical procedures. There are two circuit breakers located under the seat section. If the receptacles' maximum load is exceeded, the circuit breakers interrupt power to the receptacles.

Maximum Load.....115 VAC, 5 amps

To reset the circuit breaker...
Push circuit breaker switch.

MA7992i

Maintenance

Calling for Service

Note

Model / serial number information is required when calling for service.

If service is required, contact your authorized Midmark dealer.

To contact Midmark directly:

1-800-MIDMARK (1-800-643-6275)

8:00 am until 5:00 pm Monday - Friday (EST) [excluding standard U.S. holidays]

www.midmark.com

Cleaning

Upholstery

Equipment Alert

The upholstery is resistant to most medicinal-type stains, but may be damaged by solvents and dyes. Immediately remove any fluids spilled on the upholstery.

Wash your upholstery weekly with a mild liquid soap and water mixture, rinse with clear water and dry completely to remove disinfection cleaner build-up.

Disinfect your upholstery with a solution of standard bleach and water mixed 1 in 10 (10%) or chlorine based cleaners. Follow this with a clear water rinse and thorough drying of material. See current CDC Guideline for Disinfection & Sterilization in Healthcare Facilities.

To minimize damage caused by disinfectant cleaner residue build-up, do not allow disinfectants to pool on the upholstery surface. Once the approved contact time has been obtained, remove and dry excess liquid remaining on the surface.

Detailed care and maintenance instructions are included with your product. This information is also available on www.midmark.com in the Technical library under the User Information tab for your product.

Painted Metal / Plastic Surfaces

Clean the painted metal and plastic surfaces weekly using a clean soft cloth, and mild cleaner. Periodic applications of common furniture wax will ease cleaning, and maintain the finished luster of the surfaces.

Preventative Maintenance

Periodically inspect the following areas:

- Power cord(s) should be free of cuts or other visible damage.
- All fasteners should be in place and tightened securely.
- All mechanical functions should operate properly.

Periodically lubricate the back hinge to maintain quiet, smooth, operation (use light machine oil).

Have an authorized service technician inspect your equipment every six months.

Specifications

Weights, Dimensions, Electrical Specifications

| Specifications Chart | |
|---|---|
| Patient Weight (<i>maximum</i>): | 650 lbs (295 kg) |
| Paper Roll (<i>maximum size</i>): | (See <i>Paper Roll & Tear Strap</i> page) |
| Range of Motion & Dimensions: | (See <i>Range of Motion & Dimensions</i> page) |
| Weight of Table: w/upholstery w/packaging & skid (<i>no uph.</i>) Uph. w/packaging (<i>shipped separately</i>) | 470 lbs (213 kg) 487 lbs (221 kg) 45 lbs (20 kg) |
| Power Cord Length: | 8 ft. (244 cm) |
| Electrical Requirements: | (See <i>Regulatory Compliance Charts</i>) |
| Foot / Hand Control Voltage: | 14 VAC, SELV (Safety Extra Low Voltage) |
| Duplex Receptacles (<i>maximum load</i>) | 115 VAC, 5A, 50/60 Hz |
| Fuses: <i>IEC inlet fuses (2)</i> <i>F1 & F2 (on main PC board):</i> <i>F3 (on main PC board):</i> <i>F4 (on main PC board)</i> <u><i>Models w/ drawer heaters</i></u> <i>distribution board fuses (2)</i> | 6.3A, 250V, Type-T, 5 x 20 mm 10A, 250V, Type-T, 5 x 20 mm 160mA, 250V, Type-T, 5 x 20 mm 6.3A, 250V, Type-T, 5 x 20 mm 800mA, 250V, Type-T, 5 x 20mm |
| Duty Cycle (Motor Run Time): | Intermittent Operation [30 seconds ON / 5 minutes OFF] <i>Continuous Operation</i> |
| Receptacles, Drawer Heater: | |
| Classifications: | Class 1, Type B Applied Part |
| Protection against ingress of fluids: | IPXO [<i>Foot Control only: IPX1</i>] |

WARNING



Equipment not suitable for use in the presence of a flammable anesthetic mixture.

Model Identification / Compliance Chart



| Model | Description | Complies To: | | | | | Electrical Ratings: | | |
|---------|--|---------------|----------------------------------|---------------|--------------------------|------------|---------------------|------|----------------|
| | | UL 60601-1 | CAN / CSA 22.2, #601.1-M90 | EN 60601-1 | EN 60601-1-2 (EMC) | NFPA 99 | VAC +/- 10% | Amps | Cycles (Hz) |
| 625-001 | Two-function table (Base / Back), w/ receptacles, pelvic tilt, & drawer heater | • | • | | • | • | 115 | 12 | 50/60 |
| 625-003 | Two-function table (Base / Back), w/ receptacles, pelvic tilt, drawer heater, & wireless controls | • | • | | • | • | 115 | 12 | 50/60 |
| 625-005 | Two-function table (Base / Back), w/ pelvic tilt | | | • | • | | 230 | 4 | 50/60 |
| 625-006 | Two-function table (Base / Back), w/ pelvic tilt, drawer heater | | | • | • | | 230 | 4.5 | 50/60 |

Fire Code Ratings: All upholstery complies with California Bureau of Home Furnishing Technical Bulletin 117 and California Code of Regulations, Sect. 93120-93120.12, Title 17. Optional upholstery is available that complies with California Bureau of Home Furnishing Technical Bulletin 133.

Range of Motion & Dimensions

Range of Motion & Dimensions (continued)

MA8369i

Warranty Information

SCOPE OF WARRANTY Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

Midmark warrants to the original retail purchaser that during the applicable warranty period it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications.

APPLICABLE WARRANTY PERIOD The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product and shall be one (1) year for all warranted products and components.

OBTAINING WARRANTY SERVICE Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at www.midmark.com; by phone at 1-800-MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380. It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS: This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Products considered to be of a consumable or sterile nature;
- (4) Accessories or parts not manufactured by Midmark;
- (5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) Costs and expenses of routine maintenance and cleaning; and
- (7) Representations and warranties made by any person or entity other than Midmark.
- (8) Matching of color, grain or texture except to commercially acceptable standards;
- (9) Changes in color caused by natural or artificial light;
- (10) Custom manufactured products;
- (11) Alterations or modifications to the product by any person or entity other than Midmark; and
- (12) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market.

SOFTWARE; WITH RESPECT TO SOFTWARE THAT IS A PRODUCT OR COMPONENT THEREOF, MIDMARK DOES NOT WARRANT THAT THE SOFTWARE: (1) IS ERROR FREE; (2) CAN BE USED WITHOUT PROBLEMS OR INTERRUPTIONS; OR (3) IS FREE FROM VULNERABILITY TO INTRUSION OR ATTACK BY VIRUSES OR OTHER METHODS.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN. WARRANTY DISCLAIMER THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

NO AUTHORIZATION No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

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