Attention!

The electrical and plumbing supply connectors inside the J-Box (Junction Box) must be installed by an electrician and plumber in compliance with local building codes. All required hardware is to be supplied by contractor.

J-Box must be installed within reach of the chair umbilical [4 or 8 foot (1.2 or 2.4 m)].

Equipment Alert

When installing 115~ (VAC) receptacle in the J-Box, position it so that all the other components needed in the box will also fit.

Grounding reliability can only be achieved if this equipment is connected to an equivalent receptacle marked Hospital Grade (NEMA 5-15R, HOSPITAL GRADE).

Electrical Supply Requirements

Have an electrician install a 1/2 in. (13 mm) conduit box with quad or equal receptacle. Top of box must be no higher than 4-1/2 in. (114 mm) above finished surface of floor.

Water Supply Requirements (optional)

Use 1/2 in. (13 mm) pipe (supplied by contractor) protruding 2 in. (51 mm) above finished surface of floor. Manual shut-off valve is to be installed by a plumber.

Water pressure must be 30-50 (max.) PSI at the shut-off valve. Flush water plumbing clean before making final connections to dental equipment.

Recommended water supply hardness level <60 ppm of calcium carbonate.

Air Supply Requirements

Use 1/2 in. (13 mm) pipe (supplied by contractor) protruding 2 in. (51 mm) from finished surface of floor. Manual shut-off valve is to be installed by a plumber.

Air pressure must be 70-100 (max.) PSI at the shut-off valve. Flush air plumbing clean before making final connections to dental equipment.

Air supply to be filtered, dry and oil-free per NFPA 99 requirements or equivalent local codes.

Drain Requirements (optional)

Not less than 1-1/2 in. (38 mm) I.D. pipe from downpipe to trap, decreasing to 3/4 in. (19 mm) pipe from trap to connection box. Terminate with 3/4 in. (19 mm) NPT pipe protruding 1 in. (26 mm) above finished surface of floor. Install adapter for 5/8 in. (16 mm) drain tubing. Place trap in line and vent to conform with local codes. [All hardware supplied by contractor].

Vacuum Requirements (optional)


Plumbing to J-Box box should be specified by central vacuum supplier and terminated in J-Box with 3/4 in. (19 mm) O.D. tube protruding 1 in. (26 mm) above finished surface of floor.

Installing the J-Box


Note

Loose tubing is provided for optional equipment and may not be required for your installation. Reference the tubing diagram in the Installation Folder for proper connections.



Step 1: Place J-box over the utility connections. Brackets are present to screw J-Box to the floor if desired.

Step 3: Remove umbilical clamp and position the umbilical in the opening as shown. Replace umbilical clamp to secure.


Step 2: Use compression fittings to connect air and water regulators to the facility utilities.



Step 4: Connect the water, air and vacuum lines for the components in this installation.

Step 5: After testing the delivery handpieces to ensure all connections, place the lid on the J-Box.

DA2387-011

J-Box Dimensions and Facilities Connections

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