
Midmark Dental Light

- [English](#)
- [Español](#)
- [Français](#)

NOTE: This user guide applies to Halogen and LED-style lights.

Compliance Information & Symbols	
Serial Number Location	2
Intended Use	3
Electromagnetic Interference	3
Transportation / Storage /	
Operating Conditions	3
Disposal of Equipment	3
Authorized CE Representative	3
Calling For Service	3
Symbols	4

Operation and Adjustments	
Operate the Light	5
Adjust the Intensity	6
Composite Setting	7
Light Focus Adjustment	8
Pivot Joint Tension Adjustments	9




Cleaning and Maintenance	
Barriers	10
Cleaning Procedures	10
Cleaning the Dental Light	11
Bulb and Parts Replacement	12

Specifications and Warranty	
Light Specifications	13
Dental Light Warnings	15
EMC - Manufacturer's Declaration and	
Guidance	16
Warranty Information	19



For installation information go to:
[Technical Library at Midmark.com](#)

Serial Number Location

Ceiling Mounted Lights

Wall/Cabinet Mounted Lights



Chair Mounted Lights

Caution

Be aware that the Chair Mounted light arm can separate from the lower arm if obstructed by a stationary object as the chair is lowered. Personal injury may result if this potentially hazardous situation is not avoided.

Compliance Information

Intended Use

This product is intended to be used as a component in a dental environment to provide lighting during dental procedures.

Electromagnetic Interference

Midmark dental operatory components are designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this operatory, remove the interfering device from room and / or plug product into an isolated circuit.

Transportation / Storage / Operating Conditions

Transportation / Storage Temperature:.....23°F to 100°F (-5°C to 38°C)
Relative Humidity..... 10% to 90% (non-condensing)
Atmospheric Pressure 7.2 PSI to 15.3 PSI (50 kPa to 106 kPa)
Operating Temperature Range: 59°F to 95°F (15°C to 35°C)

Disposal of Equipment

At the end of this product's life, the unit, accessories and other consumable goods may be contaminated from normal use. Consult local codes and ordinances for proper disposal of this equipment and other consumable goods.

Calling For Service

Direct all service inquiries to your authorized Midmark dealer. When calling for service, you must provide the following information:

- Model / serial number
- Date of purchase
- Symptom(s) of malfunction

Symbols

These symbols may appear on your equipment and/or in the manuals.

WARNING

Indicates a potentially hazardous situation which could result in serious injury if not avoided.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury if not avoided. It may also be used to alert against unsafe practices

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage if not avoided.

Refer to instruction manual/booklet

**Type B,
Applied Part**

Catalogue Number

Serial Number

Manufacturer


Fragile

**Proper Shipping
Orientation**

Relative Humidity

**Transportation &
Storage
Temperature**

**Atmospheric
pressure**

Keep Dry

**Maximum stacking
height (Refer to “n”
number on package.)**



**Transportation
Orientation**

Operate the Light

Press the Power Button on the touch pad to advance thru the operation modes *AUTO-On* / *Manual ON* / *OFF*.

To use the Auto-On function...

- A) Press the Power button until the Auto-On indicator light illuminates.
- B) Pull the lighthouse down (below horizontal) to turn the light ON.
- C) Move the light up (above horizontal) to turn the light OFF.


DA3455i

Adjust the Intensity


To Adjust Light Intensity...

Press (and release) the Intensity Button to advance through the intensity settings: High / Medium / Low.

Note: Indicator light shows current setting.

View Light Specifications


DA3457i

Composite Setting



This setting provides illumination for a procedure without the normal intensity that can prematurely cure composite dental materials.

To activate the Composite Setting...

Press (and release) the Composite Button to produce an extra, low-intensity light beam.

View Light Specifications

DA3460i

Note

The Procenter Delivery Unit does not have a Composite Button.

Note


The Elevance Chair wireless control also has a Composite Button.

Light Focus Adjustment (Halogen ONLY)

This adjustment controls the sharpness of the light pattern at a specified distance. The light's focus is set at the factory for a distance of 27 inches (69 cm). This distance is measured from the lighthead to the patient.

Step 1: Remove the back cover
Pry the back cover off light.


Note: May use a flat head screwdriver to remove.

Step 2: Set the target distance

Position the lighthead at desired distance from target (ex. wall) then turn light ON.

Note: Distance must be between 21-39 in. (53-99 cm).

Step 3: Focus Light

Use a Phillips head screw driver to adjust the Focus Screw.

Note: Suggested light focus is on oval shape approximately 3 x 8 in (8 x 20 cm).

Note

Light intensity decreases as the target focus area grows.


DA3461i

Pivot Joint Tension Adjustments

To adjust the Arm Pivot or Yoke Pivot Joint Tension...

- A. Remove screws and Pivot Cover using a 3/32" Allen wrench.
- B. Tighten or loosen the Rub Pad socket head screws at the Arm Pivot Joint or the Yoke Pivot Joint.

Note: Use a 3/32" Allen wrench to remove cover and make these adjustments.


DA3503i

To adjust the Arm Pivot Joint Tension...

- A. Adjust the set screw with 9/64" Allen wrench.

Halogen Shown

DA3462i

To adjust the Lighthouse Pivot Joint Tension...

- A. Remove Screw Cover.
- B. Loosen set screw.
- C. Adjust Tension Screw as necessary.
- D. Tighten set screw.
- E. Install Screw Cover.

Note: Use 5/64" Allen wrench for set screw.

Cleaning and Maintenance

Barriers

Single-use barriers and disposable items significantly reduce the need for chemical cleaners, thus prolonging the life of the equipment. Barrier material must be impervious to moisture / fluids.

Examples of protective barriers:

- Plastic covers (available from your dealer or equipment manufacturer)
- Clear plastic wrap
- Plastic bags
- Plastic sheets
- Plastic tubing
- Plastic-backed paper
- Materials similar to those listed here

Caution

Midmark assumes no responsibility or liability for any result, expressed or implied. These are suggested practices, based on the best information available at this time this is written.

Cleaning Procedures

Use cleaners that are appropriate for the situation, such as warm water and mild detergents, or a 10% solution of bleach with water.

NOTE

Every dental practice is different, and no single cleaner is the best choice for every facility. Listed below are some organizations to assist you in choosing the best cleaners available for your practice.

- **Organization for Safety & Asepsis Procedures:**


<http://www.osap.org>

- **American Dental Association:**

<http://www.ada.org>

- **Dept. of Health & Human Resources
Centers for Disease Control & Prevention (CDC):**

<http://www.cdc.gov>

DA148100i

**Read all labels
carefully!**


DA177400i

Cleaning the Dental Light

At the beginning of each work day...
Use barrier protection to cover surfaces that must be protected from contamination.

Note: Replace barriers for each patient.


To clean the outside of the light shield... (Halogen ONLY)
Gently wipe with a soft cloth using mild soap and water.

Caution

The bulb and surrounding components may be HOT! Allow the bulb to cool for at least five minutes before removing the light shield. Never operate the light with light shield removed.

To clean LED light ...
Use damp soft cloth with mild soap gently wiping exterior surface.

To clean the inside of the light shield and reflector... (Halogen ONLY)

- A) Turn light off and allow it to cool for at least 5 minutes.
- B) Unlock the two levers
- C) Gently wipe with a soft cloth using a mild soap and water..
- D) Replace the light shield and lock the two release levers

Bulb and Parts Replacement (Halogen ONLY)

Equipment Alert

Never touch a new bulb with bare hands. Leave plastic sleeve on bulb to protect it from bodily oils. Bodily oils cause "hot spots" on the bulb that will reduce bulb life.

Caution

The bulb and surrounding components may be HOT! Allow the bulb to cool for at least five minutes before removing the light shield. Never operate the light with light shield removed.


To replace the bulb...

- Allow bulb to cool for at least 5 minutes.
- Unlock the two release levers, remove light shield, then remove the bulb.
- Leave plastic sleeve on new bulb as you insert it into the light bulb socket, then remove the plastic sleeve.
- Install the light shield and lock both release levers.

Caution

Never operate the light when the Light Shield is removed.

DA3463i

Refer to Parts Manuals in the Technical Library on www.midmark.com for all Replacement Parts

Light Specifications

Midmark Dental Light Models:		
	Halogen	LED
Bulb (Light Assembly):	17 VAC, 95W, GY 6.35 size	No User Serviceable Parts
Fuse Rating:	Type 3AG, T2.5AL, 250 VAC 1/4 x 1-1/4 in.	Type 3AG, T2.5AL, 250 VAC 1/4 x 1-1/4 in.
Fuse Rating (TLM):	Type T, 1.6A, 250 VAC, 5 mm x 20 mm, High Breaking Capacity (HBC)	Type T, 1.6A, 250 VAC, 5 mm x 20 mm, High Breaking Capacity (HBC)
Color Temperature	Bulb...3000° K Light Assembly... 5000° K (+/- 500° K) (Off Reflector)	5700 Kelvin (Operation Mode)
Light Intensity:		
Low Setting	<1900 ft.-Candles (<20,500 lux)	<1275 ft.-Candles (<13,700 lux)
Medium Setting	<2800 ft.-Candles (<30,100 lux)	<2100 ft.-Candles (<22,600 lux)
High Setting	<3500 ft.-Candles (<37,700 lux)	<2375 ft.-Candles (<25,500 lux)
Composite Setting	< 750 ft.-Candles (<8,100 lux)	<2150 ft.-Candles (<23,100 lux)
Pattern Size	3-3/4" x 7-7/8" (9.5 cm x 20 cm)	3-1/3" x 7-1/2" (8.5 cm x 19.1 cm)
Weight of Light	27 pounds (12.2 Kilograms)	27 pounds (12.2 Kilograms)
Weight of Light and Arm:		
L/R Version:	40.1 lbs (18.2 kg)	40.1 lbs (18.2 kg)
Console Version:	28.5 lbs (12.9 kg)	28.5 lbs (12.9 kg)
Track Version:	19.5 lbs (8.8 kg)	19.5 lbs (8.8 kg)

Light Specifications...continued

Midmark Dental Light Model: (Halogen)		
Mounting Style	Electrical Rating:	Complies to the applicable requirements of:
Elevance and Ultra Chairs: 153818 (-001 thru -006) 153914 (-002)	See power supply output (153941 ...)	Class I, Type B Applied Part, Ordinary Equipment (IPX0), Continuous Operation ES/IEC/EN 60601-1, CAN/CSA C22.2 60601-1 (Safety Standards), EN/IEC 60601-1-2 (EMC Standards)
Knight Chair: 153937 (-001 thru -004)	Power Supply Model No.:153935-001, -002	Class I, Type B Applied Part, Ordinary Equipment (IPX0), Continuous Operation
Ceiling: 153637-001 Wall: 153640-001 Cabinet: 153641-001 Universal: 153639-001	115/230VAC, 2A/1A, 50/60 Hz	ES/IEC/EN 60601-1, CAN/CSA C22.2 60601-1 (Safety Standards), EN/IEC 60601-1-2 (EMC Standards)
Track Light/Monitor: 153829 (-000 thru -008) / 153963 (-000 thru -008)	115VAC, 2 A, 50/60 Hz	ES/IEC 60601-1, CAN/CSA C22.2 60601-1 (Safety Standards), EN/IEC 60601-1-2 (EMC Standards)
Track: 153656-001	115/230VAC, 1.2/0.6 A, 50/60 Hz	UL 60601-1:2003 R6.03 CAN/CSA C22.2 No. 601.1-M90

Midmark Dental Light Model: (LED)		
Mounting Style	Electrical Rating:	Complies to the applicable requirements of:
Elevance and Ultra Chairs: 153901-(001 thru -008) 153915 (-002, -004)	See power supply output (153941 ...)	Class I, Type B Applied Part, Ordinary Equipment (IPX0), Continuous Operation ES/IEC/EN 60601-1, CAN/CSA C22.2 60601-1 (Safety Standards), EN/IEC 60601-1-2 (EMC Standards)
Knight Chair: 153938 (-001 thru -004)	Power Supply Model No.:153935-001, -002	Class I, Type B Applied Part, Ordinary Equipment (IPX0), Continuous Operation
Ceiling: 153905-002 Wall 153908-002 Cabinet: 153909-002 Universal: 153904-002	115/230VAC, 190mA/330mA, 50/60HZ	ES/IEC/EN 60601-1, CAN/CSA C22.2 60601-1 (Safety Standards), EN/IEC 60601-1-2 (EMC Standards)
Track Light/Monitor: 153906 (-010 thru -018) / 153964 (-010 thru -018)	115VAC, 1.3 A, 50/60 Hz	ES/IEC 60601-1, CAN/CSA C22.2 60601-1 (Safety Standards), EN/IEC 60601-1-2 (EMC Standards)
Track: 153910-002	115/230VAC, 300 mA, 50/60 Hz	UL 60601-1:2003 R6.03 CAN/CSA C22.2 No. 601.1-M90

Dental Light Warnings

WARNING

No modification of this equipment is allowed.

WARNING

Should damage or loose parts be identified, the equipment should be removed from use and serviced. Failure to do so may result in injury.

WARNING

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of the Dental Light, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen, air, or nitrous oxide.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

Caution

To completely isolate the light or light/monitor devices from electrical mains supply:

For lights provided with a power cord, the power cord must be unplugged.

For lights or light/monitor devices wired directly to mains supply, the wiring must be connected to an appropriately sized circuit breaker and the circuit breaker must be turned off.

EMC - Manufacturer's Declaration and Guidance

Guidance and manufacturer's declaration – electromagnetic emissions

The Midmark Dental Light is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Dental Light should assure that it is used in such an environment.

<i>Emissions Test</i>	<i>Compliance</i>	<i>Electromagnetic environment – guidance</i>
RF Emissions CISPR 11	Group 1	The Midmark Dental Light uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The Midmark Dental Light is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communications equipment and the Midmark Dental Light

The Midmark Dental Light is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Midmark Dental Light can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the Midmark Dental Light as recommended below, according to the maximum output power of the communications equipment.

<i>Radiated maximum output power of transmitter W</i>	<i>Separation distance according to frequency of transmitter m</i>		
	150 kHz to 80 MHz $d = [3.5/V1] \sqrt{P}$	80 MHz to 800 MHz $d = [3.5/E1] \sqrt{P}$	800 MHz to 2.5 GHz $d = [7/E1] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity

The Midmark Dental Light is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Dental Light should assure that it is used in such an environment.

<i>Immunity Test</i>	<i>IEC 60601 Test Level</i>	<i>Compliance Level</i>	<i>Electromagnetic environmental - guidance</i>
			Portable and mobile RF communications equipment should be used no closer to any part of the Midmark Dental Light including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3V rms 150kHz to 80MHz	3V	$d = [3.5/\sqrt{P}] \sqrt{P}$
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	$d = [3.5/E1] \sqrt{P}$ 80 MHz to 800 MHz $d = [7/E1] \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation in meters (m). Field strength from fixed RF transmitters, as determined by the electromagnetic site survey, ^a Should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Midmark Dental Light is used exceeds the applicable RF compliance level above, the Midmark Dental Light should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Midmark Ultra Chair.

^b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity

The Midmark Dental Light is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Dental Light should assure that it is used in such an environment.

<i>Immunity Test</i>	<i>IEC 60601 Test Level</i>	<i>Compliance Level</i>	<i>Electromagnetic environmental - guidance</i>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 2, 4, 6 kV contact ± 2, 4, 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast Transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line(s) to line(s) ±2kV lines(s) to earth	± 0.5, 1 kV ± 0.5, 1, 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T for 0.5 cycle) 40% U_T (60% dip in U_T for 5 cycles) 70% U_T (30% dip in U_T for 25 cycles) <5% U_T (>95% dip in U_T for 5 s)	0 % UT (100 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Midmark Dental Light requires continued operation during power mains interruptions, it is recommended that the Midmark Dental Light be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Warranty Information

Midmark Limited Warranty - Dental Products

SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will at Midmark's option repair or replace components of the dental products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this limited warranty is limited to the repair or replacement of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the first retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows: Effective March 1, 2018 these applicable warranty periods are measured from the date of invoice to the original user, shall be as follows:

1. OPERATORY PRODUCTS

- a. Five (5) years for all products (except for the items in (b) through (e)).
- b. Two (2) years for upholstery (chairs and stools).
- c. "KINK-VALVE" module carries a ten (10) year warranty.
- d. The original light bulb on a new light carries a one (1) year warranty.
- e. Accessories not manufactured by Midmark are excluded including but not limited to Bien Air handpiece systems, Dentsply Cavitron scaler, Satelec scaler and curing light, and Sopro cameras.

2. ORAL SURGERY PRODUCTS are warranted for a period of one (1) year.

3. STERILIZER PRODUCTS are warranted for a period of one (1) year.

4. ULTRASONIC CLEANER PRODUCTS are warranted for a period of two (2) years.

5. AIR AND VACUUM PRODUCTS

- a. PowerAir® oil-less compressors – Five (5) years or 3,500 hours of use, whichever occurs first.
- b. PowerVac® and PowerVac® G dry vacuums – Five (5) years or 10,000 hours of use, whichever occurs first (except that the vacuum pump warranty term is ten (10) years or 20,000 hours of use, whichever occurs first).
- c. Classic Series® wet-ring vacuums – Five (5) years or 10,000 hours of use, whichever occurs first.
- d. PowerMax surgical suction – Two (2) years.
- e. Hg5 Series Amalgam Separator - One (1) year. (f) Midmark manufactured accessories - One (1) year.

6. SYNTHESIS™ DENTAL CASEWORK AND ARTIZAN® EXPRESSIONS PRODUCT

- a. Five (5) years for all products and components including door and drawer fronts, casters and slides, except for the items in (b), (c) and (d).
- b. Three (3) years for electrical components such as task lights/LED lights, cords, controls and accessories.
- c. Two (2) years for sliding track monitor mount and components and upholstery. (d) One (1) year for countertops and resin, including accessories.

7. IMAGING PRODUCTS are warranted for a period of two (2) years except for the ClearVision CR reader which is warranted for a period of one (1) year.
8. MIDMARK Replacement Parts and Accessories carry a ninety (90) day warranty

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for the following;

1. defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
2. products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environment conditions and electrical power requirements;
3. products considered to be of a consumable or sterile nature;
4. accessories or parts not manufactured by Midmark;
5. charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark
6. costs and expenses of routine maintenance and cleaning;
7. representations and warranties made by any person or entity other than Midmark;
8. matching of color, grain or texture except to commercially acceptable standards;
9. changes in color caused by natural or artificial light;
10. custom manufactured products;
11. alterations or modifications to the product by any person or entity other than Midmark; and
12. Products that would otherwise be covered under Sections 1 and 2 of this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical, animal health or dental market, as the case may be, in which purchaser intends to use the product.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

WARRANTY DISCLAIMER

THIS LIMITED WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, or implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

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