

# Elevance<sup>®</sup> Standard Delivery System

[Includes: Delivery Unit and Assistant's Unit]

[English](#)  
[Español](#)  
[Français](#)

|                                                           |   |
|-----------------------------------------------------------|---|
| <b>Compliance Information and Symbols</b>                 |   |
| <a href="#">Indications for Use, Calling Service etc.</a> | 2 |
| <a href="#">Symbols</a>                                   | 3 |
| <b>Product Features</b>                                   |   |
| <a href="#">Standard Water Bottler System</a>             | 4 |
| <a href="#">Optional Water Selector Switch System</a>     | 4 |
| <a href="#">Master Switch and USB Port</a>                | 5 |
| <a href="#">Work Surface Weight Limit</a>                 | 5 |
| <a href="#">Manual Shut-Offs - Water &amp; Air Supply</a> | 6 |
| <a href="#">Pressure Regulator Valves</a>                 | 6 |
| <a href="#">Flex Arm (Brake)</a>                          | 7 |
| <a href="#">Quick-Connects/ HVE/Saliva Ejectors</a>       | 8 |

**Integration Control Module<sup>®</sup> (ICM)**  
[Link to Software Specific ICM Screens](#)

|                                                                                         |    |
|-----------------------------------------------------------------------------------------|----|
| <b>Operation and Adjustments</b>                                                        |    |
| <a href="#">Operating Handpieces</a>                                                    | 9  |
| <a href="#">Syringe Adjustment Knobs</a>                                                | 10 |
| <a href="#">Coolant Air Control Knob</a>                                                | 10 |
| <a href="#">Bypass Toggles</a>                                                          | 10 |
| <a href="#">Syringe Operation</a>                                                       | 11 |
| <a href="#">Adjustment Knobs - Assistant's Unit</a>                                     | 11 |
| <b>Cleaning and Maintenance</b>                                                         |    |
| <a href="#">Barriers</a>                                                                | 12 |
| <a href="#">Cleaning, Disinfecting and Maintenance</a>                                  | 12 |
| <a href="#">Handpiece Accessories</a>                                                   | 12 |
| <a href="#">General Cleaning / Visual Inspection</a>                                    | 12 |
| <a href="#">Waterline Maintenance</a>                                                   | 12 |
| <a href="#">Cleaning &amp; Disinfecting Assistance</a>                                  | 12 |
| <a href="#">Beginning of Each Day...</a>                                                | 13 |
| <a href="#">Flushing System- Between Patients</a>                                       | 14 |
| <a href="#">Flushing System- End of Day</a>                                             | 15 |
| <a href="#">Air/Oil Separator - Clean / Maintain</a>                                    | 16 |
| <a href="#">Maintain / Replace Regulator Filters</a>                                    | 16 |
| <a href="#">Specs / Replacement Parts</a>                                               | 17 |
| <a href="#">EMC - Manufacturer's Declaration and Guidance</a>                           | 18 |
| <a href="#">Warranty Information</a>                                                    | 21 |
| For installation information go to:<br><a href="#">Technical Library at Midmark.com</a> |    |

# **Compliance Information and Symbols**

## ***Indications for Use***

Midmark instrument delivery systems are intended to provide dental professionals with air, water, and suction along with low-voltage electricity to operate dental handpieces, syringes, and accessories during dental examinations and treatments.

## ***Electromagnetic Interference***

Midmark dental operatory components are designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this operatory, remove the interfering device from room and / or plug product into an isolated circuit.

## ***Transportation / Storage / Operating Conditions***

Transportation / Storage Temperature:.....23°F to 100°F (-5°C to 38°C)

Relative Humidity:..... 10% to 90% (non-condensing)

Atmospheric Pressure: ..... 7.2 PSI to 15.3 PSI (50 kPa to 106 kPa)

Operating Temperature Range: .....59°F to 95°F (15°C to 35°C)

## ***Disposal of Equipment***

At the end of this product's life, the unit, accessories and other consumable goods may be contaminated from normal use. Consult local codes and ordinances for proper disposal of this equipment and other consumable goods.

## ***Calling For Service***

Direct all service inquiries to your authorized Midmark dealer. When calling for service, you must provide the following information:

Model / serial number  
Date of purchase  
Symptom(s) of malfunction

### ***WARNING***

*Do not modify this equipment without authorization of the manufacturer.*

### ***Caution***

*Federal law restricts this device to sale by or on the order of a licensed dental practitioner.*

# Symbols

These symbols may appear on your equipment and/or in the manuals. Warnings and cautions are provided in the manuals where applicable.

## WARNING

*Indicates a potentially hazardous situation which could result in serious injury if not avoided.*

## Caution

*Indicates a potentially hazardous situation which may result in minor or moderate injury if not avoided. It may also be used to alert against unsafe practices*

## Equipment Alert

*Indicates a potentially hazardous situation which could result in equipment damage if not avoided.*

## Note

*Amplifies a procedure, practice, or condition.*

Type B,  
Applied Part

Type BF,  
Applied Part

Protective  
Earth Ground

Proper Shipping  
Orientation

Fragile

USB

Catalogue Number


Serial Number


Manufacturer

Refer to instruction  
manual/booklet

Pressure  
Limit

Temperature  
Limit

Humidity  
Limit

Maximum stacking  
height (Refer to “n”  
number on package.)

Keep Dry

# Product Features


## Standard Water Bottle System

The standard bottle water system features an air supply switch to shut off air supplied to the bottle system.

### Equipment Alert

The Air Supply Switch and the Water Selector Switch are mounted in the same location. Operation of the switch depends on which type of water system your delivery has.

### Before filling bottle on a water selector switch system...

- A) Turn master switch OFF.
- B) Press and hold syringe air button until all residual air pressure is released (you can hear it stop hissing).


## Optional Water Selector Switch System

The optional water selector switch system uses a toggle switch to select either city  or bottle  water. Before removing the water bottle from this system, turn the master switch OFF and hold down the syringe air button to release residual air in the system.

### Note

Use water with a hardness level <60 ppm of calcium carbonate.

DA3653i

### To fill the Water Bottle...


- A) Verify system is depressurized.
- B) Unscrew water bottle from the quick disconnect mounting cap.
- C) Fill bottle with water.
- D) Screw water bottle into quick disconnect mounting cap.

Note: You will hear a hissing noise when air is supplied to the system if the bottle is improperly installed.

## Master Switch and USB Port

The master switch must be ON to operate the delivery system.

**To turn the delivery system ON or OFF...**  
Move the master switch to the desired position.

### NOTE

Turning the system OFF may affect recorded data (e.g. pictures recorded using a camera). Be sure to retain any needed data before turning system OFF.

### WARNING

**Risk of Electrical Shock Hazard:**  
Accessory equipment connected to the USB port must be certified according to IEC 60601-1 safety standards for electromedical equipment.


## Work Surface Weight Limit

### Caution

Exceeding the work surface weight limit may result in personal injury or damage the equipment.

Maximum load on flex arm mounted units is 10 lb (4.5 kg)


## Manual Shut-Off Valves for Water and Air Supply

Use manual shut-off valves to stop the facility supplied air and water.

### Recommended Settings:

|                                      |        |
|--------------------------------------|--------|
| City Water Regulator Gauge Setting   | 30 psi |
| Air Pressure Regulator Gauge Setting | 80 psi |


### Equipment Alert

Delivery components were designed to operate at the recommended settings. Poor performance or damage to equipment may result if recommended settings are not maintained.

## Pressure Regulator Valves

Use pressure regulator valves to control the air and water pressure supplied to delivery instruments.

**Shut-off Valves**  
Rotate clockwise 1/4 turn to Shut off.


### Pressure regulator adjustment...

- A) Pull up knob and turn to adjust.
- B) Watch regulator gauge as you turn knob to achieve desired setting
- C) Push knob DOWN to lock.

Turn Clockwise  
for more pressure

Turn Counter-Clockwise  
for less pressure


DA3027I

## ***Flex Arm (Brake)***

The flex arm contains a brake mechanism that prevents the unit from lowering until you release it by touching and holding a brake release touch sensor.

### ***To position Elevance delivery unit...***


- A) Press and hold any brake release touch sensor.*
- B) Move to desired position.*

## Quick-Connect Ports

Quick-Connect Ports (air [yellow] and water [blue]) allow you to quickly connect accessories to the delivery system.


Each Quick-Connect Port features a hinged clip that secures the connecting line to a raised ring found on each port opening.

## HVE and Saliva Ejectors

### To use HVE and Saliva Ejectors...

- A) Verify that facility vacuum system is ON.
- B) Move lever to towards the OPEN position to activate and moderate suction for the HVE / Saliva Ejector.

# Operation and Adjustments

## Operating Handpieces

Handpieces deliver air and water (if water selector switch is ON) when the footpedal is pressed. The water selector switch on the foot control toggles the water selection ON or OFF.


Water can also be turned ON or OFF using a water button, located on the delivery unit and on the assistant's unit control pad. This water button **overrides** the water selector switch on the foot control. The most convenient setting is to have the foot control selector switch in the water OFF position and then use one of the water buttons on the control pads to turn the water OFF and ON. The assistant's unit has a blue dot that illuminates when the water is ON. The bottom bar of the display screen informs: *Water: On* or *Water: Normal* (meaning the water override is OFF and the **selector switch is the active control**).

### **Handpiece Operation.**



- A) Turn master switch ON.
- B) Move water selector switch to desired setting.
- C) Pull handpiece from holder.
- D) Press footpedal to deliver drive air / water to the handpiece.

### **Note**

Water control is also impacted by the control software. Always read the display screens and indicators to be aware of the active water setting.

*Illuminates when  
Water Button is ON*

## Syringe Adjustment Knobs

Water and Air is manually controlled by the buttons on the syringe. The volume of water and air supplied to those buttons is manually controlled by the syringe air and syringe water knobs shown here.

**Adjust the Syringe Air**  **and Water**   
Rotate syringe knobs to modify the volume of air and water supplied to the syringe. Operation of the syringe is manual. See **Syringe Operation**.

## Coolant Air Control Knob

The Coolant Air control knob controls the volume of coolant air supplied to all handpieces.

**Adjust the Coolant Air**   
Rotate the coolant air knob until the desired coolant air flow to the handpieces is achieved.

## Bypass Toggles for Handpiece Air and Water


Bypass controls provide a method to manually supply air and water to all the air-driven handpieces if the ICM loses power. Manual control is available when the water and drive air bypass toggles are flipped into the BYPASS positions.

**Adjust the Bypass Coolant Water**   
Rotate the coolant water adjustment knob until the desired coolant water flow to the handpieces is achieved.

*Note:* When water toggle in BYPASS position adjustment knob adjusts coolant water for all handpieces.



**Adjust the Bypass Drive Air**   
Rotate the drive air adjustment knob until the desired drive air flow to the handpieces is achieved.

*Note:* When drive air toggle in BYPASS position, the adjustment knob adjusts drive air for all handpieces.

## Syringe Operation

This page shows operation and adjustment for a standard three-way syringe.


**To adjust the water stream pressure and water spray pattern...**  
Rotate the syringe water adjustment knob as desired.

**To adjust the air spray pressure and water spray pattern...**  
Rotate the syringe air adjustment knob as desired.

## Adjustment Knobs - Assistant's Units

Elevance assistant's units have screw-type control knobs under the assistant's head to adjust air and water flow to the syringe.

# **Cleaning, Disinfecting and Maintenance**

## **Barriers**


Midmark recommends the use of disposable barriers on all clinician controls that may be in contact with clinician hands and fingers during dental procedures. The use of barriers significantly reduces the need for chemical cleaners, thus prolonging the life of the equipment.

Only use barrier material that is intended for use with dental equipment. Midmark recommends the use of an FDA market-cleared barrier such as Pinnacle Cover-all™. Follow barrier manufacturer instructions for proper use of these products.

## **Cleaning and Disinfecting**

In addition to the use of barriers Midmark recommends the use of an EPA registered and FDA market-cleared cleaner/disinfectant such as Cavicide™ to be used on all clinician controls or surfaces that may come in contact with dental instruments during dental procedures.

Follow cleaner/disinfectant manufacturer instructions for proper use of the product. Care should be taken to avoid excessive application and pooling of liquids.

## **Handpiece Accessories**

Only use dental handpiece accessories with the delivery system that are FDA market-cleared and refer to manufacturer's instructions for proper cleaning and disinfecting. Either an autoclavable syringe tip or a single use disposable syringe tip may be used.

### **Equipment Alert**

#### **AUTOCCLAVABLE SYRINGE TIP STERILIZATION**

*The autoclavable syringe tips supplied with the delivery system must be sterilized prior to use with each patient, including initial use. Be sure to thoroughly rinse and clean syringe tips prior to sterilization, any debris may reduce the effectiveness of the sterilization. Recommended sterilization process is steam autoclave. Recommended parameters are 125°C (250°F) and 106 kPa (15 PSI) for 40 minutes at temperature and pressure.*

## **General Purpose Cleaning**

For general purpose cleaning, use cleaners that are appropriate for the situation, such as warm water and mild detergents, or a 10% solution of bleach with water.

## **Visual Inspection**

After cleaning, visually inspect the product for deterioration of covers and touch pads. Do not use the delivery system if excessive discoloration, cracking, or other signs of wear are noticeable (see Calling for Service instructions).

## **Waterline Maintenance**

Waterline maintenance is necessary to keep the count of heterotrophic bacteria from rising higher than desired levels. The desired level for a specific location should be determined by any local or regional guidelines. For example, The United States Centers for Disease Control and Prevention (CDC) guideline for heterotrophic bacteria is less than or equal to 500 CFU/mL (colony forming units per milliliter). Midmark recommends keeping this level under 200 CFU/mL.

Treatment can come in many forms. The most popular methods on the market currently are tablets and straw/cartridge based systems. Midmark recommends the use of a straw/cartridge based system that keeps the bacteria levels in check.

Regular monitoring should also take place to ensure that heterotrophic bacteria is not exceeding the desired limit. If the level is higher than desired, a shock treatment of the waterlines should be performed. When performing a shock treatment, be sure to check with the manufacturer of the regular treatment regimen being used to ensure chemical compatibility. Monitoring frequency should be established by your practice. As a suggestion, Midmark would recommend that you begin by monitoring on a monthly basis, and make adjustments to the frequency based on test results.

Per the CDC, routine flushing of the waterlines should be performed between every patient. Extra flushing maybe needed within Midmark equipment when tablets are used. Undissolved tablet particles can gather over time in places within the waterlines, obstructing the line and causing water flow to slow. By flushing the waterlines, water flow is maximized and should push any undissolved particles through.

## Cleaning and Disinfecting Assistance

For assistance with cleaning and disinfecting instructions contact the Midmark Technical Service Department at 1-800-Midmark; it is helpful to provide the delivery system model number and serial number when asking for assistance.

Additional information is available from the organizations listed below:

|                                                                                                                   |                                                                                                                                                          |
|-------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Organization for Safety &amp; Asepsis Procedures:</b><br><a href="http://www.osap.org">http://www.osap.org</a> | <b>Dept. of Health &amp; Human Resources, Centers for Disease Control &amp; Prevention (CDC):</b><br><a href="http://www.cdc.gov">http://www.cdc.gov</a> |
| <b>American Dental Association:</b><br><a href="http://www.ada.org">http://www.ada.org</a>                        | <b>European Dental Association:</b><br><a href="http://www.eda-eu.org">http://www.eda-eu.org</a>                                                         |




## Cleaning and Disinfecting the Delivery System - At the Beginning of Each Day...

### **At the beginning of each day...**

- A) Thoroughly clean and visually inspect the delivery system (See General Purpose Cleaning and Visible Inspection).
- B) Fill the water bottle with fresh water and perform a flush operation (Refer to ICM Screens guide for more details.)

*Note: Use water with a hardness level <60 ppm of calcium carbonate.*



- C) Install barriers on all delivery system clinician controls.
- D) Install syringe tips (sterilized or disposable).

## Flushing the Delivery System - Between Patients...

### **Between patients...**





- A) Remove all disposable barriers, syringe tips, and other dental instruments.
- B) Clean and disinfect the delivery system (see *Cleaning and Disinfecting instructions*).
- C) Perform a flush operation (Refer to *ICM Screens guide* for more details.)
- D) Install barriers on all delivery system clinician controls.
- E) Install syringe tips (sterilized or disposable).

## Flushing the Delivery System - At the End of Each Day...

### At the end of each day...

- A) Remove all disposable barriers, syringe tips, and other dental instruments.
- B) Clean and disinfect the delivery system (see *Cleaning and Disinfecting instructions*).
- C) Fill water bottle with fresh water and perform a flush operation (Refer to *ICM Screens guide* for more details.)
- D) Turn master switch OFF. Press and hold the footpedal until all pressure is released.


DA3663I

## Air/Oil Separator - Cleaning and Maintaining

### Periodically check the air/oil separator...

- A) Turn master switch OFF.
- B) Remove container by squeezing sides in while pulling out.
- C) Remove saturated gauze.
- D) Clean container and mounting cap.
- E) Install clean gauze and reinstall the container.


## Maintaining /Replacing Regulator Filters (in J-box)

### Replacing Regulator Filters...







- A) Turn off water and air supplies using manual shut off valves.
- B) Unscrew filter cap. (Use 9/16" wrench.)
- C) Unscrew retainer nut and remove filter and o-ring.
- D) Install new filter and o-ring. Secure with retainer nut.

Note: Install the filter and retainer nut with the ridged side up (as shown).

- E) Reinstall filter cap.

| <b>Specifications</b>                                                                                                         |                                                                                                                               |
|-------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Electrical Rating for<br>Elevance Chair:                                                                                      | See power supply output (153941 ...)                                                                                          |
| Complies to the applicable<br>requirements of:                                                                                | ES/IEC/EN 60601-1, CAN/CSA C22.2<br>No. 60101-1 (Safety Standards), EN/IEC<br>60601-1-2 (EMC Standards)                       |
| Air Supply Requirements:                                                                                                      | Maximum 80-100 PSI (5.52 - 6.9 bar)                                                                                           |
| Water Supply<br>Requirements:                                                                                                 | Maximum 30 - 50 PSI (2.1 -3.4 bar)<br>Use water with a hardness level <60 ppm<br>of calcium carbonate.                        |
| Certifications:                                                                                                               | Midmark Corporation.....ISO-9001                                                                                              |
| Classifications:                                                                                                              | Class I, Type B Applied Part, except as<br>noted for optional accessories, Ordinary<br>Equipment (IPX0), Continuous Operation |
| Optional Accessories:                                                                                                         |                                                                                                                               |
| Acteon Satelec Scaler                                                                                                         | Type BF Applied Part                                                                                                          |
| Acteon Satelec Curing<br>Light                                                                                                | Type B Applied Part                                                                                                           |
| Acteon Sopro Intra-oral<br>Camera                                                                                             | Type BF Applied Part                                                                                                          |
| Bien Air Electric Motor                                                                                                       | Type BF Applied Part                                                                                                          |
| Dentsply Cavitron Scaler                                                                                                      | Type B Applied Part                                                                                                           |
| Handpiece tubing and<br>connectors intended to<br>be used with ISO 7785-1<br>or ISO 7785-2 compliant<br>air-driven handpieces | Type B Applied Part                                                                                                           |
| Sota Claris Intra-oral<br>Camera                                                                                              | Type B Applied Part                                                                                                           |

| <b>Replacement Parts</b>                                                              |                                                                                     |                                                             |                                                                                                |
|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------|------------------------------------------------------------------------------------------------|
| <b>Saliva Ejector Valve</b><br>Part Number<br>029-3975-00                             |  | <b>HVE Handpiece Valve</b><br>Part Number<br>029-3974-00    |             |
| <b>Solids Collector Screen</b><br>Part Number<br>(1) 053-1839-00<br>(50) 002-10260-00 |  | <b>Syringe Tip</b><br>Part Number<br>120514                 |             |
| <b>Saliva Ejector Tip</b><br>Part Number<br>016-1116-00                               |  | <b>Regulator Valve Filter</b><br>Part Number<br>053-1166-00 |             |
| <b>Saliva Ejector Lever</b><br>Part Number<br>016-1113-00                             |  | <b>HVE Lever</b><br>Part Number<br>016-1109-00              | <br>DA2989I |

### **WARNING**

**Equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen, air, or nitrous oxide.**

**Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.**

## EMC - Manufacturer's Declaration and Guidance

### Guidance and manufacturer's declaration – electromagnetic emissions

The Midmark Elevance Unit dental device is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Elevance Unit dental device should assure that it is used in such an environment.

| <b>Emissions Test</b>                                | <b>Compliance</b> | <b>Electromagnetic environment – guidance</b>                                                                                                                                                                                                                                                                                                                                                                                                                            |
|------------------------------------------------------|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF Emissions CISPR 11                                | Group 1           | The Midmark Elevance Unit dental device uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.<br><br>The Midmark Elevance Unit dental device is suitable for use in all establishments, including domestic establishment and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF Emissions CISPR 11                                | Class A           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Harmonic emissions IEC 61000-3-2                     | Class A           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Complies          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |

### Recommended separation distances between portable and mobile RF communications equipment and the Midmark Elevance Unit dental device

The Midmark Elevance Unit dental device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Midmark Elevance Unit dental device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the Midmark Elevance Unit dental device as recommended below, according to the maximum output power of the communications equipment.

| <b>Radiated maximum output power of transmitter<br/>W</b> | <b>Separation distance according to frequency of transmitter<br/>m</b> |                                                |                                                 |
|-----------------------------------------------------------|------------------------------------------------------------------------|------------------------------------------------|-------------------------------------------------|
|                                                           | 150 kHz to 80 MHz<br>$d = 1.2 \times \sqrt{P}$                         | 80 MHz to 800 MHz<br>$d = 1.2 \times \sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2.3 \times \sqrt{P}$ |
| 0.01                                                      | 0.12                                                                   | 0.12                                           | 0.23                                            |
| 0.1                                                       | 0.37                                                                   | 0.37                                           | 0.74                                            |
| 1                                                         | 1.17                                                                   | 1.17                                           | 2.34                                            |
| 10                                                        | 3.69                                                                   | 3.69                                           | 7.38                                            |
| 100                                                       | 11.67                                                                  | 11.67                                          | 23.34                                           |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Guidance and manufacturer's declaration – electromagnetic immunity

The Midmark Elevance Unit dental device is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Elevance Unit dental device should assure that it is used in such an environment.

| <i>Immunity Test</i>                                                                                   | <i>IEC 60601 Test Level</i>                                                                                                                                                                             | <i>Compliance Level</i>                                                                                   | <i>Electromagnetic environmental - guidance</i>                                                                                                                                                                                                                                                                                        |
|--------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                         | ± 6 kV contact<br>± 8 kV air                                                                                                                                                                            | ± 6 kV contact<br>± 8 kV air                                                                              | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                                                                                          |
| Electrical fast Transient / burst<br>IEC 61000-4-4                                                     | ± 2 kV for power supply lines<br><br>± 1 kV for input/output lines                                                                                                                                      | ± 2 kV for AC and DC power lines<br><br>I/O lines not tested, all less than 3 meters                      | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                                    |
| Surge<br>IEC 61000-4-5                                                                                 | ±1kV line(s) to line(s)<br>±2kV lines(s) to earth                                                                                                                                                       | ±1kV line(s) to line(s)<br>±2kV lines(s) to earth                                                         | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                                    |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | <5% $U_T$<br>(>95% dip in $U_T$ for 0.5 cycle)<br><br>40% $U_T$<br>(60% dip in $U_T$ for 5 cycles)<br><br>70% $U_T$<br>(30% dip in $U_T$ for 25 cycles)<br><br><5% $U_T$<br>(>95% dip in $U_T$ for 5 s) | V Dip >30% of $U_T$ for 10ms<br><br>V Dip < 60% of $U_T$ for 100ms<br><br>V Dip > 95% of $U_T$ for 5000ms | Mains power quality should be that of a typical commercial or hospital environment. If the user of the Midmark Elevance Unit dental device requires continued operation during power mains interruptions, it is recommended that the Midmark Elevance Unit dental device be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8                                             | 3 A/m                                                                                                                                                                                                   | 3 A/m                                                                                                     | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                                                              |

**NOTE:**  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration – electromagnetic immunity

| <i>Immunity Test</i>          | <i>IEC 60601 Test Level</i> | <i>Compliance Level</i> | <i>Electromagnetic environmental - guidance</i>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-------------------------------|-----------------------------|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                               |                             |                         | <p>Portable and mobile RF communications equipment should be used no closer to any part of the Midmark Elevance Unit dental device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p>                                                                                                                                                                                                                                   |
| Conducted RF<br>IEC 61000-4-6 | 3V rms<br>150kHz to 80MHz   | 3V                      | $d = 1.2 \times \sqrt{P}$ 80 MHz to 800 MHz<br>$d = 2.3 \times \sqrt{P}$ 800 MHz to 2.5 GHz                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Radiated RF<br>IEC 61000-4-3  | 3V/m<br>80MHz to 2.5GHz     | 3V/m                    | <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation in meters (m).</p> <p>Field strength from fixed RF transmitters, as determined by the electromagnetic site survey, <sup>a</sup> Should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location is which the Midmark Elevance Unit dental device is used exceeds the applicable RF compliance level above, the Midmark Elevance Unit dental device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Midmark Elevance Unit dental device.

<sup>b</sup> Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.

# Warranty Information

## Midmark Limited Warranty - Dental Products

### SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will at Midmark's option repair or replace components of the dental products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this limited warranty is limited to the repair or replacement of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the first retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

### APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows: Effective March 1, 2018 these applicable warranty periods are measured from the date of invoice to the original user, shall be as follows:

#### 1. OPERATORY PRODUCTS

- a. Five (5) years for all products (except for the items in (b) through (e)).
- b. Two (2) years for upholstery (chairs and stools).
- c. "KINK-VALVE" module carries a ten (10) year warranty.
- d. The original light bulb on a new light carries a one (1) year warranty.
- e. Accessories not manufactured by Midmark are excluded including but not limited to Bien Air handpiece systems, Dentsply Cavitron scaler, Satelec scaler and curing light, and Sopro cameras.

#### 2. ORAL SURGERY PRODUCTS are warranted for a period of one (1) year.

#### 3. STERILIZER PRODUCTS are warranted for a period of one (1) year.

#### 4. ULTRASONIC CLEANER PRODUCTS are warranted for a period of two (2) years.

#### 5. AIR AND VACUUM PRODUCTS

- a. PowerAir® oil-less compressors – Five (5) years or 3,500 hours of use, whichever occurs first.
- b. PowerVac® and PowerVac® G dry vacuums – Five (5) years or 10,000 hours of use, whichever occurs first (except that the vacuum pump warranty term is ten (10) years or 20,000 hours of use, whichever occurs first).
- c. Classic Series® wet-ring vacuums – Five (5) years or 10,000 hours of use, whichever occurs first.
- d. PowerMax surgical suction – Two (2) years.
- e. Hg5 Series Amalgam Separator - One (1) year. (f) Midmark manufactured accessories - One (1) year.

#### 6. SYNTHESIS™ DENTAL CASEWORK AND ARTIZAN® EXPRESSIONS PRODUCT

- a. Five (5) years for all products and components including door and drawer fronts, casters and slides, except for the items in (b), (c) and (d).
- b. Three (3) years for electrical components such as task lights/LED lights, cords, controls and accessories.
- c. Two (2) years for sliding track monitor mount and components and upholstery. (d) One (1) year for countertops and resin, including accessories.

7. IMAGING PRODUCTS are warranted for a period of two (2) years except for the ClearVision CR reader which is warranted for a period of one (1) year.

8. MIDMARK Replacement Parts and Accessories carry a ninety (90) day warranty

### EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for the following:

1. defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
2. products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environment conditions and electrical power requirements;
3. products considered to be of a consumable or sterile nature;
4. accessories or parts not manufactured by Midmark;
5. charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark
6. costs and expenses of routine maintenance and cleaning;
7. representations and warranties made by any person or entity other than Midmark;
8. matching of color, grain or texture except to commercially acceptable standards;
9. changes in color caused by natural or artificial light;
10. custom manufactured products;
11. alterations or modifications to the product by any person or entity other than Midmark; and
12. Products that would otherwise be covered under Sections 1 and 2 of this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical, animal health or dental market, as the case may be, in which purchaser intends to use the product.

### EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

### WARRANTY DISCLAIMER

THIS LIMITED WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

### STATUE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, or implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

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