354 / 355 Mobile Light Installation

Applies to Models:

354 (-022)


355 (-022 / -024 / -037)

Special Tools:

none

Language of origin: English

Step 1: Inspect components.


MA9926I

Step 2: Install transformer housing.

- A) Position caster base insert on base.
- B) Secure transformer housing to base with spacer disk and nut*.


*Note: Tighten nut to 35 - 45 ft-lbs. (47.5 - 61 N•m)

Step 3: Install vertical tube.

- A) Remove tape securing wires, then feed wires thru bottom of vertical tube.
- B) Slide sleeve onto vertical tube.
- C) Connect wires from vertical tube to corresponding wires from transformer housing. (ex. brown to brown)
- D) Slide vertical tube onto mounting tube of transformer housing.
- E) Secure vertical tube with four #10-24 x 3/8" hex screws*.

Note: Use 1/8" hex wrench (supplied).

15A22 271


Equipment Alert

Align three screw holes on the vertical tube with holes in ball pivot housing of the counterbalance assembly before installing any screws. Failure to comply may result in damage to screws.

Step 4: Install counterbalance assembly.

- A) Align three screw holes, then insert counterbalance into vertical tube.
- B) Install two #10-24 x 3/8" hex screws* into the holes marked "A".
- C) Pull two wires from counterbalance thru bottom access hole.
- D) Connect wires from counterbalance to corresponding wires from vertical tube, then push wires inside tube.
- E) Slide sleeve into position, then secure with #10-24 x 3/8" phillips screw.

Note: Use 1/8" hex wrench (supplied)

Install Bulb: 355 only!

WARNING

The maximum bulb wattage for this light is 100 Watts.
There is a risk of fire if this limit is exceeded.

Use Midmark part number: 002-0640-00.

Caution

- To prevent electrical shock, the light must be in **STANDBY** mode (⏻).
- Allow the light to cool before removing bulb.
- Move the lighthead away from the exam site to prevent contaminants from falling into the area.


Equipment Alert

Halogen bulbs are sensitive to body oils. Do not touch the glass portion of the bulb.
If the glass portion of the bulb is handled, wipe with a clean soft lint-free cloth.
Wipe with alcohol, then pat dry.

To install bulb...

- Turn light switch to **STANDBY** (⏻).
- Loosen** three screws on top of lighthead.
(Do not attempt to remove these screws).
- Separate top cap from lighthead housing.
- Pull bulb out of bulb socket, then discard.
- Grasp new bulb with a clean cloth, then insert into bulb socket.
- Align three top cap screws with holes, then tighten screws.

Electrical Input Setting: 230 V models only!

Equipment Alert


The voltage setting that appears in the window of the fuse drawer, must match the facility voltage. Determine facility voltage before adjusting Electrical Input Setting.

Equipment Alert

Locking tabs are fragile.
Pull gently to prevent damage.


To adjust Electrical Input Setting...

- A) Turn light switch to STANDBY (⏻).
- B) Disconnect power cord from light.
- C) Using a flat-bladed screwdriver, push in on two locking tabs, then pull fuse drawer out.
- D) Pull fuses out of fuse drawer.
- E) Pull out on two selector lock tabs, then pull voltage selector out of fuse drawer.
- F) Rotate voltage selector until the correct setting aligns with the voltage window.
- G) Insert voltage selector into fuse drawer.

Electrical Information

Step 5: Connect power cord.

Electrical Ratings:	
120 VAC models	354: 120 VAC, 60 Hz, 1.5 amps, single phase 355: 120 VAC, 60 Hz, 1.25 amps, single phase
230 VAC models	355: 220 / 230 / 240 VAC, 50/60 Hz, 0.7 amps, single phase
Fuse Ratings:	
120 VAC models	354 / 355: 1.6A, 250V, 5mm x 20mm, Type T, low breaking capacity, IEC 127-2/3
230 VAC models	355: 0.8A, 250 VAC, 5mm x 20mm, Type FST, Time Lag, IEC 127-2/3
Grounding reliability can only be achieved if this equipment is connected to an equivalent receptacle marked Hospital Grade (NEMA 5-15R, HOSPITAL GRADE)	

Step 6: Check for proper operation.

A) Turn ON / STANDBY switch to ON "I".


(Light should illuminate).

(If light does not illuminate, check wire connections).

B) Rotate the lighthead, lamp tube, and cross tube thru range of motion. Release the lighthead at different positions.


(Light should move freely, and should not drift when released).

(If drifting occurs, perform appropriate tension adjustment procedure).

To adjust ball pivot tension...


- A) Remove screw, then slide sleeve down.
- A) Insert screwdriver into adjustment hole.
- B) Move the screwdriver LEFT / RIGHT to adjust tension setting (LEFT: **decreases** tension / RIGHT: **increases** tension)
- C) Remove screwdriver, and check for desired tension.
Repeat until desired tension is achieved.
- D) Slide fixed arm sleeve up, then install screw.

To adjust cross tube counterbalance...

- A) Remove housing cap from counterbalance.
- B) Loosen / tighten the 3/16" adjustment screw until the cross tube balances in horizontal position.
- C) Install housing cap onto counterbalance.


Note: Turning the adjustment screw clockwise will lower the lighthead.
Turning the adjustment screw counterclockwise will raise the lighthead.

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