Elevance® Continental-Style Delivery Unit

[English](#)

[Español](#)

[Français](#)

For Model:
153800-003

DA40371

User Guide

Style P

© Midmark Corporation 2016

TP202 20-42-FO-00014 Rev A1 C2169 003-2631-99 Rev. HA3 (12/12/18)

Product Information

Dealer :

Date of Purchase:

Model / Serial Number:

**Midmark Authorized Service
Company:**


Table of Contents

Integration Control Module®

For information regarding software / interactive screens, refer to applicable software version of ICM Screens. www.midmark.com

Important Information

<i>Symbols</i>	4
<i>Product Registration</i>	5
<i>Indications for Use</i>	5
<i>Electromagnetic Interference</i>	5
<i>Transportation / Storage / Operating</i>	5
<i>Disposal of Equipment</i>	6
<i>Safety Instructions</i>	6

Product Features

<i>Master Switch</i>	7
<i>Water Control</i>	8
<i>Water Bottle Systems:</i>	
<i>Standard Water Bottle System</i>	9
<i>Water Selector System (optional)</i>	10
<i>Filling the Water Bottle (all systems)</i>	11
<i>Tray (Optional)</i>	12
<i>USB Port</i>	13
<i>Flex Arm Brake</i>	14
<i>Quick Connect Ports</i>	15
<i>Manual Shut-off Valves and Regulators</i>	16

Delivery Operation

<i>Operating the Syringe</i>	17
<i>Syringe Air / Water Adjustments</i>	18
<i>Normal Operation of Smart Handpieces</i>	19
<i>Adjusting the Volume of Coolant Air</i>	20
<i>Bypass Mode Operation</i>	21

Maintenance

Daily Maintenance:

<i>At the beginning of each day</i>	22
<i>Between patients</i>	23
<i>At the end of each day</i>	24

Periodic Maintenance:

<i>Air / Oil Separator</i>	25
<i>Replace Regulator Filters</i>	26

Cleaning / Infection Control

<i>Barriers</i>	27
<i>Cleaning and Disinfecting</i>	27
<i>General Purpose Cleaning</i>	27
<i>Visual Inspection</i>	27
<i>Waterline Maintenance</i>	28
<i>Handpiece Accessories</i>	28

Servicing Equipment

<i>Calling for Service</i>	29
<i>Replacement Parts</i>	29

Specifications / Compliance

<i>Specifications Chart</i>	30
-----------------------------------	----

EMC - Manufacturer's Declaration

<i>and Guidance</i>	31
---------------------------	----

Warranty Information

<i>Limited Warranty</i>	34
-------------------------------	----

Important Information - Safety Symbols

WARNING
Indicates a potentially hazardous situation which could result in serious injury.

Caution
Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices


Equipment Alert
Indicates a situation which could result in equipment damage.

Note
Amplifies a procedure, practice, or condition.

Symbol Glossary

These symbols may appear on your equipment and/or in the manuals.

	Manufacturer		Serial Number		Catalogue Number
	Type B, Applied Part		Keep Dry		Pressure Limit
	Type BF, Applied Part		USB		Temperature Limit
	Protective Earth Ground		Refer to instruction manual/booklet		Humidity Limit
	Proper Shipping Orientation		Fragile		Maximum stacking height (Refer to “n” number on package.)

Product Registration

To register your product online go to: www.midmark.com

Indications for Use

Midmark instrument delivery systems are intended to provide dental professionals with air, water, and suction along with low-voltage electricity to operate dental handpieces, syringes, and accessories during dental examinations and treatments.

Electromagnetic Interference

Midmark dental operatory components are designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this operatory, remove the interfering device from room and / or plug product into an isolated circuit.

Transportation / Storage / Operating Conditions

Transportation / Storage Temperature:..... 23°F to 100°F (-5°C to 38°C)
Relative Humidity:..... 10% to 90% (non-condensing)
Atmospheric Pressure: 7.2 PSI to 15.3 PSI (50 kPa to 106 kPa)
Operating Temperature Range:59°F to 95°F (15°C to 35°C)

Disposal of Equipment

At the end of this product's life, the unit, accessories and other consumable goods may be contaminated from normal use. Consult local codes and ordinances for proper disposal of this equipment and other consumable goods.

Safety Instructions

WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen, air, or nitrous oxide.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

WARNING

Do not modify this equipment without authorization of the manufacturer.

Caution

Federal law restricts this device to sale by or on the order of a licensed dental practitioner.


Product Features

Master Switch

The master switch must be ON to operate the delivery system. The master switch controls air and electric to the delivery system.

To turn the delivery system ON or OFF...
Move the master switch to the desired position.

Note: Turning the system OFF may affect recorded data (e.g. pictures recorded using a camera). Be sure to retain any needed data before turning system OFF.


Water Control

Water delivered to the handpieces can be turned ON/OFF using the two controls shown here.

To turn water ON using the selector switch on the foot control...

- A) Turn master switch ON.
- B) Move water selector switch on foot control to water ON position.


Note: ICM screen displays Water: Normal.

To turn water ON using water button on the ICM...

- A) Turn master switch ON.
- B) Move water selector switch on foot control to water OFF position.
- C) Toggle the water button on the ICM to display Water: On.

Note: Some assistant's units have a control pad water button that functions the same as the water button on the ICM.

To turn water OFF...

- A) Turn master switch ON.
- B) Move water selector switch on foot control to water OFF position.
- C) Toggle water button to display Water: Normal on the ICM screen.

Water Bottle Systems

The Elevance Continental-Style Delivery Unit may be ordered with either the Standard Water Bottle System, or the optional Water Selector System.

Standard Water Bottle System


The Standard Water Bottle System supplies water to the syringe and handpieces from the water bottle. Use water with a hardness level <60 ppm of calcium carbonate.

The water bottle must be pressurized to operate properly.


Standard Water Bottle System

To pressurize the water bottle...

Move the air supply switch to ON.

Air Supply Switch Location


DA3907i

Water Bottle Systems - continued

Water Selector System (optional)

The Water Selector System allows the operator to use either the water bottle, or city water in the delivery system.

Water Selector System
To select the water source...
Move the water selector switch to the desired setting.

Water Bottle Systems - continued

Filling the Water Bottle (all systems)

Equipment Alert

The water bottle must be depressurized before it is removed.

Standard Water Bottle System **Before filling the water bottle...**

- A) Turn air supply switch OFF.

Note: Refer to 'Standard Water Bottle System' for location.

Water Selector System **Before filling the water bottle...**



- A) Turn master switch OFF.

Note: Refer to 'Master Switch' section for location.

All Systems

To fill the water bottle...

- A) Press and hold syringe air button until all residual air pressure is released ('hissing' noise will stop).
B) Remove water bottle.
C) Fill bottle with water with a hardness level <60 ppm of calcium carbonate.
D) Install water bottle.

DA3910i


Tray (Optional)

Maximum load on the tray area should not exceed 10 pounds (4.5 kg).

Caution

Exceeding the recommended weight limit may result in personal injury or damage the equipment.

USB Port


The delivery unit has one standard USB port located on the syringe side of the unit.

WARNING

Risk of Electrical Shock Hazard:

Accessory equipment connected to the USB port must be certified according to IEC 60601-1 safety standards for electromedical equipment.


Flex Arm Brake

The flex arm brake prevents unwanted, downward movement of the unit.


To release the flex arm brake...

- A) *Touch and hold any brake release sensor.*
- B) *Move unit to desired position.*

Quick Connect Ports

Quick connect ports [air (yellow) and water (blue)] allow you to quickly connect accessories to the delivery system. Each quick connect port features a hinged clip that secures the connecting line to a raised ring found on each port opening.

Manual Shut-off Valves and Regulators

Use manual shut off valves to turn facility supplied air and water ON and OFF.

Use regulators to control the pressure of the air and water supplied to the handpieces.


Equipment Alert

Delivery components are designed to operate at the recommended settings. Poor performance or damage to equipment may result if recommended settings are not maintained.

Recommended Settings:	PSI	kPa
City Water Gauge	30	207
Air Pressure Gauge	80	552

To adjust regulator pressure...

- Pull up on regulator knob.
- Turn clockwise to increase pressure, counter-clockwise to decrease pressure.
- Push knob down to lock.


DA3912i

Delivery Operation

Operating the Syringe

To deliver air and water to syringe...

- A) Turn master switch ON.
- B) Turn the water ON
- C) Pull the syringe.
- D) Press the appropriate button(s) on the syringe to deliver the desired fluid(s) from the syringe tip.


DA3918i

Syringe Air / Water Adjustments

The volume of air and water delivered by the syringe can be set using the syringe control knobs located under the cover on the back of the delivery unit.

To adjust the syringe air / water volume setting...

- A) Rotate syringe control knobs to adjust the volume of air and water supplied to syringe.
- B) Press and hold a button on the syringe as you rotate the corresponding control knob to observe the effect of the adjustment.

DA3853i

Normal Operation of Smart Handpieces



Smart handpieces have a control screen that appears on the ICM when the handpiece is pulled. Users can adjust settings for the drive air and water for each smart handpiece using the control screen. Reference the ICM Screens guide for more detailed information about using the ICM.

Smart handpieces include Bien Air electric motors, air turbines and scalers.

To deliver drive air and water to handpieces...

- A) Turn master switch ON.
- B) Turn the water ON.
- C) Pull a handpiece and observe the ICM screen that appears.
- D) Press the footpedal (fully) to deliver the air and water at the settings shown on the ICM screen.

Note: Deliver less than the displayed air/water settings by pressing the footpedal partially.

Equipment Alert

Normal operation occurs when the toggle controls located under the cover on the back of the delivery unit are in Normal position.


Adjusting the Volume of Coolant Air

Coolant Air


Coolant air is delivered thru the handpiece to the work surface. Users can adjust the volume of coolant air delivered using the coolant air control knob.

Air turbines and electric motors are equipped with coolant air; scalers are not.

The coolant air control knob can be used in both Normal and Bypass modes of operation.


Coolant Air Control Knob

To adjust the coolant air to handpieces...

- A) Turn master switch ON.
- B) Pull a handpiece and hold your hand in front of it.
- C) Rotate the coolant air control knob as you press the footpedal to feel the effects of the adjustment.

Note The coolant air control knob regulates the coolant air for all handpieces equipped with coolant air.


DA3853i

Bypass Mode Operation

Bypass toggles provide a method to manually operate air driven handpieces if Normal operation becomes impossible due to an ICM power loss or malfunction.

Bypass mode operation is activated by flipping the air and water bypass toggles to the BYPASS positions. The bypass control knobs for drive air and water become active.


Maintenance

Daily Maintenance

At the beginning of each day...

- A) Thoroughly clean and visually inspect the delivery system*.
- B) Fill water bottle with water with a hardness level <60 ppm of calcium carbonate, then perform a flush operation.
(Refer to: ICM Screens guide for details).
- C) Install barriers on all delivery system clinician controls*.
- D) Install syringe tips (sterilized or disposable).

*Note: Refer to "Cleaning / Disinfecting" section in this document.


Daily Maintenance - continued

Between patients...

- A) Remove all disposable barriers, syringe tips, and other dental instruments.
- B) Clean and disinfect the delivery system*.
- C) Perform a flush operation. (Refer to: ICM Screens guide for details).
- D) Install barriers on all delivery system clinician controls.
- E) Install syringe tips (sterilized or disposable).

*Note: Refer to "Cleaning / Disinfecting" section in this document.


Daily Maintenance - continued

At the end of each day...

- A) Remove all disposable barriers, syringe tips, and all dental instruments.
- B) Clean and disinfect the delivery system*.
- C) Fill the water bottle with fresh water and perform a flush operation. (Refer to: ICM Screens guide for details).
- D) Turn master switch OFF. Press and hold the footpedal until all pressure is released.


*Note: Refer to "Cleaning / Disinfecting" section in this document.

Periodic Maintenance

Inspect / clean the air / oil separator


- A) Turn master switch OFF.
- B) Remove container by squeezing sides while pulling out.
- C) Remove saturated gauze.
- D) Clean container and mounting cap.
- E) Install clean gauze, then reinstall the container.

Periodic Maintenance - continued

Replace the regulator filters.

- A) Turn off water and air supplies using manual shut off valves.
- B) Unscrew filter cap. (Use 9/16" wrench.)
- C) Unscrew retainer nut and remove filter and o-ring.
- D) Install new filter and o-ring. Secure with retainer nut.
- E) Reinstall filter cap.

DA3924i

Cleaning / Infection Control

Barriers

Midmark recommends the use of disposable barriers on all clinician controls that may be in contact with clinician hands and fingers during dental procedures. The use of barriers significantly reduces the need for chemical cleaners, thus prolonging the life of the equipment.

Only use barrier material that is intended for use with dental equipment. Midmark recommends the use of an FDA market-cleared barrier such as Pinnacle Cover-all™. Follow barrier manufacturer instructions for proper use of these products.

Cleaning and Disinfecting

In addition to the use of barriers Midmark recommends the use of an EPA registered and FDA market-cleared cleaner/disinfectant such as Cavicide™ to be used on all clinician controls or surfaces that may come in contact with dental instruments during dental procedures.

Follow cleaner/disinfectant manufacturer instructions for proper use of the product. Care should be taken to avoid excessive application and pooling of liquids.

Cleaning and Disinfecting Assistance

For assistance with cleaning and disinfecting instructions contact the Midmark Technical Service Department at 1-800-Midmark, it is helpful to provide the delivery system model number and serial number when asking for assistance.

Additional information is available from the organizations listed below:	
Organization for Safety & Asepsis Procedures: http://www.osap.org	Dept. of Health & Human Resources, Centers for Disease Control & Prevention (CDC): http://www.cdc.gov
American Dental Association: http://www.ada.org	European Dental Association: http://www.eda-eu.org

General Purpose Cleaning

For general purpose cleaning, use cleaners that are appropriate for the situation, such as warm water and mild detergents, or a 10% solution of bleach with water.

Visual Inspection

After cleaning, visually inspect the product for deterioration of covers and touch pads. Do not use the delivery system if excessive discoloration, cracking, or other signs of wear are noticeable (see Calling for Service instructions).

Waterline Maintenance

Waterline maintenance is necessary to keep the count of heterotrophic bacteria from rising higher than desired levels. The desired level for a specific location should be determined by any local or regional guidelines. For example, The United States Centers for Disease Control and Prevention (CDC) guideline for heterotrophic bacteria is less than or equal to 500 CFU/mL (colony forming units per milliliter). Midmark recommends keeping this level under 200 CFU/mL.

Treatment can come in many forms. The most popular methods on the market currently are tablets and straw/cartridge based systems. Midmark recommends the use of a straw/cartridge based system that keeps the bacteria levels in check.

Regular monitoring should also take place to ensure that heterotrophic bacteria is not exceeding the desired limit. If the level is higher than desired, a shock treatment of the waterlines should be performed. When performing a shock treatment, be sure to check with the manufacturer of the regular treatment regimen being used to ensure chemical compatibility. Monitoring frequency should be established by your practice. As a suggestion, Midmark would recommend that you begin by monitoring on a monthly basis, and make adjustments to the frequency based on test results.

Per the CDC, routine flushing of the waterlines should be performed between every patient. Extra flushing maybe needed within Midmark equipment when tablets are used. Undissolved tablet particles can gather over time in places within the waterlines, obstructing the line and causing water flow to slow. By flushing the waterlines, water flow is maximized and should push any undissolved particles through.

Handpiece Accessories

Only use dental handpiece accessories with the delivery system that are FDA market-cleared and refer to manufacturer's instructions for proper cleaning and disinfecting. Either an autoclavable syringe tip or a single use disposable syringe tip may be used.

Equipment Alert

AUTOCCLAVABLE SYRINGE TIP STERILIZATION

The autoclavable syringe tips supplied with the delivery system must be sterilized prior to use with each patient, including initial use. Be sure to thoroughly rinse and clean syringe tips prior to sterilization, any debris may reduce the effectiveness of the sterilization. Recommended sterilization process is steam autoclave. Recommended parameters are 125°C (250°F) and 106 kPa (15 PSI) for 40 minutes at temperature and pressure.

Servicing Equipment

Calling for Service

If service is required, contact your authorized Midmark dealer.
[Model / serial number information will be required when calling for service]

To contact Midmark directly:

1-800-MIDMARK (1-800-643-6275)
8:00 am until 5:00 pm Monday - Friday (EST)
www.midmark.com

Replacement Parts

Syringe Tip Part number: 12051	
Regulator Valve Filter Part number: 053-1166-00	

Specifications / Compliance

Specifications Chart	
Electrical Rating for Elevance Chair:	See power supply output (153941 ...)
Complies to the applicable requirements of:	ES/IEC/EN 60601-1, CAN/CSA C22.2 No. 60101-1 (Safety Standards), EN/IEC 60601-1-2 (EMC Standards)
Air Supply Requirements:	Maximum 80-100 PSI (5.52 - 6.9 bar)
Water Supply Requirements:	Maximum 30 - 50 PSI (2.1 -3.4 bar) Use water with a hardness level <60 ppm of calcium carbonate.
Certifications:	Midmark Corporation.....ISO-9001
Classifications:	Class I, Type B Applied Part, except as noted for optional accessories, Ordinary Equipment (IPX0), Continuous Operation
Optional Accessories:	
Acteon Satelec Scaler	Type BF Applied Part
Acteon Satelec Curing Light	Type B Applied Part
Acteon Sopro Intra-oral Camera	Type BF Applied Part
Bien Air Electric Motor	Type BF Applied Part
Dentsply Cavitron Scaler	Type B Applied Part
Handpiece tubing and connectors intended to be used with ISO 7785-1 or ISO 7785-2 compliant air-driven handpieces	Type B Applied Part
Sota Claris Intra-oral Camera	Type B Applied Part

EMC - Manufacturer's Declaration and Guidance

Guidance and manufacturer's declaration – electromagnetic emissions		
The Midmark Elevance Unit dental device is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Elevance Unit dental device should assure that it is used in such an environment.		
<i>Emissions Test</i>	<i>Compliance</i>	<i>Electromagnetic environment – guidance</i>
RF Emissions CISPR 11	Group 1	The Midmark Elevance Unit dental device uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The Midmark Elevance Unit dental device is suitable for use in all establishments, including domestic establishment and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communications equipment and the Midmark Elevance Unit dental device			
The Midmark Elevance Unit dental device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Midmark Elevance Unit dental device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the Midmark Elevance Unit dental device as recommended below, according to the maximum output power of the communications equipment.			
<i>Radiated maximum output power of transmitter W</i>	<i>Separation distance according to frequency of transmitter m</i>		
	150 kHz to 80 MHz $d = 1.2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \times \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.69	3.69	7.38
100	11.67	11.67	23.34
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guidance and manufacturer's declaration – electromagnetic immunity

The Midmark Elevance Unit dental device is intended for use in the electromagnetic environment specified below. The customer or the user of the Midmark Elevance Unit dental device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environmental - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for AC and DC power lines I/O lines not tested, all less than 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line(s) to line(s) ±2kV lines(s) to earth	±1kV line(s) to line(s) ±2kV lines(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T for 0.5 cycle) 40% U_T (60% dip in U_T for 5 cycles) 70% U_T (30% dip in U_T for 25 cycles) <5% U_T (>95% dip in U_T for 5 s)	V Dip >30% of U_T for 10ms V Dip < 60% of U_T for 100ms V Dip > 95% of U_T for 5000ms	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Midmark Elevance Unit dental device requires continued operation during power mains interruptions, it is recommended that the Midmark Elevance Unit dental device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

<i>Immunity Test</i>	<i>IEC 60601 Test Level</i>	<i>Compliance Level</i>	<i>Electromagnetic environmental - guidance</i>
			Portable and mobile RF communications equipment should be used no closer to any part of the Midmark Elevance Unit dental device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3V rms 150kHz to 80MHz	3V	$d = 1.2 \times \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \times \sqrt{P}$ 800 MHz to 2.5 GHz
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation in meters (m). Field strength from fixed RF transmitters, as determined by the electromagnetic site survey, ^a Should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location is which the Midmark Elevance Unit dental device is used exceeds the applicable RF compliance level above, the Midmark Elevance Unit dental device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Midmark Elevance Unit dental device.

^b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.

Limited Warranty

Midmark Limited Warranty - Dental Products

SCOPE OF WARRANTY

Midmark Corporation (“ Midmark”) warrants to the original retail purchaser that it will at Midmark’s option repair or replace components of the dental products manufactured by Midmark (except for components not warranted under “Exclusions”) that are defective in material or workmanship under normal use and service. Midmark’s obligation under this limited warranty is limited to the repair or replacement of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the first retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows: Effective March 1, 2018 these applicable warranty periods are measured from the date of invoice to the original user, shall be as follows:

1. OPERATORY PRODUCTS

- a. Five (5) years for all products (except for the items in (b) through (e)).
- b. Two (2) years for upholstery (chairs and stools).
- c. “KINK-VALVE” module carries a ten (10) year warranty.
- d. The original light bulb on a new light carries a one (1) year warranty.
- e. Accessories not manufactured by Midmark are excluded including but not limited to Bien Air handpiece systems, Dentsply Cavitron scaler, Satelec scaler and curing light, and Sopro cameras.

2. ORAL SURGERY PRODUCTS are warranted for a period of one (1) year.

3. STERILIZER PRODUCTS are warranted for a period of one (1) year.

4. ULTRASONIC CLEANER PRODUCTS are warranted for a period of two (2) years.

5. AIR AND VACUUM PRODUCTS

- a. PowerAir® oil-less compressors – Five (5) years or 3,500 hours of use, whichever occurs first.
- b. PowerVac® and PowerVac® G dry vacuums – Five (5) years or 10,000 hours of use, whichever occurs first (except that the vacuum pump warranty term is ten (10) years or 20,000 hours of use, whichever occurs first).
- c. Classic Series® wet-ring vacuums – Five (5) years or 10,000 hours of use, whichever occurs first.
- d. PowerMax surgical suction – Two (2) years.
- e. Hg5 Series Amalgam Separator - One (1) year. (f) Midmark manufactured accessories - One (1) year.

6. SYNTHESIS™ DENTAL CASEWORK AND ARTIZAN® EXPRESSIONS PRODUCT

- a. Five (5) years for all products and components including door and drawer fronts, casters and slides, except for the items in (b), (c) and (d).
- b. Three (3) years for electrical components such as task lights/LED lights, cords, controls and accessories.
- c. Two (2) years for sliding track monitor mount and components and upholstery. (d) One (1) year for countertops and resin, including accessories.

7. IMAGING PRODUCTS are warranted for a period of two (2) years except for the ClearVision CR reader which is warranted for a period of one (1) year.

8. MIDMARK Replacement Parts and Accessories carry a ninety (90) day warranty

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for the following;

1. defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
2. products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environment conditions and electrical power requirements;
3. products considered to be of a consumable or sterile nature;
4. accessories or parts not manufactured by Midmark;
5. charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark
6. costs and expenses of routine maintenance and cleaning;
7. representations and warranties made by any person or entity other than Midmark;
8. matching of color, grain or texture except to commercially acceptable standards;
9. changes in color caused by natural or artificial light;
10. custom manufactured products;
11. alterations or modifications to the product by any person or entity other than Midmark; and
12. Products that would otherwise be covered under Sections 1 and 2 of this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical, animal health or dental market, as the case may be, in which purchaser intends to use the product.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

WARRANTY DISCLAIMER

THIS LIMITED WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, or implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

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