

# Folding Exam Table Backer Kit [002-1692-00]

## Applies to Models:

8012-00-001

8012-00-XXX

8012\*

Language of origin: English



### Equipment Alert

The wall structure must be capable of supporting the table load.  
Weight: 90 lb. (41 kg) Torque (max) 1660 ft. lb. (2250 N-m)




### Attention!

Steps 1 and 2 must be performed from the side of the wall **that the table will be mounted to**.

### Step 1: Position the template on the wall.

- Draw a **level** line on the wall at 24 1/2" (158 cm) from the floor.  
(Be sure the line spans two wall studs)
- Locate / mark the center of the two wall studs.
- Align the bottom of the template with the line drawn (B).
- Align the left or right edge of the template with the centerline of the appropriate stud (the left edge / stud is shown in the illustration).




### Equipment Alert

For interior wall installation only.

### Step 2: Drill holes for mounting plate.

- Use a long drill bit to drill six pilot holes at the locations marked on the template.  
The holes **must** go thru both sides of the wall. **Keep the drill level.**  
**Do not drill thru the studs!** If the template is aligned correctly, the hole location marks will be offset from the studs.
- Expand each pilot hole using a 1" drill bit, or a hole saw from **both** sides of the wall.







### Attention!

Steps 3 and 4 must be performed from the **back of the wall** that the table will be mounted to. This step requires two wooden supports (4" x 4" x 11") not supplied.


### Step 3: Cut access hole for wall supports.

- A) Draw a **level** line on the wall 2-3/8" **below** the center of the bottom holes drilled in Step 3.
- B) Locate / mark the center of the two wall studs.
- C) Draw a **level** line on the wall 5-1/2" from the centerline of each stud as shown.
- D) Use a saw to cut out the 5" x 12" section of wallboard marked in Step 4 (A / B / C).



### Step 4: Install two wall supports.

- A) On the side where the stud is between the access hole and the 1" holes:  
Position the wood support against the stud.  
Secure with two screws (not supplied).
- B) On the side where the 1" holes are between the access hole and the stud:  
Position the wood support so that it does **not** obstruct the 1" holes.  
Secure with two screws (not supplied).




VA1856i

**Attention!**

Steps 5 and 6 must be performed from the **back** of the wall that the table will be mounted to.


**Step 5: Install standoffs onto backer plate.**

- A) Thread six carriage bolts into six standoffs.
- B) Insert bolt / standoff into keyed hole in backer plate.
- C) Use pliers to tighten bolt / standoff.
- D) Repeat for remaining bolts / standoffs.

**Step 6: Install backer plate.**

- A) Align standoffs with 1" holes, then push backer plate in until it is against the wall.
- B) Secure plate to steel studs with two screws.
- C) Install four hole plugs.

Hole Plug

**Equipment Alert**

The backer plate can be utilized for installation with other wall construction types (i.e. concrete, block, brick, etc.), however, the contractor/installer must supply appropriate hardware in those instances.

**Step 7: Refer back to installation manual 003-2602-00 to complete Vet Exam Folding table installation.**