Veterinary Dental Suction Accessory [9A560001]

Applies to Models:
8000 (-001 thru -008)

Special Tools:
1/8", 7/32" and 3/8" Drill Bits

Intended Use:

This product is intended to provide suction for dental procedures being conducted by qualified veterinary professionals in both the office and hospital environments. This product is for veterinary use only.

WARNING

Verify all tubing and electrical wiring has been moved from path of drilling. Check tubing and wiring after drilling for accidental damage. Must be installed on left side of Midmark 1000 (as shown).



Step 1: Drill mounting holes.

- Turn off Midmark 1000 and unplug from receptacle.
- Open back door of Midmark 1000. Release pressure from tank by turning valve until it aligns with tubing (open position). After air is out of tank, close valve.
- Align top edge of template with bottom edge of top cover and wrap right edge of template around to front cover of Midmark 1000. Tape template in place, mark hole locations and remove template.
- Use a straight edge to verify holes are aligned and perpendicular with top cover.
- Move wires and tubing from inside the wall behind marked holes.
- Drill three pilot holes with 1/8" drill bit.
- Re-drill top and bottom holes through with 7/32" drill bit.
- Re-drill center hole through with 3/8" drill bit.
- Check wires and tubing for any damage from drilling. Replace if damaged.

Note


Tubing color may vary.

Wire & Tube Clearance

Step 2: Install backer plate.

- A) Reach inside the back of Midmark 1000, position backer plate locating studs through drilled holes.
- B) Install nut on top stud.
- C) Hang suction assembly on top stud and tighten nut.
- D) Install nut onto bottom stud.
- E) Route tubing from on/off switch through center drilled hole.
- F) Install collection container.


Caution

Verify disconnected tube is from **BACK** side of regulator gauge for 65 psi.
If not system will be operating at 100 psi. .

Step 3: Connect tubing inside Midmark 1000.

- A) Connect tube from suction assembly to Tee fitting provided.
- B) Disconnect tubing from back side of air regulator gauge and insert into Tee fitting.
- C) Connect short tubing from Tee fitting to back side of air regulator (where tube was removed in step B).
- D) Shut back door of Midmark 1000.


Note: In step B, there are two tubes connected to back of air regulator. Either tube may be disconnected.

To disconnect/connect tubing from tee..

- A) Pull down release ring on tee to disconnect tube.
- B) Push tube into ring to connect.

Completed Connections Shown

Tube Connection
on Back Side of Air
Regulator Gauge

Short Tube Connected to Tee
Connect to Back Side of Air Regulator Gauge


Tube From
Suction Assembly
*Shown in Front of Tank for Clarity

Removed From Back of
Regulator - Now Connected to Tee

VA1870i


Step 4: Install cover.

- A) Route HVE hose through strain relief on side of cover and secure onto accessory barb fitting.
- B) Place cover onto assembly, guiding on/off switch through grommet on top of cover.
- C) Install two painted screws.
- D) Place handpiece in holder.
- E) Perform functional check.

To Perform functional check...


- A) Plug in and turn on Midmark 1000.
- B) Allow tank to fill.
- C) Turn vacuum on.
- D) Verify suction at handpiece.

Equipment Alert

HVE has a wire to block large solids from getting stuck in system. If HVE loses suction, solids might be lodged on wire. To clean out large solids, pull top of HVE off and remove solids lodged in wire.

Wire

VA1871i

Cleaning

Note

Only autoclaveable suction tips or single use, disposable tips may be used. Refer to manufacturer's manual for recommended cleaning procedures.

Note

Dental unit waterline treatment is important for the health and safety of the patient and all clinic personnel. Due to the long, narrow waterline tubing in dental units, standing water can lead to microbial colonization and biofilm growth. Biofilm is a community of bacterial cells and other microbes that adhere to surfaces and form a protective slime layer. If left untreated these microbes can break off, transferring bacteria to the patient, or even become aerosolized into the clinic environment.

Waterline treatment is necessary to keep the count of heterotrophic bacteria from rising higher than desired levels. The desired level for a specific location should be determined by any local or regional guidelines. Regular monitoring should also take place to ensure that heterotrophic bacterial is not exceeding the desired limit. Monitoring frequency should be established by your facility. As a suggestion, Midmark would recommend that you begin by monitoring on a monthly basis, and make adjustments to the frequency based on test results.

Treatment can come in many forms. The most popular methods on the market currently are tablets and straw/cartridge based systems. Midmark recommends the use of DentaPure by Crosstex, part number DP365B. Refer to the manufacturer's instructions for use.

Extra flushing may be needed within Midmark equipment when tablets are used. Undissolved tablet particles can gather over time in places within the waterlines, obstructing the line and causing water flow to slow. By flushing the waterlines, water flow is maximized and should push any undissolved particles through.


Every veterinary practice setting is different and no one disinfectant is the best choice for every facility. There are several organizations educated in disinfection procedures. Refer to web sites of accredited sources that can assist veterinary personnel in choosing what is best for their practice. Listed below are a few of the sites:


Organization for Safety & Asepsis Procedures
<http://www.osap.org>

Dept. of Health & Human Resources Centers for Disease Control & Prevention (CDC)
<http://www.cdc.gov>

To clean vacuum tubing and HVE instrument...

- A) Empty the collection container and reinstall.
- B) Mix chosen treatment system, per manufacturer's instructions, into a clean bucket.
- C) Turn unit ON and suction the treatment solution through the HVE instrument.
- D) Empty solution from the collection container and reinstall.
- E) Fill bucket with 16 oz. of water.
- F) Turn unit ON and suction water from bucket through the HVE instrument.
- G) Empty the collection container and reinstall.

9A560001 Veterinary Dental Suction Accessory (includes items 1 thru 21)		
Item	Description	Qty.
1	Grommet	1
2	• Painted Screw	2
3	Cover	1
4	Nut, Keps (#8-32)	1
5	Assembly	1
6	Vacuum Generator Assembly	1
7	Collection Container Assembly, 1000 ML	1
8	Foam Pad, Reticulated	1
9	Nut, (3/4-16)	1
10	Tubing	5 ft.
11	Mounting Bracket	1
12	• Backer Plate Assembly	1
13	Pan HD, PHL Screw (#10-24 x 3/8")	2
14	• Nut Keps #10-32)	2
15	HVE Holder	1
16	Flat HD, PHL Screw (#8-32 x 3/4")	1
17	HVE Handpiece	1
18	Wire, Obstruction	1
19	Snap Bushing	1
20	• 1/4" O.D. Tubing	.5
21	• Tee Union, 1/4"	1
22	Dental Aspirator Suction Tips (qty 25)	1
23	HVE Suction Tips (qty 25)	1

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