M9[®] / M11[®] Self -Contained Steam Sterilizer

[English](#)

[Español](#)

[Français](#)

For Models:

M9 - 043

M11 - 043

SA1986p

User Guide

Product Information

Attention Canadian Users - Action Required!

Thank you for purchasing an M9®/M11® Steam Sterilizer.

To comply with changes in Canadian regulatory requirements you must now provide documentation, a Manufacturer's Data Report (MDR), to your appropriate provincial regulatory agency when purchasing and using a steam sterilizer.

To obtain the MDR from Midmark®, record the sterilizer Serial Number and the ASME National Board Number of your sterilizer. See the illustration below for the location of these numbers. Then, contact Midmark customer service at 1.800.MIDMARK and provide the recorded numbers.

Midmark will then complete the MDR and send it to you. Once you have the MDR, you must submit it to the appropriate provincial regulatory agency.

We apologize for any inconvenience and appreciate your assistance in helping maintain compliance to Canadian regulations for sterilizers. If you have any questions or concerns, please contact 1.800.MIDMARK.


Table of Contents

Important Information

Safety Instructions	4
Intended Use	4
Sterilization Technology	4
Electromagnetic Interference	4
Safety Symbols	4
Shipping Symbols	5
Transportation / Storage Conditions	5
Accessories, Tools and Service Parts	5
Included with Sterilizer	6
Component Location	7
Controls and Indicators	8
Sterilization Monitoring Guidelines	10

Installation

Operating Environment	11
Location Requirements	11
Electrical Requirements	12
Connecting the Power Cord	13
User Settings	14

Operation

Before Operating the Sterilizer	15
Filling the Reservoir	15
Qualification Testing	16
Guidelines for Loading	16
Standard Cycle Parameters	24
Cycle Operation	25
Post-Sterilization Processing	28
Programmable Cycle Buttons	29

Maintenance

Maintenance Messages	31
Daily Maintenance	32
Weekly Maintenance	33
Monthly Maintenance	35
Extended Use Maintenance	37

Thermal Printer (Optional)

Operating the Printer	38
Printer Tape Description	38
Example of Typical Printout of a Program Cycle	39
Paper Roll Removal / Installation	40

Tray and Cassette Tool (optional)

Using the Optional Tray / Cassette Tool	41
---	----

Troubleshooting

Troubleshooting Chart	42
Informational Messages	43
Accessing the last 5 Error Codes	43
Error Messages	43
Calling for Service	45

Specifications / Compliance

Specifications Chart:	
M9®	46
M11®	47
Water Purity Specifications	48

Warranty Information

Scope of Warranty	49
-------------------------	----

Important Information

Safety Instructions

The primary concern of Midmark® is that this equipment is operated and maintained with the safety of the patient and staff in mind. To assure safe and reliable operation:

- Read and understand this manual before attempting to install or operate the sterilizer.
- Assure that the appropriate personnel are informed on the contents of this manual.
(This is the responsibility of the purchaser).
- Assure that this manual is located near the sterilizer, or if possible, permanently affixed to the sterilizer.

Intended Use

The Midmark and Ritter M9® and M11® Steam Sterilizers can be used in veterinary offices, hospitals, clinics, laboratories and other facilities to sterilize heat and moisture stable reusable items (including dental handpieces) that are compatible with steam sterilization. Refer to Guidelines for loading and Standard Cycle Parameters in this manual for detailed information.

Sterilization Technology

The Midmark M9 and M11 utilize a dynamic-air-removal method to effectively remove air in the chamber and load while generating steam. There are two commonly used dynamic-air-removal systems; one type is pre-vacuum which utilizes a vacuum pump, and the other type is steam-flush pressure-pulse (SFPP) which utilizes a series of steam flushes and pressure pulses without the need for additional vacuum components. The Midmark M9 and M11 sterilization cycles utilize the SFPP type of air removal.

Electromagnetic Interference

This sterilizer is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this product:

- Remove interfering device from room
- Plug sterilizer into isolated circuit
- Increase separation between sterilizer and interfering device
- Contact Midmark if interference persists

Authorized Representatives

WARNING

In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.

Customers in Mexico should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Sofía Zavala | Asuntos Regulatorios

Carlos Nafarrate SA de CV

Avenida Las Américas 601,

Ladrón de Guevara, 44600,

Guadalajara, Jalisco, México | Tel (33) 3669 9800

Safety Symbols

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage.

Shipping Symbols

Proper shipping orientation

Keep dry

Temperature Limit

Fragile

Max. stacking height
(Refer to number on package)

Humidity Limit

2 Person Lift Caution

>75 lbs (34 kg)

Transportation / Storage Conditions

Equipment Alert

Water must be drained from the unit's reservoir before transporting / storing below +32°F (0°C).

Ambient Temperature Range:

-22°F to 140°F (-30°C to +60°C)

Relative Humidity:

10% to 90% (non-condensing)

Atmospheric Pressure:

7.2 psia to 15.4 psia (49.6 kPa to 106.4 kPa)

Accessories, Tools and Service Parts

Unless noted, items can be used on the M9® and M11®.

Accessories	Order Number
Speed-Clean, 1 (16oz. [.47 liter]) bottle	002-0396-00
Speed-Clean, 1 case (12 - 16oz. [.47 liter]) bottles	002-0396-05
Thermal Printer	9A599001
Data Logger	9A682001
Cassette Rack (Horizontal)	9A215001 (M11 only)
Cassette Rack (Vertical)	9A215002 (M11 only)


Common Service Parts / Tools	Order Number
Printer Thermal Paper Roll (Refill)	060-0016-00
Printer Stick-able Thermal Paper (Refill)	060-0016-01
Door and Dam Gasket Kit (M11 only)	002-0504-00
Door and Dam Gasket Kit (M9 only)	002-0361-01
Pouch Rack Kit (M11 only)	002-2108-00
Pouch Rack Kit (M9 only)	002-2108-01
Tray/Cassette Tool	9A307001
VistaCool direct-to-drain thermal reduction system	9A586002 (Dual)

Included with Sterilizer

Speed-Clean


Care / Operation Card


Speed-Clean SDS Sheet

Mikmark Corporation	07/24/2018 Rev P Date Created: June 2007 Revised: February 2015 Page 1 of 6
SAFETY DATA SHEET	
Section 1: Identification	
1.1. Product Identifier	
Special Case	
1.2. Recommended Use and Restrictions on Use Surface Cleaner	
1.3. Supplier's Details	
Developed for	Mikmark Corporation 60 Vista Drive, P.O. Box 286, Versailles, OH 45386-0286, USA Telephone: 614-937-0246 E-Mail: info@Mikmark.com
1.4. Emergency Telephone Number:	
(011) 725-666-08	
Section 2: Hazard Identification	
2.1. Hazard Classification	
Skin Irritation, Category 2	
Eye Irritation, Category 2A	
2.2. Label Elements	
Signal Word	
Warning	
Symbols:	
Exclamation Mark	
	
GHS07	
Hazard Statements:	
P261+P353+P361+P531: Pre-wash or pre-rinse skin contact areas with water before using. After use, thoroughly wash the lower and upper extremities.	
Precautionary Statements:	
Prevention:	
Wear eye protection.	
Wear protective gloves.	
Wash hands after handling this material.	
Response:	
IF IN EYES: Immediately flush eyes with water for 15 minutes, intermittently lifting the lower and upper eyelids.	
IF ON SKIN: Flush skin with running water for 15 minutes.	
2.3. Hazard Not Otherwise Classified	
None	


Power Cord

Pouch Rack
(M11® Rack Shown)


2 Small Trays

2 Large Trays

MA9010-1i


Component Location

SA1943i


Controls and Indicators

The sterilizer controls and indicators are shown in the illustrations on this and the following page. The accompanying tables describe the function of each control / indicator.

Control		Function
Display	See illustration	Indicates cycle selected, cycle temperature and exposure time for the selected cycle. During the cycle, display shows messages describing status of cycle. When cycle enters sterilization mode, remaining cycle time is displayed as well as temperature and pressure. Display also shows error message if a malfunction occurs. Refer to the <i>Troubleshooting</i> section of this manual for a detailed explanation of <i>Informational / Error Messages</i> .
Unwrapped button		A program cycle designed to process unwrapped instruments at: 270°F (132°C) for 3:00 minutes / 30 minute drying cycle.
Pouches button		A program cycle designed to process instruments in combination paper / plastic sterilization pouches or wrapped instruments at: 270°F (132°C) for 4 minutes / 30 minute drying cycle.
Packs button		A Program cycle designed to process packs of instruments or textiles at: 250°F (121°C) for 30 minutes / 30 minute drying cycle.
Handpieces button		A program cycle for dental handpieces which runs at: 270°F (132°C) for 4 minutes / 30 minute drying cycle.
Start button		Initiates selected program or, when SELECT CYCLE is displayed, pressing Start will activate heater for 10 minutes.
Stop button		Terminates selected program or function.
1 or 2 buttons NOTE: All material run in these cycles must be validated for sterilization by the user.		Programmable cycle buttons that allow an operator to create two different programmed cycles for special applications. Sterilization time and temperature, along with drying time and venting procedure can be adjusted or changed. Caution <i>These cycles are not FDA cleared and validation of sterility of items processed using them is the responsibility of the user.</i>

Controls and Indicators - continued...

Control		Function
P button		Programming mode button that allows operator to change temperature, time, dry time and/or venting procedure. Used in conjunction with buttons 1 or 2. (Refer to Programming Mode).
+ (plus) button		Allows temperature, time, vent mode or dry time to be increased or changed when in location 1 or 2 and the P (programming) mode is activated.
- (minus) button		Allows temperature, time, vent mode or dry time to be decreased or changed when in location 1 or 2 and the P (programming) mode is activated.
Door Handle	refer to illustration	For latching / opening door.
Water Level Indicator / Reservoir Drain Tube	refer to illustration	Shows amount of water in reservoir. Tube also used for drainage of reservoir into suitable container.
Fill Opening	refer to illustration	Access for filling reservoir with water.
Pressure Relief Valve Test Lever	refer to illustration	Allows operator to check pressure relief valve.
Thermal Printer (Optional)	refer to illustration	The printer (optional equipment) can be used to provide a permanent record of time, temperature and pressure during a cycle.

Sterilization Monitoring Guidelines

Note

This information below is provided for reference only. Contact appropriate state/local agencies for specific sterilization guidelines for your office. Additional information on infection control is available from the Centers for Disease Control and Prevention (CDC), Organization for Safety and Asepsis Procedures (OSAP) , Association for the Advancement of Medical Instrumentation (AAMI) and Association for Professionals in Infection Control and Epidemiology (APIC).

Physical Monitors

Temperature and pressure measuring devices can help detect sterilizer malfunctions. The sterilizer's control system aborts the cycle and displays a message if physical conditions go outside established limits. The optional thermal printer can be used to create a record of each load's actual cycle time, temperature and pressure.

Note

Use only FDA cleared chemical and biological indicators designed for steam sterilization that are compatible with the particular sterilization cycle temperature and exposure time being monitored. Use sterility monitors with each sterilization load. If a sterilizing cycle is terminated prematurely, reprocess instruments to ensure sterility of the load. Follow manufacturer's instructions for proper disposal of used indicators.

Chemical Indicators

Chemical indicators are designed to verify that conditions in the sterilizer chamber were adequate to achieve sterilization. They do not validate that a processed item is sterile. If a chemical indicator shows a failure, items in that load are considered non-sterile. Potential causes for sterilization failure include: improper cleaning, packing, loading, or a sterilizer malfunction. Determine the cause of any sterilization failure and remedy the situation before running the next cycle. Only FDA cleared chemical indicators labeled for use with the steam sterilization cycle parameters. e.g. temperature and exposure time, of the M9® / M11® Sterilizers should be used for monitoring the cycles. Follow the chemical indicator's instructions for proper storage, use, interpretation and disposal.

Biological Indicators


Biological indicators are microbiological devices designed to accompany items being sterilized to monitor adequacy of the sterilization process. If a biological indicator shows a failure, items in that load are considered non-sterile. Potential causes for sterilization failure include: improper cleaning, packing, loading, or a sterilizer malfunction. Determine the cause of any sterilization failure and remedy the situation before running the next cycle. Only FDA cleared biological indicators labeled for use with the steam sterilization cycle parameters. e.g. temperature and exposure time, of the M9 / M11 Sterilizers should be used for monitoring the cycles. Follow the biological indicators instructions for proper storage, use, interpretation and disposal.

Installation

Operating Environment

Ambient Temperature Range:	68°F to 104°F (+20°C to +40°C)
Relative Humidity:	< 80% (non-condensing) (Pollution Degree 2, in accordance to IEC664)
Normal Operating Altitude:	< 9842 ft. (3000 m) above sea level
Device approved for INDOOR USE ONLY .	
Device to be operated in a relatively dust-free environment. (Pollution Degree 2, in accordance to IEC664)	
Device should be connected to a power source with over-voltage limits less than 1500 watts from mains to ground. (Installation Category II in accordance to IEC664)	
The M9® and M11® will emit 5000 BTU / HR during operation.	

Location Requirements

Location Requirements - continued...

Support Surface

- Material should be water-resistant material. (Ex. laminate, stainless steel, stone, etc.)
- Surface must be level to ensure chamber fills with correct water level.
Improper water level in the chamber could cause a sterilizer malfunction.
- Surface should meet minimum dimensions listed below:

Dimensions

Depth (front to back) **M11®** - 60.5 cm (23.8") **M9®** - 51.8 cm (20.4")

Clearance Requirements

To ensure proper air circulation and to allow access to the reservoir fill port and drain coupling, adhere to the minimum clearance requirements listed below. If the sterilizer will be operated in continuous cycles, locate sterilizer where steam will not damage materials or equipment in the surrounding area.

Back of Unit - Back Wall..... 5.1 cm (2")
Front Support Surface - Front Sterilizer 2.5 cm (1")
Sides of Unit - Side Wall 5.1 cm (2")
Distance above Unit for Printer Access..... 12.7 cm (5")
Maximum Upper Cabinet Shelf Overhang **M11** - 46 cm (18") **M9** - 38 cm (15")
Under Cabinet or Shelf..... **M11** - 58 cm (23") **M9** - 56 cm (22")

Relocation Requirements for Sterilizer

- Disconnect power cord from electrical outlet and allow sterilizer to cool.
Drain water from reservoir or do not tip sterilizer, allowing water to spill.

Electrical Requirements

WARNING

For 115 VAC models: Use 104 - 127 VAC, 50/60 Hz alternating current only.

For 230 VAC models: Use 207 - 253 VAC, 50/60 Hz alternating current only.

Failure to do so may result in electric shock to personnel and / or damage to sterilizer.

Note

For safety, the unit must be connected to a properly polarized and grounded receptacle. Always use a power cord with grounding connections that match the receptacles in your location.

115 VAC Unit:

115 VAC, 50/60 Hz, 12 amp
Dedicated Supply Circuit*: 120 VAC, 50/60 Hz, 15 amp
Max. Power Consumption: 1425 Watts

230 VAC Unit:

230 VAC, 50/60 Hz, 6.4 amp
Dedicated Supply Circuit*: 230 VAC, 50/60 Hz, 10 amp
Max. Power Consumption: 1500 Watts

*** Power source must have over-voltage limits less than 1500 watts from mains to ground.**
(Installation Category II in accordance to IEC664)

Connecting the Power Cord

WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen, air, or nitrous oxide.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

WARNING

Check the serial number label on back panel of sterilizer to verify voltage rating for the unit. Failure to connect sterilizer to an appropriate power supply could result in damage to the unit and electrical shock to personnel.

Equipment Alert



For optimal performance, allow sterilizer to reach room temperature before operating.

To connect the power cord...

- A) Plug power cord into receptacle on back of sterilizer.
- B) Plug power cord into a properly polarized and grounded receptacle rated for a minimum of 15 amps. A dedicated circuit only used for the sterilizer is recommended.
- C) M9® and M11® are not equipped with an on/off switch, the display operates off very low power. (example: microwave oven display)

Note: When power is connected, the messages shown below will appear on the display.


Display:





* These screens will display the total number of cycles run on the unit, the model number (M9 or M11), the software version number, serial number date and time.

User Settings

User setting mode will enable the user to set the desired unit of measure and adjust the time clock. To enter the user settings mode, unplug and replug the power cord while pressing and holding the  button until the User Settings screen below appears.

Action	Description	Display
Press 	Press start to initialize the User Settings menus.	
Adjust Units of Measure 	The "+" and "-" buttons alternate the setting between English or Metric	
	When the desired Units of Measure appear on the display...	
Press 	To store the desired Units of Measure This brings up the Time Clock programming display	
Adjust Time Clock 	The "+" and "-" buttons will adjust the values of the following settings - Year / Month / Day / Hour / AM or PM / Minute.	
	When the desired values appear on the display...	
Press 	To store setting and bring up the next setting in the list. Repeat as necessary for all settings. When all settings are completed pressing the "P" will initiate normal power up.	

NOTE: User Settings for Date and Clock must be updated manually for day light savings.

Operation

Before Operating the Sterilizer

WARNING

Do not use this sterilizer for sterilizing volatile substances or for any purpose other than its intended design. Burns and toxic or explosive conditions could result.

Do not force door handle at any time. Chamber pressure may cause door to open with extreme force. If door handle does not move freely, allow unit to cool and depressurize for 40 minutes before opening door. Failure to adhere could result in serious personal injury.

Do not run the sterilizer without the tray plate in place. If the sterilizer malfunctions, immediately unplug sterilizer and call for service; do not attempt to repair the sterilizer yourself. Doing so could result in serious injury.

Caution

Programmable cycles 1 and 2 are provided for those applications requiring sterilization parameters different than the preset cycles. These cycles are not FDA cleared for medical use. All material processed in these cycles must be validated by the user to ensure sterility of the processed load.

Equipment Alert

For optimal performance, allow the sterilizer to reach room temperature before operating.

Filling the Reservoir


To fill reservoir...

- A) Open door to unit.
- B) Pour distilled water into fill opening until water level reaches the top of the fill level label on the water level indicator tube.

Equipment Alert

Use distilled water or water that meets the referenced water purity specifications. Failure to comply may result in sterilizer malfunction and/or premature failure due to excessive corrosion.

Qualification Testing

Your sterilizer should be tested after sterilizer installation, malfunctions, relocation, major repairs, being turned off or shut down for more than three days and after sterilization process failure. Qualification testing should be performed prior to placing the sterilizer in service. If multiple cycle types are used, e.g. “Pouches” and “Packs” each cycle type should be qualified. Qualification testing should include at least one Biological Indicator (BI) (sometimes referred to as Spore Tests) and one Chemical Indicator (CI). The test pack should be placed on the bottom tray near the chamber door and performed with items routinely processed and considered to be the most difficult to sterilize. Additional items should be placed in the chamber along with the Biological Indicator and Chemical Indicator so that chamber is fully loaded (don't exceed the maximum capacities listed in the tables under “Guidelines for Loading” in this manual). Three consecutive test runs, for each cycle type tested, with negative results from the BIs and the appropriate readings from all physical monitors and chemical indicators demonstrating complete sterilization, provide verification that the sterilizer has been properly installed (or reinstalled after relocation) or repaired to the manufacturer's specifications and that it will function effectively in the facility in which it is installed. All items processed during qualification testing should be quarantined until the results of the biological testing for all three test runs are available.

Guidelines for Loading

All items must be processed in accordance with Centers for Disease Control and Prevention (CDC), “Guidelines for Infection Control in Dental Healthcare Settings” – 2003, MMWR; 52 (no.RR-17) and “Guidelines for Disinfection and Sterilization in Healthcare Facilities” – 2008, which states:

“Items to be sterilized should be arranged to permit free circulation of the sterilizing agent (e.g., steam, chemical vapor, or dry heat); manufacturer's instructions for loading the sterilizer should be followed.”

Equipment Alert

Loads must be placed on trays at all times – unless the optional cassette racks are used – otherwise, serious instrument or equipment damage may occur.

Types of Items that can be processed in the M9® and M11®

Before placing any instrument in the M9 or M11, check with the instrument manufacturer to be sure the materials are compatible with steam sterilization and to verify the acceptability of sterilization parameters. The M9 and M11 are designed to sterilize the following:

- Metal instruments
- Rubber / plastic devices (*ex. suction cannulas, impression trays, etc.*)
- Wrapping / bundling materials (*ex. CSR wrap, instrument pouches, etc.*)
- Cassettes (which fit in the sterilizer trays or the cassette rack accessories)
- High / low speed handpieces
- Surgical instruments (*ex. ophthalmologic instruments*)

Types of Items that can be processed in the M9® and M11® - Continued...

Equipment Alert

Do not sterilize items composed of any of the following materials in the M9 or M11.

- Corrosion sensitive metal (ex. carbon steel, iron, etc.)
- Fragile items susceptible to breaking under pressure / high temperature
- Biomedical waste
- Plastics that may break down or produce residue when exposed to steam / high temperatures.

Examples

Polyethylene, Styrene, Cellulosics, ABS, PVC, Acrylic (Plexiglass™), PPO (Noryl™), Latex and Neoprene

Preparing Items for Sterilization

WARNING

Clean and dry instruments thoroughly before placing them into tray. Improper cleaning may result in non-sterile instruments or damage to the unit. Follow instrument manufacturer's guidelines and CDC recommendations for handling and cleaning instruments prior to sterilization.

Instruments must be thoroughly cleaned to remove all residual matter, such as debris, disinfectant residuals, blood, organic tissue, etc. General cleaning guidelines are listed below but the device manufacturer's instructions for proper cleaning and preparation of the device for sterilization should always be followed:

- Clean instruments immediately after use to avoid drying of residual matter.
- The use of automated cleaning equipment (e.g. ultrasonic cleaner or washer-disinfector) is recommended over manual cleaning for clinician safety and cleaning effectiveness.
- After cleaning, thoroughly rinse instruments with tap water to remove any loosened debris or residual cleaning fluid. The purity of tap water varies significantly thus, it's recommended the final rinse be done with water of adequate quality to avoid instrument staining. After rinsing instruments should be inspected for damage, debris, detergent residue and then dried before packaging.
- If the instrument manufacturer's instructions require lubrication of the instruments after cleaning, wipe off excess lubricant before packing for sterilization or loading into the sterilizer.

Immediate Use Sterilization

The M9 and M11 are capable of *Immediate Use* sterilization - sterilizing unwrapped instruments for immediate use. Place a surgical cotton towel, paper tray liner or CSR wrap folded to fit on the tray bottom before putting unwrapped items in the tray. Arrange unwrapped items on the towel so they do not touch one another (See Photo 1).

Immediate Use Sterilization - Continued...

Photo 1

Please consider the following when choosing whether or not to sterilize your instruments unwrapped:

- The sterility of unwrapped instruments is compromised upon exposure to a non-sterile environment. Follow CDC guidelines for using unwrapped, sterilized instruments.
- Due to the sensitive nature of some types of surgery (including, but not limited to ophthalmological), instruments used in such procedures must be wrapped or pouched in order to reduce their exposure to sterilization process residues. The water reservoir should also be drained and refilled with fresh distilled water on a daily basis when processing instruments for these procedures on a routine basis

WARNING

Do not overload the chamber! Adequate space is required around items in trays for steam circulation and drying. Failure to allow adequate space will compromise sterilization and drying. Items and packaging should be completely dry when removed from the sterilizer to minimize the potential for recontamination.

General Guidelines

- Use only M9® and M11® trays in their appropriate sterilizer. Using other trays could restrict air / steam flow to items resulting in inadequate sterilization and drying.
- All items must fit within the tray and not extend over the lip of the sterilizer tray. Instruments must not scrape the chamber walls when sliding the tray into the chamber.
- Jointed items must be sterilized in an open position so all surfaces are exposed to the steam (See Photo 2).

General Guidelines - Continued...


Photo 2

- Handpieces and instruments must be arranged in a single layer on the trays (not piled or stacked), to permit proper steam flow and drying.
- Glassware or utensils capable of holding water, e.g. bottles, basins, beakers, etc., should be positioned on the tray with the open side down so any water condensate drains from the container (See Photo 3). When sterilizing glassware check with the manufacturer to make sure it is compatible with steam sterilization.

Photo 3

- Rinse tubing with distilled water or water that meets the referenced water purity specifications and do not dry prior to sterilizing. Arrange tubing on the tray so there are no sharp bends and the tubing ends are open and unobstructed (See Photo 4).

Photo 4

General Guidelines - Continued...


- Follow the device manufacturer's instructions for disassembly of multi-part instruments prior to packaging / sterilization to assure all parts are adequately exposed to the steam.
- If items are being sterilized and stored for later use they must be packaged, e.g. pouched, wrapped, etc. and completely dry when removed from the sterilizer chamber to avoid potential recontamination.
- Variations in load configuration, size, wrapping materials and the environment may require the operator to increase the default drying time to assure all packaging and instruments are completely dry. See "Cycle Operation" section of this manual for instructions on adjusting drying time.
- When sterilizing a load that contains one or more handpieces, utilize the hand piece cycle, not Pouches or Unwrapped Cycle.

Pouching and Wrapping Items


The M9® and M11® are capable of sterilizing pouched or wrapped items to preserve sterility after processing.

- When pouching or wrapping items, use only sterilization pouches and wraps that have been cleared by the FDA and labeled for use with the steam sterilization cycle being used. Follow the manufacturer's instructions for use.
- Instruments made from different materials (stainless steel, carbon steel, plastic, etc.) should not be mixed in the same pouch or wrapped pack to avoid potential instrument damage.
- Pouches, wraps and the included instruments should not touch the chamber wall to allow proper steam circulation and avoid potential instrument damage.
- For proper steam circulation and drying the preferred orientation of pouches is resting on their edge, best accomplished using the Midmark Pouch Rack supplied with the M9/M11 sterilizer. If additional pouch racks are needed order kit P.N. 002-2108-00 (M11- 6 slot) or 002-2108-01 (M9- 5 slot). If more than two (2) pouch racks are used in a single load additional dry time may be required.

M11 Pouch Rack
002-2108-00

M9 Pouch Rack
002-2108-01

- The use of the Pouch Rack requires the removal of some of the sterilizer trays. When using the Pouch Rack load a single pouch per rack slot. When using paper / plastic sterilization pouches the pouches should be oriented in the rack so the plastic side of one pouch faces the paper side of the adjacent pouch.
- Pouches loaded directly in the trays will dry best if loaded with plastic side up.
- Do not layer pouches in the trays. Pouches should be loosely packed with handpieces and instruments single height loaded (not piled or stacked), to permit proper steam flow and penetration to the items (See Photo 2).

Pouching and Wrapping Items - Continued...

- If pouches or wrapped packs are labeled using a marking pen the ink should be nontoxic. On paper-plastic pouches only mark on the plastic side of the pouch. On wrapped packs mark the labeling information on the indicator tape or apply a separate label.

Shown below are some recommended pouch load configurations for the M9® (Photo 5 and 6) and M11® (Photo 7 and 8):

Caution

When sterilizing pouches, always use pouch racks on the bottom tray. Do not place pouches flat on the bottom tray.


M9 Pouches with Pouch Rack

Photo 5

M9 Pouches Mixed

Photo 6

M11 Pouches with Pouch Rack

Photo 7

M11 Pouches Mixed

Photo 8

- When using cassettes in the M9 / M11, follow the cassette manufacturer's instructions for use.
- Do not wrap items too tightly. Sterilization can be compromised if an item is excessively wrapped and the wrap is more likely to tear if wrapped too tightly.
- Cassettes can be loaded on the trays but they must fit within the boundaries of the tray and they should not touch each other, the tray above, or the chamber wall to allow proper steam circulation. The total instrument load, including the cassettes, should not exceed the limits listed in maximum capacities tables.
- When using the optional Horizontal and Vertical cassette racks (9A215001 and 9A215002) in the M11 the maximum load including the cassettes can be increased to 12 lbs. (5.44 kg). Additional dry time may be required.

Pouching and Wrapping Items - Continued...

Shown below are some recommended cassette load configurations for the M9® (Photos 9 and 10) and M11® (Photos 11, 12, 13 and 14):

Caution

Always use pouch racks on the bottom tray. Stand cassettes in the bottom tray on edge. Do not place cassettes flat on the bottom tray.

M9 Large Cassette
Photo 9

M9 Mixed Load with Small Cassette and Pouches
Photo 10

M11 Mixed Load with Cassette and Pouches
Photo 11

M11 with Vertical Cassette Rack
Photo 12

M11 with Horizontal Cassette Rack
Photo 13

M11 with Horizontal Cassette Rack and Top Tray
Photo 14


Textile Loads

- Clean textiles recommended for steam sterilization can be processed in the Packs cycle. Verify the acceptability of the Packs cycle sterilization parameters with the textile manufacturer to assure the textiles are compatible.
- All textile packs must not exceed the thickness limits specified in the maximum capacity tables below and must fit within the boundaries of the sterilizer trays.
- If multiple packs are placed on a single tray maintain a minimum of 1/4" (6.4 mm) between packs for proper steam circulation and drying. (Photo 15)

Equipment Alert

Do not use towels or packaging containing chlorine bleach residue. Failure to comply may cause rusting or discoloration of the chamber / trays and significantly shorten the life of the sterilizer.

Textile load

Photo 15

Load Size

Successful sterilization is dependent on correct loading of the sterilizer. Do not overload the sterilizer chamber! Adequate space must be maintained around all items placed in the chamber to assure proper steam circulation and adequate drying. The charts below are provided as a reference regarding the maximum loads that can be processed in the sterilizers but maintaining proper spacing between all items processed to assure good steam circulation and drying should be the guiding factor in determining the maximum load that can be processed.

M9® / M9D® Maximum Capacities **			
Load Type	Large Tray	Small Tray	Total
Solid Instruments	2.4 lbs. (1089 grams) or	1.6 lbs. (726 grams) or	8.0 lbs. (3629 grams) or
Handpieces	2 handpieces with other instruments 2.4 lbs. (1089 grams) or	2 handpieces with other instruments 1.6 lbs. (726 grams) or	8 handpieces with other instruments 8.0 lbs. (3629 grams) or
Packs*	1475 cm³ ≤ 3.2 cm thick (90 in³ ≤ 1.25 in. thick)	901 cm³ ≤ 1.9 cm thick (55 in³ ≤ .75 in. thick)	4752 cm³ (290 in³)
M11® Maximum Capacities **			
Load Type	Large Tray	Small Tray	Total
Solid Instruments	2.7 lbs., (1225 grams) or	1.8 lbs., (816 grams) or	9.0 lbs., (4082 grams) or
Handpieces	2 handpieces with other instruments 2.7 lbs. (1225 grams) or	2 handpieces with other instruments 1.8 lbs. (816 grams) or	8 handpieces with other instruments 9.0 lbs. (4082 grams) or
Packs*	2376 cm³ ≤ 3.8 cm thick (145 in³ ≤ 1.5 in. thick)	108 in³ ≤ 1.5 in. thick (1770 cm³ ≤ 3.8 cm thick)	8275 cm³ (505 in³)
* Allow a minimum of 1/4" (6.4 mm) space between each pack and from the chamber wall.			
** The default dry time may need increased due to variations in load configuration, wrapping materials and the environment to completely dry the chamber contents at these capacities.			

Standard Cycle Parameters

Cycle Type	Cycle Parameters			Drying Time ²	Items to be Sterilized (Always consult the item manufacturer's recommendations for sterilization.)
	Temp. (Min.)	Time	Press. ¹ (Ref.)		
 Unwrapped	270°F (132°C)	3 min.	27.1 psi (186 kPa)	30 min.	• Instruments loose on a tray. • Open glass or metal canisters. • Tubing not used in surgical procedures. (Max. length - 40" and Min. inside diameter - .187") • Loose items manufacturers recommend for exposure at 270°F (132°C) for 3 minutes. <i>Note: The sterility of unwrapped items is compromised on exposure to a non-sterile environment.</i>
 Pouches	270°F (132°C)	4 min.	27.1 psi (186 kPa)	30 min.	• Pouched or loosely wrapped instruments. • Multiple layers of instruments separated by fabric. • Wrapped trays of loose instruments. • Wrapped cassettes. • Wrapped items manufactures recommend for exposure at 270°F (132°C) for 4 minutes.
 Packs	250°F (121°C)	30 min.	15 psi (104 kPa)	30 min.	• Textiles and surgical packs wrapped for sterilization. • Items, except liquids, manufacturers recommend for exposure at 250°F (121°C) for 30 minutes.
 Handpieces	270°F (132°C)	4 min.	27.1 psi (186 kPa)	30 min.	• Dental handpieces (wrapped or unwrapped). <i>Note: Verify acceptability of sterilization parameters with handpiece manufacturer.</i>
1. The pressure shown in this table is for reference only. It's the ideal pressure of saturated steam at the sterilization temperature. The pressure shown on the sterilizer display may be higher. 2. Dry time can be changed from 0 to 60 minutes. Refer to Cycle Operation.					

Cycle Operation

Refer to the following steps for a detailed description of cycle operation:


Equipment Alert

The sterilizer will not operate unless the door is closed and latched properly.

Step 1: Close and latch the door.

- A) Lift the door handle, then push the door closed.
- B) While pushing in on the door, slide the door handle down to engage the latch.


Cycle Buttons
(Programmable)

Cycle Buttons
(Pre-set)

Note

Pressing the **Start** button when 'SELECT CYCLE' is displayed, at beginning or end of a cycle, activates the heater for 10 minutes. The display flashes 'ADDITIONAL HEAT'. This allows the operator to preheat the chamber before starting a cycle or to add additional time to the Dry mode at the end of a cycle. Pressing the **Stop** button will end the 'ADDITIONAL HEAT' time.

Equipment Alert

Using an incorrect sterilization program could result in non-sterile goods and may damage instruments. Consult instrument manufacturer for specific sterilization instructions.

Step 2: Select the desired cycle.

- A) Press the appropriate cycle button on the display panel.
(Refer to "Cycle Parameters" chart for time / temperature specifications).

Cycle Operation - continued...

After the cycle button is pressed, the parameters for that cycle will appear on the display.



Note

On units using the metric display, °F will display as °C and PSI will display as KPA.

Pressing **P** enables operator to change DRY time from 0 to 60 minutes in 1 minute increments on a pre-set cycle.

Pressing **−** decreases time. Pressing **+** increases time.

If the new Drying Time is desired on a regular basis pressing "P" after adjusting the dry time will store the new time as the default dry time for the selected cycle.

WARNING

STOP button may be pressed at any time to stop or interrupt a cycle. Goods must not be considered sterile if this occurs before the Dry Cycle. Sterilizer will return to SELECT CYCLE mode.

Step 3: Press 'Start' button to initiate cycle.

You will hear a "beep" for two seconds, indicating the cycle has started.


Cycle Operation - continued...

After the 'Start' button is pressed, the stage / status of the current cycle will appear on the display. The chart below illustrates the display messages that will appear during each stage of the cycle.

Stage of Cycle	Description	Display
Filling	Chamber begins filling with water.	FILLING CHAMBER
	When water reaches the proper level...	CHAMBER IS FULL
Heating	Display changes as temperature and pressure in chamber changes. <i>(Metric units may be displayed if desired).</i>	HEATING – UNWRAPPED XXX°F XX.X PSI
Sterilizing	Sterilizing begins when correct temperature / pressure is reached. Time remaining in cycle counts down while current temperature / pressure in chamber is continuously updated.	STERILIZING MM:SS XXX°F XX.X PSI
Ready to Vent	'READY TO VENT' is displayed when 10 seconds remain in sterilization cycle.	READY TO VENT XXX°F XX.X PSI
Fast Vent	When time runs out in sterilizing mode, the vent valve opens. Steam / water are released back into the reservoir. The display changes as temperature and pressure in chamber changes.	FAST VENT XXX°F XX.X PSI
	Caution <i>Keep clear when M9®/M11® door is ready to open! Failure to do so could result in severe burns from steam being released.</i>	
Door To Open (automatic)	An audible signal is emitted to indicate that the door is about to open. When pressure in chamber reaches zero, the door actuates to partially open (drying mode) position.	DOOR TO OPEN XXX°F XX.X PSI
This chart continues on the following page...		

Cycle Operation - continued...

Stage of Cycle	Description	Display
	Caution <i>The processed loads may still be wet if the Dry Cycle is aborted prior to completion. Always allow processed loads to dry in the sterilizer before they are handled to avoid re-contamination.</i>	
Drying	Time of Dry Cycle is counted down. If desired, the Dry Cycle can be aborted by pressing the STOP button.	
	When Drying time reaches 0:00...	
Cycle Complete	An audible signal is emitted for 10 seconds.	

Caution

- The processed load and inner surfaces will be hot. Avoid contact with hot surfaces. Failure to do so could result in serious burns.
- Do not apply downward pressure to the open sterilizer door while unloading the chamber. Doing so could cause the sterilizer to tip resulting in serious burns or injury from the trays or instruments sliding out of the chamber.

Step 4: Remove processed load from the chamber.

- A) Refer to “Using the Optional Tray/Cassette Tool” later in this manual.
- B) Sterilizer is now ready for another cycle.


Post-Sterilization Processing

After sterilization is complete, all items must be handled in accordance with accepted and documented standards, such as the Centers for Disease Control and Prevention (CDC) documents, “Guidelines for Infection Control in Dental Healthcare Settings” - 2003, MMWR; 52 (no. RR-17) and “Guidelines for Disinfection and Sterilization in Healthcare Facilities” - 2008, as well as any local requirements that may apply.

Qualified personnel responsible for infection control should prepare a protocol for handling sterilized items. This protocol should be followed by all personnel responsible for handling sterilized items.

Programmable Cycle Buttons

Cycle buttons 1 2 can be used for custom applications that are not covered by standard cycle programs. These cycles have not been FDA cleared for medical use. The programmable cycles allow the user to adjust cycle parameters in order to sterilize items that cannot be sterilized in any of the standard pre-programmed cycles. Once a custom program has been stored, it can be used just by pushing the 1 or 2 button. Use the instructions in the chart below to set time / temperature parameters for these buttons. (If you wish to change the settings, these buttons may be reprogrammed at any time).

Note:
The programmable cycles can be programmed for the monthly cleaning cycles with no dry time to expedite monthly maintenance.

NOTE:
Pressing the STOP button during this procedure will abort the changes and revert to the original settings.






Action	Description	Display
Press desired button 1 2	This selects the button that will be programmed.	PROGRAM 1 XXX°F, XX MINUTES
		FAST VENT XXX MINUTE DRY
Press P	This brings up the sterilization temperature programing display.	STERILIZATION TEMP: XXX.X°F
	Caution Sterilization temperature can be adjusted from a minimum of 250°F (121°C) to a maximum of 275°F (135°C). At high altitudes sterilization temperatures above 270° F (132° C) may cause the pressure relief valve to leak.	
Adjust Sterilization Temperature + -	The “+” and “-” buttons adjust the temperature by 1°increments.	STERILIZATION TEMP: XXX.X°F
	When the desired temperature appears on the display...	<P> FOR NEXT <STOP> TO CANCEL

This continues on the following page...

Programmable Cycle Buttons - continued...

NOTE:

Pressing the STOP button during this procedure will abort the changes and revert to the original settings.

Action	Description	Display
Press 	To store the temperature. This brings up the sterilization time programing display showing the current value.	
	Caution <i>Sterilization time can be adjusted from a minimum of 3 minutes to a maximum of 45 minutes. Because it is important to properly coordinate the cycle time with the sterilization temperature the minimum selectable time will be determined by the temperature setting.</i>	
Adjust Sterilization Time 	The "+" and "-" buttons adjust the time by 1 minute increments.	
	When the desired time appears on the display...	
Press 	To store the sterilization time. This brings up the vent speed programing display and the current value.	
Adjust Vent Speed 	In FAST vent, the valve fully opens and vents the chamber.	
	In SLOW vent, the valve opens for a fraction of a second (once per minute) to slowly vent the chamber.	
	Pressing "+" sets it to: FAST Pressing "-" sets it to: SLOW	
	When the desired setting appears on the display...	
Press 	To store the setting and bring up the dry time programing display and the current value.	

This chart continues on the following page...

Programmable Cycle Buttons - continued...

NOTE:
Pressing the STOP button during this procedure will abort the changes and revert to the original settings.

Action	Description	Display
Adjust Dry Time 	Dry time can be adjusted from 0 to 60 minutes.	
	The “+” and “-” buttons adjust the time by 1 minute increments.	
	When the desired time appears on the display...	
Press 	To store the setting and complete the programming process. The display will show the new cycle parameters.	
		
	Note <i>The programmed settings are retained under Program # button <1> or <2>. Even if power is interrupted, or the unit is unplugged the setting will be retained.</i>	
	Caution <i>All materials processed using these cycles must be validated by the user. Failure to do so could result in incomplete sterilization.</i>	

Maintenance

Maintenance Messages

To assure proper operation and maximum sterilizer life, carefully follow all recommendations for periodic maintenance. Recommended maintenance is easy to do and takes very little time.

One of the MOST important steps you can take to prevent problems with your sterilizer, is to ensure that ONLY distilled water or water that meets the referenced water purity specifications - NOT TAP WATER - is used in the sterilizer. Since the sterilizer operates with high water temperatures, any minerals dissolved in the water will form mineral deposits. This can prevent valves from opening or closing properly and can also lead to corrosion in the chamber and tubing.

Maintenance reminders will be displayed on the screen at the appropriate intervals to assist the operator. These reminders are removed from the display screen once a cycle is started.

The user is responsible for establishing a periodic maintenance procedure to assure correct operation of equipment and reliable sterilization of loads. Contact your local Midmark® distributor or service representative to develop a program for planned maintenance.

Daily Maintenance

Equipment Alert

If the sterilizer is used frequently to process dental handpieces that have been lubricated or dipped in dental milks, drain the water from the reservoir daily. Refill the reservoir with distilled water or water that meets the referenced water purity specifications.

Clean External Surfaces

- A) Wash the exterior of the sterilizer each day according to your facility's procedure for clinical contact surfaces, noting the following: (Use only quaternary disinfectants to disinfect unit. Staining, pitting, discoloration, or softening could occur if phenolic, iodophor, or glutaraldehyde-based disinfectant is used on plastic surfaces of the unit. Also, use of alcohol or aerosol spray cleaner / disinfectant containing substantial amounts of alcohol in the formula can damage the faceplate).
- B) Wring excess solution from the cloth.
- C) Using soft cloth, wipe all external surfaces.
- D) Follow the instructions provided with the cleaner / disinfectant used regarding rinsing and drying of the external surfaces.


Clean sterilizer door / dam gaskets

Caution

To prevent burns, allow unit to cool before cleaning gaskets and mating surfaces.

- A) Examine gaskets for possible damage.
- B) Clean gaskets and mating surfaces with a damp cloth.

Weekly Maintenance

Equipment Alert

Failure to change water may result in sterilizer malfunction. Do not use bleaching agents or any abrasive materials / substances in chamber (i.e. bleach, steel wool, wire brush, scouring powder, etc.). Failure to comply may result in damage to the chamber and/or other components.

Note

Every seven days, the autoclave will automatically display the PERFORM WEEKLY MAINTENANCE message. If power is disconnected, the cycle of weekly messages will be reset.


Caution

To prevent burns, allow unit to cool before draining reservoir.

Clean Chamber / Trays (including Rack and Plate)

- A) Disconnect the upper portion of the reservoir drain tube from the panel clips, bend it downward, and drain the reservoir water into a suitable container, e.g. a bucket and dispose of the water.
- B) Remove the trays, tray rack and tray plate from the sterilizer.
(Refer to the following page for instructions on removing / installing the tray rack and tray plate).
- C) Wash trays, rack, plate and inside of chamber with mild soap or Speed-Clean and distilled water or water that meets the referenced water purity specifications.
- D) Refill reservoir with distilled water or water that meets the referenced water purity specifications.

Weekly Maintenance - continued...


WARNING

*Allow unit to cool before removing or installing tray rack and plate.
Handle metal tray rack carefully to avoid injury. Do not run sterilizer without tray plate in place.*

To remove tray rack / plate...

- A) Remove trays.
- B) Using a screwdriver, pry plate up while pulling tray rack / plate out of chamber.

Equipment Alert

Install tray rack / plate with angled end of plate toward the back of the chamber.

Do not allow plate to contact the water level sensor.

To install tray rack / plate...

- A) Insert the tray rack into the tray plate.
- B) Place back of tray plate in chamber.
- C) Press down on tray rack, while sliding into chamber.

Monthly Maintenance

WARNING

Do not process instruments while flushing system.

Equipment Alert

Use only **Speed-Clean** to flush system. Failure to flush system with Speed-Clean may result in the premature failure of sterilizer components.

Note

Every 28 days, the sterilizer will automatically display the **PERFORM MONTHLY MAINTENANCE** message. If power is disconnected, the cycle of monthly messages will be reset.


Caution

To prevent burns, allow unit to cool before draining reservoir.

Clean Chamber / Plumbing

- A) With a cooled chamber, drain the sterilizer's reservoir and refill with clean distilled water or water that meets the referenced water purity specifications. Add one ounce of Speed-Clean sterilizer cleaner directly to the bottom of chamber.
- B) Run one **Pouches** cycle. 
- C) Press **Stop** button  when Dry Cycle begins. (Dry Cycle is not needed during maintenance.)
- D) Drain reservoir and refill a second time with clean distilled water or water that meets the referenced water purity specifications. .
- E) Rinse by running one **Unwrapped** cycle . Push the "Stop" button  when the drying cycle begins.
- F) Drain and refill reservoir with clean distilled water or water that meets the referenced water purity specifications, then allow sterilizer to cool.
- G) Remove tray plate, rack and trays. Wipe off with a damp cloth.

Equipment Alert

Use care when wiping the inside of the chamber. Failure to comply may result in damage to the heating element, steam temperature probe and/or level sensor probe.


Equipment Alert

Do not operate sterilizer without filters in place.

- H) *Remove and clean filters. The filters are intended to prevent debris from causing valve failures. Between regular monthly cleanings if the fill or vent times become too long or items will not dry the filters should be cleaned. (Refer to the illustration for location of filter screens.)*
- I) *Grasp filter and gently pull away from chamber wall while twisting slightly. (If necessary, pliers may be used to remove filters)*
- J) *Clean filters with Speed-Clean and distilled water. A small stiff bristled brush or ultrasonic cleaner may be helpful. Rinse filters with distilled water. Replace filter(s) if debris cannot be removed by cleaning.*
- K) *Wipe out the inside of chamber.*
- L) *Install filters. (Press inward, toward chamber wall while twisting slightly).*
- M) *Install tray plate, rack and trays.*

Remove / Clean Door and Dam Gaskets

- A) *Remove door and dam gaskets from chamber door, then remove the gasket ring from the door gasket.*
- B) *Clean gaskets and ring with Speed-Clean, distilled water and a soft brush.*
- C) *Inspect gaskets for damage / shrinking / swelling. Replace gaskets if damage is apparent.*
- D) *Press gasket ring into the channel in the door gasket and reinstall the gasket in the door.*
- E) *Install dam gasket.*

Check Pressure Relief Valve

(must be checked each month to assure it functions properly)


- A) Press **Unwrapped** button 
- B) Press **Start** button 

Caution

During the pressure relief valve check, steam will be vented from under the sterilizer. To keep from being burned, place a steam barrier (a rolled up towel) around the bottom of the sterilizer.

- C) Wait until pressure in chamber reaches 20 PSI (138 kPa).


- D) Pull upward firmly on the pressure relief lever for approximately 3 seconds, then release. (Steam should discharge freely from beneath rear of unit when lever is pulled. If the valve does not close completely when lever is released, pull lever again and release quickly so that it snaps closed. Repeat this until valve seats properly).

- E) Press **Stop** button 
(This aborts the cycle to prevent overheating).

Equipment Alert

If excessive force is required to open the pressure relief valve, or if the valve will not reseal properly, the valve must be replaced. (Refer to "Calling for Service" in this manual).

Extended Use Maintenance

The M9® and M11® are designed and tested to provide exceptional reliability throughout their service life. However, like all electro-mechanical devices they are subject to wear and degradation with use.

To ensure the integrity, performance and safety of all major components it is the responsibility of the user to have the sterilizer performance / operation verified by a Midmark® Authorized Service Provider at least every 10 years or 10,000 cycles, whichever comes first. After 10 years or 10,000 cycles of use, an annual inspection by a Midmark Authorized Service Provider is recommended.

Thermal Printer (Optional)

Operating the Printer

After the optional printer is installed and the sterilizer is plugged in, the printer is automatically powered up and initialized. No user intervention or setup is required.

Printer Tape Description

The printer will print the following information for each sterilizer cycle:

- Model Number: indicates sterilizer model and software version.
- Cycle Number: reflects actual cycle count of sterilizer.
- Sterilizer ID: Serial Number of sterilizer.
- Operator: a line is printed so the operator's signature can be recorded on the printer tape.
- Date and Time: of start of each cycle
- "BEGIN selected CYCLE": to indicate the beginning of the cycle selected by the operator.
- Summary of selected cycle set points.

Once the cycle starts, the printer will print the words "FILLING CHAMBER" to show that the sterilizer is filling with water.

Once the sterilizer begins the Heating Phase of the sterilization cycle, the word "HEATING" is printed and the printer will print the chamber temperature, pressure and elapsed time in 2 minute increments until the heating phase is completed.

When the sterilizer enters the sterilization phase, the word "STERILIZING" is printed and the printer will print the chamber temperature, pressure and elapsed time in 1 minute increments until the sterilization phase is completed.


When the sterilizer has completed the sterilization phase of the program cycle, the printer will print the words "VENTING CHAMBER" to show that the steam pressure is being exhausted from the chamber.

When the sterilizer enters the drying phase, the word "DRYING" is printed and the printer will print the words "DRYING START" and show temperature, pressure and elapsed time in 5 minute increments starting with 5:00. The printer continues to print the elapsed time in 5 minute increments until the drying phase is completed. The final record for the drying phase will include the words "DRYING COMPLETE". In the event the drying time is programmed to a time that isn't divisible by 5, the final printed record for the drying phase will reflect the actual programmed drying time in 1 minute increments, e.g. a programmed dry time of 12 minutes will have 5,10, and 12 minutes printed on the printer tape.

When the sterilizer has completed the drying phase of the sterilization cycle, the printer will print a summary of the sterilization cycle with the duration of each phase of the cycle and the Total Cycle Time. Following the summary the printer will print "CYCLE COMPLETE".

[NOTE: If drying cycle is aborted, "DRYING COMPLETE" and "CYCLE COMPLETE" will not print].

Example of Typical Printout of a Program Cycle

Paper Roll Removal / Installation

Equipment Alert

Use only thermal paper rolls with a maximum diameter of 48 mm (1.89") and width of 58 mm (2.28").

To remove paper...

- A) Insert finger into groove at the back of the small cover and lift upward until lid is released from locked position. To avoid damage do not use excessive force.
- B) The printer lid will hinge up towards the back of the printer.
- C) Remove any remaining paper.


To install the paper roll...

- A) Unroll 7.6 to 10.2 cm (3" to 4") of paper.
- B) Place the paper roll into cradle, unrolled side down.
- C) Hold unrolled edge of paper out the front of the printer and shut the lid by applying equal amounts of pressure on each side of lid until it latches.
- D) Press the paper feed button for additional paper.

Note

A pink stripe that gets progressively darker on the paper indicates when paper is low.

Printer Paper Refill	Order Number
Printer Thermal Paper Roll	060-0016-00
Printer Stick-able Thermal Paper	060-0016-01

Using the Optional Tray / Cassette Tool

Equipment Alert

Only use the 9A307001 Tray / Cassette Tool with Midmark® manufactured trays. Slide the proper plate into the handle as shown.

Note: There are three plates included. Each plate is designed for a specific purpose. (M9® trays, M11® trays, cassettes)

To remove M9 / M11 trays and cassettes...


- A) Hook the top, saw-toothed tab of tool to top center of tray lip / cassette.
- B) Rotate tool downward until the tool plate is completely under tray / cassette.
- C) Check to ensure tray / cassette is being held securely and then remove from chamber.

Note: Cassette tool can handle cassettes up to 3.8 cm (1 1/2") thick.

Caution

Trays / cassettes may be HOT - use care when removing or transporting. Hold the tray level and slightly elevated to prevent it from shifting and becoming dislodged. Failure to comply may result in personal injury due to burns.

Troubleshooting

Troubleshooting Chart

Use the following chart to assist in correcting minor problems with the sterilizer.


Problem	Possible Cause	Solution
Sterilizer does not operate (no display).	Sterilizer power cord came loose from supply outlet or back of sterilizer.	Assure Sterilizer power cord is plugged into outlet and sterilizer.
	No power at Sterilizer supply outlet.	Check circuit breaker for supply outlet. If problem recurs, unplug unit power cord and contact an authorized service technician. (see "Calling for Service")
	Fuse open on main P.C. Board.	Unplug unit power cord and contact an authorized service technician. (see "Calling for Service")
Steam escaping from pressure relief valve.	Pressure relief valve not properly seated.	Reseat pressure relief valve. (see "Monthly Maintenance")
Sterilization failure evidence from process monitor (chemical indicator, biological indicator, etc.).	Sterilization conditions were not present at location of the indicator.	Reload sterilizer in accordance with guidelines for loading the M9® or M11®. If problem recurs, unplug unit power cord and contact an authorized service technician. (see "Calling for Service")
	Indicator is out-of-date, is inappropriate for sterilization cycle, or has otherwise malfunctioned.	Use an indicator, appropriate for the load and cycle selected, that has been stored properly. Contact the indicator manufacturer for additional information on proper selection, use, storage and potential misapplication or malfunction.
Door gasket leaks.	Door gasket is damaged or dirty.	Clean or replace door gasket. (see "Weekly Maintenance")
Items are not dry at end of Dry Cycle.	Sterilizer is improperly loaded.	Reload sterilizer in accordance with guidelines for loading the M9 or M11. If the particular load configuration being used is consistently not dry at the end of the dry cycle increase the cycle drying time (see note on page 26) or reduce the load.
	Filter screen(s) clogged.	Clean or replace filter screen(s). (see "Monthly Maintenance")
Printer not printing.	Printer cables are not connected properly.	Ensure that printer cables are connected to printer.
	Printer malfunction. (status light flashes multiple times per second)	Unplug sterilizer power cord, wait 15 seconds and then plug sterilizer power cord back in.
	Printer is out of paper. (status light flashes once per second)	Insert a new paper roll. (see paper roll removal and installation section)
	Wrong paper	Assure printer paper is thermal paper
Display goes blank / black / is non legible	Electromagnetic Interference	Press any key except Stop on the keypad to reset the display without aborting the cycle. Pressing stop will also reset the display but will also abort the cycle in process.

Informational Messages

The chart below lists the informational messages that may appear during operation.

Message	Possible Cause	Solution
INITIALIZING SYSTEM	Unit power cord was just plugged in - standard informational message.	Normal operation will occur after a 4 second pause.
TOTAL CYCLES XXX M9, VX.X.X	Unit power cord was just plugged in - standard informational message. (Displays # of cycles, model number (M9® or M11®) and software version)	Normal operation will occur after this message.
MM/DD/YYYY 8:07 Am VXXXXXXX	Unit power cord was just plugged in - standard informational message. (Displays date, time and serial number.)	Normal operation will occur after this message.
PERFORM WEEKLY MAINTENANCE	This message is displayed every 7, 14 and 21 days after the unit is plugged into a power source to prompt the operator to perform weekly maintenance described in this manual.	Perform weekly maintenance. The message will automatically clear after the next cycle is completed.
PERFORM MONTHLY MAINTENANCE	This message is displayed every 28 days to prompt operator to perform monthly maintenance described in this manual.	Perform monthly maintenance. The message will automatically clear after the next cycle is completed.

Accessing the last 5 Error Codes

Action	Description	Display
Press and hold  during power up	This menu shows the last five error codes displayed on the unit. <i>NOTE: “1” is the most recent error code; “5” is the oldest.</i>	
Press the “stop” button 	To exit the error code menu press the stop button.	

Error Messages

The chart below and on the following pages, lists the error messages that may appear during operation.

Caution

If an error occurs more than once, stop using the sterilizer. Record the message or error code, unplug unit, then call an authorized service representative. (see “Calling for Service”)
If an error message includes: “Items Not Sterile”, the items in sterilizer shall not be considered sterile; they must be run through a successful sterilization cycle.

Error Messages - continued...

Message	Possible Cause	Solution
C010: POWER UP MODE SYSTEM PWR LOSS ITEMS NOT STERILE PUSH STOP TO RESTART	Unit had loss of power during cycle.	Press STOP button to restart.
CO60: POWER UP MODE SYSTEM HARDWARE	Power was interrupted briefly or an internal glitch.	Unplug unit power cord for 1 minute and then plug back in. If problem persists, contact an authorized service representative. (see "Calling For Service")
C102: FILL MODE STOP PRESSED	STOP button was pressed during cycle.	Press STOP button to restart.
C103: THROUGH C105 HEATUP, STERILIZE OR VENT MODE STOP PRESSED	STOP button was pressed during cycle.	Wait briefly (up to one minute) while chamber pressure / temperature dissipates. Press STOP button to return to select cycle mode where a new cycle may be initiated.
C106: DOOR MODE STOP PRESSED	STOP button was pressed during cycle.	Press STOP button to restart.
C232: FILL MODE WATER LOW	Water level in reservoir is too low.	Refill water reservoir with distilled water, or water that meets the referenced water purity specifications. Wait briefly (up to one minute). Press STOP button to return to select cycle mode where a new cycle may be initiated.
	Fill / Vent filter (in bottom of chamber) is clogged.	Clean Fill / Vent filter. (see "Monthly Maintenance")
C326: DOOR MODE DOOR CLOSED	Door latch safety switch is still making contact after door motor operated.	Open door.
C382: FILL DOOR OPEN	Sterilizer detects that door safety switch contacts opened.	Close the sterilizer door. (Cycle will continue where left off.)
C383: HEATUP MODE DOOR OPEN	Sterilizer detects that door safety switch contacts opened.	Wait briefly (1 minute) while chamber pressure / temperature dissipates. Press STOP button to return to select cycle mode. Initiate a new cycle.
C384: STERILIZE MODE DOOR OPEN	Sterilizer detects that door safety switch contacts opened.	Unplug unit power cord for 1 minute and then plug back in. If problem persists, contact an authorized service representative. (see "Calling For Service")
C533 through C633: STEAM TEMP LOW or STEAM TEMP HARDWARE or PRESSURE LOW	Sterilizer detects that temperature and / or pressure is outside the limits for normal operation.	Unplug unit power cord for 1 minute and then plug back in. If problem persists, contact an authorized service representative. (see "Calling For Service")

This chart continues on the following page...

Error Messages - continued...

Message	Possible Cause	Solution
C642 THROUGH C647 PRESSURE HIGH	Chamber pressure is outside the limits for normal operation.	Wait briefly (up to one minute) while chamber pressure / temperature dissipates. Press STOP button to return to select cycle mode. Initiate a new cycle.
C660 THROUGH C677 PRESSURE HARDWARE or PRESSURE OVERLIMIT	Chamber pressure is outside the limits for normal operation.	Unplug unit power cord for 1 minute and then plug back in. If problem persists, contact an authorized service representative. (see "Calling For Service")
C980 THROUGH C987 HI-LIMIT OPEN	High limit switch has opened for at least 1/4 second during specific operational mode.	Unplug unit power cord for 30 minute and then plug back in. If problem persists, contact an authorized service representative. (see "Calling For Service")

Calling for Service

Note

Please mark down any displayed code(s) and be sure to relay this information to the service technician.

Contact your **Midmark® Authorized Dealer**, or log on to www.midmark.com/technical-library.
Model and serial number information will be required when calling for service.

If service is required, contact Midmark directly:
US 844.856.1232 / Canada 937.526.8585
Monday-Thursday 8:00 am until 6:00 pm
Friday 8:00 am until 5:00 pm (ET)

Specifications Chart: M9®

Physical Dimensions:

Overall Length with Plug	51.8 cm (20.38")
Overall Width	38.9 cm (15.3")
Overall Height with Printer	40.9 cm (16.1")
Counter Area	49.1cm x 59.4 cm (19.3" x 23.4")
Chamber	22.9 cm x 38 cm (9" dia. x 15" deep)
Standard Tray, Large	18.6 cm x 30.5 cm x 2.2 cm (7 5/16" x 12" x 7/8")
Standard Tray, Small	14.3 cm x 30.5 cm x 2.2 cm (5 5/8" x 12" x 7/8")

Weights:

Weight with Empty Reservoir	73 lb. (33.1 kg)
Weight with Shipping Carton	81 lb. (36.7 kg)
Water Reservoir Capacity	1.1 Gallons (4.1 liters) to Full Mark usable volume is 0.5 gallons (1.9 liters)

Electrical Rating:

Note: A separate (dedicated) circuit is required for this sterilizer. Sterilizer should not be connected into an electrical circuit with other appliances or equipment unless circuit is rated for the additional load.

115 V models	115 VAC (+/- 10%), 12 Amp, 50/60 Hz
230 V models	230 VAC (+/- 10%), 6.4 Amp, 50/60 Hz

Fuse Ratings:

115 VAC	F1.....0.25 Amp, 250 V, Slo-blo, 1/4" x 1 1/4"
	F2.....15 Amp, 250 V, Fast Acting, 1/4" x 1 1/4"
230 VAC	F1.....0.125 Amp, 250 V, Slo-blo, 5 x 20 mm
	F2.....8 Amp, 250 V, Fast Acting, 5 x 20 mm
Safety Valve Setting	40 psi (275.8 kPa)
Heat Emission	5000 BTU / hr during operation

Certifications

ASME Boiler and Pressure Vessel Code,
Section VIII, Division 1.
Canadian Registration Number available

Midmark is an ISO 13485 and ISO 9001 Certified
Company

UL61010-1, 2nd Edition
IEC 61010-2-040:2005, 1st Edition
CAN/CSA C22.2, No. 61010-1, 2nd Edition
CSA C22.2, No. 61010-2-040-07 Part 2-040, -
1st Edition
IEC 60601-1-2 Edition 4.0 2014-02
CISPR II:2009, A1:2010 Class A, Group 1
IEC 61000-3-2:2014, Part 3-2
IEC 61000-3-3:2014, Part 3-3
ANSI/AAMI ST55:2010(R)2014
ANSI/AAMI/IEC62304

Specifications Chart: M11®

Physical Dimensions:

Overall Length with Plug	60.5 cm (23.8")
Overall Width	45.2 cm (17.8")
Overall Height with Printer	50.0 cm (18.1")
Counter Area	55.4 cm x 68.1 cm (21.8" x 26.8")
Chamber	28 cm x 46 cm (11" dia. x 18" deep)
Standard Tray, Large	22.9 cm x 38 cm x 2.9 cm (9" x 15" x 1 1/8")
Standard Tray, Small	16.8 cm x 38 cm x 2.9 cm (6 5/8" x 15" x 1 1/8")

Weights:

Weight with Empty Reservoir	99 lb. (44.9 kg)
Weight with Shipping Carton	131 lb. (59.4 kg)
Water Reservoir Capacity	1.4 Gallons (5.3 liters) to Full Mark usable volume is 1.0 gallons (3.8 liters)

Electrical Rating:

Note: A separate (dedicated) circuit is required for this sterilizer. Sterilizer should not be connected into an electrical circuit with other appliances or equipment unless circuit is rated for the additional load.

115 V models	115 VAC (+/- 10%), 12 Amp, 50/60 Hz
230 V models	230 VAC (+/- 10%), 6.4 Amp, 50/60 Hz

Fuse Ratings:

115 VAC	F1.....0.25 Amp, 250 V, Slo-blo, 1/4" x 1 1/4"
	F2.....15 Amp, 250 V, Fast Acting, 1/4" x 1 1/4"
230 VAC	F1.....0.125 Amp, 250 V, Slo-blo, 5 x 20 mm
	F2.....8 Amp, 250 V, Fast Acting, 5 x 20 mm
Safety Valve Setting	40 psi (275.8 kPa)
Heat Emission	5000 BTU / hr during operation

Certifications

ASME Boiler and Pressure Vessel Code,
Section VIII, Division 1.
Canadian Registration Number available

Midmark is an ISO 13485 and ISO 9001 Certified
Company

UL61010-1, 2nd Edition
IEC 61010-2-040:2005, 1st Edition
CAN/CSA C22.2, No. 61010-1, 2nd Edition
CSA C22.2, No. 61010-2-040-07 Part 2-040, -
1st Edition
IEC 60601-1-2 Edition4.0 2014-02
CISPR II:2009, A1:2010 Class A, Group 1
IEC 61000-3-2:2014, Part 3-2
IEC 61000-3-3:2014, Part 3-3
ANSI/AAMI ST55:2010(R)2014
ANSI/AAMI/IEC62304

Water Purity Specifications

Water Purity Specification Chart		
AAMI ST-46 (ref.)		
Evaporate Residue	$\leq 15 \text{ mg/l}$	$\leq 15 \text{ ppm}$
Silica	$\leq 2 \text{ mg/l}$	$\leq 2 \text{ ppm}$
Iron	$\leq 0.2 \text{ mg/l}$	$\leq 0.2 \text{ ppm}$
Cadmium	$\leq 0.005 \text{ mg/l}$	$\leq 0.005 \text{ ppm}$
Lead	$\leq 0.05 \text{ mg/l}$	$\leq 0.05 \text{ ppm}$
Rest of heavy metals, excluding Iron, Cadmium and Lead ¹	$\leq 0.1 \text{ mg/l}$	$\leq 0.1 \text{ ppm}$
Chloride	$\leq 3 \text{ mg/l}$	$\leq 3 \text{ ppm}$
Phosphate	$\leq 0.5 \text{ mg/l}$	$\leq 0.5 \text{ ppm}$
Conductivity (at 20°C) ²	$\leq 50 \text{ }\mu\text{s/cm}$	TDS $\leq 25 \text{ ppm}$
pH	6.5 to 8	
Appearance	Colorless, clean, without sediment	
Hardness	$\leq 0.1 \text{ mmol/l}$ (.585 gr/gal)	$\leq 10 \text{ ppm}$

Midmark Sterilizer Warranty

1. SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will at Midmark's option repair or replace components of the sterilizer products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this limited warranty is limited to the repair or replacement, at Midmark's option of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the first retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

2. APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of invoice to the original user, shall be as follows:

- (1) STERILIZER PRODUCTS are warranted for a period of one (1) year.
- (2) MIDMARK Replacement Parts and Accessories carry a ninety (90) day warranty.

3. OBTAINING WARRANTY SERVICE

Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via our website at www.midmark.com, by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380, or by phone at: 1.800.MIDMARK.

It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

4. EXCLUSIONS

This limited warranty does not cover and Midmark shall not be liable for the following;

- (1) defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) products which are not installed, used and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environment conditions and electrical power requirements;
- (3) products considered to be of a consumable or sterile nature;
- (4) accessories or parts not manufactured by Midmark;
- (5) charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) costs and expenses of routine maintenance and cleaning;
- (7) representations and warranties made by any person or entity other than Midmark;

Warranty Information (continued)

(8) matching of color, grain or texture except to commercially acceptable standards;

(9) changes in color caused by natural or artificial light;

(10) custom manufactured products;

(11) alterations or modifications to the product by any person or entity other than Midmark;

(12) Products that would otherwise be covered under Sections 1 and 2 of this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue outside of their market.

5. EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF DATA, LOSS OF USE, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

6. WARRANTY DISCLAIMER

THIS LIMITED WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

7. STATUE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period. In the event multiple warranty periods exist with respect to a product, the ninety (90) day period provided for herein shall begin to run from expiration of the warranty period for the component to which the claim relates.

8. NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

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

Caring for your M9 / M11 sterilizer...

Daily/Weekly Maintenance

Caution
To prevent burns, allow unit to cool before cleaning gaskets and internal surfaces.

Daily:

- **Clean External Surfaces and Gaskets**
- A) Wash exterior of sterilizer according to your facility's procedure for clinical contact surfaces. Use gels, aqueous or disinfectants to disinfect unit. See User Guide for more information.
- B) Cleanse door gaskets for damage. Clean gaskets and making them flexible with a damp cloth.


Weekly:

Caution
To prevent burns, allow unit to cool before draining reservoir.

Clean reservoir (surface)

- Drain water from the reservoir using drain tube on front of unit.
- Remove tray, tray rack, and tray plate.
- Wash reservoir items and make chamber with mild soap or Speed-Clean and water.
- Rinse reservoir with clean water.

Monthly Maintenance

WARNING
Do not process instruments while flushing system.

Equipment Alert
Use gels, Speed-Clean in flush system. Failure to flush system with Speed-Clean may result in the premature failure of sterilizer components.

Flush System

- With a closed chamber, drain reservoir and fill with clean water.
- Add one ounce of Speed-Clean sterilizer cleaner directly to the bottom of chamber.
- Rinse one **Pressure** cycle.
- Press **Stop** button **1** when Dry Cycle begins.
- Drain reservoir and fill with clean water.
- Repeat by running one **Unchanged** cycle **2**, press the **Stop** button **2** when the drying cycle begins.
- Drain and refill reservoir with clean water. Then allow sterilizer to cool.

Clean Filter Screens (fill / vent, and air)

Equipment Alert
Do not remove screens without filters in place.


- Allow sterilizer to cool. Then remove trays / rack / tray plate.
- Slide in the **Pressure** filter for position of the screen.
- Push / Pull filter to remove. (Use pliers if necessary).
- Clean filter with Speed-Clean and water. (Replace filter if screen cannot be removed by cleaning).
- Install filter in place, rack, and trays.
- Check **Pressure Filter** label: **Replace / Clean / Inspect Gasket**.

*Water when referenced is defined as distilled water or water that meets the water purity specifications found in the User Guide.
* Refer to the User Guide for complete list of maintenance procedures

Hoja SDS de Speed-Clean

Midmark Corporation		077-601540 Rev P Data Control Issue 2007 Revised February 2015 Page 1 of 5
SAFETY DATA SHEET		
Section 1: Identification		
1.1. Product Identifier		
Speed-Clean		
1.2. Recommended Use and Restrictions on Use		
Sterilizer Cleaner		
1.3. Supplier Details		
Revised Date: Midmark Corporation Address: 60 Veda Drive, P.O. Box 286, Versailles, OH 43084-0286, USA Telephone: 1-800-MIDMARK		
1.4. Emergency Telephone Number:		
(977) 525-6108		
Section 2: Hazard Identification		
2.1. Hazard Classification		
Skin Irritation: Category 2 Eye Irritation: Category 2A		
2.2. Label Elements		
Signal Word: Warning		
Symbols		
Exclamation Mark		
Hazard Statements		
P303+P361+P353: In case of contact with skin, wash thoroughly with plenty of water.		
P305+P351+P338: In case of contact with eyes, rinse thoroughly with plenty of water.		
P501: Dispose of contents and container in accordance with local, state, and federal regulations.		
Precautionary Statements		
P201: Hazardous for the Environment		
P202: Hazardous for the Environment		
P203: Hazardous for the Environment		
P204: Hazardous for the Environment		
P205: Hazardous for the Environment		
P206: Hazardous for the Environment		
P207: Hazardous for the Environment		
P208: Hazardous for the Environment		
P209: Hazardous for the Environment		
P210: Hazardous for the Environment		
P211: Hazardous for the Environment		
P212: Hazardous for the Environment		
P213: Hazardous for the Environment		
P214: Hazardous for the Environment		
P215: Hazardous for the Environment		
P216: Hazardous for the Environment		
P217: Hazardous for the Environment		
P218: Hazardous for the Environment		
P219: Hazardous for the Environment		
P220: Hazardous for the Environment		
P221: Hazardous for the Environment		
P222: Hazardous for the Environment		
P223: Hazardous for the Environment		
P224: Hazardous for the Environment		
P225: Hazardous for the Environment		
P226: Hazardous for the Environment		
P227: Hazardous for the Environment		
P228: Hazardous for the Environment		
P229: Hazardous for the Environment		
P230: Hazardous for the Environment		
P231: Hazardous for the Environment		
P232: Hazardous for the Environment		
P233: Hazardous for the Environment		
P234: Hazardous for the Environment		
P235: Hazardous for the Environment		
P236: Hazardous for the Environment		
P237: Hazardous for the Environment		
P238: Hazardous for the Environment		
P239: Hazardous for the Environment		
P240: Hazardous for the Environment		
P241: Hazardous for the Environment		
P242: Hazardous for the Environment		
P243: Hazardous for the Environment		
P244: Hazardous for the Environment		
P245: Hazardous for the Environment		
P246: Hazardous for the Environment		
P247: Hazardous for the Environment		
P248: Hazardous for the Environment		
P249: Hazardous for the Environment		
P250: Hazardous for the Environment		
P251: Hazardous for the Environment		
P252: Hazardous for the Environment		
P253: Hazardous for the Environment		
P254: Hazardous for the Environment		
P255: Hazardous for the Environment		
P256: Hazardous for the Environment		
P257: Hazardous for the Environment		
P258: Hazardous for the Environment		
P259: Hazardous for the Environment		
P260: Hazardous for the Environment		
P261: Hazardous for the Environment		
P262: Hazardous for the Environment		
P263: Hazardous for the Environment		
P264: Hazardous for the Environment		
P265: Hazardous for the Environment		
P266: Hazardous for the Environment		
P267: Hazardous for the Environment		
P268: Hazardous for the Environment		
P269: Hazardous for the Environment		
P270: Hazardous for the Environment		
P271: Hazardous for the Environment		
P272: Hazardous for the Environment		
P273: Hazardous for the Environment		
P274: Hazardous for the Environment		
P275: Hazardous for the Environment		
P276: Hazardous for the Environment		
P277: Hazardous for the Environment		
P278: Hazardous for the Environment		
P279: Hazardous for the Environment		
P280: Hazardous for the Environment		
P281: Hazardous for the Environment		
P282: Hazardous for the Environment		
P283: Hazardous for the Environment		
P284: Hazardous for the Environment		
P285: Hazardous for the Environment		
P286: Hazardous for the Environment		
P287: Hazardous for the Environment		
P288: Hazardous for the Environment		
P289: Hazardous for the Environment		
P290: Hazardous for the Environment		
P291: Hazardous for the Environment		
P292: Hazardous for the Environment		
P293: Hazardous for the Environment		
P294: Hazardous for the Environment		
P295: Hazardous for the Environment		
P296: Hazardous for the Environment		
P297: Hazardous for the Environment		
P298: Hazardous for the Environment		
P299: Hazardous for the Environment		
P300: Hazardous for the Environment		
P301: Hazardous for the Environment		
P302: Hazardous for the Environment		
P303: Hazardous for the Environment		
P304: Hazardous for the Environment		
P305: Hazardous for the Environment		
P306: Hazardous for the Environment		
P307: Hazardous for the Environment		
P308: Hazardous for the Environment		
P309: Hazardous for the Environment		
P310: Hazardous for the Environment		
P311: Hazardous for the Environment		
P312: Hazardous for the Environment		
P313: Hazardous for the Environment		
P314: Hazardous for the Environment		
P315: Hazardous for the Environment		
P316: Hazardous for the Environment		
P317: Hazardous for the Environment		
P318: Hazardous for the Environment		
P319: Hazardous for the Environment		
P320: Hazardous for the Environment		
P321: Hazardous for the Environment		
P322: Hazardous for the Environment		
P323: Hazardous for the Environment		
P324: Hazardous for the Environment		
P325: Hazardous for the Environment		
P326: Hazardous for the Environment		
P327: Hazardous for the Environment		
P328: Hazardous for the Environment		
P329: Hazardous for the Environment		
P330: Hazardous for the Environment		
P331: Hazardous for the Environment		
P332: Hazardous for the Environment		
P333: Hazardous for the Environment		
P334: Hazardous for the Environment		
P335: Hazardous for the Environment		
P336: Hazardous for the Environment		
P337: Hazardous for the Environment		
P338: Hazardous for the Environment		
P339: Hazardous for the Environment		
P340: Hazardous for the Environment		
P341: Hazardous for the Environment		
P342: Hazardous for the Environment		
P343: Hazardous for the Environment		
P344: Hazardous for the Environment		
P345: Hazardous for the Environment		
P346: Hazardous for the Environment		
P347: Hazardous for the Environment		
P348: Hazardous for the Environment		
P349: Hazardous for the Environment		
P350: Hazardous for the Environment		
P351: Hazardous for the Environment		
P352: Hazardous for the Environment		
P353: Hazardous for the Environment		
P354: Hazardous for the Environment		
P355: Hazardous for the Environment		
P356: Hazardous for the Environment		
P357: Hazardous for the Environment		
P358: Hazardous for the Environment		
P359: Hazardous for the Environment		
P360: Hazardous for the Environment		
P361: Hazardous for the Environment		
P362: Hazardous for the Environment		
P363: Hazardous for the Environment		
P364: Hazardous for the Environment		
P365: Hazardous for the Environment		
P366: Hazardous for the Environment		
P367: Hazardous for the Environment		
P368: Hazardous for the Environment		
P369: Hazardous for the Environment		
P370: Hazardous for the Environment		
P371: Hazardous for the Environment		
P372: Hazardous for the Environment		
P373: Hazardous for the Environment		
P374: Hazardous for the Environment		
P375: Hazardous for the Environment		
P376: Hazardous for the Environment		
P377: Hazardous for the Environment		
P378: Hazardous for the Environment		
P379: Hazardous for the Environment		
P380: Hazardous for the Environment		
P381: Hazardous for the Environment		
P382: Hazardous for the Environment		
P383: Hazardous for the Environment		
P384: Hazardous for the Environment		
P385: Hazardous for the Environment		
P386: Hazardous for the Environment		
P387: Hazardous for the Environment		
P388: Hazardous for the Environment		
P389: Hazardous for the Environment		
P390: Hazardous for the Environment		
P391: Hazardous for the Environment		
P392: Hazardous for the Environment		
P393: Hazardous for the Environment		
P394: Hazardous for the Environment		
P395: Hazardous for the Environment		
P396: Hazardous for the Environment		
P397: Hazardous for the Environment		
P398: Hazardous for the Environment		
P399: Hazardous for the Environment		
P400: Hazardous for the Environment		
P401: Hazardous for the Environment		
P402: Hazardous for the Environment		
P403: Hazardous for the Environment		
P404: Hazardous for the Environment		
P405: Hazardous for the Environment		
P406: Hazardous for the Environment		
P407: Hazardous for the Environment		
P408: Hazardous for the Environment		
P409: Hazardous for the Environment		
P410: Hazardous for the Environment		
P411: Hazardous for the Environment		
P412: Hazardous for the Environment		
P413: Hazardous for the Environment		
P414: Hazardous for the Environment		
P415: Hazardous for the Environment		
P416: Hazardous for the Environment		
P417: Hazardous for the Environment		
P418: Hazardous for the Environment		
P419: Hazardous for the Environment		
P420: Hazardous for the Environment		
P421: Hazardous for the Environment		

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Daily:

- **Clean External Surfaces and Gaskets**
- A) Wash exterior of sterilizer according to your facility's procedure for clinical contact surfaces. Use gels, aqueous or disinfectants to disinfect unit. See User Guide for more information.
- B) Remove door gaskets for storage. Clean gaskets and making sure with a damp cloth.

Weekly:

- **Clean Internal Surfaces and Gaskets**
- A) Wash interior of sterilizer according to your facility's procedure for clinical contact surfaces. Use gels, aqueous or disinfectants to disinfect unit. See User Guide for more information.
- B) Remove door gaskets for storage. Clean gaskets and making sure with a damp cloth.

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Clean reservoir (surface):

- A) Drain water from the reservoir using drain tube on front of unit.
- B) Remove tray, tray rack, and tray plate.
- C) Wash removed items and make chamber with mild soap or Speed-Clean and water.
- D) Refill reservoir with clean water.

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Equipment Alert
Use gels, Speed-Clean to flush system. Failure to flush system with Speed-Clean may result in the premature failure of sterilizer components.

Flush System:

- A) With a closed chamber, drain reservoir and fill with clean water.
- B) Add one ounce of Speed-Clean sterilizer cleaner directly to the bottom of chamber.
- C) Press Stop button (red) when Dry Cycle begins.
- D) Drain reservoir and fill with clean water.
- E) Press to start one (one) cycle (red) press the 'Stop' button (red) when the drying cycle begins.
- F) Drain and refill reservoir with clean water. Then allow sterilizer to cool.

Clean Filter Screens (fill, vent, and air):

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

- A) Allow sterilizer to cool. Then remove trays / rack / tray plate.
- B) Pull / front filter to remove. (Use gels if necessary).
- C) Pull / front filter to remove. (Use gels if necessary).
- D) Clean filter with Speed-Clean and water. (Replace filter if gels cannot be removed by cleaning).
- E) Install filter with water.
- F) Install filter with water.
- G) Check Pressure Relief Valve.
- H) Return / Clean / Repeat Cycle.


Fiche SDS Speed-Clean



Midmark Corporation		077-601540 Rev P Data Control Issue 2007 Revised February 2015 Page 1 of 5
SAFETY DATA SHEET		
Section 1: Identification		
1.1 Product Identifier		
Speed-Clean		
1.2 Recommended Use and Restrictions on Use		
Sterilizer Cleaner		
1.3 Supplier Details		
Revised Date: Midmark Corporation Address: 60 Veda Drive, P.O. Box 286, Versailles, OH 43140-0286, USA Telephone: 1-800-843-0000		
1.4 Emergency Telephone Number:		
(977) 525-6108		
Section 2: Hazard Identification		
2.1 Hazard Classification		
Skin Irritation Category 2 Eye Irritation Category 2A		
2.2 Label Elements		
Signal Word Warning		
Symbols		
Exclamation Mark		
Hazard Statements		
P303+P361+P353: In case of contact with skin, wash thoroughly with plenty of water.		
P305+P351+P338: In case of contact with eyes, rinse thoroughly with plenty of water.		
P308+P313: For all other effects, consult a physician.		
Precautionary Statements		
S2: Keep away from children.		
S3: Wear protective gloves.		
S4: Wash hands after handling this material.		
Response		
P303+P361+P353: In case of contact with skin, wash thoroughly with plenty of water.		
P305+P351+P338: In case of contact with eyes, rinse thoroughly with plenty of water.		
P308+P313: For all other effects, consult a physician.		
2.3 Hazard Not Otherwise Classified		
None		

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© Midmark Corporation 2018