


Cable Management for the Sidestream LoFlo Sensor

The following illustrations demonstrate proper cable management for the sidestream LoFlo sensor.



Always coil the cable in the direction of the natural cable coil, greater than 3 inches [7.6 cm] in diameter. The use of cable ties is acceptable, and will aid in maintaining the diameter of the cable during storage. Position the cable tie to avoid sharp kinks at the connector and sensor strain reliefs.



DO NOT kink the sensor cable, or bend the sensor cable less than 3 inches [7.6 cm] in diameter.



DO NOT coil and band the cable, as this will kink the cable and bend the cable smaller than 3 inches [7.6 cm] in diameter.



DO NOT fold the cable over on itself.



DO NOT wrap the cable around the sensor in any manner.



DO NOT wrap the cable around the sensor in any manner.




DO NOT wrap the cable around the optional pole mount in any manner.



DO NOT wrap the cable around the optional pole mount in any manner.



DO NOT coil the sensor cable around a ventilator circuit support arm, monitor stand, or any object less than 3 inches [7.6 cm] in diameter.



DO NOT pinch and straighten the cable during cleaning, handling, or use.