
Ritter®

228/229 Barrier-Free® Exam Chair

User Guide


For Models: 228 (-001/-002/-003), 229 (-002/-003)

English
Français
Español

Product Information

Model/Serial Number Label
Located on the Back Mounting Frame
(exact location may vary)

Warranty Registration

Scan the following QR code for access to the Examination Chair warranty registration.

Table of Contents

Important Information

Symbol Glossary	4
Intended Use	5
Electrical Requirements	5
Electromagnetic Interference	5
Cybersecurity	6-7
EMC - Manufacturer's	
Declaration and Guidance	8
Transportation/Storage Conditions	9
Operating Conditions	9
Disposal of Equipment	9
Authorized Representatives	10
Safety Instructions	11-12

Operation

Active Sensing Technology®	13
Base Up/Down Function	14
Back Up/Down Function (228)	15
Back Up/Down Function (229)	16
Paper Roll + Tear Strap	17
Stirrups	18
Footrest + Debris Tray	19
Drawer w/ Heater	20
Pelvic Tilt	21
Duplex Receptacles	22
Clean Assist Roller System	23

Maintenance

Calling for Service	24
Cleaning	24-25

Specifications/Compliance

Specifications Chart	26
Model Identification/Compliance Chart	27
Range of Motion	28
Dimensions	29

Accessories + Common Parts

Accessory/Common Part	30
-----------------------------	----

Warranty Information

Limited Warranty	31
Warranty Registration	31

Important Information

Symbol Glossary

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage.

Note

Amplifies an examination, practice or condition.

Duty Cycle (motor run time)
30 seconds ON, 5 minutes OFF

Serial Number

Proper shipping orientation

Keep dry

Fragile

Maximum stacking height
(palletted units)

Type B, Applied Part
(Upholstery + Stirrups)

Fuse rating specification

Protective earth ground

Dangerous voltage/shock hazard

Patient Weight Limit

Refer to manual

Manufacturer

Catalogue Number

Operating Instructions

Intended Use

This product is intended to be used as a chair to provide positioning and support of patients during general examinations conducted by qualified medical professionals in a standard medical office environment.

Standard Electrical Requirements

Before connecting the Midmark device power cord into an electrical power source outlet, it is recommended to contact a local licensed electrician or a local authority having jurisdiction over the internal building power distribution to verify that the device when powered up will comply with all local electrical codes. Use the electrical power rating marked on the device when determining the appropriate electrical branch circuit and outlet requirements. For 115V rated device models with a current rating at or above 12 Amps, a dedicated branch circuit may be required if the branch circuit protection is provided by a 15 Amp breaker.

WARNING

To completely isolate the chair from electrical mains supply, power cord must be unplugged.

When using high frequency surgical devices or endocardial catheters, use non-conductive material to insulate patient and clinician/operator from metal portions of the chair.

Failure to comply may result in electric shock or burns to patient.

Electromagnetic Interference

This product is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this product:

- Remove interfering device from room
- Plug chair into isolated circuit
- Increase separation between chair and interfering device
- Contact Midmark if interference persists

Note

The 228/229 Barrier-Free® Exam Chair is MR Unsafe and is not recommended for use in a Magnetic Resonance (MR) or MRI environment.

Cybersecurity

The 228/229 Barrier-Free® Exam Chair is designed to align with the cybersecurity labeling recommendations outlined in the FDA guidance, Cybersecurity in Medical Devices: Quality System Considerations and Content of Premarket Submissions. Compliance is subject to ongoing assessment and may evolve with regulatory updates.

Device Security Overview

The 228/229 Barrier-Free Exam Chair is designed with integrated security features intended to support safe operation in clinical environments. The device does not support wireless communication and is intended for use in controlled clinical settings. Users should follow best practices for access control and system updates to maintain cybersecurity protections.

Access Control

Access to the 228/229 Barrier-Free Exam Chair should be restricted to authorized clinical personnel. Care should be taken within the clinical environment to ensure that only authorized users are performing chair operations.

Physical and Port Security

The 228/229 Barrier-Free Exam Chair includes one RJ45 port which allows a foot controller to be connected to the chair. On the PCB assembly itself are a JTAG and UART connector. These are not accessible to end users in the clinical environment. These connectors have been assessed for risks within our Cybersecurity Risk Assessment and Threat Model. To prevent unauthorized access, it is recommended that any unused ports be restricted or monitored. Only Midmark-approved accessories should be used to maintain system integrity and cybersecurity compliance.

Software Bill of Materials (SBOM) Availability

A Software Bill of Materials (SBOM) is available upon request to provide transparency regarding the software components included in the 228/229 Barrier-Free Exam Chair. This SBOM helps users understand the device's software dependencies and potential security considerations within the software supply chain. Users may obtain the SBOM by contacting Midmark technical support.

Software Updates and Cybersecurity Maintenance

Midmark periodically releases firmware and security updates to maintain device performance and protect against known vulnerabilities. Updates are deployed as necessary, with users receiving notifications through official Midmark communication channels. To ensure compliance with cybersecurity best practices, users should contact Midmark Technical Services to install any applicable updates in accordance with institutional IT policies.

FAILURE TO INSTALL UPDATES IN A TIMELY MANNER MAY INCREASE CYBERSECURITY RISKS AND IMPACT COMPLIANCE WITH IT SECURITY POLICIES.

Cybersecurity - continued

Security Event Logging and Incident Response

The 228/229 Barrier-Free® Exam Chair is designed to automatically record system events, such as configuration changes and detected anomalies. These logs are intended to be securely stored within the device and can be reviewed by Midmark authorized personnel as needed.

- In the event of a security issue, users should:
- Report the incident to Midmark Support following institutional security procedures.
- Document details of the observed issue, including error messages and device behavior. Follow the recommended remediation steps provided by Midmark.

Midmark monitors cybersecurity risks and endeavors to provide mitigation guidance as necessary if a significant security vulnerability is identified. However, Midmark does not guarantee that all vulnerabilities will be detected or mitigated.

Backup and Recovery

The 228/229 Barrier-Free Exam Chair does not store user data and does not support user-configurable system settings. As a non-networked medical procedure table, no backup or recovery actions are required by the user. In the event of device malfunction, service and restoration are performed through Midmark Technical Services.

End of Support and Risk Communication

Midmark aims to provide software and security updates for a minimum of 7 years after device deployment. After this period, Midmark may continue to make security updates available, but general cybersecurity risks could increase. Users will generally receive advance notification regarding end-of-support timeliness and recommended actions. It is strongly recommended that users consider transitioning to an updated device before the end-of-support date to maintain cybersecurity protection.

Secure Device Decommissioning

The 228/229 Barrier-Free Exam Chair does not store sensitive user data and does not support user-configurable system settings. As a non-connected medical procedure table, no data removal or cybersecurity-specific decommissioning actions are required by the user. Device decommissioning should be performed in accordance with standard equipment disposal procedures and applicable local regulations. For disposal or end-of-life servicing, contact Midmark Technical Services.

Manufacturer Disclosure Statement for Medical Device Security (MDS2)

For additional cybersecurity details, including network security measures, communication mechanisms, and software update management, users should refer to the 228/229 Barrier-Free Exam Chair Manufacturer Disclosure Statement for Medical Device Security (MDS2). This document provides structured cybersecurity information relevant to IT and security teams managing Midmark devices. The MDS2 form is available upon request through Midmark technical support.

EMC - Manufacturer's Declaration and Guidance

Note

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Standard	Test	Level
Emissions		
CISPR 11 edition 6.2 (2019)	Conducted and Radiated RF	CISPR 11 Group 1 Class A
Immunity		
IEC 61000-4-2 edition 2.0 (2008)	Electrostatic Discharge	±8 kV Contact, ±15 kV Air
IEC 61000-4-3 edition 3.2 (2010)	Radiated RF EM Fields	3V/m 80MHz-2.7 GHz 80% AM at 1kHz
IEC 61000-4-3 edition 3.2 (2010)	Proximity Fields from RF Wireless Communications Equipment	IEC 60601-1-2, Table 9
IEC 61000-4-4 edition 3.0 (2012)	Electrical Fast Transients	±2 kV at 100kHz repetition frequency
IEC 61000-4-5 edition 3.1 (2017)	Surges	Line to Line: ±1 kV
IEC 61000-4-5 edition 3.1 (2017)	Surges	Line to Ground: ±2 kV
IEC 61000-4-6 edition 4.0 (2013)	Conducted RF	3V: 0.15MHz -80MHz 6V in ISM Bands between 0.15MHz - 80MHz 80% AM at 1kHz
IEC 61000-4-8 edition 2.0 (2009)	Rated Power Frequency Magnetic Fields	30 A/m
IEC 61000-4-11 edition 2.1 (2017)	Voltage Dips	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle
IEC 61000-4-11 edition 2.1 (2017)	Voltage Dips	70 % UT; 25/30 cycles Single phase: at 0°
IEC 61000-4-11 edition 2.1 (2017)	Voltage Interruptions	0 % UT; 250/300 cycles
IEC 61000-4-39 edition 1.0 (2017)	Proximity Magnetic Fields	IEC 60601-1-2, Table 11
ANSI C63.18 edition 2.0 (2014)	Radio Frequency	Separation Distance Cellular 5G, FDD ≥ 30 cm Wireless Power Transfer (WPT) ≥ 30 cm

Transportation/Storage Conditions

Ambient Temperature Range: -30°C to +60°C (-22°F to +140°F)
Relative Humidity 10% to 90% (non-condensing)
Atmospheric Pressure 500hPa to 1060hPa (0.49atm to 1.05atm)

Operating Conditions

Ambient Temperature Range: +10°C to +40°C (+50°F to +104°F)
Relative Humidity 30% to 75% (non-condensing)
Atmospheric Pressure 700hPa to 1060hPa (0.69atm to 1.05atm)

Disposal of Equipment

At the end of product life, the chair, accessories, and other consumable goods may have become contaminated due to normal medical use. Consult local codes and ordinances for proper disposal of equipment, accessories and other consumable goods.

Authorized Representatives

WARNING

In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.

Customers in Mexico should direct all questions, incidents, and complaints to the Midmark Authorized Representative:

Carlos Nafarrate SA de CV
Avenida Las Américas 601,
Ladrón de Guevara,
44600, Guadalajara, Jalisco, México
Tel (33) 3669 9800
<https://www.nafarrate.com/> [nafarrate.com]

Customers in Panama should direct all questions, incidents, and complaints to the Midmark Authorized Representative:

TECNIDENTAL PANAMA S.A.
VILLA MAGNA, VILLA ZAITA
OFFICE B8 PANAMA
Panama
Email: tecnidentalpanama@gmail.com
Tel: 507 6416-2235
+507 6492-0884
<https://tecnidentalpa.com/> [tecnidentalpa.com]

Customers in Costa Rica should direct all questions, incidents, and complaints to the Midmark Authorized Representative:

Emergo Costa Rica S.R.L.
Edificio Ribleh 2 Piso, Diagonal a la Biblioteca Municipal
Heredia, Belén
La Ribera, Costa Rica
Tel: +(506) 2201-14 53
Email: LAvigilance@ul.com

Customers in Colombia should direct all questions, incidents, and complaints to the Midmark Authorized Representative:

Emergo Colombia
Calle 76 No. 4-28 of. 803
Bogotá - Colombia 110221
Tel: (57) 313 4419785
CC: LST.AUS.ProgramCommColombia@ul.com

Safety Instructions

Note

Maintenance and service should not be performed on the chair while in use with a patient.

WARNING

No modification of this equipment is allowed.

Do not use outlet strips or extension cords with the duplex receptacles.

Connecting equipment to the duplex receptacles effectively leads to creating a Medical Electrical System and can result in a reduced level of safety. Refer to IEC 60601-1 clause 16 for requirements related to Medical Electrical Systems.

To avoid risk of electric shock, this equipment must only be connected to a supply with earth ground.

Equipment is not suitable for use in the presence of a flammable anesthetic mixture.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

This product can expose you to chemicals including Antimony Oxide (Antimony Trioxide), which is known to the state of California to cause cancer. For more information go to www.p65warnings.ca.gov

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the 228/229 Barrier-Free® Examination Chair including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Safety Instructions - continued

WARNING

Be sure that all personnel and equipment are clear of the chair before activating any function. Keep a minimum 6 in. (15.25 cm) clearance surrounding the chair. Failure to do so could result in personal injury.

At the chair's lowest height [16.5 in. (42 cm)], the clearance between the bottom of the drawer and the floor is approximately 4 in. (10.16 cm). Make sure patient (esp. feet) and all objects are clear of this area before lowering chair.


Equipment Alert

Chair must be positioned in normal use such that the power cords can be easily accessed to unplug.

If the chair malfunctions, immediately release controller, unplug the chair power cord, and assist the patient from the chair. If chair continues to malfunction, call for service.

Failure to comply with these instructions may result in injury to patient and/or damage to equipment.

Operation

Active Sensing Technology® (not available on all models)

Note

The chair will sound a series of "beeps" when powered on.


WARNING

Be sure that all personnel and equipment are clear of the chair before activating any function. Failure to do so could result in personal injury.

The Active Sensing Technology prevents the chair from moving Down if the bottom of the drawer contacts another object. If the Active Sensing Technology is activated, it will sound a series of "Beeps," disable the Base Down function and move chair in the Base Up direction. To resume Base Down operation, remove object from under the bottom of the drawer and press the Base Up/Down buttons on controller.

Note

The Active Sensing Technology is not available on models without drawer.

Base Up/Down Function

Attention!

The Active Sensing Technology® will stop the Base Down function and retract Up if the bottom of the drawer contacts another object.

(You will hear a series of “Beeps” if the Active Sensing Technology is activated.)

To resume Base Down operation:

Remove object from under the bottom of the drawer.

Press Base Down button.



Note

This control port is intended for a wired Ritter® foot or hand control only.

Note

The chair automatically stops at its maximum and minimum height.


*To move the chair Up or Down...
Press and hold the appropriate button.*

Back Up/Down Function - 228

To position the back section...

- A) While holding the back section in place, squeeze the Back Release Handle (either side).
- B) Move the back section.
- C) Release handle to lock back section in position.

Back Up/Down Function - 229

Note

The back section automatically stops at its maximum and minimum height.

To position the back section...

Press and hold the appropriate button.


Paper Roll + Tear Strap

Paper Roll Size (maximum):..... 61 cm wide x 9.2 cm diameter
(24 in. x 3-5/8 in.)


Note

One additional paper roll may be stored under the back section.

WARNING

To avoid injury, do not put head or extremities under the paper tear strap.

To install paper tear strap...
Attach strap to snaps on each side of chair.

Stirrups

Caution

In office environments where the room temperature exceeds 35°C (95°F), the surface of the stirrup foot hold may reach 46°C (114.8°F) while in its stowed position. To prevent potential discomfort when in use, stirrup covers can be installed or socks should be worn by the patient.


Caution

Be sure the stirrups are locked in place before using.
The stirrups will not support the patient's entire weight.

To position the stirrups...

- A) Pull the stirrup out, then unfold.
- B) Lift the stirrup slightly, then move it left or right as desired.
- C) Release stirrup to lock in position.

Footrest + Debris Tray (not available on all models)


WARNING

Do not use footrest to support patient's entire weight.

To extend footrest...

Pull footrest out to desired position.

To access the debris tray...


Pull debris tray out to desired position.

Drawer w/ Heater (not available on all models)

The drawer heater warms to 33°C (92°F) to 44°C (112°F) measured at the center of the drawer at room temperature and nominal input voltage.

To turn the drawer heater ON/OFF...
Press and release the drawer heater switch.

Note: Switch will illuminate when heater is ON.

Pelvic Tilt (not available on all models)


WARNING

To avoid risk of injury do not put extremities under the seat.


To elevate seat section...

Lift seat until pelvic support drops into place.

To return seat section to flat position...

- A) Lift seat section slightly.
- B) Pivot pelvic support handle upward.
- C) Lower seat section slowly.


Duplex Receptacles (not available on all models)

The duplex receptacles on each side of chair provide power for accessories used during medical examinations. There are two circuit breakers located under the seat section. If the receptacles' maximum load is exceeded, the circuit breakers interrupt power to the receptacles.

Maximum Load.....115 VAC, 4 amps

To reset the circuit breaker...
Push circuit breaker switch.

Clean Assist Roller System (not available on all models)

Control levers located under the back section on each side, activate and deactivate the Clean Assist Roller System. The rollers are located on each corner of the chair to allow positioning of chair.

WARNING

To prevent injury, never activate/deactivate Clean Assist Roller System or move chair with patient aboard. Use caution when moving chair.

Equipment Alert

Be sure to activate both control levers.

WARNING

To avoid risk of fall, the Clean Assist Roller System must be disengaged while a patient ingresses or egresses the chair.

To activate Clean Assist Roller System...


- Raise back section to 30 degrees or greater.
- Ensure chair is at least two inches higher than its lowest position.
- Move control levers on each side of chair to the Clean Assist Roller System engagement position.
- Move chair to the lowest position to activate rollers.

Equipment Alert

The intended use of the Clean Assist Roller System is for chair positioning within the exam room only.


Clean Assist Roller
System Engagement
Position

Equipment Alert

After Clean Assist Roller System has been activated, disconnect power cord and foot control and stow away from chair to prevent damage to cords.

Clean Assist Roller
System Disengagement
Position

To deactivate Clean Assist Roller System...

- Move chair to desired location.
- Reconnect power cord and foot control.
- Raise chair at least 2 inches to deactivate Clean Assist Roller System.

Note: An audible click will occur to confirm deactivation.

Note

Control levers will automatically reset to the Clean Assist Roller System disengagement position.

Maintenance

Calling for Service

Note

Model/serial number information is required when calling for service.

Note

Detailed service information is available on [midmark.com](https://www.midmark.com) in the technical library for your product.

If service is required, contact Midmark directly: US 844.856.1230/Canada 937.526.8585
Monday-Thursday 8:00AM until 6:00PM Friday 8:00AM until 5:00PM ET.

Cleaning

WARNING

The upholstery should be inspected before each use and replaced if punctured or otherwise damaged.

Equipment Alert

The upholstery is resistant to most medicinal-type stains, but may be damaged by solvents and dyes. Immediately remove any fluids spilled on the upholstery.

Upholstery

Whenever upholstery has been soiled, follow a two-step process of cleaning and disinfection as follows. Wash it with a mild liquid soap and water mixture, rinse with clean water and dry completely. After it is dry, apply an EPA registered disinfectant per the manufacturer's instructions. Alternatively, there are EPA registered products that are both cleaner and disinfectant in one that can be used in a single step per manufacturer instructions. To minimize damage caused by chemical disinfectants, do not allow them to pool on surfaces. Once the recommended contact time has elapsed, rinse with clean water and dry the disinfected surface. Complete a visual inspection to confirm that no additional contamination remains and repeat as necessary.

For additional cleaning and disinfection recommendations, follow the current [Guideline for Disinfection + Sterilization in Healthcare Facilities from the CDC](https://www.cdc.gov/infection-prevention/disinfection/).

Detailed care and maintenance instructions are included with your product. Additionally, the QR code below will take you to a list of suggested disinfectants for use on your upholstery. Another disinfectant option is to mix standard bleach and water at a ratio of 1:10, respectively.

Cleaning - continued

High-Touch Metal and Plastic Chair Surfaces

Whenever a chair surface has been soiled, follow a two-step process of cleaning and disinfection as follows. Wash it with a mild liquid soap and water mixture, rinse with clean water and dry completely. After it is dry, apply an EPA registered disinfectant per the manufacturer’s instructions. Alternatively, there are EPA registered products that are both cleaner and disinfectant in one that can be used in a single step per manufacturer instructions. To minimize damage caused by chemical disinfectants, do not allow them to pool on surfaces. Once the recommended contact time has elapsed, rinse with clean water and dry the disinfected surface. Failure to do so could result in premature damage to these surfaces. Complete a visual inspection to confirm that no additional contamination remains and repeat as necessary.

For additional cleaning and disinfection recommendations, follow the current [Guideline for Disinfection + Sterilization in Healthcare Facilities from the CDC](#). The table below lists suggested types of disinfectants at acceptable levels of active ingredient concentration.

<i>Disinfectant Type</i>	<i>Maximum Active Ingredient Concentration</i>
Bleach	0.53%
Hydrogen Peroxide	1.40%
Quaternary Ammonium	Total Alcohol: 18% Total Ammonium Compounds: 0.30%

Another disinfectant option is to mix standard bleach and water at a ratio of 1:10, respectively.

Periodic Maintenance

Periodically inspect the following areas:

- Power cord(s) should be free of cuts or other visible damage.
- All fasteners should be in place and tightened securely.
- All mechanical functions should operate properly.

Periodically lubricate the back hinge to maintain quiet, smooth, operation (use light machine oil).

Have a qualified service technician inspect your equipment every twelve months.

WARNING
Should damage and or loose parts be identified, the examination chair should be removed from use and serviced. Failure to do so may result in injury.

Specifications/Compliance


<i>Specifications Chart</i>	
Patient Weight (max):	227 kg (500 lbs)
Paper Roll (maximum size):	61 cm long x 9.2 cm diameter (24 in. x 3-5/8 in.)
Range of Motion + Dimensions:	[See Range of Motion + Dimensions pages]
Weight of Chair: w/upholstery w/packaging & skid (no uph.) Uph. with packaging (shipped separately)	197.3 kg (435 lbs) 213.2 kg (470 lbs) 18.1 kg (40 lbs)
Stirrup Weight Capacity:	56.7 kg (125 lbs) each/113.4 kg (250 lbs) combined
Power Cord Length:	244 cm (8 ft.) long
Electrical Requirements:	[See Model Identification/Compliance Chart]
Foot Control Voltage:	3.3 VDC
Duplex Receptacles Maximum Load: (not available on all models)	115 VAC, 4A, 60 Hz
Fuse Ratings: F1 (on Main PC board) F2 (on Main PC board)	T6.3AH, 250V, 5 x 20 mm T250mAH, 250V, 5 x 20 mm
Duty Cycle (Motor Run Time):	Intermittent Operation [30 seconds ON/5 minutes OFF]
Receptacles, Drawer Heaters:	Continuous Operation
Classifications:	Class I, Type B Applied Part, Intermittent Operation [30 seconds ON/5 minutes OFF]
Certification:	ES 60601-1:2012, CSA C22.2 No. 60601-1:14; (Safety Standards)
Compliance:	IEC 60601-1-2:2020 Ed. 4.1 (EMC Standards)
Protection against ingress of fluids:	IPX1 Foot Control

Model Identification/Compliance Chart


Model	Description	Complies To:			Electrical Ratings:		
		ANSI/AAMI ES60601-1	CAN/CSA-C22.2 No. 60601-1:14	IEC 60601-1-2	VAC	Amps	Cycles (Hz)
228-001	Ritter® Basic Chair	•	•	•	115	9	60
228-002	Ritter® Basic Chair w/ Receptacles, Pelvic Tilt + Drawer Heater	•	•	•	115	11.5	60
228-003	Ritter® Basic Chair w/ Receptacles, Recessed Rollers, Pelvic Tilt + Drawer Heater	•	•	•	115	11.5	60
229-002	Ritter® Power Back Chair w/ Receptacles, Pelvic Tilt + Drawer Heater	•	•	•	115	11.5	60
229-003	Ritter® Power Back Chair w/ Receptacles, Recessed Rollers, Pelvic Tilt + Drawer Heater	•	•	•	115	11.5	60

Fire Code Ratings: All upholstery complies with California Bureau of Home Furnishing Technical Bulletin 117 and California Code of Regulations, Sect. 93120-93120.12, Title 17. Optional upholstery is available that complies with California Bureau of Home Furnishing Technical Bulletin 133.

Range of Motion:

Dimensions:

Accessories + Common Parts

<i>Accessory/Common Part</i>		<i>Part Number</i>
Patient Support Rail	(Type B Applied Part)	9A600001
Patient Support Rail Plus	(Type B Applied Part)	9A600002
Stainless Steel Debris Tray (not available for 228-001)		9A610001
Assist Arms	(Type B Applied Part)	9A615001
253 LED Light Spline Mounting Hardware		9A619001
Knee Crutch Mounting Brackets		9A641001
Articulating Knee Crutch w/Mounting Brackets (9A641001)	(Type B Applied Part)	9A645001
229 (-002/-003) Hand Control Upgrade		9A681001
Flat Premium 28" Upholstery Top	(Type B Applied Part)	002-2008-XXX*
Flat UltraFree 28" Upholstery Top	(Type B Applied Part)	002-10361-XXX*
Premium Thick Footrest Pad	(Type B Applied Part)	002-10494-XXX*
Premium UltraFree Thick Footrest Pad	(Type B Applied Part)	002-10495-XXX*

*XXX indicates your three-digit color code. Call 1.800.MIDMARK to contact Customer Service to obtain information regarding Midmark's current available color selection. To maintain reliability and performance use only genuine Midmark replacement parts.


Consult the documentation provided with the accessories for their connection, disconnection, intended use and operation.

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