



Technical Advisory Notice	003-10167-00	Rev. A03	Page 1 of 2
Midmark Imaging Service Kit Part Number Updates			Release Date: April 2018



PURPOSE: This notice provides additional information for the following products:

- *Vantage Panoramic X-Ray System, Vantage Cephalometric X-Ray System*
- *VetProDC IO Vet X-Ray System, Preva IO Dental X-Ray System*

Specifically, this addendum provides additional information for service kits with part numbers that have changed.

Attention

The equipment must only be installed and operated in accordance with the safety procedures and operating instructions given in the Installation Manual and in the User Manual for these products, for the purposes and applications for which it was designed. Modifications and/or additions to the equipment may only be carried out by Midmark Corporation or by third parties expressly authorized by Midmark Corporation to do so. Such changes must comply with legal requirements as well as with the generally-accepted technical rules. It is the responsibility of the user to ensure that existing legal regulations regarding installation of the equipment with respect to the building are observed.

Midmark Imaging Service Kits

Vantage – Panoramic and Cephalometric X-Ray Systems

FORMER Part Number	CURRENT Part Number	Vantage Service Kit Description
60-A2083	002-10277-00	ASSY, SECONDARY COLLIMATOR, CEPH
60-A1004	002-10278-00	SENSOR ASM. CEPHALOMETRIC,
60-A2009	60-A2092	PRIMARY COLLIMATOR ASM
60-A1010	002-10294-00	SENSOR ASM. PANORAMIC
60-A1014	60-A2064	TUBE HEAD WITH OIL, PANORAMIC

Preva / VetProDC – Intraoral X-Ray Systems

FORMER Part Number	CURRENT Part Number	Preva / VetProDC Service Kit Description
30-A2195	002-10279-00	MODULAR BLD, 70 MM CONE,
30-A2196	002-10280-00	MODULAR BLD, 60 MM CONE, WHITE,
30-A2198	002-10281-00	MODULAR BLD, 30X40MM CONE,
30-A2199	002-10282-00	MODULAR BLD, 20X30MM CONE,
30-A2200	002-10283-00	KIT, CONE, 70MM CONE, MODULAR, LONG, GRAY
30-A2201	002-10284-00	KIT, CONE, 60MM, MODULAR LONG, WHITE
30-A2203	002-10285-00	KIT, CONE, 30X40MM, MODULAR, LONG, WHITE
30-A2221	002-10286-00	MODULAR BLD, 35x45MM CONE, GRAY
30-A2222	002-10287-00	MODULAR BLD, 35x45MM CONE, WHITE
30-A2223	002-10288-00	MODULAR BLD, 35x45MM CONE LONG, GRAY
30-A2224	002-10289-00	MODULAR BLD, 35x45MM CONE LONG, WHITE
30-A2228	002-10290-00	CONE, 60MM, MODULAR BLD, GRAY,
30-A2229	002-10291-00	KIT, CONE, 60MM, MODULAR LONG, GRAY
30-A1027	002-10292-00	TUBE HEAD WITH OIL, PREVA,
30-A1027-w	002-10293-00	TUBE HEAD WITH OIL, PREVA RAL, 9010



Technical Advisory Notice	003-10167-00	Rev. A03	Page 2 of 2
Midmark Imaging Service Kit Part Number Updates			Release Date: April 2018

This page intentionally blank.

**Manufacturer**

Midmark Corporation
675 Heathrow Drive
Lincolnshire, IL 60069 USA
Phone: 847-415-9800
Fax: 847.415-9801
www.midmark.com

**Contacting Technical Support**

Phone: 1-800-MIDMARK (1-800-643-6275)

Email: imagingtechsupport@midmark.com

Technical Library

www.midmark.com/technical-library