


Specifications: Shor-Line 8-Foot Wide Stainless Steel Cage Bank Assembly



902-0109-29

This 8-foot wide stainless steel cage bank specification includes dimensions and material pertaining to options included in this particular cage bank configuration.

| Stainless Steel Cage Bank | |
|--|---|
| Physical Dimensions | |
| Overall Width (A) | 244.32 cm (96.19 in) |
| Overall Height (B) | 174.07 cm (68.53 in) |
| Overall Depth (C) | 71.75 cm (28.25 in) |
| Floor Surface Area | 3.27 sq m (35.20 sq ft) |
| Load Bearing | 66.98 cm (26.37 in) – Center line of the front and back tube framing |
| Weight | |
| Weight | 185 kg (408 lbs) |
| Material | |
| Cages/Base Framing | Series 304 stainless steel* |
| Sound Dampening Pads | 76.2 x 203.2 x 1.5 mm (approx 1/8 in) |
| Hinges | Stainless steel reinforced polymer hinge assembly (US Patent 7805808) |
| Casters | Reinforced nylon caster with polyurethane locking wheels |
| * Series 304 is referred to as medical grade due to its composition. | |



| Cage Bank Configuration | |
|-------------------------|----------------------------------|
| Top Row | 2 @ 30"W x 30"H, 1 @ 36"W x 30"H |
| Bottom Row | 2 @ 48"W x 30"H |
| Platform | Standard mobile platform |



Equipped Options

- Mobile platform includes two 48" long stainless steel frames with nylon casters
- Cage bank stainless steel top and end panels
- Cage bank ships fully assembled

Note: Specifications include the standard wire cage door for weight calculations.

⚠ This product can expose you to chemicals, including nickel (metallic), which are known to the state of California to cause cancer. For more information, go to: [P65Warnings.ca.gov](https://www.P65Warnings.ca.gov)