


Specifications: Shor-Line 3-Foot Wide Stainless Steel Cage Bank Assembly



902-0328-50

This 3-foot wide stainless steel cage bank specification includes dimensions and material pertaining to options included in this particular cage bank configuration.

| Stainless Steel Cage Bank | |
|----------------------------------------------------------------------|-----------------------------------------------------------------------|
| Physical Dimensions | |
| Overall Width (A) | 91.77 cm (36.13 in) |
| Overall Height (B) | 177.95 cm (70.06 in) |
| Overall Depth (C) | 71.75 cm (28.25 in) |
| Floor Surface Area | 1.77 sq m (19.05 sq ft) |
| Load Bearing | 66.98 cm (26.37 in) – Center line of the front and back tube framing |
| Weight | |
| Weight | 96 kg (212 lbs) |
| Material | |
| Cages/Base Framing | Series 304 stainless steel* |
| Sound Dampening Pads | 76.2 x 203.2 x 1.5 mm (approx 1/8 in) |
| Hinges | Stainless steel reinforced polymer hinge assembly (US Patent 7805808) |
| * Series 304 is referred to as medical grade due to its composition. | |



| Cage Bank Configuration | |
|-------------------------|-----------------|
| Top Row | 2 @ 18"W x 18"H |
| Middle Row | 2 @ 18"W x 24"H |
| Bottom Row | 1 @ 36"W x 28"H |
| Platform | No platform |



Equipped Options

- Cage bank ships fully assembled

Notes:

- Specifications include the standard wire cage door for weight calculations.
- Surface should be water-resistant material (e.g., laminate, stainless steel, solid surface materials, stone, concrete).
- Do NOT use trim moulding—the front edge of the platform must be clear of obstruction.
- Ensure base moulding does not interfere with the cage front drain lip, alignment plate.

 This product can expose you to chemicals, including nickel (metallic), which are known to the state of California to cause cancer. For more information, go to: [P65Warnings.ca.gov](https://www.P65Warnings.ca.gov)