


Specifications: Shor-Line 4-Foot Wide Stainless Steel Cage Bank Assembly



902-0420-52

This 4-foot wide stainless steel cage bank specification includes dimensions and material pertaining to options included in this particular cage bank configuration.

Stainless Steel Cage Bank	
Physical Dimensions	
Overall Width (A)	122.25 cm (48.13 in)
Overall Height (B)	185.50 cm (73.03 in)
Overall Depth (C)	71.75 cm (28.25 in)
Floor Surface Area	1.63 sq m (17.55 sq ft)
Load Bearing	66.98 cm (26.37 in) – Center line of the front and back tube framing
Weight	
Weight	105 kg (232 lbs)
Material	
Cages/Base Framing	Series 304 stainless steel*
Sound Dampening Pads	76.2 x 203.2 x 1.5 mm (approx 1/8 in)
Hinges	Stainless steel reinforced polymer hinge assembly (US Patent 7805808)
Casters	Reinforced nylon caster with polyurethane locking wheels
PVC Floor	14-gauge perforated sheet steel with 3/8" diameter openings coated with 1/8" thick non-toxic PVC
* Series 304 is referred to as medical grade due to its composition.	



Cage Bank Configuration	
Top Row	2 @ 24"W X 30"H
Bottom Row	1 @ 48"W X 36"H
Platform	Low-profile mobile platform



Equipped Options

- Stainless steel low-profile frame mobile platform with nylon casters
- Cage bank stainless steel left end panel
- PVC-coated cage floors (PVC floor height – 2 1/4" above cage floor)
- Cage bank ships fully assembled

Note: Specifications include the standard wire cage door for weight calculations.