




This 5-foot wide stainless steel cage bank specification includes dimensions and material pertaining to options included in this particular cage bank configuration.

| Stainless Steel Cage Bank | |
|--|--|
| Physical Dimensions | |
| Overall Width (A) | 152.73 cm (60.13 in) |
| Overall Height (B) | 185.50 cm (73.03 in) |
| Overall Depth (C) | 71.75 cm (28.25 in) |
| Floor Surface Area | 2.02 sq m (21.74 sq ft) |
| Load Bearing | 66.98 cm (26.37 in) – Center line of the front and back tube framing |
| Weight | |
| Weight | 128 kg (283 lbs) |
| Material | |
| Cages/Base Framing | Series 304 stainless steel* |
| Sound Dampening Pads | 76.2 x 203.2 x 1.5 mm (approx 1/8 in) |
| Hinges | Stainless steel reinforced polymer hinge assembly (US Patent 7805808) |
| Casters | Reinforced nylon caster with polyurethane locking wheels |
| PVC Floor | 14-gauge perforated sheet steel with 3/8" diameter openings coated with 1/8" thick non-toxic PVC |
| * Series 304 is referred to as medical grade due to its composition. | |



| Cage Bank Configuration | |
|-------------------------|-----------------------------|
| Top Row | 2 @ 30"W x 30"H |
| Bottom Row | 1 @ 60"W x 36"H |
| Platform | Low-profile mobile platform |



Equipped Options

- Stainless steel low-profile frame mobile platform with nylon casters
- PVC-coated cage floors (PVC floor height – 2 1/4" above cage floor)
- Cage bank ships fully assembled

Note: Specifications include the standard wire cage door for weight calculations.

⚠ This product can expose you to chemicals, including nickel (metallic), which are known to the state of California to cause cancer. For more information, go to: P65Warnings.ca.gov