


# Specifications: Shor-Line 4.5-Foot Wide Stainless Steel Cage Bank Assembly



902-1800-01

This 4.5-foot wide stainless steel cage bank specification includes dimensions and material pertaining to options included in this particular cage bank configuration.

Stainless Steel Cage Bank	
Physical Dimensions	
Overall Width (A)	137.64 cm (54.19 in)
Overall Height (B)	204.70 cm (80.59 in)
Overall Depth (C)	71.75 cm (28.25 in)
Floor Surface Area	3.48 sq m (37.46 sq ft)
Load Bearing	66.98 cm (26.37 in) – Center line of the front and back tube framing
Weight	
Weight	186 kg (411 lbs)
Material	
Cages/Base Framing	Series 304 stainless steel*
Sound Dampening Pads	76.2 x 203.2 x 1.5 mm (approx 1/8 in)
Hinges	Stainless steel reinforced polymer hinge assembly (US Patent 7805808)
Casters	Reinforced nylon caster with polyurethane locking wheels
* Series 304 is referred to as medical grade due to its composition.	



Cage Bank Configuration	
Top Row	3 @ 18"W x 18"H
Middle Row 1	3 @ 18"W x 18"H
Middle Row 2	3 @ 18"W x 18"H
Bottom Row	3 @ 18"W x 18"H
Platform	Standard mobile platform



### Equipped Options

- Mobile platform includes two 48" long stainless steel frames with nylon casters
- Cage bank components ship separately and require on-site assembly

*Note: Specifications include the standard wire cage door for weight calculations.*

⚠ This product can expose you to chemicals, including nickel (metallic), which are known to the state of California to cause cancer. For more information, go to: [P65Warnings.ca.gov](https://www.P65Warnings.ca.gov)