


Specifications: Shor-Line 3-Foot Wide Stainless Steel Cage Bank Assembly



902-0102-12

This 3-foot wide stainless steel cage bank specification includes dimensions and material pertaining to options included in this particular cage bank configuration.

Stainless Steel Cage Bank	
Physical Dimensions	
Overall Width (A)	91.77 cm (36.13 in)
Overall Height (B)	204.62 cm (80.56 in)
Overall Depth (C)	71.75 cm (28.25 in)
Floor Surface Area	1.77 sq m (19.05 sq ft)
Load Bearing	66.98 cm (26.37 in) – Center line of the front and back tube framing
Weight	
Weight	109 kg (241 lbs)
Material	
Cages/Base Framing	Series 304 stainless steel*
Sound Dampening Pads	76.2 x 203.2 x 1.5 mm (approx 1/8 in)
Hinges	Stainless steel reinforced polymer hinge assembly (US Patent 7805808)
Casters	Reinforced nylon caster with polyurethane locking wheels
* Series 304 is referred to as medical grade due to its composition.	



Cage Bank Configuration	
Top Row	2 @ 18"W x 24"H
Middle Row	2 @ 18"W x 24"H
Bottom Row	1 @ 36"W x 24"H
Platform	Standard mobile platform



Equipped Options

- Mobile platform includes two 48" long stainless steel frames with nylon casters
- Cage bank stainless steel top and end panels
- Cage bank ships fully assembled

Note: Specifications include the standard wire cage door for weight calculations.