




This 8-foot wide stainless steel cage bank specification includes dimensions and material pertaining to options included in this particular cage bank configuration.

Stainless Steel Cage Bank	
Physical Dimensions	
Overall Width (A)	244.17 cm (96.13 in)
Overall Height (B)	174.07 cm (68.53 in)
Overall Depth (C)	71.75 cm (28.25 in)
Floor Surface Area	3.32 sq m (35.74 sq ft)
Load Bearing	66.98 cm (26.37 in) – Center line of the front and back tube framing
Weight	
Weight	186 kg (410 lbs)
Material	
Cages/Base Framing	Series 304 stainless steel*
Sound Dampening Pads	76.2 x 203.2 x 1.5 mm (approx 1/8 in)
Hinges	Stainless steel reinforced polymer hinge assembly (US Patent 7805808)
Casters	Reinforced nylon caster with polyurethane locking wheels
Back Panel	Vented
* Series 304 is referred to as medical grade due to its composition.	



Cage Bank Configuration	
Top Row	2 @ 48"W x 30"H
Bottom Row	2 @ 48"W x 30"H
Platform	Standard mobile platform



Equipped Options

- Mobile platform includes two 48" long stainless steel tube frames with nylon casters
- Vented back panel
- Cage bank stainless steel end panels
- Quiet close latch
- Cage bank ships fully assembled

Note: Specifications include the standard wire cage door for weight calculations.

⚠ This product can expose you to chemicals, including nickel (metallic), which are known to the state of California to cause cancer. For more information, go to: [P65Warnings.ca.gov](https://www.P65Warnings.ca.gov)