


Suites are available in several sizes. This specification includes dimensions and materials pertaining to options included in this particular model. Units can be combined to create various configurations.

| Specifications                                                                            |                                                          |
|-------------------------------------------------------------------------------------------|----------------------------------------------------------|
| Physical Dimensions                                                                       |                                                          |
| Overall Height (A)                                                                        | See "Dimensions"                                         |
| Overall Width (B)                                                                         | See "Dimensions"                                         |
| Overall Depth (C)                                                                         | See "Dimensions"                                         |
| Unit Interior Dimensions                                                                  | 22.5" W x 28.5" H x 28" D                                |
| Litter Interior Dimensions                                                                | 16.5" W x 28.5" H x 28" D                                |
| Weight                                                                                    |                                                          |
| Weight                                                                                    | See "Dimensions"                                         |
| Materials                                                                                 |                                                          |
| Unit Panels                                                                               | Polyvinyl Chloride                                       |
| Door Options:<br>Tri-Door<br>Tri-Door w/Viewing Dome<br>Full Poly Door<br>Full Grill Door | Stainless Steel Grill                                    |
|                                                                                           | Polycarbonate                                            |
|                                                                                           |                                                          |
|                                                                                           |                                                          |
| Legs                                                                                      | Stainless Steel                                          |
| Casters                                                                                   | Reinforced Nylon Caster with Polyurethane Locking Wheels |
| Hinges                                                                                    | Injection-Molded Fiberglass-Nylon Thermoplastic Polymer  |




| Configurations |            |                                    |
|----------------|------------|------------------------------------|
| Option A       | Top Row    | 1 @ 18"W x 30"H<br>1 @ 24"W x 30"H |
|                | Bottom Row | 1 @ 18"W x 30"H<br>1 @ 24"W x 30"H |
| Option B       | Top Row    | 1 @ 24"W x 30"H<br>1 @ 18"W x 30"H |
|                | Bottom Row | 1 @ 24"W x 30"H<br>1 @ 18"W x 30"H |


| Dimensions                                                                                                                                 |             |                        |                    |                     |                                          |                                 |
|--------------------------------------------------------------------------------------------------------------------------------------------|-------------|------------------------|--------------------|---------------------|------------------------------------------|---------------------------------|
| Base Option                                                                                                                                |             | A - Exterior Height    | B - Exterior Width | C - Exterior Depth  | Floor Surface Area                       | Weight*                         |
| Option A                                                                                                                                   | No Platform | 152.40 cm (60 in)      | 106.68 cm (42 in)  | 73.03 cm (28.75 in) | 5.34 sq m (57.50 sq ft)                  | Configurations vary per options |
|                                                                                                                                            | Mobile      | 173.99 cm (68.50 in)   |                    |                     |                                          |                                 |
|                                                                                                                                            | Stationary  | 175.90 cm (69.25 in)** |                    |                     |                                          |                                 |
| Option B                                                                                                                                   | No Platform | 152.40 cm (60 in)      | 106.68 cm (42 in)  | 73.03 cm (28.75 in) | 5.34 sq m (57.50 sq ft)<br>(Living Unit) | Configurations vary per options |
|                                                                                                                                            | Mobile      | 173.99 cm (68.50 in)   |                    |                     | 4.01 sq m (43.13 sq ft)<br>(Litter Unit) |                                 |
|                                                                                                                                            | Stationary  | 175.90 cm (69.25 in)** |                    |                     |                                          |                                 |
| * Weights vary depending on base platform, back panel, and door type options.<br>** Heights may vary depending on adjustable feet setting. |             |                        |                    |                     |                                          |                                 |

\* Weights vary depending on base platform, back panel, and door type options.

\*\* Heights may vary depending on adjustable feet setting.



NO PLATFORM



Casters (TYP)

MOBILE PLATFORM



Adjustable  
Legs (TYP)

STATIONARY PLATFORM

NOTE: Caster height is 21.60 cm (8.5 in), adjustable feet height is 23.50 cm (9.25 in). Pass-through included.