



Rear Storage Pod [9A313]

Contents:


Installation 1
Exploded View 2

Intended Use:

To provide additional storage space at exam table.

Special Tools:

none




Step 1: Raise back section; then slide storage pod onto top cover.



Equipment Alert

Do not lower back section with Rear Storage Pod lid open.
Doing so may result in damage to lid and/or storage unit.

MA607700i



MA612900i

9A313001 Storage Pod		
Item	Description	Qty.
1	• Storage Container Lid	1
2	• Storage Divider	1
3	• Storage Container	1

Serial Numbers / Dates:

all