

Work Surface (18" x 13") 152396

Intended Use:

Work Surface for Dentist, Hygienist, and / or Assistant

Requirements / Restrictions:

Work Surface maximum load is 40 lbs. (18.14 kgs.)


Components:

Screws:

#10 x 5/8"
wood screw 4
#10 - 32 x 1/2" 8
Spacer 1
Plate 1
Work Surface 1

Special Tools:

none




DA178300i

Work Surface, 152396-01

Cart (18" x 13")

ACCESSORY PARTS LIST



DA178301i

DA178301i

Item	Part No.	Description	Qty	Item	Part No.	Description	Qty
	152396-01	Work Surface (<i>pebble</i>)		3	• 121369-50	• Plate, Mounting (<i>pebble</i>)	1
1	• 063-9943-00	• Shelf (<i>pebble</i>)	1	4	• 040-0010-48	• Screw (#10 x 5/8", wood)	4
2	• P12577	• Screw (#10-32 x 1/2")	8	5	• 121370-50	• Spacer (<i>pebble</i>)	1

Always Specify Model & Serial Number