







Notes:



Leviton 42080-2xS or equivalent
 Leviton 5G108-Rx5 or equivalent RJ45
 Leviton 41106-Rx6 or equivalent RJ11



Leviton 42080-2xS or equivalent
 Leviton 41106-Rx6 or equivalent RJ11
 Leviton 5G108-Rx5 or equivalent RJ45



Leviton 42080-2xS or equivalent
 2x Leviton 5G108-Rx5 or equivalent RJ45

* Total length of wiring between devices,
 including cabling outside the wall.

Midmark
Digital Vital Signs Device/625-004
Pre-wire Diagram
004-10042-00 REV. AA2 1/9/20

MA10177i