







Notes:



Leviton 42080-1xS or equivalent
Leviton 41106-Rx6 or equivalent RJ11



Leviton 42080-2xS or equivalent
Leviton 41106-Rx6 or equivalent RJ11
Leviton 5G108-Rx5 or equivalent RJ45



Leviton 42079-1xS or equivalent
Leviton 5G108-Rx5 or equivalent RJ45

* Total length of wiring between devices,
including cabling outside the wall.

Midmark
IQvitals Zone®/626
Pre-wire Diagram
004-10045-00 REV. AA2 7/22/20

MA10223i