




Notes:




Table / Midmark Digital Vital Signs Device (RJ11)

Leviton 42080-1xS or equivalent
Leviton 41106-Rx6 or equivalent RJ11





Table / Midmark Digital Vital Signs Device (RJ11)

Midmark Digital Vital Signs Device / PC (RJ45)

Leviton 42080-2xS or equivalent
Leviton 41106-Rx6 or equivalent RJ11
Leviton 5G108-Rx5 or equivalent RJ45



Midmark Digital Vital Signs Device / PC (RJ45)

Leviton 42079-1xS or equivalent
Leviton 5G108-Rx5 or equivalent RJ45

* Total length of wiring between devices, including cabling outside the wall.

Midmark
Digital Vital Signs Device/626
Mobile PC Pre-wire Diagram
004-10287-00 REV. AA1 7/22/20

MA11757i