


Dual Unit Install with Air/Water Separator

Note
Waste Water/Vacuum Hookup Kit must be purchased and installed in conjunction with either the CV3 or CV5 Dual Hookup Kit.



Note
When the exhaust is routed to and air water separator the larger twin air water separator for twin vacuums should be installed.

- Instructions are for:**
- Dual Hookup Kit:
77000034 - CV3(R) -or- 77000950 - CV5(R)
 - Waste/Water/Vacuum Hookup Kit:
77000145
 - Twin Air Water Separator (recommended)
77000246

ArtFull